

Spray Dryer

GB211C

The First Edition

- Thank you for purchasing GB series product of Yamato Scientific Chongqing Co., Ltd.
- To use this unit properly, read this "Instruction Manual" thoroughly before using this unit. Keep this instruction manual around this unit for referring at anytime.



WARNING!

Carefully read and thoroughly understand the important warning items described in this manual before using this unit.

Yamato Scientific Co.,Ltd.

1. Safety precautions.....	1
Explanation of symbols	1
List of symbols	2
Warning · Cautions.....	3
2. Before using this unit	5
Precautions when installing the unit.....	5
Service socket capacity.....	8
Temperature output terminal	8
3. Names and functions of the parts	9
Main unit + GF301C set (GB211C-A)	9
Main unit + GF200 set (GB211C-B)	11
GB211C+GF301C (GB211C-A) operation interface overview.....	13
GB211C+GF200 (GB211C-B) operation interface overview	14
Description of switch and indicator lamp in the interface	15
Value display and input description in the interface	16
4. Operating procedures	17
PUSH mode installation	17
Preparation before operation	20
Preparation before operation (GB211C+GF301C installation).....	22
Preparation before operation (GB211C+GF200 installation)	27
GB211C+GF301C operating method.....	30
GB211C+GF200 operating method	35
Use of automatic needle spray nozzle (for GB211C+GF301C)	41
Granulation automatic liquid sending function (for GB211C+GF200)	42
If want to interrupt the sample processing, or when nozzle blockage occurs	43
The relation between rotate speed of liquid sending pump and liquid amount/between blower power and dry air amount (reference).....	44
Calibration of temperature sensor.....	46
5. Handling Precautions	47
Drying Method under Appropriate Condition	48
Caution during operation.....	49
6. Maintenance Method	50
Daily Inspection and Maintenance	50
7. Long storage and disposal	56
When not using this unit for long term / When disposing.....	56
Matters to consider when disposing of the unit.....	56
8. When a trouble occurs	57
Safety device and error indications	57
Trouble Shooting.....	60
9. After Service and Warranty	62
When requesting a repair.....	62

10. Specification	63
11. Wiring diagram.....	66
GB211C wiring diagram	66
12. System diagram	67
Standard mode system diagram	67
PUSH mode system diagram.....	68
Granulation mode system diagram (GB211C+GF200)	69
13. Operation principle	70
14. Replace parts list	72
15. List of Dangerous Substances.....	76
16. Standard installation manual.....	77


1. Safety precautions


Explanation of symbols

About symbols

A variety of pictograms are indicated in this operating instruction and on products for safe operation. Possible results from improper operation ignoring them are as follows.

Be sure to fully understand the descriptions below before proceeding to the text.

 **Warning** Indicates a situation which may result in death or serious injury (Note 1.)

 **Caution** Indicates a situation which may result in minor injury (Note 2) and property damages (Note 3.)

(Note 1) Serious injury means a wound, an electrical shock, a bone fracture or intoxication that may leave after effects or require hospitalization or outpatient visits for a long time.

(Note 2) Minor injury means a wound or an electrical shock that does not require hospitalization or outpatient visits for a long time.

(Note 3) Property damage means damage to facilities, devices and buildings or other properties.

Meanings of symbols



This pictogram indicates a matter that encourages the user to adhere to warning ("caution" included).

Specific description of warning is indicated near this pictogram.



This pictogram indicates prohibitions

Specific prohibition is indicated near this pictogram.



This pictogram indicates matters that the user must perform

Specific instruction is indicated near this pictogram.

1. Safety precautions

List of symbols

Warning



General warnings



Danger!: High voltage



Danger!: High temperature



Danger!: Moving part



Danger!: Hazard of explosion

Caution



General cautions



Electrical shock!



Burning!



Caution for no liquid heating!



Caution for water leak!



For water only



Poisonous material

Prohibitions



General bans



Fire ban



Do not disassemble



Do not touch

Compulsions



General compulsions



Connect ground wire



Install levelly



Pull out the power plug



Regular inspection

1. Safety precautions

Warning · Cautions

Warning



Do not use this unit in an area where there is flammable or explosive gas

Never use this unit in an area where there is flammable or explosive gas. This unit is not explosion-proof. An arc may be generated when the power switch is turned on or off, and fire/explosion may result. (Refer to page 76 “15. List of Dangerous Substances”.)



Always ground this unit

Always ground this unit on the power equipment side in order to avoid electrical shock due to a power surge.



Apply the source of rated power

Be sure to apply the source of rated power or more. Applying non-rated voltage or non-rated power supply may cause the fire or electric shock.



Prohibition of use for error

If a smoke or abnormal smell may be occurred, turn off the power switch of the main unit immediately, and turn off the original power source, and finally contact to either the dealer you purchased this unit or our sales office. Leaving the failure may cause the fire or electric shock. Since the repairing of this unit is dangerous for non-specified service person, never repair the unit by the customer himself.



Do not use the power cord if it is bundled or tangled

Do not use the power cord if it is bundled or tangled. If it is used in this manner, it can overheat and fire may be caused.



Do not damage power cord

Do not damage power cord by bending, pulling, or twisting forcibly. It may cause the fire or electric shock. Besides, operating the unit with the something put on the cord may cause overheat, and result in fire.



Never use explosive or flammable material with this unit.

Never use an explosive material, a flammable material or a material containing them. An explosion or an electrical shock may result. See section “15. List of Dangerous Substances” on page 76. **Connect GB211C+GF301C (GB211C-A) with the optional GAS series product to form an enclosed and low-oxygen circulation system, which is able to use the organic solvent sprays without the risk of explosion. When using the organic solvents, pay special attention to their explosion conditions, especially the mixture of multiple organic solvents. Please read the GAS series product instruction manual for operations.**



Never try to touch a hot part.

Some parts of the unit are hot during and immediately after operation. Take special care for possible scald.



Never try to disassemble or modify the unit.

Never try to disassemble or alter the unit. A malfunction, a fire or an electric shock may result.

1. Safety precautions

Warning · Cautions

Caution

During a thunder storm

During a thunderstorm, turn off the power key immediately, then turn off the circuit breaker and the main power. If this procedure is not followed, fire or electric shock may be caused.

After power failure

When power is shut off during operation (while the blower is operating or the liquid is being sent) due to turning the ELB to "OFF" or a power failure, all operation modes will reset to the initial states after recovery. When the temperature inside the chamber is higher, keep operating the blower until it cools down to 45°C or below after recovery from a power failure.


Do not perform unattended operation

Do not perform unattended operation.

After running out of the samples, it will cause idling and nozzle blocking, the outlet temperature will rise, the sample hose will be disconnected from the nozzle and the remaining sample will flow out, which may result in unexpected accident.

About countermeasures against static electricity

The cyclone may be charged with static electricity depending on the samples, operating environment or conditions. Implement countermeasures against static electricity such as installing the accessorial earth clips on the clamp (4 positions) at the connection of the cyclone, or installing an antistatic brush (optional) to the body of the cyclone.



2. Before using this unit


Precautions when installing the unit

Warning

1. Always ground this unit



- The power supply is single-phase 200-230V \sim 50/60Hz. Please entrust the professional personnel of the nearest electrical construction shop or agent to connect the power supply.
- The protection impedance of the machine is 0.5 Ω or less. Carry on the grounding construction according to the technical requirements of the electrical equipment at the location of the customer. If the technical requirements are not clear, the construction and acceptance of grounding works according to the grounding resistance of less than 4 Ω .
- Be sure to connect the earth wire (the green cable of power cord) to the grounding conductor or ground terminal to prevent accidents caused by electric leakage.
- Do not connect the earth wire to gas or water pipes. If not, the fire disaster may be caused.
- Do not connect the earth wire to the grounding of telephone wire or lightning conductor. If not, the fire disaster may be caused.



Core Wire Color	In-house Wiring
Brown	Voltage Side
Blue	Voltage Side
Yellow green	Ground Side

- Before connecting the power cord, turn off the protection switch on the power supply device.
- This device does not have accessorial plug, please select the plug and terminal that match the power supply capacity according to the connected power supply device.
- Please note the color of each core.

2. Please use a dedicated power supply



Use a power supply that matches the power supply capacity.

Electrical capacity: single phase 200-230V \sim 50/60Hz 17-21A (Protection circuit breaker operating current is 25A)

If the power supply is ON, but the equipment does not start, check whether the main power supply voltage is lower or share a power cord with other machines. Please use the power cord separately from other machines. In consideration of the safety of the machine power connection, please entrust the seller, agent or electrical construction shop to carry on.

3. Choose a proper place for installation




Do not install this unit in a place where:

- Rough or dirty surface.
- Flammable gas or corrosive gas is generated.
- Ambient temperature below 5°C or above 30°C.
- Ambient temperature fluctuates violently.
- There is direct sunlight.
- There is excessive humidity and dust.
- There is a constant vibration.
- Place where the water is easy-to-be splashed.



Install this unit on a stable place with the space as shown below.



2. Before using this unit

Precautions when installing the unit

Warning


4. Do not use this unit in an area where there is flammable or explosive gas



Never use this unit in an area where there is flammable or explosive gas. This unit is not explosion-proof. An arc may be generated when the power switch is turned ON or OFF, and fire/explosion may result.



Refer to page 76 "15. List of Dangerous Substances".




5. Do not modify



Modification of this unit is strictly prohibited. This could cause a failure.


Modification



6. Do not topple or tilt this unit



Set this unit to the flattest place. Setting this unit on rough or slope place could cause the vibration or noise, or cause the unexpected trouble or malfunction.



7. Place the unit



Due to sudden earthquake, impact, etc., the product may collapse or move, and then be damaged.

It is best to avoid places where there are many people and take safety precautions.

2. Before using this unit


Precautions when installing the unit



Warning

8. Do not use explosive or flammable substances

- ⊘ Never use explosive substances, flammable substances and substances that include explosive or flammable ingredients in this unit. Explosion or fire may occur. **GB211C+GF301C (GB211C-A) supports organic solvents by connecting it to the optional GAS series product. Carefully read the operation manual of GAS series product and take special care for handling of organic solvents.** Refer to page 76 “15. List of Dangerous Substances”.



9. Absolutely prohibit the use of toxic or biohazardous substances

- ⊘ This product and GAS are not developed for biosafety purpose, do not have the ability to treat toxic or biohazardous substances, and the use of toxic or biohazardous substances is absolutely prohibited, for example: polychlorinated biphenyls, cyanide, virus or bacteria.

10. Absolutely prohibit the use of substances containing unknown ingredients

- ⊘ The thermal decomposition of ingredients with unknown properties may cause explosion, fire, poisoning or other accidents.

11. Handling of power cord

- ⊘ Do not entangle the power cord. This will cause overheating and possibly a fire.
Do not bend or twist the power cord, or apply excessive tension to it. This may cause a fire and electrical shock.
Do not lay the power cord under a desk or chair, and do not allow it to be pinched in order to prevent it from being damaged and to avoid a fire or electrical shock.
Keep the power cord away from any heating equipment such as a room heater. The cord's insulation may melt and cause a fire or electrical shock.
- ! If the power cord becomes damaged (wiring exposed, breakage, etc.), immediately turn off the power at the rear of this unit and shut off the main supply power. Then contact your nearest dealer for replacement of the power cord. Leaving it may cause a fire or electrical shock. Connect the power plug to the receptacle which is supplied appropriate power and voltage.

2. Before using this unit

Service socket capacity

Service socket capacity



Apply the 200-230V~ service socket for this unit.

Connecting the service socket with its capacity over 1A will blow out the fuse, and the power source to the service socket is shut down.


For resetting, replace the fuse in the fuse holder on the right side of the back of the unit.

If the rated current exceeds 1A, please use another power supply.

Temperature output terminal

The temperature output signals of the Outlet (outlet temperature) and the Inlet (inlet temperature), aiming at the measured temperature 0-300°C, the current output is 4-20mA.


[Current output 4-20mA: Measured temperature 0-300°C]



Conversion formula: Current output I (mA) = Measured temperature T (°C) \div 18.75 + 4

Measured temperature T (°C) = 18.75 \times [Current output I (mA) - 4]

When you connecting to the voltage input of the recorder, connect a fixed resistor (shunt resistor) of 300Ω or lower to the voltage input of the recorder.




+ Inlet-temp
- 4-20mA:0-300°C

+ Outlet-temp
- 4-20mA:0-300°C

3. Names and functions of the parts

Main unit + GF301C set (GB211C-A)


Front side



3. Names and functions of the parts

Main unit + GF301C set (GB211C-A)


Back side



3. Names and functions of the parts

Main unit + GF200 set (GB211C-B)


Front side



3. Names and functions of the parts






Main unit + GF200 set (GB211C-B)

Back side




3. Names and functions of the parts

GB211C+GF301C (GB211C-A) operation interface overview



3. Names and functions of the parts

GB211C+GF200 (GB211C-B) operation interface overview



※ The not shown GB211C+GF200 interfaces are consistent with GB211C+GF301C, which will not be detailed here.

3. Names and functions of the parts




Description of switch and indicator lamp in the interface

In the operation interface on the touch screen, the action state of each switch button can be confirmed by checking if the indicator lamp is on. The appearance of the switch button is characterized by a button frame inside which is an effective area for operation.

Type 1: Manual operation button

This kind of button is not only a functional switch but also an indicator lamp for the current state of the function. It has only two states of ON/OFF. The operator will switch the state once every time he operates, and the indicator will switch accordingly.


This kind of button has no power-off memory function, and it will be in OFF state after power outage recovery. Its corresponding function can only be activated by the operator.

	Initial state: OFF The temperature controller and heater are not working, and the indicator lamp is in standby state (yellow).
	Click once: the state switches from OFF to ON The temperature controller and heater are working, and the indicator lamp is in operating state (green).
	Click again: the state switches from ON to OFF The temperature controller and heater are not working, and the indicator lamp is in standby state (yellow). The subsequent operations repeat the above actions.

Type 2: Function select button and function enable indicator lamp


This kind of button is only the switch of function, and the usage status of the function is shown by a separate indicator lamp.

This kind of button has power-off memory function, after power outage recovery, will keep the state before power outage. It is characterized by the ability to maintain the last setting of operating parameters.




Inlet temp.
controller
enable button

Status indicator
lamp of inlet
temp. controller
(green means
enabled)



Outlet temp.
controller
enable button

Status indicator
lamp of outlet
temp. controller
(gray means
not enabled)



Automatic
cleanout
needle
start/stop
button

Status
indicator
lamp of
cleanout
needle

Type 3: Status indicator lamp

An indicator lamp indicating the current status of each functional unit. This indicator lamp is internal automatic function and has no corresponding switch.

If the indicator lamp is gray or invisible, it indicates that the function is not running. If the indicator lamp is green or visible, it indicates that the function is running. If the indicator lamp is red, it indicates that the function is in the alarm status.

3. Names and functions of the parts

Value display and input description in the interface

In the operation interface of touch screen, the values can be divided into two types: only display but no input; display and input.

Type 1: Value display

There is no input box for the value display. Its background and the bottom color are the same, and the display value is white.

The value display will display the data in real time. When an alarm occurs, some specific numerical values will change their colors to prompt the operator.


For example, if the temperature exceeds the upper limit allowed by the equipment, the value will change to red; and the danger level of oxygen concentration is displayed in green, yellow and red colors.

Type 2: Value input

There is input box for the value input. Its input box is white, and the display value is black.


Click the white input box to pop up the value input popup. Input the required value and press ENTER button to complete the input.

The value input will display the operator's previous input data in real time. Except for the values used in system setting, the inputs of other values are power-off memory type. The operator only needs to operate once in the initial setting, then no operation is required later as long as the value is not changed. (If it is not used for more than 14 days, data may be lost due to PLC internal power supply exhaustion, and may need to be reset.) ※ In order to keep data for a long time, a button battery can be installed on PLC to provide continuous power for PLC. Please refer to P. 54 "About the use of PLC batteries".



- ※ Special note: when the data setting in the system parameter interface is completed, click the Save Parameter button in the lower left, the parameters will be saved only after the button turns green. If the modified data is not saved or the equipment is powered off, the data will be restored to the previous data before modification.

Set the system information correctly at the first startup.



4. Operating procedures

PUSH mode installation

The spray drying mode of this unit is divided into standard mode and PUSH mode.

For the differences between the two modes, see P. 67 "12. System diagram". The biggest differences in their uses are as below:


In standard mode of spray drying, the air in all pipes and glassware is under negative pressure, and the generated hot exhaust gas will pass through the blower.

In PUSH mode of spray drying, the air in all pipes and glassware is under positive pressure, and the generated hot exhaust gas will not pass through the blower.

The factory setting of this unit is standard mode, if need to use PUSH mode, please perform the following operations to switch it to PUSH mode during installation.

- ※ It is not recommended for customer to switch back and forth between the two modes, which may lead to mixed samples.
- ※ After the standard mode is used, if you must switch to PUSH mode, clean the removed pipes and blower filter first, and then check whether there is any foreign matter attached to the blower. If there are foreign matters in the fan blades or air duct of blower, replace the blower components. Otherwise, the attachments of all pipes and fan blades will be blown into the drying chamber and contaminate the samples.


- (1) Please open the doors on the front and left side of the base and remove the pipe according to the following figure.



4. Operating procedures

PUSH mode installation


- (2) Please remove the 4 screws, open the right side door, remove the clamp of the air intake filter, and then remove the air intake filter, as shown in the following figure.



4. Operating procedures

PUSH mode installation

- (3) Please reassemble the exhaust port connecting pipe and exhaust pipe according to the following figure.




- (4) Please keep the spare parts properly for later use.

4. Operating procedures

Preparation before operation

(1) Connection of power cord

The power cord of this unit is earthed 3-core including the earth wire, and the yellow green wire must be earthed.



(2) Connection of the exhaust pipe


If concerned about the hot air and fine powder discharged from the blower, please connect the exhaust pipe, attached to the exhaust port, to the air pipe and exhaust the gas outdoors.

In consideration of safety and environmental protection, please choose GWS411C cleaning unit to clean and filter the exhaust gas before discharging it outdoors.

If organic solvents or flammable and explosive substances are used, please choose GAS series organic solvent recovery unit to carry on closed condensation recovery.

(3) The attached pressure-proof hose is used to connect the port ($\phi 7$) on the back of the upper frame with the pressurized air devices such as compressor. Please tighten by using the hose clamp. When using compressor, etc. to provide the pressure air source, use a pressure reducing valve to keep the pressure constant.

※ Please adjust the output pressure of the reducing valve to 0.3-0.6MPa and keep it constant. Using pressure greater than 0.6MPa will damage the pipe. If the pressure is less than 0.3MPa, the automatic nozzle cannot be opened. The non-constant pressure will cause the spraying to form uneven droplets.




(4) Unpack the glass components and make sure that there is no broken glass or missing items.

4. Operating procedures

Preparation before operation

- (5) Set the system parameters before using the system for the first time or when it has not been used for a long time.



- ① Turn ON the ELB on the right side of the main unit.


ATOMIZING AIR



POWER




- ② Turn ON the **Power** switch on the operation panel of the main unit, touch screen display.



- ③ Press the time display area for more than 5 secs to jump to the system parameter screen.


※ If the equipment is transported, stored or not used for a long time, leading to the exhaustion of PLC internal power supply, the previous time setting will become unreliable (will also stop updating). Reset the system time of the equipment based on the local time. Otherwise the running curve will show the wrong point in time.



※ Select the power supply frequency based on the local power supply. Otherwise, the blower cannot operate normally, which may lead to abnormal temperature or even fire. **No operating parameters are set before delivery.**

※ Please adjust the recording time interval of the running curve according to the needs of the test. The interval time from 1 to 60 secs can be set, **and the factory parameter is 5 secs.**

※ It is not necessary to set the inlet temp. and outlet temp. correction at the first use. If the sensor has deviation after a long-time use, **please refer to P. 46 "Calibration of temperature sensor"**.



- ④ Click the **Save Parameters** button until it turns green. Finally click the Return button to return to the initial screen.

4. Operating procedures


Preparation before operation (GB211C+GF301C installation)

GB211C main unit + GF301C mini spray installation procedures (GB211C-A)


(1) Install the distributor and aluminum honeycomb assembly onto the top of the unit.

Install the O-ring P16 into the $\phi 20$ groove on the top of the distributor.

(install using three M6 x 20 hex bolts, spring washers, flat washers each)




(2) Insert the pipe in the center of the distributor and tighten it clockwise.






4. Operating procedures

Preparation before operation (GB211C+GF301C installation)

- (3) Install the O-ring P110 in the groove of the lifter.
Install the tray for product collecting container to the upper row of four screw holes at the front of the main unit with knurled screws.




- (4) Install the drying chamber taking care to align the groove with the lifter positioning pin. Turn on the power switch, hold the drying chamber with your hand, keep pressing the **Raise** button on the touch screen to raise the lifter. Stop raising the lifter when the upper part of the drying chamber is close to the flange of the distributor. Then press the **Raise** button again and again to make the upper part of the drying chamber enter the inside of the distributor flange. Finally, keep pressing the **Raise** button until the lifter stops and the **BLOWER** button is displayed. Install the outlet temperature sensor at the pipe of the glassware connecting port and insert the plug into the connector at the main unit.



4. Operating procedures

Preparation before operation (GB211C+GF301C installation)


(5) Connect the cyclone following the step numbers below.



4. Operating procedures

Preparation before operation (GB211C+GF301C installation)

- (6) Insert the automatic needle spray nozzle vertically into the mounting hole on the top as shown in the figure below. Then according to the labels on the nozzle, the attached connecting pipes are used to connect the corresponding air pipe ports at the back of the unit body.



- (7) When the characteristics of the samples change greatly due to the influence of temperature, or the samples are easy to be evaporated, lead to the blockage of the pipe inside the nozzle, it is recommended to use the cooling function of the spray nozzle.

Please connect the coolant inlet and coolant outlet of the automatic needle spray nozzle with the cooling water circulator (such as CF312C) which is sold separately to circulate the coolant through the inside of the nozzle.

The whole outer wall of the automatic needle spray nozzle will be filled with flowing coolant, and the liquid and air pipes inserted into the drying chamber through the center of the nozzle will be completely wrapped, ensuring that the temperature of the samples before entering the drying chamber is basically unchanged.


4. Operating procedures

Preparation before operation (GB211C+GF301C installation)

- (8) When the setting is changed from GB211C-B to GB211C-A, please confirm the suction port connecting position of the main unit and switch over.

Remove the 4 screws, open the right side door, and assemble as shown in the figure below.

※ The default setting of this unit is GB211C-A. Therefore, you can skip this step when installing GB211C-A.




4. Operating procedures


Preparation before operation (GB211C+GF200 installation)

GB211C main unit + GF200 granulation set installation procedures (GB211C-B)

- (1) Please remove the 4 screws, open the right side door, and switch the position of the suction port. Assemble as shown in the figure below. (Make sure the filter case presses the microswitch before switching to the granulation control system)




- (2) Please insert the nozzle bushing in the nozzle insertion port on the main unit. The direction is shown as below.



4. Operating procedures

Preparation before operation (GB211C+GF200 installation)

- (3) Unpack the attachment GF200 and make sure that there is no broken glass or missing items.
- (4) Please install each part according to the serial number in the figure below.




- (5) Turn on the power switch, keep pressing the **Raise** button to raise the lifter. Stop raising the lifter when the upper part of the flow layer chamber is close to the flange of the filter chamber. Then press the **Raise** button again and again to make the upper part flange of the flow layer chamber fit to the packing and the filter chamber flange. Keep pressing the **Raise** button until the lifter stops and the **BLOWER** button is displayed.
- (6) Please install the outlet temperature sensor at the connector of the central sleeve, and insert the plug into the connector on the side of the main unit.
- (7) Please use the suction hose to connect the pipe connector of the central sleeve to the connection port at the lifter, and fix it with hose clamps.

4. Operating procedures

Preparation before operation (GB211C+GF200 installation)

(8) Insert the granulating spray nozzle vertically into the top mounting hole as shown below. Insert the connecting pipe with **NOZZLE B** label into the **NOZZLE A** port, and the other end of **NOZZLE B** connecting pipe is connected to the nozzle quick connector.




4. Operating procedures

GB211C+GF301C operating method


The following will use the standard sample setting method as an example for reference.
Sodium chloride water solution, 100g, solid concentration 5wt %

(1) Mini spray accessories are installed in the order described above (P.22 to P.26).




(2) Turn ON the ELB on the right side of the main unit.

ATOMIZING AIR




(3) Turn ON the **Power** switch on the operation panel of the main unit, touch screen display.



(4) Please switch the **Blower** ON to jump to the running screen for blower output power setting.

e.g.: blower output power 20.0% (the air volume is about 0.4m³/min)

(Please refer to P. 44 "The corresponding table as below is for blower output power and average dry air amount")




※ Click the **Earth** button at the upper right of the touch screen to enter the selection screen of language switching, and you can select the display language (Chinese, Japanese, English). After selecting, click the **x** button at the upper right of the screen to close the screen.



4. Operating procedures

GB211C+GF301C operating method



(5) There are inlet temp. controller and outlet temp. controller respectively on the operation screen, which are used for display and temp. setting.

By clicking the icon of **inlet temp.** or **outlet temp.**, you can select the inlet temp. or outlet temp. at will. After selecting, the inlet control or outlet control indicator lamp in the inlet temp. or outlet temp. controller will light up.


PV value displays the real-time temp. of the temp. sensor, SP value is black characters on white background, click to set the operating temp.

※The setting range of each temp. controller is different.

Inlet temp. setting range: 0 — 240°C


Outlet temp. setting range: 0 — 100°C

e.g.: select the inlet control, set the inlet temp. as 150°C.




(6) Turn ON the **Heater** switch, and the heater begins to work.

※ After heating, the function switch between inlet temp. control and outlet temp. control will become invalid to prevent the mistake contact in the experiment. If need to switch the control function, firstly turn OFF the **Heater** switch, switch to inlet temp. control or outlet temp. control, and then turn ON the **Heater** switch.




(7) Set the liquid-sending hose as shown on the left, turn the pull rod of the liquid-sending pump CCW to open the pump head, put the liquid-sending hose in it, and then turn the pull rod CW to make the liquid-sending hose stuck. Insert the other end of the liquid-sending hose into the liquid-sending interface of the spray nozzle. Please use the distilled water as the sample.




4. Operating procedures

GB211C+GF301C operating method

ATOMIZING AIR



ON



(8) After the inlet and outlet temp. reach the desired temp., set the spray flow and liquid-sending speed, and turn ON the **liquid-sending** switch to transport the distilled water.

e.g.: When the outlet temp. reaches about 80°C, the spray flow is set as 10L/min and the liquid-sending speed is set as 10rpm (about 2.5mL/min). (Please refer to P.44 “The corresponding table as below is for rotate speed of liquid sending pump and average liquid sending amount”) Adjust the liquid-sending speed to make the outlet temp. be slightly lower than 75°C.

(9) In order to stabilize the outlet temp. and inlet temp. at the desired temp., please adjust the dry air volume, spray flow and liquid-sending speed again.

e.g.: Adjust the liquid-sending speed to make the outlet temp. be slightly lower than 75°C.

— Hint —


- When the inlet temp. is constant, the influences of each setting on the outlet temp. are as follows.
Sample liquid-sending volume → small: outlet temp. → high
Dry air volume → large: outlet temp. → high
Sample concentration (external factor) → high: outlet temp. → high
- If increase the spray flow, the spray droplets will become micronized.
- The volume of spray flow is in direct proportion to the diameter of nozzle orifice.
- When the samples are replaced from the distilled water to the actually used samples, the outlet temp. will become slightly higher due to the non-evaporative part (solid part).

(10) When the outlet temp. is stable, replace the samples with the actually used samples. At this point, the outlet temp. will change more or less, if necessary, please adjust the liquid-sending speed again.

e.g.: Replace the samples with 100g sodium chloride 5% solution

4. Operating procedures

GB211C+GF301C operating method




- (11) During normal spraying, when the sample cannot be sprayed, the orifice of the spray nozzle may be blocked. Operate the **needle** button to squeeze out the blockage, or set the automatic needle to prevent the orifice of the spray nozzle from being blocked. Please refer to P. 41 "Use of automatic needle spray nozzle (for GB211C+GF301C)".

When the blockage at the orifice of the spray nozzle cannot be cleaned out by needle, please click the **Liquid Sending** button to stop the liquid sending, and then long press the **Liquid back** button to make the sample in the liquid-sending hose return to the sample container.

According to P.34 "End process", stop the unit.


According to P.50 "About cleaning after use", thoroughly clean the spray nozzle. After drying and assembling, continue to test.



- (12) During normal spraying, if the conical misty samples sprayed from the nozzle becomes irregular, it may be due to the attachment of samples near the orifice of the spray nozzle. Please press the **Pulse jet** button on the nozzle to see whether the attachment can be squeezed out.

If it still cannot be cleaned out, please click the **Liquid Sending** button to stop the liquid sending, and then long press the **Liquid back** button to make the sample in the liquid-sending hose return to the sample container.

Turn the **Heater** button OFF, wait until the inlet temperature is lower than 60 ° C and the outlet temperature is lower than 50 ° C, take out the nozzle and scrape off the attachment near the nozzle with a knife.




- ※ The **Liquid back** button is a manual/auto button. Press the button, the liquid-sending pump reverses, and release the button, the pump stops reversing. Long press the button for 5 secs, the pump will automatically reverse, even if release the button, it will not stop. At this time, click the **Liquid back** button again to stop the reverse.

- ※ Liquid Sending and Liquid back cannot operate at the same time. When the liquid-sending pump is sending liquid, the operation of **Liquid back** button is invalid. Similarly, when the liquid-sending pump is making liquid back, the operation of **Liquid Sending** button is also invalid.

4. Operating procedures

GB211C+GF301C operating method

ATOMIZING AIR




OFF


~ End process ~

(13) When the sample liquid sending is finished, replace the samples with the distilled water again to clean the nozzle. Clean for about 5mins, turn OFF the **liquid-sending** switch, and then adjust the spray flow to 0.

e.g.: After about 15mins, when the process of 100g sending liquid is finished, please replace the samples with the distilled water.



(14) Turn the **Heater** button OFF, when the inlet temp. is lower than 60°C and the outlet temp. is lower than 50°C, please turn the **Blower** button OFF. Otherwise, you will be prompted:




※ When the inlet temp. is above 60°C or the outlet temp. is above 50°C, do not stop the operation of the blower by forcibly cutting off the power. Otherwise, the malfunction may occur.

ATOMIZING AIR



(15) Turn OFF the **Power** switch.



(16) Remove the container fixing clamp and take out the product collecting container. At this point, please note that the back of the cyclone cover also has powder attached.

e.g.: Amount of collected powder is about 3-3.5g.

(17) Wash the containers according to P.50 "About cleaning after use".


※ When use a sample such as sodium chloride that corrodes metals, break down the spray nozzle and wash thoroughly.

4. Operating procedures

GB211C+GF200 operating method


The following will use the standard sample setting method as an example for reference.

(1) The granulation set is installed in the order described above (P.27 to P.29).




(2) Turn ON the ELB on the right side of the main unit.

ATOMIZING AIR




(3) Turn ON the **Power** switch on the operation panel of the main unit, touch screen display.



(4) Please switch the **blower** ON to jump to the running screen for blower output power setting.

e.g.: blower output power 20.0% (the air volume is about 0.4m³/min)

(Please refer to P. 44 "The corresponding table as below is for blower output power and average dry air amount")




※ Click the **Earth** button at the upper right of the touch screen to enter the selection screen of language switching, and you can select the display language (Chinese, Japanese, English). After selecting, click the **x** button at the upper right of the screen to close the screen.



4. Operating procedures

GB211C+GF200 operating method



(5) There are inlet temp. controller and outlet temp. controller respectively on the operation screen, which are used for display and temp. setting.

By clicking the icon of **inlet temp.** or **outlet temp.**, you can select the inlet temp. or outlet temp. at will. After selecting, the inlet control or outlet control indicator lamp in the inlet temp. or outlet temp. controller will light up.


PV value displays the real-time temp. of the temp. sensor, SP value is black characters on white background, click to set the operating temp.

※The setting range of each temp. controller is different.

Inlet temp. setting range: 0 — 240°C


Outlet temp. setting range: 0 — 100°C

e.g.: select the inlet control, set the inlet temp. as 150°C.




(6) Turn ON the **Heater** switch, and the heater begins to work.

※ After heating, the function switch between inlet temp. control and outlet temp. control will become invalid to prevent the mistake contact in the experiment. If need to switch the control function, firstly turn OFF the **Heater** switch, switch to inlet temp. control or outlet temp. control, and then turn ON the **Heater** switch.




(7) Set the liquid-sending hose as shown on the left, turn the pull rod of the liquid-sending pump CCW to open the pump head, put the liquid-sending hose in it, and then turn the pull rod CW to make the liquid-sending hose stuck.

※ When the sample cannot be sprayed, the orifice of the spray nozzle may be blocked. At this time, press the needle of the nozzle (P.50 “About cleaning after use” Spray nozzle exploded drawing). The needle (P.50 “About cleaning after use” Spray nozzle exploded drawing) will squeeze out the blockage at the orifice.



4. Operating procedures

GB211C+GF200 operating method




- (8) Place the liquid sending hose in a container that contains binder, turn the **liquid-sending** switch ON, and when the binder is close to the nozzle inlet, turn the **liquid-sending** switch OFF. At this time, adjust the liquid-sending speed of the pump to an appropriate setting.

※ Please refer to the table below for the use of binder.

e.g.: Binder: Polyvinyl carbazole 50g
(Actual usage is about 20g)
Adjust the liquid-sending speed to 12mL/min


Types and characteristics of binders (reference material)

Types	Characteristics
Gelatin	Due to the low concentration and weak viscosity, the high concentration solution needs to be heated before spraying.
Dextrin	Weak viscosity, very helpful in forming tablets
Potato starch	This is superior in granule property and inexpensive. Utilized in medical and food fields.
Gum Arabic	Spray after heating. This requires a lot of binder.
CMC (carboxymethyl cellulose)	This presents a higher viscosity at a lower temperature. This tends to leave more powder residues.
MC (methyl cellulose)	Strong viscosity, suitable for coarse particles granulation



- (9) After the outlet temperature is stable, please turn off the **Heater** switch to stop heating and then turn off the **Blower** switch to stop the air supply.


e.g.: The outlet temperature is stable at about 60°C



- (10) Click the **Lower** button to lower the flow layer chamber, and place the sample uniformly on the microporous plate of the flow layer chamber.

e.g.: Sample Sintered alumina 300g

※ As the flow layer chamber is in a high temperature state, please wear the heat resistant gloves for operation.




- (11) Click the **Raise** button to raise the flow layer chamber, and make the flow layer chamber be in close contact with the filter.

※ Perform the operations (9) to (11) quickly. If the inlet temperature is above 60°C or the outlet temperature is above 50°C for a long time without air blowing, a fault may occur.

4. Operating procedures


GB211C+GF200 operating method



(12) Turn on the **Blower** and **Heater** switches, start the operation of blower and heater, and make the sample flow. Adjust the air volume so that the height of the flow layer is consistent with the position of the silicone rubber cap of the flow layer chamber.

— Operational hint —

Too much powder attached to the filter may decrease the air volume. Press the **Pulse jet** switch, the pressurized air will be blown into the filter to remove the powder. Please use the **Pulse jet** switch regularly to remove the powder attached to the filter.



ATOMIZING AIR



ON

(13) After the outlet temperature is stable, set the spray flow and turn on the **liquid-sending** switch.

e.g.: Spray flow 10L/min Liquid-sending speed 12mL/min

When the flow of the sample slows down, turn off the **Liquid-sending** switch and reduce the spray flow to a minimum.

The nozzle does not work after the liquid-sending is completely stopped. Therefore, to avoid clogging the nozzle, use the **Liquid back** switch to return the binder to a position close to the nozzle hose connection port. The **Liquid back** switch is activated when it is pressed.


e.g.: After spraying for 30secs, turn the **liquid-sending** switch OFF.

Reduce the spray flow to 1L/min.

※ The **Liquid back** button is a manual/auto button.

Press the button, the liquid-sending pump reverses, and release the button, the pump stops reversing. Long press the button for 5 secs, the pump will automatically reverse, even if release the button, it will not stop. At this time, click the **Liquid back** button again to stop the reverse.

※ Liquid Sending and Liquid back cannot operate at the same time. When the liquid-sending pump is sending liquid, the operation of **Liquid back** button is invalid. Similarly, when the liquid-sending pump is making liquid back, the operation of **Liquid Sending** button is also invalid.



4. Operating procedures


GB211C+GF200 operating method

(14) Repeat spraying and drying operations in step (13) until the desired particle size is obtained.

The particle size will become gradually larger in the second session of spraying and drying and after, and thus you need to gradually shorten the time of binder spraying and gradually increase air amount.


In order to simplify the above operations, the granulation automatic liquid sending function can be used. Please refer to P.42 "Granulation automatic liquid sending function (for GB211C+GF200)".

e.g.: Repeat spraying and drying of polyvinyl carbazole for 5 times (approx. 20g is used in total)




Operational hint

When flow condition of samples has degraded, increase the blower air amount and press the **Mixer** switch in order to disperse the samples evenly to always keep a good flow conditions. And when the spraying is unstable, press the upper tip of the nozzle to remove the blockage at the nozzle tip. The **Mixer** switch is activated while it is pressed. Press the **Mixer** switch for 3secs or more if you want to operate continuously. Press it again to cancel the continuous operation. The mixer will also stop when the blower is turned OFF.




(15) When the sample reaches the desired particle size, turn off the **liquid-sending** switch, reduce the spray flow to a minimum, and then through the **liquid back** switch to return the binder to a position close to the nozzle hose connection port.

※ **Note that the heating cannot be stopped immediately at this time. The sample should be fully dried before stop heating.**



(16) After the sample is fully dried, please turn off the **heater** switch and stop heating. Then reduce the output power of the blower, and use the **pulse jet** switch to remove the attachment on the filter.




(17) Turn off the **blower** switch to stop the air supply. Then click the **Lower** button to lower the flow layer chamber. Lower the flow layer chamber to the lowest position, take out the chamber and collect the granulated samples.

(18) If a large number of samples need to be processed in multiple batches, repeat the steps (9) to (17).


4. Operating procedures

GB211C+GF200 operating method

～ End process ～




(19) When the sample granulation is finished, put the flow layer chamber back on the lifter, click the **Raise** button to raise the flow layer chamber, make the flow layer chamber and filter in close contact. (Follow the previous sequence [P.28] to operate)



(20) Turn on the **Blower** switch. Long press the **liquid back** switch, return the remaining binder in the liquid-sending hose to the container, and then turn off the **liquid back** switch.


(21) When the inlet temperature is lower than 60°C and the outlet temperature is lower than 50°C, turn off the **blower** switch to stop the air supply.

※ When the inlet temperature is above 60°C or the outlet temperature is above 50°C, if stop the operation, a malfunction may occur.



(22) Click the **Lower** button to lower the flow layer chamber. Lower the flow layer chamber to the lowest position, take out the chamber.

ATOMIZING AIR




(23) Please turn OFF the **Power** switch.

(24) Please clean the containers according to P.50 "About cleaning after use".


4. Operating procedures

Use of automatic needle spray nozzle (for GB211C+GF301C)




During normal spraying, if the sample cannot be sprayed, possibly because the orifice of the spray nozzle is blocked.

When using the automatic needle spray nozzle, press the **Needle** button on the operation screen, the automatic needle will continue to move at the frequency of poking once per second until the **Needle** button is released.



When it is estimated that the blocking will occur in the test, click the icon of needle controller to turn on the automatic needle function while switching from the distilled water to the samples. The manual mode indicator lamp of the needle controller goes off and the automatic mode indicator lamp lights on. When the automatic needle function is not used, just need to click the icon of needle controller to turn off the function.




When the automatic needle function is enabled, the countdown of the needle timer starts. When the timing (PV value) is 0, the needle activates. Later the timer resets and starts the next timing.

The time of the timer can be set by clicking the SP value of the needle controller.

When the timing is reached, the action times of automatic needle can be set by clicking the action times setting button. Its setting method is to click once, then the needle action times increases by 1, successively increasing to the maximum value of 3, and then click again to return to the value of 1, and cycle again.


※ The maximum action times of automatic needle can be set for three times. The test data show that if it acts for three times, the blockage still cannot be cleared, then there is no effect to act for more times. At this time, you can consider reducing the interval time of the timer, and you can activate the needle when the blockage has not accumulated to be unable to clear it, so as to prevent blockage.



When using the automatic needle function, if the blockage is found between the two actions of automatic needle, the manual **needle** button can also be clicked, and its action mode is the same as manual needle. The action of manual needle does not affect the timing and action of automatic needle. Manual and automatic modes can be carried out in parallel without cross influence.

4. Operating procedures

Granulation automatic liquid sending function (for GB211C+GF200)





In the process of granulation, repeatedly transport the binder to achieve the particle size requirements. The granulation automatic liquid sending function can simplify the operations and improve the operation precision.

In the granulation operation screen, click the icon of Liquid pump.




In the pop-up screen, set the time of auto liquid sending and auto liquid back of binder. After setting, click the **AUTO liquid feeding** button. When the button turns green, it means that the granulation automatic liquid sending function is started (see the figure below). Click the upper right corner of the Setting window to close the Setting window.



According to the above picture, when in auto liquid sending, liquid sending for 15secs and liquid back for 5secs. Click the Liquid sending button on the operation panel, and the liquid pump will firstly send liquid for 15secs, then reverse for 5secs after liquid sending, and return the binder to a position close to the nozzle hose connection port. After liquid back, the liquid pump stops and waits for the next operation.

If need to stop this function, return to the Setting window again and click the **AUTO liquid feeding** button. When the button turns yellow, it means that the granulation automatic liquid sending function is stopped (see the figure below).



- ✘ Although this quick function can reduce the difficulty of granulation operation and improve the precision, the particle size will become gradually larger in the second session of spraying and drying and after, and thus you need to gradually shorten the time of binder sending and gradually increase air amount.
- ✘ This quick function can only be used in granulation mode.

4. Operating procedures

If want to interrupt the sample processing, or when nozzle blockage occurs

If want to interrupt the sample processing, or when nozzle blockage occurs, please stop sending liquid according to P.34 "End process".

In addition, if want to process another samples, please firstly recycle the products in the collecting container, clean it according to P.50 "About cleaning after use", and then change to another sample for spray test.

4. Operating procedures

The relation between rotate speed of liquid sending pump and liquid amount/between blower power and dry air amount (reference)

The corresponding table as below is for rotate speed of liquid sending pump and average liquid sending amount (the calibration liquid is water at 23°C).

Please refer to it during operation. Please pay attention that the density and viscosity of solution have a great influence on the liquid sending amount. If the liquid sending amount is less, the liquid sending pump head may be not pressed tightly, or there is block, deform or damage for the liquid sending hose.

Rotate speed of liquid sending pump	Avg. liquid sending amount (ml/min)
0.0	0
10.0	2.5
20.0	5.1
30.0	7.5
40.0	9.9
50.0	12.5
60.0	14.7
70.0	17.5
80.0	20.3
90.0	23.3
100.0	26.4

The corresponding table as below is for blower output power and average dry air amount. (Measured value of blower operating independently at 220V ~ standard voltage)

Please refer to it during operation. If the air amount is too low, the blower filter or the suction filter may be blocked. Clean the filter according to the maintenance method (see P. 51 "Filter Cleaning"). In addition, please note that when the power frequency is different at 50Hz or 60Hz, the dry air amount corresponding to each power is also different. For details, please refer to the following table.

In the case of 230V/50Hz power source	
Blower output power (%)	Avg. dry air amount (m ³ /min)
5	0.18
10	0.26
20	0.41
30	0.53
40	0.62
50	0.70
60	0.77
70	0.83
80	0.87
90	0.89
100	0.91

In the case of 230V/60Hz power source	
Blower output power (%)	Avg. dry air amount (m ³ /min)
5	0.17
10	0.26
20	0.41
30	0.53
40	0.63
50	0.71
60	0.77
70	0.83
80	0.87
90	0.90
100	0.91


4. Operating procedures

The relationship between the output power of the blower and the maximum operating temperature of the inlet (Reference)

The following table is the corresponding curve between the output power of the blower and the maximum operating temperature of the inlet. (Under the condition of 220V supply voltage)

Please use it as reference during operation.

- ※ When the air volume is too large, the wind speed is too fast, resulting in a reduction of the heat transfer time of the air in the heating sleeve, and the air cannot raise enough temperature in a short time, and eventually the maximum operating temperature of the inlet will drop.



4. Operating procedures

Calibration of temperature sensor

After long-term use, the temperature sensor will have temperature drift, resulting in the deviation of operating temperature.

In this case, you can calibrate the temperature sensor to restore the control temperature to the normal state.

If need to calibrate the temperature sensor, please contact the agent or Yamato Scientific. Calibration of temperature sensor is a paid service.


After calibration, the correction parameters of inlet temperature or outlet temperature (in the red box below) will be modified according to the calibration report of the sensor. Please do not change them before the next calibration.

The screenshot displays the Yamato HMI interface for the GB211C-A stand-alone operation. At the top, the Yamato logo is on the left, and the date/time '2023/03/16 13:35:29' is on the right. Below the logo, it shows 'HMI Ver:1.01 PLC Ver:1.01'. The main interface is divided into several sections: a 'Time setting' section with fields for Year (2000), Month (01), Day (01), Hour (00), Minute (00), and Second (00); an 'Inlet temp' section with 'CAL:' (0.0 °C) and 'SC:' (1.000) fields; an 'Outlet temp' section with 'CAL:' (0.0 °C) and 'SC:' (1.000) fields; a 'Frequency of power supply' section with '50Hz' and '60Hz' buttons; and a '2 second recording period' section. A red box highlights the 'CAL:' and 'SC:' fields for both Inlet and Outlet temperatures. At the bottom left is a 'Save parameters' button, and at the bottom right is a home icon.

5. Handling Precautions


Warning

1. Substances that cannot be used


-  Never use explosive, flammable or substance that contains them. Otherwise, an explosion or a fire may occur. See P. 76 "15. List of Dangerous Substances".

Connect GB211C+GF301C (GB211C-A) with the optional GAS series product to form an enclosed and low-oxygen circulation system, which is able to use the organic solvent sprays without the risk of explosion. When using the organic solvents, pay special attention to their explosion conditions, especially the mixture of multiple organic solvents. Please read the GAS series product instruction manual for operations.

2. If a problem occurs


-  If smoke or strange odor comes out of this unit, turn off the main power supply right away, and pull out the plug. Immediately contact the sales agent or our business office for maintenance. If continue to operate, fire or electric shock may result. Never perform repair work yourself.

3. Do not touch the part with high temperature


-  During or just after operation, the temperature of drying chamber, cyclone and surrounding area is higher. Do not touch these parts to avoid scalding.

Caution


1. Do not put anything on this unit

-  Do not put anything on this unit. It will cause injury if fall.


2. During a thunder storm

-  During a thunderstorm, turn off the power key immediately, then turn off the circuit breaker and the main power. If this procedure is not followed, fire or electrical shock may be caused.


3. Do not use corrosive sample

-  Stainless steel SUS304 is used for the interior; however, it may be corroded by strong acid etc. In addition, the sealing strip and silicon rubber may be corroded by some kind of solvent like acid, alkali, oil, halogen, etc. Do not use the sample containing those substances.


4. Recovery after power outage

-  During operation, the machine stops due to power outage. When the power is supplied again, it will be restored to the initial state.

5. Take measures against toppling and falling

-  It may cause injure to a person if this unit falls down or moves by a sudden earthquake, impact, etc. For the sake of safety, please take measures against toppling and falling.

6. Do not disassemble glassware and pipes when the outlet temperature is above 50℃.

-  Do not disassemble glassware and pipes when the outlet temperature is above 50℃. Otherwise, there is a risk of scalding.
The heating pipe inside the machine expands in size at high temperatures. At this time, if disassemble the glassware and pipes for cleaning, the glassware and pipes will shrink due to cold, and the port size will become smaller. The size mismatch will occur when they are re-installed, and the installation by force will cause damage.

5. Handling Precautions

Drying Method under Appropriate Condition

- (1) The best appropriate drying condition is differed depending on the sample to be dried. There are some data for reference. Please consult the agent for details.
- (2) The attachment on the drying chamber is particularly significant because of high sample concentration, low inlet temperature, too high or too low spray air pressure, or too much liquid sending amount of the sample. When there is an abnormality, the above reasons can be taken into account, and please adjust appropriately.
- (3) During operation, when the spray direction is changed due to the sample attachment at the tip of spray nozzle, turn ON the pulse jet switch, and blow off the attachment from the tip of the nozzle using the pressurizing air. If the attachment cannot be blown off, disassemble the spray nozzle and use ultrasound, etc. to clean it.
- (4) As for the reasons of the sample attachment on the cyclone part, it can consider about whether the solvent (distilled water or ion-exchange water) is not fully evaporated, or the unique characteristics of the sample (low melting point, absorbability, etc.). In order to make the moisture content of the powder as little as possible, the heat of the sample is the more the better, it's able to increase the inlet temperature and dry air flow, or reduce the liquid sending amount of the sample. That is to reduce the difference between the inlet temperature and the outlet temperature. When the sample has its unique characteristics, please add additives to adjust the sample.

- (5) In the case that the hygroscopicity is high, the product may become the moist powder in the container. Change the drying condition following the method in (4), or, if required, heat up the container for product before operation.
- (6) The orifice of the spray nozzle is 460 μ . If the sample is blocked with suspension at orifice part impetuously, use the 508 μ and 711 μ nozzles prepared for the orifice as optional (Nozzle main body P.错误!未定义书签。 “错误!未找到引用源。”, the nozzle main body, the needle, and the ring in the exploded view of the spray nozzle are common with the 406 μ nozzle) These 508 μ and 711 μ nozzles are differed on the point of the size of the spray pattern and particle diameter of the drop slightly compared to the 406 μ one, and these differences may affect the interference status. Refer to the Graph 1 for the relation between spray air pressure and spray airflow rate (atmospheric conversion).
When the gas source pressure is different, the maximum adjustable spray flow is different. Take the standard orifice 711 μ as an example, the maximum adjustable spray flow under different pressures is shown in the right figure.




Fig.1

MPa

- (7) The too small powder (few μ or less) among dried ones is impossible to be collected, and exhausted to the outside through the blower. If this exhausted amount of the too small powder becomes more, decrease either spray airflow rate or spray air pressure. Also, since the particle diameter becomes smaller as the density of the sample is lower, adjust the density of the sample if required.

5. Handling Precautions

Caution during operation

- (1) When connecting the power supply, be sure to ground it.
- (2) Pressurized air should be controlled at a stable pressure of 0.3-0.6MPa. Using pressure greater than 0.6MPa will damage the pipe.
- (3) The outlet temperature will cause deterioration of the material of the suction/exhaust hose, the material of the filter and the performance of the blower. Please do not use it over 130℃ for a long time. The heater will stop automatically when the temperature exceeds 140℃.
- (4) Check the glass chambers are fixed to the specified position with no gap, and then turn on the switches of blower and heater.
- (5) The unit is not explosion-proof. Do not use any solvent that contains flammable organic solvents as the samples. ※ **When you use an organic solvent for GB211C+GF301C (GB211C-A) , connect the optional organic solvent recovery unit (GAS).**
- (6) When the heater is ON, do not expose the end cap of drying chamber and mounting port of the spray nozzle to the non-guard status, and do supply the air to the heater part for at least 0.1-0.2m³/min.
- (7) During normal spraying, when the sample cannot be sprayed, the orifice of the spray nozzle may be blocked. Operate the needle button to squeeze out the blockage, or set the automatic needle to prevent the orifice of the spray nozzle from being blocked. Please refer to P. 41 "Use of automatic needle spray nozzle (for GB211C+GF301C)".
- (8) When the liquid sending pump cannot send samples, please check if the sample hose is crushed at the roller of the pump, the inner wall of the hose is adhered tightly without restoration, or the inner of the nozzle is blocked. Run again after all the above conditions return to normal.
- (9) Do not perform unattended operation. Because idling after the sample being used up and nozzle blockage will cause the outlet temperature to rise, or the sample hose to fall off from the nozzle resulting in sample outflow, unexpected accidents may occur.
- (10) If it is a silicone hose, easy to be eroded by halogen solvents and acids (concentrated), then the expansion or fracture may occur, need to pay attention during the operation.
- (11) When the inlet temperature is set at high temperature, if the air flow of the blower is excessive, it may not reach the set temperature due to the capacity of the heater. At this time, either turn down the air flow, or increase the set temperature for operation. The set value is inconsistent with the actual inlet temperature when increasing the set temperature. The heater will stop automatically when the inlet temperature exceeds 260℃ or when the outlet temperature exceeds 140℃.
If this unit is not operated, turn "OFF" the earth leakage breaker on the right side of the unit.
- (12) Depending on the sample used, operating environment and conditions, the cyclone may be prone to generate the static electricity. Either install auxiliary grounding terminals at 3 positions of the clamp of the cyclone connection, or install the anti-static brush (optional) on the main body of the cyclone.
- (13) If there is a leakage between the product collecting container and the metal cover at the lower of the cyclone, the dried powder will accumulate at the lower of the cyclone and may not fall into the product container. Therefore, pay special attention when installing the product container.
- (14) Because the capacity of the product container is about 750mL, when the collected powder is about 200-250g, it accounts for about 80% of its capacity. If continue to operate, it will cause the reduction of the collecting efficiency of the powder. Stop the operation for a while, and take out the collected powder.
- (15) Due to the different samples, the cyclone part may be prone to generate the static electricity, please use appropriate methods to remove it. Although the grounding effect of wrapping the metal wire on the glass part is very obvious, using the anti-static brush (optional) to vertically contact with the cyclone part is more convenient.

6. Maintenance Method

Daily Inspection and Maintenance

Warning



- Disconnect the power cable from the power source when doing an inspection or maintenance unless needed.
- Perform the daily inspection and maintenance after the machine is restored to normal temperature.
- Do not disassemble this unit.

⚠ Caution

- Use a well-drained soft cloth to wipe dirt on this unit. Do not use benzene, thinner or cleanser for wiping. Deformation, deterioration or color change may result in.

About cleaning after use

- (1) After completing the operation, remove the attachments following the process “Preparation before operation” on P. 22 in reverse order.
- (2) Clean the portion of attachment to which the powder is adhered.
- (3) Flow the distilled water into the sample tube by pressing the pump switch, and remove the contaminant attached to the inner of the part.
- (4) Remove the spray air hose and sample liquid sending hose from the spray nozzle, and disassemble the nozzle as per the following figure. After disassembling, clean it using the supersonic cleaner. Remaining the contaminant to the inner of the part may cause the insufficient spray. Therefore, clean it completely.



6. Maintenance Method


Daily Inspection and Maintenance

Filter Cleaning


● The filter in blower

Clean up the filter in blower periodically.

1. Please open the door in the front of the main body and remove the connecting pipe of the blower.



2. Please turn on the fixed fittings at the front cover of the blower and take out the filter.



3. The followings are the cleaning procedures of the filter.
 - (1) Wash the filter pressing in the water repeatedly, and air-dry it.
 - (2) Compressed air blowing.
 - (3) Vacuum cleaning with a cleaner.
 - (4) Press washing the filter after being immersed into the solvent that hot water (approx. 40°C) and neutral detergent are mixed at a rate of 5:95 one whole day and night, then rinse it with water and air-dry it.
4. Please follow the reverse order of disassembly to install.
Turn the soft side of the filter to windward when installing the filter.



6. Maintenance Method

Daily Inspection and Maintenance

● Suction filter

Regularly clean the suction filter.

1. Remove the 4 screws, open the right side door, and remove the knurled screws on the left and right sides of the suction filter.



2. Pull the filter out of the suction filter case. (Refer to the exploded drawing above)
3. The followings are the cleaning procedures of the filter.
 - (1) Wash the filter pressing in the water repeatedly, and air-dry it.
 - (2) Compressed air blowing.
 - (3) Vacuum cleaning with a cleaner.
 - (4) Press washing the filter after being immersed into the solvent that warm water (approx. 40°C) and neutral detergent are mixed at a rate of 5:95 one whole day and night, then rinse it with water and air-dry it.
4. Please follow the reverse order of disassembly to install.
Turn the soft side of the filter upward when installing the filter.


6. Maintenance Method

Daily Inspection and Maintenance

Monthly maintenance

Check the earth leakage breaker function.

- Connect the power cord and power on before test.
- Turn the breaker on.
- Push the red test switch by a ballpoint pen etc. If there is no problem, the earth leakage breaker will be turned off.



6. Maintenance Method

Daily Inspection and Maintenance


About the use of PLC batteries

● Selection of battery

When need to use a battery, please use a button battery with a production date less than two years.

The model of button battery is CR2032.

The reading of production date




The left character: the rightmost digit of Year

The right character: month (0: October, Y: November, Z: December)

Example: "8Z" indicates the production of December, 2018

● Installation of battery

Turn off the ELB, remove the 4 screws from the PDC door on the right side of the unit body by using a M4 cross screwdriver, and then open the PDC with the right side of the door as the shaft. See as below:




6. Maintenance Method

Daily Inspection and Maintenance

Open the battery holder of the CPU unit, place the battery into the battery holder and close the battery holder. The installation direction of the battery is shown in the following figure:

- ※ When replacing with a new battery, take out the old battery and then put in the new battery. The battery replacement must be completed within 5 minutes after powering off the CPU unit to ensure that the clock data is not lost. If this step is not completed within 5 minutes, the clock will stop and the time will be reset to "2001-01-01 01:01:01 Sunday".




After the installation is completed, perform the reverse operations to restore the device to its original state.

● Battery life

The maximum battery life after installation is 3 years at 25°C, regardless of whether the device is powered on or not. If used at higher temperatures, the battery life will be shortened.

● Low battery power indicator

When the battery power is about to be exhausted, the ERR/ALM indicator lamp in the front of the CPU unit will blink.



7. Long storage and disposal

When not using this unit for long term / When disposing



Caution

When not using this unit for long term...

- Turn off the earth leakage breaker and original power source for safe without fail. Also, store the glass unit after removing it from the main unit. When the glass unit is contacted to the external, it may cause the breakage.



Warning

When disposing...

- Keep out of reach of children.
- Remove the power cord.

Matters to consider when disposing of the unit

Environmental protection should be considered

We request you to disassemble this unit as possible and recycle the reusable parts considering to the environmental protection. The feature components of this unit and materials used are listed below.

Component Name	Material
Parts of Main Unit	
Exterior	Cold rolled steel plate with surface coating
Insulating material	Ceramic fiber
Sample tray	Stainless steel
Label	Polyethylene (PET) resin film
Hose	Silicon rubber, Teflon
Electrical Parts	
Heater	Stainless steel and others
Motor	Iron, Aluminum, Copper wire and others
Circuit boards	Board, capacitor, resistor, transformer, etc.
Power cord & wiring materials and others	Synthetic rubber, resins
Sensor	Stainless steel and others

8. When a trouble occurs

Safety device and error indications

The table shows possible causes of activation of the safety unit and solutions.

[Error indication]

When an abnormality occurs to the inlet temperature controller or the outlet temperature controller, the touch screen at the operation panel displays the error screen. When an abnormality occurs, confirm the error content and implement appropriate solutions.

Display	Reasons	Solutions
Er01.PLC analog module failure	① The wire connection of the PLC analog module is loose ② The PLC analog module is damaged	① Power off and restart ② If it cannot reset after power off and restart, please contact our service department or agent.
Er02.Inlet temperature transmitter disconnection	① The wire connection of inlet temperature transmitter is loose ② The inlet temperature transmitter is damaged	① Power off and restart ② If it cannot reset after power off and restart, please contact our service department or agent.
Er03.Inlet temperature sensor disconnection	① The wire connection of inlet temperature sensor is loose ② The inlet temperature sensor is damaged	① Power off and restart ② If it cannot reset after power off and restart, please contact our service department or agent.
Er04. Granulation inlet temperature sensor disconnection	① The wire connection of granulation inlet temperature sensor is loose ② The granulation inlet temperature sensor is damaged	① Power off and restart ② If it cannot reset after power off and restart, please contact our service department or agent.
Er05.Outlet temperature transmitter disconnection	① The wire connection of outlet temperature transmitter is loose ② The outlet temperature transmitter is damaged	① Power off and restart ② If it cannot reset after power off and restart, please contact our service department or agent.
Er06.Outlet temperature sensor disconnection	① The outlet temperature sensor is not installed ② The wire connection of outlet temperature sensor is loose ③ The outlet temperature sensor is damaged	① Install the outlet temperature sensor, click the alarm reset ② If it cannot reset, please contact our service department or agent.
Er07.Blower current detector malfunction	① The blower current detector is disconnected ② The blower current detector is damaged	Please contact our service department or agent.
Er08.Blower does not alarm	① The output power of the blower is set too low ② The power supply voltage is too low ③ The blower main control relay is disconnected ④ The blower main control relay is damaged ⑤ The blower speed controller is disconnected ⑥ The blower speed controller is damaged	Please contact our service department or agent.

8. When a trouble occurs

Safety device and error indications

Display	Reasons	Solutions
Er09.Blower overload alarm	① The blower is stuck ② The blower is damaged ③ The blower current detector is damaged	Please contact our service department or agent.
Er12.SSR short circuit alarm	SSR short circuit	Please contact our service department or agent.
Er13.Heater disconnection alarm	① Heater disconnection ② Heater damage	Please contact our service department or agent.
Er14.Heater overheat alarm	Heater overheat (Overheat protector activates)	Please contact our service department or agent.
Er15.Liquid sending pump is overloaded	The pump head is stuck	① Clean the pump head ② Please contact our service department or agent.
Er16.Liquid sending pump is not running	① The wire connection of liquid sending pump is disconnected ② Liquid sending pump damage	Please contact our service department or agent.
Er17.Inlet temperature controller alarm	① The inlet temperature controller reports an error ② The parameters of the inlet temperature controller are abnormal	① Power off and restart ② If it cannot reset after power off and restart, please contact our service department or agent.
Er18.Outlet temperature controller alarm	① The outlet temperature controller reports an error ② The parameters of the outlet temperature controller are abnormal	① Power off and restart ② If it cannot reset after power off and restart, please contact our service department or agent.
Er19.Inlet temperature overheat alarm	① The inlet temperature exceeds 260℃	① Stop heating, blow to drop the inlet temperature, and then click the alarm reset button in the alarm screen. ② If it cannot reset, please contact our service department or agent.
Er20.Outlet temperature overheat alarm	① The outlet temperature exceeds 140℃	① Stop heating, blow to drop the outlet temperature, and then click the alarm reset button in the alarm screen. ② If it cannot reset, please contact our service department or agent.

※ When the temperature is abnormal, the blower ON, the heater OFF and the liquid sending pump stops. The same is true when the temperature sensor is disconnected. After troubleshooting, press the "alarm reset" button to release the alarm, still hold the mode that the blower ON, the heater OFF and the liquid sending pump stops. The abnormal display can be set to Chinese, Japanese or English by language switch.

8. When a trouble occurs

Safety device and error indications

When used in connection with GAS, the following table shows the reasons and solutions when the safety device activates:

Display	Reasons	Solutions
Er30. Communication with GAS is interrupted	① The wire connecting to GAS is loose ② The wire connecting to GAS is damaged	① Reconnect the connecting wire. ② Replace the connecting wire. ③ After power off and restart, if it cannot reset, please contact our service department or agent.

※ When GAS alarm occurs, GB211C will display the error prompt of GAS, but will not display the error code of GAS. Regarding GAS troubleshooting, please refer to GAS instruction manual.

8. When a trouble occurs

Trouble Shooting

In case of the following conditions

Symptoms	Possible causes	Countermeasures
The POWER does not turn ON	<ul style="list-style-type: none"> ● The ELB is OFF ● The power switch is OFF ● Malfunction of the power supply ● The power cord is disconnected ● Malfunction of power switch 	<ul style="list-style-type: none"> ● Turn the ELB ON ● Turn the power switch ON ● Check the power supply circuit ● Replace the cord ● Replace the power switch
The blower does not activate	<ul style="list-style-type: none"> ● The power supply voltage is too low ● The blower power is too low ● The blower connector is not correctly connected ● The blower input line is disconnected ● Blower switch failure ● Blower motor failure ● Blower motor brush failure ● Blower circuit and wiring failure 	<ul style="list-style-type: none"> ● Select the appropriate power supply ● Increase the output power of the blower ● Connect correctly ● Replace the input line ● Replace the touch screen, PLC or temperature controller ● Replace the motor or motor board ● Replace the brush ● Maintain or replace the part
The temperature cannot rise	<ul style="list-style-type: none"> ● The heater button is not ON ● The heater connector is not correctly connected ● Failure of other parts causes the protection circuit activation (error display) ● The blower switch is not ON ● Protection circuit activates ● Heater disconnection ● Heater switch failure ● Heater circuit and wiring failure 	<ul style="list-style-type: none"> ● Turn on the heater button ● Connect correctly ● Solve the problem and turn ON the switch ● Turn ON the blower switch, and then turn ON the heater switch ● Check if there is an alarm ● Replace the heater ● Replace the touch screen or PLC ● Maintain the part or replace the temperature controller
The liquid sending pump does not activate	<ul style="list-style-type: none"> ● Pump speed is set to 0 ● Pump switch failure ● Pump motor failure ● Pump circuit and wiring failure ● Nozzle installation failure 	<ul style="list-style-type: none"> ● Set the pump rotate speed ● Replace the touch screen or PLC ● Replace the motor or driver ● Maintain the part ● Confirm installation status of nozzle and correct
No air flow for spray	<ul style="list-style-type: none"> ● Flowmeter is not open ● Flowmeter switch failure ● Pressurized air source failure ● Hose connection failure ● Solenoid valve failure ● Pulse jet circuit and wiring failure 	<ul style="list-style-type: none"> ● Turn on the flowmeter and adjust ● Maintain or replace the part ● Replace the solenoid valve ● Replace the touch screen or PLC ● Maintain the part

8. When a trouble occurs

Trouble Shooting

Symptoms	Possible causes	Countermeasures
Temperature controller failure	<ul style="list-style-type: none">● The outlet temperature sensor is not installed● Defective display function● Sensor failure● Overheat prevention function failure	<ul style="list-style-type: none">● Install the outlet temperature sensor correctly● Maintain the part or replace the PLC● Replace the sensor● Lower the set temperature
Cannot reach the set temperature	<ul style="list-style-type: none">● Heater capacity is insufficient due to excessive dry air flow● Regulating circuit and wiring failure	<ul style="list-style-type: none">● No abnormality. During the high-temperature operation, either reduce the dry air flow, or increase the set value of operating temperature● Maintain the part or replace the PLC

- ◆ In case if the error other than listed above occurs, please immediately cut off the power supply, pull out the power cord, and contact the sales store or our company's business, customer service center.

9. After Service and Warranty

When requesting a repair

When requesting a repair

If any trouble occurs, immediately stop operation, turn the power switch off, pull out the power plug and contact your dealer, our sales office or our customer service center.

Information necessary for requesting a repair

- Model name of the product
 - Serial number
 - Date (y/m/d) of purchase
 - Description of trouble (as in detail as possible)
- } See the warranty card or the nameplate on the unit.
See the section “3. Names and functions of the parts” on page 9.

Be sure to indicate the warranty card to our service representative.

Warranty card (attached separately)

- Warranty card is given by your dealer or one of our sales offices and please fill in your dealer, date of purchase and other information, and then store it securely.
- Warranty period is one full year from the date of purchase. Repair service for free is available according to the conditions written on the warranty card.
- For repairs after the warranty period consult your dealer, one of our sales offices or our customer service center.
Paid repair service is available on your request when the product's functionality can be maintained by repair.

Minimum holding period of repair parts

The minimum holding period of repair parts for this product is seven years after end of production. Repair parts here refer to parts necessary for maintaining performance of the product.

10. Specification

Specifications of main unit

Model		GB211C	
Suitable for solvent recovery unit※1		GAS411C/510C (optional)	
Temperature		5°C to 35°C (Indoor use only)	
Altitude		Up to 2,000 meters	
Relative humidity		≤75%RH	
Drying object samples	※1	Solution, suspension, emulsion (for flammable and explosive substances, please use when connected with GAS)	
Spray mode		Two-fluid nozzle (orifice diameter about Φ0.7)	
Spray and hot air contact mode		Vertical downward spray and parallel flow	
Water evaporation		Max. about 1500mL / Hr	
Structure	Temperature controller	PID temperature controller	
	Heater	3.2kw (220V~) Stainless steel pipe heater	
	Blower	Series excited blower	
	Liquid sending pump	Peristaltic pump	
	Blowout mechanism for pressurized air	Continuous spraying	
	Automatic cleanout needle	Use the needle in the nozzle to realize the automatic orifice cleaning (The air cylinder in the nozzle is driven by solenoid valve, electronic timer and pressure air)	
	Nozzle blower	Blow off the powder attached to the nozzle tip (using of solenoid valve, pressure air)	
Control part	Temperature adjustment range	Inlet temperature: 0-240°C, Outlet temperature: 0-100°C	
	Temperature adjustment accuracy ※2	±1°C	
	Temperature display	Digital display of inlet temperature and outlet temperature (display accuracy 0.1°C)	
	Spray air flowmeter	Measurement range: 0-30L/min (the max. flow rate is related to the pressure of the spray gas source)	
	Adjustable range of liquid sending volume	0 – 26 ml/min variable (the max. drying capacity 25 ml/min)	
	Adjustable range of drying air volume ※2	0.2-0.9m³/min (when a single blower runs independently)	
Spec.	External dimension (mm)(WxDxH) ※3	760×420×1350	
	Power supply	Single phase 200-230V~ 50/60Hz 17-21A	
	Weight	Approx. 110kg	
Attached accessories	Outlet temperature sensor	1	Sample hose: silicone, I.D.2mm×O.D.4mm×1m 2
	Fuse 250V 2A	1	Polyvinyl chloride, I.D.2mm×O.D.4mm×1m 3
	Earth wire	1	Exhaust hose: made from vinyl chloride, I.D.50mm×2m 1
	Exhaust connector	1	Hose tie: #64 1
	Clamp	1	Hose: 5m (for connecting pressurized air) 2
	Warranty card	1	Hose clamp 1
	Instruction manual	1	Air tube A Φ 6×0.6 m (NOZZLE A) 1
			Air tube B Φ 6×0.6 m (NOZZLE B) 1

※1 GB211C+GF301C (GB211C-A) is special for water soluble solvent. When using the organic solvent, it needs to be connected with the optional GAS series product. **Also, please note that GAS series product cannot be connected to GB211C+GF200 (GB211C-B) .**

※2 Performance under 230V~ power supply conditions, with an ambient temperature of 5°C-35°C and a relative humidity of ≤75% RH.

※3 The outer dimension does not include the projection part.

10. Specification

Mini granulation accessories [GF200]	Model	GF200	
	Processing capacity	50-300g	
	Spray nozzle	Two-fluid nozzle 1A	
	Flow layer chamber capacity	3L	
	Glass parts	Made of super hard glass	
	Stirring shaft	Inside the flow layer chamber (using GB210 stirring structure)	
	Dust removal of nozzle	Pulse jet type (use the pressurizing air blowout mechanism for GB210 model)	
	Weight	Approx. 13kg	
	Parts list	Flow layer chamber	1
		Central sleeve (O-ring P16 and P135 included)	1
		Butterfly nut M6	3
		Flat washer M6	6
		Hexagon bolt M6×30	3
		Spring washer M6	3
		Cap silicone	1
		Hose $\phi 38.1 \times \phi 41.7 \times 800$ L silicone	1
		Hose clamp	2
		O-ring P145	1
		pipe	1
		filter (O-ring P30 included)	1
		Filter chamber	1
		Flange	1
		Packing C t5mm silicone	1
		Spray nozzle	1
		Filter fixing nut	1
		Polyvinyl carbazole 5% solution 50g container	1
		Sintered alumina 300g container	1
		Wrench	1
		Warranty card	1


10. Specification

Mini Spray accessories [GF301C]	Model	GF301C	
	Amount of water evaporation	Max. Approx. 1500ml/h	
	Spray nozzle	Two-fluid nozzle 1A	
	Drying chamber	Made of super hard glass	
	Cyclone	Made of super hard glass	
	Product collecting container	Made of super hard glass	
	Dust removal of nozzle tip	Pulse jet type (use the pressurizing air blowout mechanism for GB211C model)	
	Weight	Approx. 11kg	
	Parts list	Cyclone	1 set
		Drying chamber	1 set
		Product collecting container	1
		Container holding band	1
		Packing 40A、50A	1 each
		Power clamp 40A、50A	1 each
		Cap	1
		Connecting ferrule (D)	1
		PFA corrugated pipe 1-1/2 3 feet (for connecting the cyclone)	1
		Hose clip	2
		Distributor (O-rings P16, P135 included)	1
		Hex bolt M6×20	3
		Flat washer M6	3
		Spring washer M6	3
		Aluminum honeycomb	1
		Pipe	1
		Spray nozzle	1
		Wrench	1
		Plastic container for 100g of 5% sodium chloride solution	1
		Warranty card	1

Please remind that this product may be changed the specification and others for revision without any announce to the user.

11. Wiring diagram


GB211C wiring diagram



Symbol	Part name	Symbol	Part name	Symbol	Part name
QR	Earth leakage breaker (30mA)	M	Liquid sending pump drive motor	AD	Analog input module
XD	AC power distribution wiring terminal block	PLC	Programmable controller	DA	Analog output module
SF	Panel power switch	HMI	Touch screen	(MAD)	Analog input and output module
2GQ	Heat dissipation fan	LAN	Network switch	X	Outlet temperature sensor quick plug
FA	Fuse	2X,3X	Network interface socket for online use	3XD	Temperature output terminal block
2XD	Service socket	KA	Intermediate relay	3BT	Granulation temperature sensor
KF	Heating control main relay	KH	Pulse jet nozzle solenoid valve	2/3KA	Lifter up/down relay
BE	Heater current inductive switch	2KH	Pneumatic cleaning nozzle	4KA	Stirring start relay
KFX	Solid state relay	3KH	GAS connection solenoid valve	5KA	Spray/granulation switching relay
EB,2EB	Heating pipe	4KH	Spray solenoid valve	2M	Lifter motor (including capacitance C)
2KF	Blower control main relay	BG	Nozzle detection switch	3M	Stirring motor (including capacitance 2C)
2BE	Blower current transmitter	BT	Inlet temperature sensor (Pt100)	2BG	Granulation function detection switch
TR	Solid state voltage regulator	2BT	Outlet temperature sensor (Pt100)	3BG	Lifter lowest limit switch
GQ	Blower	3BE	Inlet temperature transmitter	4BG	Lifter highest limit switch
TB	DC power supply (DC24V)	4BE	Outlet temperature transmitter	5~7BG	Lifter in place limit switch
QCS	Liquid sending pump driver	BB	Overheat protector		

12. System diagram


Standard mode system diagram



No.	Part name	No.	Part name
①	Heater	⑪	Three-way solenoid valve for spray air switching
②	Inlet temperature sensor	⑫	Spray flowmeter
③	Distributor	⑬	Solenoid valve for automatic nozzle
④	Drying chamber	⑭	Spray nozzle
⑤	Cap	⑮	Liquid sending pump
⑥	Outlet temperature sensor	⑯	Cooling water connection port
⑦	Cyclone	⑰	Three-way solenoid valve for pulse jet gas source switching
⑧	Product collecting container	⑱	Suction filter
⑨	Blower		Suction blind plate
⑩	Solenoid valve for pulse jet		

12. System diagram


PUSH mode system diagram



No.	Part name	No.	Part name
①	Heater	⑪	Three-way solenoid valve
②	Inlet temperature sensor	⑫	Flowmeter
③	Distributor	⑬	Solenoid valve for automatic cleanout needle
④	Drying chamber	⑭	Spray nozzle
⑤	Cap	⑮	Liquid sending pump
⑥	Outlet temperature sensor	⑯	Cooling water connection port
⑦	Cyclone	⑰	Three-way solenoid valve for pulse jet gas source switching
⑧	Product collecting container	⑱	Suction filter
⑨	Blower		
⑩	Solenoid valve		

12. System diagram

Granulation mode system diagram (GB211C+GF200)



No.	Part name	No.	Part name
①	Heater	⑪	Three-way solenoid valve for spray air switching
②	Microporous plate	⑫	Spray flowmeter
③	Flow layer chamber	⑬	Solenoid valve for automatic nozzle
④	Filter chamber	⑭	Stirring blade
⑤	Nozzle	⑮	Stirring motor
⑥	Outlet temperature sensor	⑯	Inlet temperature sensor
⑦	Liquid sending pump	⑰	Solenoid valve for pulse jet gas source switching
⑧	Blower	⑱	Suction filter
⑨	Connecting hose	⑲	Blind plate
⑩	Solenoid valve for pulse jet	⑳	Filter

13. Operation principle

Operation principle of standard mode

Refer to P.67 "Standard mode system diagram"

The sample is sent from the appropriate container to ⑭ spray nozzle with ⑮ liquid sending pump. Moreover, the compressed air from the air compressor is regulated by ⑫ flow meter, and sent to ⑭ spray nozzle. At the tip of the nozzle, the compressed air is mixed with the sample, and the mixed sample is sprayed into the ④ drying chamber. This sample becomes drop shape that the particle diameter is approx. $20\ \mu$ and the surface area is $3,000\ \text{cm}^2$ per 1 liter of sample.

On the other side, the air is sucked into the unit by ⑨ blower, and heated up by ① heater to the set temperature. Since the contact area of the heated air and the sample is very large, the approx. 90% or more of the moisture will be evaporated in the drying chamber momentarily.

The dried sample that became fine particles is sent to ⑦ cyclone after further drying, and separated from the water vapor here, and then sent to ⑧ product collecting container. The time after the sample is sprayed with the nozzle till it is collected into this container does not take 0.5 seconds. Moreover, since the sample particles are always surrounded by the solvent vapor (water vapor), the temperature does not rise extremely around the fine particles due to the vaporization heat. However, the heat-sensitive material, such as enzyme, may be damaged even if the outlet temperature is about 80°C .

The evaporated moisture is evacuated to outside via the blower.

The temperature conditions in the experiment are displayed on the display panel through the inlet temperature sensor and the outlet temperature sensor. Moreover, the airflow rate that dries the sample is regulated by the blower.

In case that the sample attachment on the nozzle tip is too much, turn on ⑩ solenoid valve to make the pressurized air blow to the nozzle tip from ③ distributor in order to remove the attachment. If necessary, remove ⑤ cap to lead the external gas into the inside of the chamber.

Operation principle of PUSH mode

Refer to P. P.68 "PUSH mode system diagram"

The sample is sent from the appropriate container to ⑭ spray nozzle with ⑮ liquid sending pump. Moreover, the compressed air from the air compressor is regulated by ⑫ flow meter, and sent to ⑭ spray nozzle. At the tip of the nozzle, the compressed air is mixed with the sample, and the mixed sample is sprayed into the ④ drying chamber. This sample becomes drop shape that the particle diameter is approx. $20\ \mu$ and the surface area is $3,000\ \text{cm}^2$ per 1 liter of sample.

On the other side, the air is pushed into the unit by ⑨ blower, and heated up by ① heater to the set temperature. Since the contact area of the heated air and the sample is very large, the approx. 90% or more of the moisture will be evaporated in the drying chamber momentarily.

The dried sample that became fine particles is sent to ⑦ cyclone after further drying, and separated from the water vapor here, and then sent to ⑧ product collecting container. The water vapor is discharged through the exhaust pipe. In PUSH mode, air is pushed into the unit by the blower. When the high-temperature solvent vapor (water vapor), which dries the sample particles, is discharged, even if there is a high-temperature air flow, it will not influence the blower, and it can be discharged to the outside through the exhaust pipe.

PUSH mode can be used for high outlet temperature test conditions, the blower will not be influenced by high temperature and thus reduce its service life.

The temperature conditions in the experiment are displayed on the display panel through the inlet temperature sensor and the outlet temperature sensor. Moreover, the airflow rate that dries the sample is regulated by the blower.

In case that the sample attachment on the nozzle tip is too much, turn on ⑩ solenoid valve to make the pressurized air blow to the nozzle tip from ③ distributor in order to remove the attachment. If necessary, remove ⑤ cap to lead the external gas into the inside of the chamber.

13. Operation principle

Operation principle of granulation mode

Refer to P. 69 "Granulation mode system diagram (GB211C+GF200)".

By switching the installation position of ⑱ suction filter, the system will automatically switch to granulation mode when the micro switch is triggered by the filter case. Please refer to P.27~P.29 to correctly install the granulation components and connecting pipes.

The ambient air enters ① heater through the suction filter, and preheats ③ flow layer chamber by using the heated air flow. After connecting the ⑨ connecting hose with ⑧ blower, the exhaust gas can be sucked out.

After preheating, be sure to adjust the output air volume of the blower to 5%, then stop the blower for a short time, use the lifter switch to lower the flow layer chamber and take out the flow layer chamber, put the sample evenly on ② microporous plate, and then use the lifter switch to raise the flow layer chamber, adjust the packing of the flow layer chamber, start the blower and heating again. Adjust to the appropriate air volume so that the sample flowing does not exceed the position of the silicone rubber cap. When too much sample is attached on the surface of ⑳ filter, it may cause the air volume being reduced, it's able to turn on ⑩ solenoid valve for pulse jet to clean the filter surface.

After the ⑥ outlet temperature is stable, adjust the ⑫ flowmeter to an appropriate spray air source flow, then start ⑦ liquid sending pump to input the binder, at the same time you can choose to turn on ⑮ stirring motor according to the actual situation. After the binder sprayed from ⑤ nozzle is fully mixed with the sample, the fluidity of the sample will slow down, and the particle size of the sample will gradually become larger (It is suggested to obtain the desired particle size by intermittent input of a small amount of binder).

When the sample reaches the desired particle size and is completely dried, the heater can be turned off. When the inlet temperature is lower than 60℃ and the outlet temperature is lower than 50℃, the blower can be turned off for a long time.

Use the lifter switch to lower and take out the flow layer chamber, and collect the granulated sample. The experiment can be finished only after the glass container, filter, connecting pipe, nozzle, liquid sending hose, etc. are thoroughly cleaned.

14. Replace parts list

Replacement parts for GB211C

	Part name	Specifications	Manufacturer	Code No.
※	O-ring	NPF22.4	Yamato Scientific	A081902016
※	O-ring	NPF38	Yamato Scientific	A081902055
※	O-ring	P110 Viton	Yamato Scientific	B081902002
※	O-ring	P145 Viton	Yamato Scientific	B081902003
※	Liquid sending hose	Φ2*Φ4 (silicone rubber)	Yamato Scientific	B080807050
※	Liquid sending hose	Φ2*Φ4 (PVC)	Yamato Scientific	B080807051
※	Exhaust pipe	φ50×2m	Yamato Scientific	B080807002
	Heat-resistant hose	MD25 Φ50*380L	Yamato Scientific	B080807041
※	Filter screen	GB210-40360	Yamato Scientific	B040300010
※	Suction filter	DL410-40540	Yamato Scientific	B080199005
※	Filter screen	PS-600 T20X140X140	Yamato Scientific	B040300035
※	Blower	119625-00 220V 1100W	AEMETEK	A080103029
※	Peristaltic pump head	YZ10 φ2*φ4	KONAP	A080400090
※	Flowmeter	LZB-10WBF 3~30L/min	Yamato Scientific	A040409027
	Solenoid valve	AG43-02-4-E-AC220V	CKD	A040403102
	Direct acting 3-way valve	3PA210-06-P-3	CKD	A040499060
	Heat-resistant pipe	I.D.φ48 L=1000mm	Yamato Scientific	A080807103
※	Clamp	40KF Center ring with O-ring	Yamato Scientific	A041500010
	Stainless steel vacuum corrugated pipe	40KF L=500	Yamato Scientific	A041500008
	Nozzle bushing	GB210E_01_04-16	Yamato Scientific	H090201031
※	Heater	ADL311SE_01_03_02	Yamato Scientific	H090101045
	Temperature sensor (inlet)	ADL311SC_03_01-03	Yamato Scientific	H090101056
	Temperature sensor (outlet)	ADL311SC_03_01-04	Yamato Scientific	H090101057
	Temperature sensor (granulation)	GB210C_03_01-02	Yamato Scientific	H090201034

14. Replace parts list

Part name	Specifications	Manufacturer	Code No.
Touch screen	NB7W-TW11B	Yamato Scientific	A020400014
PLC	CP2E-N40DT-A	Yamato Scientific	A020300078
Earth leakage breaker	BV-DN 1P+N 25A 30mA	Yamato Scientific	A010410002
SSR	KS15/D-38Z40-L with protective cover RPC-1	Yamato Scientific	A011006024
Power relay (KF, 2KF)	HF116F-3/024DF2HTFW	Yamato Scientific	A011001013
Relay (KA ~ 4KA)	G2R-2-SN DC24V	Yamato Scientific	A011001015
Relay (5KA)	MY4NJ DC24V/MY4N-GS DC24 BY OMZ/C	Yamato Scientific	A011001006
Switching power	LRS-100-24	Yamato Scientific	A010801045
Solid state voltage regulator	Single-phase SZTV-1 60A	Yamato Scientific	A010799007
Liquid sending hose	Φ2*Φ4	Yamato Scientific	B080807050
Fuse	250V 2A	Yamato Scientific	A010301005
Micro switch	SS-01GL2	Omron	B011505003
Solenoid valve	VX3334Q-02-1G1-B	SMC	B040403001
Lifter motor	4RK25GN-CW2L2 for lifter	ORIENTAL MOTOR	A011603026
Reducer	4LB10N-1 for lifter	ORIENTAL MOTOR	B011601007
Micro switch	V-104-1A4	Omron	B011505006
Slide rail	C115-1148A CA=1876	Japan Accuride	B080400005
Motor	2IK6A-CW2L2 for stirring	ORIENTAL MOTOR	A011603025
Overheat protector	350℃	Yamato Scientific	B020103001

Note: Parts marked with ※ are consumable parts.

14. Replace parts list

Replacement parts for GF301C

	Part name	Specifications	Manufacturer	Code No.
	Drying chamber system	GF300-30000 super hard glass	Yamato Scientific	B080699008
	Cyclone set	GF300-30060 super hard glass	Yamato Scientific	B080699005
	Container holding band	GF300-40000 Stainless steel	Yamato Scientific	LT00027540
	Nozzle set	S00360-00-316L-ASSM	Yamato Scientific	Q110901009
※	O-ring	P16 4 types D Viton	Yamato Scientific	4210026021
※	Aluminum honeycomb	91.5mm* 45mm*27mm	Yamato Scientific	A080199062
※	O-ring	P135 4 types D Viton	Yamato Scientific	F0020073
※	Cap	GF300-40100 Silicone	Yamato Scientific	B082002008
	Connecting ferrule (D)	GF300-40080		LT00027543
※	PFA corrugated pipe	1·1/2 3 feet (915mm)	IIDA GOMU	LT00027545
	Hose clip	JCS-Win-2A φ35~50	Okada Industry	LT00027550
※	Packing	40A Silicone	OSAME INDUSTRIES	B081903072
※	Packing	50A Silicone	OSAME INDUSTRIES	B081903070
	Clamp	40A	OSAME INDUSTRIES	B080913014
	Clamp	50A	OSAME INDUSTRIES	B080913013
	Product collecting container	GF300-30090	Yamato Scientific	B080604007

Note: Parts marked with ※ are consumable parts.

14. Replace parts list

Replacement parts for GF200

	Part name	Specifications	Manufacturer	Code No.
	Spray nozzle set	GF200-30000	Yamato Scientific	LT00028787
※	O-ring	P135 4 types D fluorine rubber	Yamato Scientific	4210026044
※	O-ring	P16 4 types D fluorine rubber	Yamato Scientific	4210026021
※	O-ring	P30 4 types D fluorine rubber	Yamato Scientific	4210026026
※	O-ring	P145 4 types D fluorine rubber	Yamato Scientific	4210026045
※	Packing C	GF200_40310 Silicone rubber	Yamato Scientific	LT00028762
	Filter chamber	GF200_30120 super hard glass	Yamato Scientific	LT00028126
	Rolling bearing	SST-1260	Yamato Scientific	4180126002
	Packing (A)	GF300_40160 Silicone rubber	Yamato Scientific	LT00024524
	Packing (B)	GF300_40150 Silicone rubber	Yamato Scientific	LT00024523
	Microporous plate	GF200_40160 SUS304	Yamato Scientific	LT00028128
	Cap	GF300_40100 Silicone rubber	Yamato Scientific	LT00027544
	Flow layer chamber	GF200_30070 super hard glass	Yamato Scientific	LT00028129
	Cleanout needle	GF200-40000 14293-6-1/16-SS	Yamato Scientific	LT00028432
※	Packing	CP-4042-2-TEF	Yamato Scientific	3280016002
※	Washer (A)	CP104369-TEF	Yamato Scientific	3280016003
※	Washer (B)	CP3612-TEF	Yamato Scientific	3280016006
※	O-ring	JASO-1017 hycar	Yamato Scientific	4210076002
	Filter	GF200-30100 polyester	Yamato Scientific	LT00028130
	Filter packing	GF200_40260 Silicone rubber	Yamato Scientific	LT00028458
	Hose	SRDH (GS type) L800mm	Yamato Scientific	3040080004
	Hose clip	JCS-Win-2A φ35~50	Yamato Scientific	LT00027550

Note: Parts marked with ※ are consumable parts.

15. List of Dangerous Substances

Never use explosive substances, flammable substances and substances that include explosive or flammable ingredients in this unit.



GB211C+GF301C (GB211C-A) supports organic solvents by connecting it to the optional GAS series product. Carefully read the operation manual of GAS series product and take special care for handling of organic solvents.

Explosive substance	Explosive substance	①Nitroglycol, glycerine trinitrate, cellulose nitrate and other explosive nitrate esters
		②Trinitrobenzen, trinitrotoluene, picric acid and other explosive nitro compounds
		③Acetyl hydroperoxide, methyl ethyl ketone peroxide, benzoyl peroxide and other organic peroxides
	Flammable substances	Explosive substances Metal "lithium", metal "potassium", metal "natrium", yellow phosphorus, phosphorus sulfide, red phosphorus, celluloids, calcium carbide (a.k.a, carbide), lime phosphide, magnesium powder, aluminum powder, metal powder other than magnesium and aluminum powder, sodium dithionous acid (a.k.a., hydrosulphite)
		Oxidizing substances ①Potassium chlorate, sodium chlorate, ammonium chlorate, and other chlorates ② Potassium perchlorate, sodium perchlorate, ammonium perchlorate, and other perchlorates ③ Potassium peroxide, sodium peroxide, barium peroxide, and other inorganic peroxides ④Potassium nitrate, sodium nitrate, ammonium nitrate, and other nitrates ⑤Sodium chlorite and other chlorites ⑥Calcium hypochlorite and other hypochlorites
		Flammable substances ①Ethyl ether, gasoline, acetaldehyde, propylene chloride, carbon disulfide, and other substances with ignition point at a degree 30 or more degrees below zero. ② n-hexane, ethylene oxide, acetone, benzene, methyl ethyl ketone and other substances with ignition point between 30 degrees below zero and less than zero. ③Methanol, ethanol, xylene, pentyl acetate, (a.k.a.amyl acetate) and other substances with ignition point between zero and less than 30 degrees. ④Kerosene, light oil, terebinth oil, isopenthyll alcohol(a.k.a. isoamyl alcohol), acetic acid and other substances with ignition point between 30 degrees and less than 65 degrees.
	Combustible gas	Hydrogen, acetylene, ethylene, methane, ethane, propane, butane and other gases combustible at 15°C at one air pressure.

16. Standard installation manual

※ Follow the items below to make installation. (Check the procedures separately for optional parts or products of special specifications.)

Model	Serial number	Date	Installation manager (company name)	Installation manager	Judgment

No	Item	Implementation method	Table of contents No. Section for reference in manual	Judgment
Specification				
1	Accessories	Check of quantity according to the accessory columns	10.Specifications	
2	Installation	<ul style="list-style-type: none"> Visual check of the environmental status Caution: Surrounding environment 	2. Before using this unit <ul style="list-style-type: none"> At the installation site... 	
Operation related matters				
1	Source voltage	<ul style="list-style-type: none"> Measure customer side voltage (ELB etc.) with a multimeter Measure voltage while the heater is operating (Shall meet the standards) Caution: Use a power supply that meets the standard when you are going to install it on a plug or an ELB. 	2. Before using this unit <ul style="list-style-type: none"> Be sure to connect the earth wire... Use the dedicated outlet for power supply 4. Operating procedures Preparations (1) & (2) 10. Specifications <ul style="list-style-type: none"> Power supply 	
2	Installation of the attachment	Preparations <ul style="list-style-type: none"> Connecting the exhaust duct Connection to the compressor Connection of the spray nozzle cooling mechanism (as necessary) Checking the contents of GF301C set Installation of the distributor Installation of the nozzle guide block Installation of the drying chamber Installation of the temperature sensor Installation of the cyclone, the product collecting container, the cap, and the hose <ul style="list-style-type: none"> Insert the spray nozzle from the ceiling of the main unit and then connect the liquid sending hose and the pressurized air hose 	4. Operating procedures, preparations <ul style="list-style-type: none"> (3) Connection of the exhaust duct (4)Rear of the upper frame... (5)Cooling the spray nozzle (6)Mini spray... (7)On the top of the main unit... (8)...In the center of the distributor... (9)Lifter positioning... (10)...The temperature sensor... (11)...The cyclone... Operating method <ul style="list-style-type: none"> Set referring to the left drawing in section(7) 	

16. Standard installation manual

No	Item	Implementation method	Table of contents No. Section for reference in manual	Judgment
3	Operation start (test run)	Perform the test run <ul style="list-style-type: none"> • ELB and the power switch ON • Set the setting select to INLET and set the INLET temperature to 150°C • Installation of the mini spray attachment • Set the BLOWER switch ON and to air amount 20.0% (0.4m³/min) • Set the liquid sending pump to be 40rpm • Turn the heater switch ON • Setting the liquid sending hose and distilled water • Spraying pure water Set the spray pressure to 10L/min when the outlet temperature has risen to around 80°C. Adjust liquid sending speed so that the outlet temperature will be slightly lower than about 75°C <ul style="list-style-type: none"> • Change from distilled water to the sample and shift to the powder collecting operation 	4. Operating procedures Operating method <ul style="list-style-type: none"> • (1)·(2) • (3) on the operation panel... • (4)Mini spray... • (5)Blower switch... • (6) the heater switch... • (7) the liquid sending tube... • (8)·(9) • (10) When the outlet temperature has become stable... 	
4	Operation stop	Stop operation <ul style="list-style-type: none"> • Change from the sample to distilled water and wash inside the spray nozzle Approx.5 min → PUMP switch OFF→ Adjust spray pressure to 0 <ul style="list-style-type: none"> • Turn the HEATER switch OFF • Turn the BLOWER switch OFF when the outlet temperature dropped to 50°C or less • Turn the POWER switch OFF • Collect powder • Clean the containers according to the maintenance method 	4. Operating procedures Operating procedures <ul style="list-style-type: none"> • (11)When sample has been... • (12)Turn the heater OFF... • (12)Turn the heater OFF... • (13)...The power switches... • (14)...The container holding band... • (15)...to the maintenance method 6. Maintenance procedures	
説明				
1	Description of operation	Description of operation of each part to the customer according to the manual	1. Safety precautions to 15. List of Dangerous Substances	
2	Error codes	Description of the error codes and countermeasures to the customer according to the manual	8. When a trouble occurs to 9. After-sales service and warranty	
3	Maintenance & inspection	Description of operation of each part to the customer according to the manual	6. Maintenance procedures <ul style="list-style-type: none"> • Daily inspection/care 	
4	Completion of installation Matters to note	<ul style="list-style-type: none"> • Indicate the installation date and the manager name on the nameplate of the main unit. • Fill in the warranty card with necessary matters and hand it over directly to the customer. • Description of after-sales service route 	9. After-sales service and warranty	

Responsibility

Please follow the instructions in this document when using this unit. Yamato Scientific has no responsibility for the accidents or breakdown of device if it is used with a failure to comply. Never conduct what this document forbids. Unexpected accidents or breakdown may result in.

Note

- ◆ The contents of this document may be changed in future without notice.
- ◆ Any books with missing pages or disorderly binding may be replaced.

Instruction Manual
Spray Dryer
GB211C
The first edition May. 24, 2024
Revision

YAMATO SCIENTIFIC CO., LTD.
Harumi Triton Square Y-36F, 1-8-11 Harumi,
Chuo-ku, Tokyo 104-6136, Japan
Tel : +81-3-5548-7122
Fax : +81-3-5548-0132
<https://www.yamato-scientific.com/>