



Ware Drying Oven

**Model: DG400C/410C/440C/450C
800C/810C/840C/850C**

- Second Edition -

- Thank you very much for purchasing this Yamato DG series Ware Drying Oven.
- Please read the “Operating Instructions” and “Warranty” before operating this unit to assure proper operation. After reading these documents, be sure to store them securely together with the “Warranty” at a handy place for future reference.

⚠Warning!

Before operating the unit, be sure to read carefully and fully understand important warnings in the operating instructions.

Yamato Scientific Co.,Ltd.

Table of Contents

1. Safety precautions	1
Explanation of pictograms	1
List of symbols	2
Warning·Cautions	3
2. Before operating the unit.....	5
Precautions when installing the unit	5
3. Names and functions of parts.....	11
Main body	11
Operation panel	15
Explanation of characters.....	16
4. Operating procedures.....	17
List of operation modes and functions.....	17
Operation mode·function setting keys and characters.....	19
Operating procedures (settings for overheat prevention device)	20
Operating procedures (fixed temperature operation)	21
Operating procedures (quick auto stop operation).....	22
Operating procedures (auto stop operation)	24
Operating procedures (auto start operation).....	26
Useful functions (calibration offset function)	28
Useful function (setting lock function).....	29
Useful function (power outage compensation function)	30
Air exhaust fan (only for DG450C/850C).....	31
Bactericidal lamp (only for DG450C/850C).....	32
5. Cautions on handling.....	33
6. Maintenance procedures.....	37
Daily inspection/maintenance.....	37
7. When the unit is not to be used for a long time or when disposing	40
When the unit is not to be used for a long time or when disposing	40
Notes about disposition.....	40
8. Troubleshooting	41
Safety device and error codes.....	41
When a malfunction is suspected.....	42
9. After sales service and warranty	43
When requesting a repair.....	43
10. Specifications	44
11. Wiring diagram.....	46
12. List of replacement parts.....	49
13. List of dangerous materials	53
14. Standard installation manual	54

1. Safety precautions

Explanation of pictograms

About pictograms

A variety of pictograms are indicated in this operating instruction and on products to assure safe operation. Possible results from improper operation ignoring them are classified as follows.

Be sure to fully understand the descriptions below before proceeding to the text.



Warning

Indicates a situation which may result in death or serious injury (Note 1)



Caution

Indicates a situation which may result in minor injury (Note 2) and property damage (Note 3).

(Note 1) Serious injury means a wound, an electrical shock, a bone fracture or intoxication that may leave after effects or require hospitalization or outpatient visits for a long time.

(Note 2) Minor injury means a wound or an electrical shock that does not require hospitalization or outpatient visits for a long time.

(Note 3) Property damage means damage to facilities, devices and buildings or other properties.

Meanings of pictograms



This pictogram indicates a matter that encourages the user to adhere to warning ("caution" included).

Specific description of warning is indicated near this pictogram.



This pictogram indicates prohibitions

Specific prohibition is indicated near this pictogram.



This pictogram indicates matters that the user must perform.

Specific instruction is indicated near this pictogram.

1.Safety precautions

List of symbols

Warning



General warnings



Danger!:
High voltage



Danger!:
High temperature



Danger!:
Moving part



Danger!:Hazard of
explosion

Caution



General cautions



Electrical shock!



Burning!



Caution for no
liquid heating!



Caution for water
leak!



For water only



Poisonous
material

Prohibitions



General bans



Fire ban



Do not
disassemble



Do not touch

Compulsions



General
compulsions



Connect ground
wire



Install levelly



Pull out the power
plug



Regular
inspection

1. Safety precautions

Warning-Cautions

Warning



Never operate the unit in an atmosphere containing flammable or explosive gas

Never operate the unit in an atmosphere containing flammable or explosive gas. Otherwise, an explosion or a fire may result since the unit is not explosion-proof. See section "Never operate the unit in an atmosphere containing flammable or explosive gas. Otherwise, an explosion or a fire may result since the unit is not explosion-proof." See section "

13. List of dangerous materials" on page 53 about flammable or explosive gas.



Be sure to connect the ground wire

Be sure to connect the ground wire correctly. Otherwise, electrical leak may result and cause an electrical shock or a fire.



Ban on operation when an abnormality occurs

When a smoke or an unusual odor is seen or sensed, immediately turn the power switch on the main unit off and pull out the power cord (plug) from the power supply. A fire or an electrical shock may result.



Never use electrical power cords bundled

When these are used bundled, they might overheat causing a fire.



Take care not to damage electrical power cords

Avoid tightly bend, pull with a strong force or twist to prevent electrical power cords from damaging. A fire or an electrical shock may result.



Never use an explosive or a flammable material with this unit

Never use an explosive material, a flammable material or a material containing them. An explosion or an electrical shock may result. See section "

13. List of dangerous materials" on page 53 about flammable or explosive gas.



Never try to touch a hot part

Some parts of the unit are hot during and immediately after operation. Take special care for possible burning.

**Never try to disassemble or alter the unit**

Some parts of the unit are hot during and immediately after operation. Take special care for possible burning.

**Caution****When a thunder is heard**

When a thunder is heard, turn the main power off immediately. A malfunction, fire or an electrical shock may result.

2. Before operating the unit

Precautions when installing the unit

1. Carefully select an installation site.




Take special care not to install the unit at a place described below:

- On uneven or dirty floor
- Where combustible gas or corrosive gas exists
- Where the ambient temperature is 35°C or more
- Where temperature fluctuates widely
- Where dust or humidity is excessive
- Where subject to direct sunlight
- Where vibration is severe




Install this unit at a place with spaces shown below.



2. Install the unit on a level surface



Install the unit on a level surface, or else it may cause unexpected malfunctions like vibration, noise, etc.



3. Installation



The unit might fall down or move by an earthquake or an impact resulting a product damage. We recommend to make safety measures when installing the unit at a busy place. Take appropriate safety measures to prevent the unit from tripping over.

2. Before operating the unit

Precautions when installing the unit

4. Ensure adequate ventilation




Do not operate if the vents at the side and back of product are blocked.

If blocked, it may cause the internal temperature rising and performance degradation of product, also result in accident, failure or fire.

5. Never operate the unit in a place containing liquid substance



Never operate the unit in a place containing liquid substance. Once liquid entering into the product, it may cause failure, electric shock or fire.



6. Never operate the unit in an atmosphere containing flammable or explosive gas





Never operate the unit in an atmosphere containing flammable or explosive gas.

Otherwise, an explosion or a fire may result since the unit is not explosion-proof.




See “

13. List of dangerous materials” on page 53 about flammable or explosive gas.



Flammable gas



2. Before operating the unit

Precautions when installing the unit

7. Use dedicated distribution board and outlet

Use a power distribution board or outlet that meets the electrical capacity of the unit.

Power capacity:	DG400C	AC115V	9A	DG800C	AC115V	12A
	DG440C	AC115V	9A	DG840C	AC115V	12A
	DG410C	AC220V	5A	DG810C	AC220V	6.5A
	DG450C	AC220V	5A	DG850C	AC220V	6.5A

- * When the unit does not start even when you turn the power ON, check if for low main voltage or if the unit is connected to the same power supply line as other devices and connect it to another line if necessary.
- * Avoid connecting too many devices using a branching outlet or extending a wire, the temperature controlling function may degrade due to voltage drop.



Do not connect the unit to any parts or lines such as a gas pipe, a water pipe or a telephone line. Otherwise, an accident or a malfunction may result.

8. Handling of a power cord



Never use electrical power cords bundled. When these are used bundled, they might overheat causing a fire.

Do not convert, forcibly bend, twist or pull the power cord. Otherwise, a fire or an electrical shock may result.

Do not place the power cord under a desk or a chair, or sand between objects to avoid it from being damaged. Otherwise, a fire or an electrical shock may result.

Do not place the power cord close to a stove or other heat generating device. Sheath of the cord may burn and result in a fire or an electrical shock.



If the power cord should be damaged (exposure of core wire or disconnection), immediately turn the main unit off, pull out the power cord (plug) out of the power supply and ask your dealer to replace the cord. Otherwise, a fire or an electrical shock may result.

Connect the power cord to an appropriate wall outlet.

9. Be sure to connect the ground wire




- When there is no ground terminal available, class D grounding work is necessary and please consult your dealer or our nearest sales office.

- Be sure to connect the ground wire to the wall outlet securely.




We recommend use a ground type outlet



When there is no ground terminal.

In this case, class D grounding work is necessary and please consult your dealer or our nearest sales office.

When a bipolar type outlet is used



Insert the ground adaptor included as an option, into a power plug confirming the polarity of the outlet. Connect the ground wire (green) of the ground adaptor to the ground terminal on the power supply equipment.



Never connect the ground wire to gas pipe, water pipe, or telephone line. Otherwise, an accident or a malfunction may result.


2. Before operating the unit

Precautions when installing the unit

(1) Lift the foot of caster adjustor

Lift the feet of 4 caster adjustors as shown in the right figure, confirm the 4 casters can turn freely, and then move the unit.

※ If move to an uneven surface, it would give excessive impact to the casters and cause damage of casters. Please uplift the unit to move.




(2) Choose a proper installation place

Flat place (4 casters touch land totally), put down the feet of adjustors as check if the unit shakes or tilts.

(3) Set the water tray


Set the water tray for water containing when drying.



(4) Install shelf board

- Install the shelf rest brackets at left and right shelf supports at proper height.
- Slide the shelf board inward.
- ※ Pay attention to the height of left and right shelf rest brackets.
- Confirm shelf board would not fall down and get loosening.


Withstand load of each shelf board is 15kg in uniform loading. Place the samples separately.



2. Before operating the unit

Precautions when installing the unit

- Too many samples may influence proper temperature control. Be sure to use shelf boards and place samples apart each other so as to make free space of 30% or more to assure proper temperature accuracy.




Make at least 30% of space

- (5) Never put samples at the bottom of inner chamber
- If directly put samples at the bottom of inner chamber and operate the unit, it may influence the temperature performance. It may cause not only corrosion, damage or rust of inner chamber, but also sample damage or fire.
 - The sensor is set at the inner wall of chamber, do not make samples contact the walls. Place the samples on the shelf boards.
- (6) Pay attention not to use the following samples
- ① Samples containing flammable or explosive substances
This product is not explosion proof, never use the unit to dry samples containing flammable or explosive substances.
 - ② Corrosive samples
Pay attention not to use corrosive samples. This product may be corroded by corrosive substances like strong acid, although it mostly made of stainless steel SUS304. Otherwise, the sealing strip may be corroded by organic solvent such as acid, alkali and oil.
- (7) Open the vents when operating
The vents are set on the top of product. Usually, operate with the 2 vents opened. When operating, adjust the opening as per the moisture discharge of samples.
- ! High-temperature steam would erupt from the vents, it may cause burn injury, do not observe at the vents or touch them.**
- (8) Operate when door is fully closed
- Confirm that the locking mechanism at right side of door is totally locked, and then operate the unit.


3. Names and functions of parts

Main body

DG400C/410C front




DG400C/410C Back




3. Names and functions of parts

Main body

DG440C/450C front




DG440C/450C Back




3. Names and functions of parts

Main body

DG800C/810C front




DG800C/810C back




3. Names and functions of parts

Main body

DG840C/850C front




DG840C/850C back



3. Names and functions of parts

Operation panel









No.	Name	Operation/action
①	RUN/STOP key	Used for starting/stopping operation.
②	▼▲ keys	Used for selecting settings.
③	TIMER key	Key for selecting timer operation settings. Quick auto stop operation, auto stop operation or auto start operation can be selected.
	SUB MENU key (Long press of the Timer key)	Key for setting calibration offset temperature, the key lock function or the power outage compensation function.
④	RUN lamp	Illuminates during fixed temperature operation and blinks during timer operation.
⑤	HEATER lamp	Illuminates while heater power is on.
⑥	Measured temperature screen	Displays measured temperature in the bath · set characters · alarm information.
⑦	Set temperature screen	Displays a set temperature, timer settings and timer remaining time.

3. Names and functions of parts

Explanation of characters

Characters on the controller are explained in this section.

Characters	Identifier	Name	Application
	AStP	Auto stop setting	Used for setting auto stop operation.
	AStr	Auto start setting	Used for setting auto start operation.
	End	Time up	Displayed when timer operation has ended. See pages 22 and 24.
	cAL	Calibration offset setting	Used for inputting a calibration offset temperature See section "Using the calibration offset function" on page 28.
	Lock	Key lock of settings	Key locks settings to prevent their alteration See section "Using the lock function" on page 29.
	Pon	Power outage compensation setting	Selects operations after recovery from power outage. See section "Using the power outage compensation function" on page 30.

*See the section "Operation mode-function setting keys and characters" on page 19 for characters of operation modes and functions.

4. Operating procedures

List of operation modes and functions

Operation modes of the unit are as shown below:

No	Name	Description	Page
1	Fixed temperature operation	Turning the ELB on to enter the operation setting mode. Proceed to temperature setting that uses ▼▲ keys. Pressing the RUN/STOP key longer to start operation, and pressing the RUN/STOP key longer again to stop operation.	P.21
2	Quick auto stop operation	Used when you want to “stop fixed temperature operation being performed automatically in several hours.” Press the TIMER key during fixed temperature operation to display “AStP.” Set a duration before stop with the ▼▲ keys. Pressing the RUN/STOP key starts quick auto stop operation and activates the timer in the middle of it to automatically stop it after the set period of time.	P.22
3	Auto stop operation	Used when you want to “set automatic stop for fixed temperature operation when making settings for it.” Press the TIMER key to display “AStP.” Set a duration before stop with the ▼▲ keys. Pressing the RUN/STOP key starts auto stop operation.	P.24
4	Auto start operation	Used when you want to “start operation automatically after several hours” after power is turned on. Press the TIMER key to display “AStr.” Set a duration before stop with the ▼▲ keys. Pressing the RUN/STOP key starts auto start operation.	P.26
* Operation mode cannot be changed while the unit is in operation. First stop operation before changing the mode.			

4. Operating procedures

List of operation modes and functions


Functions of the unit are as shown below:

No	Name	Description	Page
1	Overheat prevention function	<p>Automatic overheat prevention function:</p> <p>This function is linked to the unit set temperature and has been set to so that it is automatically activated (returned automatically) at a temperature 12°C higher than the set temperature in the bath.</p> <p>Standalone overheat prevention device:</p> <p>When the temperature in the bath reaches the set temperature of the overheat prevention device, its heater circuit trips to shut off controller operation.</p> <p>The temperature can be set with the manual dial on the hydraulic overheat prevention device installed at the right side of the unit.</p>	P.20
2	Calibration Offset function	<p>Calibration offset function compensates any differences between the target temperature in the bath and the control temperature of the controller (sensor temperature.)</p> <p>The function can compensate to either plus or minus side for the whole temperature band of the unit.</p> <p>This compensation can be set with the SUB MENU keys.</p>	P.28
3	Setting lock function	<p>This function locks the set operation status.</p> <p>The lock can be set or released with the SUB MENU key.</p>	P.29
4	Power outage compensation function	<p>This function returns the main unit operation to the resume status after recovery from power outage, or keeps the current stop status.</p> <p>This compensation can be set with the SUB MENU keys.</p>	P.30

4. Operating procedures

Operation mode-function setting keys and characters

Key operations and characters in the diagram below are used for operation mode and function settings.



4. Operating procedures

Operating procedures (settings for overheat prevention device)

As a safety measure for preventing overheat, a hydraulic overheat prevention device (manual return) is installed.


Temperature setting range and functions

The temperature setting range for the standalone overheat prevention device is “50°C ~ 120°C.”

When the temperature in the bath keeps rising beyond the controller set temperature and reaches the set temperature of the overheat prevention device, the heater circuit trips and the controller operation is shut off.

When the overheat prevention device is activated, it will not be released until the ELB is turned on.

How to set temperature



Set the temperature scale
to the arrow

Setting the overheat prevention temperature

- Set the temperature scale on the hydraulic overheat prevention device installed on the right side of the unit to the arrow in the diagram shown left.
- Turn the ELB to “OFF” and wait for a while without opening the door.
- After a while, turn the ELB “ON.” (Turn the ELB “ON.”)


Caution

- ① Set temperature as “set temperature +20°C” as a rough standard and add 5°C to the setting if the device functions improperly.
- ② The temperature setting range for the standalone overheat prevention device is “50°C ~ 120°C.” Be sure to set the overheat prevention activation temperature correctly otherwise the device may not start, the overheat prevention device is activated before temperature in the bath increases completely, or a fire or other unexpected accidents may result.
The temperature is set at 120°C on shipping from the factory.
- ③ If the temperature set of standalone overheat protector is at room temperature or below, open the door at this time, the standalone overheat protector may work.
- ④ The overheat prevention device has been designed to prevent overheating of devices not to protect samples. The device does not prevent accidents caused from use of explosive or flammable substances.

4. Operating procedures

Operating procedures (fixed temperature operation)

How to start fixed temperature operation




1. Turn the ELB ON. (Turn the ELB to “ON.”)

When the ELB is turned ON, the initial values will be displayed for about four seconds, then the initial screen will appear and the current bath temperature and the previous set temperature are displayed on each of the indicators.


Measured temperature screen: Displays the current bath temperature

Set temperature screen: Displays the previous set temperature



2. Setting the temperature

Set a temperature using the ▼▲ keys.



3. Starting operation

Press the **RUN/STOP** key longer.

Fixed value operation will start and the RUN lamp and the HEATER lamp come on.

4. Stopping operation

Press the **RUN/STOP** key longer.

Operation stops, the RUN lamp goes off and the screen switches to the initial setting screen.

When you want to correct setting errors or change settings

When you want to change settings, press the ▼▲ keys on the current screen to enter the setting mode where you can change settings. Blink stops three seconds after three seconds after change and setting is completed.



Caution


- ① When you want to lower the set temperature during fixed temperature operation, note that it takes some time to reach the reset temperature since the unit has no cooling capacity.
- ② Immediately after operation has been stopped, the temperature in the bath is around the set temperature. Operation stop refers only to machine stop and time needed for decreasing the emperature in the bath is not considered.

4. Operating procedures

Operating procedures (quick auto stop operation)

Used when you want to “stop fixed temperature operation being performed automatically in several hours. Quick auto stop operation is a function to enable auto stop timer setting during operation.

Procedures for quick auto stop operation



1. Setting time period before stop during fixed temperature operation

- ① Make sure that the RUN lamp is illuminated to indicate the unit is in operation.

Press the **TIMER** key.

Characters **AStP** **85tP** are indicated on the measured temperature screen to indicate the auto stop operation mode and set duration blinks on the set temperature screen.

- ② Set a duration you want using the **▼▲** keys.


About the timer function

The maximum time that can be set for the timer is 999 hours 50 minutes.

Up to 99 hours 59 minutes, time can be set in minutes.

One hundred hours and over are set only in 10 minutes.

Keep the **▼▲** keys pressed to continuously change set time and you can quickly reach the time you want. Press the **▼▲** keys once at a time for fine adjustment.




2. Starting timer operation

When the time you want is set, press the **RUN/STOP** key while the set temperature screen is blinking.

The RUN lamp blinks and timer operation is started.

Timer starts counting when the temperature in the bath reaches the set temperature.

Once timer counting is started, the set temperature screen changes to the remaining time display.



3. Stopping and ending timer operation

Operation stops automatically when the set temperature has elapsed.

Characters **End** **End** blink on the set temperature screen to indicate operation has ended.

Press the **RUN/STOP** key for approx. one second to end the timer operation mode. The screen switches to the initial setting screen.

4. Operating procedures

Operating procedures (quick auto stop operation)

When you want to correct set temperature or set time, or change settings


When you want to change settings, press the ▼▲ keys on the current screen to enter the setting mode where you can change settings. Blinking stops three seconds after three seconds after change and setting is completed. Note, however, that temperature changes after timer activation are counted also while temperature is changing.

When you want to change settings before timer activation, press the **TIMER** key on the current screen to enter the setting mode where you can change settings. Enter a time duration from when the set temperature is reached to the time the device shall be stopped.

When you want to change settings after timer activation, press the **TIMER** key on the current screen to enter the setting mode where you can change settings. Note, however, you need to set a time calculated by adding the time already passed to the time to be added.

After change has been made, press the **RUN/STOP** key to complete the process.

When you want to stop quick auto stop operation in the middle of it, press the **RUN/STOP** key long once to stop device control once, then make settings again in the appropriate mode.


In terms of the remaining time display  a blinking dot indicates count down and an illuminating dot indicates a wait status (while temperature is increasing or decreasing to the set temperature) during which the timer has stopped counting.

4. Operating procedures

Operating procedures (auto stop operation)

This mode automatically stops fixed temperature operation after a certain time from its start set with the timer.

Procedures for auto stop operation



1. Setting a stop time

- ① After confirming the temperature you want is set, Press the **TIMER** key to display characters AStP **88.88** on the measured temperature screen that indicate auto stop operation.
The set time is displayed on the set temperature screen.
- ② Set a time you want using the **▼▲** keys.
Pressing the **▼▲** keys makes the set time blink. The time is determined when blinking stops.


About the timer function

The maximum time that can be set for the timer is 999 hours 50 minutes.

Up to 99 hours 59 minutes, time can be set in minutes.

One hundred hours and over are set only in 10 minutes.

Keep the **▼▲** keys pressed to continuously change set time and you can quickly reach the time you want. Press the **▼▲** keys once at a time for fine adjustment.



2. Starting timer operation

When the time you want is set, press the **RUN/STOP** key for about one second while characters AStP **88.88** that indicate auto stop operation are displayed on the measured temperature screen and the set time on the set temperature screen.

The **RUN** lamp blinks and timer operation is started.

Timer starts counting when the temperature in the bath reaches the set temperature.


Once timer counting is started, the set temperature screen changes to the remaining time display.

3. Stopping and ending timer operation

Operation stops automatically when the set temperature has elapsed.

Characters End **End** blink on the set temperature screen to indicate operation has ended.

Press the **RUN/STOP** key for approx. one second to end the timer operation mode. The screen switches to the initial setting screen.



4. Operating procedures

Operating procedures (auto stop operation)

When you want to correct set temperature or set time, or change settings

When you want to change settings, press the ▼▲ keys on the current screen to enter the setting mode where you can change settings. Blinking stops three seconds after three seconds after change and setting is completed. Note, however, that temperature changes after timer activation are counted also while temperature is changing.


When you want to change settings before timer activation, press the **TIMER** key on the current screen to enter the setting mode where you can change settings. Enter a time duration from when the set temperature is reached to the time the device shall be stopped.

When you want to change settings after timer activation, press the **TIMER** key on the current screen to enter the setting mode where you can change settings. Note, however, you need to set a time calculated by adding the time already passed to the time to be added.

After change has been made, press the **RUN/STOP** key to complete the process.

Auto stop operation is not available together with auto start operation.

When you want to stop auto stop operation in the middle of it, press the **RUN/STOP** key long once to stop device control once, then make settings again in the appropriate mode.

In terms of the remaining time display  a blinking dot indicates count down and an illuminating dot indicates a wait status (while temperature is increasing or decreasing to the set temperature) during which the timer has stopped counting.


4. Operating procedures

Operating procedures (auto start operation)

This mode automatically starts fixed value operation after a certain time from its start set with the timer.

However, operation does not stop automatically but needs to be stopped manually.

Procedures for auto start operation



1. Setting an operation start time

- ① After confirming the temperature you want is set, Press the **TIMER** key to display characters AStr^{85Er} on the measured temperature screen that indicate auto start operation. The set time is displayed blinking on the set temperature screen.
- ② Set a time you want using the **▼▲** keys. Pressing the **▼▲** keys makes the set time blink. The time is determined when blinking stops.


About the timer function

The maximum time that can be set for the timer is 999 hours 50 minutes.

Up to 99 hours 59 minutes, time can be set in minutes.

One hundred hours and over are set only in 10 minutes.

Keep the **▼▲** keys pressed to continuously change set time and you can quickly reach the time you want. Press the **▼▲** keys once at a time for fine adjustment.




2. Starting timer operation

When the time you want is set, press the **RUN/STOP** key for about one second while characters AStr^{85Er} that indicate auto start operation are displayed on the measured temperature screen and the set time on the set temperature screen.

Timer starts counting when the **RUN/STOP** key is pressed and RUN lamp blinks.

Display on the measured temperature screen switches from set time display to remaining time display.



3. Stopping and ending timer operation

Operation automatically starts at the set time and the **RUN** lamp comes on.

To stop operation, press the **RUN/STOP** key for approx. one second to end the timer operation mode. The screen switches to the initial setting screen.

4. Operating procedures

Operating procedures (auto start operation)

When you want to correct set temperature or set time, or change settings

When you want to change the set temperature during timer counting, press the ▼▲ keys during that status to switch the set temperature screen to the set temperature input mode, which blinks to enable change of the set temperature with the ▼▲ keys.

When you want to change the set time during timer counting, press the TIMER key during that status to switch the set temperature screen to the set time input mode, which blinks to enable change of the set time with the ▼▲ keys.

In either case, the set temperature screen will stop blinking after a while and switch to the timer count mode and the change made is determined. Note, however, when you change the set time you need to set a time calculated by adding the time already passed to the time to be added.

When operation has started after the auto start time, you cannot change the set time.

When you want to stop auto start operation in the middle of it, press the RUN/STOP key long to stop device control once, then make settings again in the appropriate mode.

In terms of the remaining time display 1.30 a blinking dot indicates count down and an illuminating dot indicates a wait status (while temperature is increasing or decreasing to the set temperature) during which the timer has stopped counting.

4. Operating procedures



Useful functions (calibration offset function)

Using the calibration offset function

Calibration offset function compensates any differences between the target temperature in the bath and the control temperature of the controller (sensor temperature.) The function can compensate in parallel to either plus or minus side for the whole temperature band of the unit.

The lock can be set or released with the **[SUB MENU]** keys.

The temperature is set at “0” on shipping from the factory.



- ① Start operation at the target set temperature and confirm the temperature in the bath with a temperature recorder after temperature has stabilized.
- ② Confirm the difference between the set temperature and that in the bath.
- ③ Press the **[TIMER]** key (**[SUB MENU]** key) long to enter the sub menu mode.
Press the **[TIMER]** key (**[SUB MENU]** key) several times to select the characters **cAL** **[cAL]** that indicates the calibration offset function.
- ④ Enter the difference between the set temperature and the temperature in the bath using the **▼▲** keys and press the **[TIMER]** key (**[SUB MENU]** key) long to exit the sub menu mode. (When you want to set the key lock function, proceed to character selection process for the key lock function without pressing the **[TIMER]** key (**[SUB MENU]** key) long.)

- * You can set either of + or – side for the offset compensation temperature.
When compensation is set for the – side, the measured temperature display decreases by the compensation temperature while the temperature in the bath increases by the same amount.
When compensation is set for the + side, the measured temperature display increases by the compensation temperature while the temperature in the bath decreases by the same amount.
- * Since too large a compensation value may result in larger difference between the actual and indicated temperatures and may present a danger, consult our nearest sales office before entering a large compensation value.
- * The device has, in addition to the calibration offset function, the two-point compensation function that adjusts offset for the lower temperature range and higher temperature range, for which adjustment temperatures have been input on shipping from the factory.
- * Consult the nearest sales office before attempting validation work for the temperature adjusting device.


4. Operating procedures

Useful function (setting lock function)

Using the lock function


This function locks the set operation status.

The temperature is set at “off” on shipping from the factory.




- ① Press the **TIMER** key (**SUB MENU** key) long to enter the sub menu mode.

Press the **TIMER** key (**SUB MENU** key) several times to select the characters Lock **LocH** that indicate the setting lock function.



- ③ “Off” is displayed on the set temperature screen. To lock settings, change to “on” using the ▲ key.

Press the **TIMER** key (**SUB MENU** key) long to exit the sub menu mode.



- (3) To release lock, press the **TIMER** key (**SUB MENU** key) long again and select the characters Lock **LocH** that indicate setting lock using the ▼▲ keys.

Lock is released when “off” is selected using the ▼ key.

* When the lock function is “on”, keys other than the **RUN/STOP** key and the **TIMER** key (**SUB MENU** key) are locked.


4. Operating procedures

Useful function (power outage compensation function)

Using the power outage compensation function


The power outage compensation function returns the main unit operation to the resume status after recovery from power outage, or keeps the current stop status.

The function is set at “on” on shipping from the factory.



- ① Press the **TIMER** key (**SUB MENU** key) long to enter the sub menu mode.

Press the **TIMER** key (**SUB MENU** key) several times to select the characters Pon **Pon** that indicate the power outage compensation function.



- ② “On” is displayed on the set temperature screen. The device keeps stop status after recovery from power outage when this setting is set to “off” using the ▼ key.


Press the **TIMER** key (**SUB MENU** key) long to exit the sub menu mode.

4. Operating procedures

Air exhaust fan (only for DG450C/850C)

- ① Use the regulating damper of air inflow at lower right to manually regulate. Loosen the fixed screws of regulating damper, it's total open if slide inward and total close if slide outward.

It's total close at factory setting.




- ② When use the air exhaust fan, please turn its switch ON. The indicator lamp of air exhaust fan lights up, and the exhaust fan activates. When not use, please turn off the regulating damper of air inflow, switch OFF.
(※ the air exhaust fan is controlled through ON/OFF of standalone switch)
- ③ When not use the air exhaust fan, please turn its switch OFF. The indicator lamp of air exhaust fan goes out, and the exhaust fan stops.

Warning

- ※ The exhaust vent of unit would emit high-temperature vapor, keep your head and hands away.
- ※ There is too much moisture contained in the gas exhausted by the unit, do not put goods easily affected with damp near the unit.

4. Operating procedures

Bactericidal lamp (only for DG450C/850C)



- ① When use the bactericidal lamp, close the door and then turn the switch of bactericidal lamp ON.
The ultraviolet light (253.7nm nearby) irradiates in the chamber.
For safety, if door opens, even the switch is ON, the bactericidal lamp would goes out.
- ② When not use the bactericidal lamp, please turn its switch OFF. The indicator lamp of sterilization goes out, and bactericidal lamp goes out.

Caution

- ※ The bactericidal lamp is controlled through ON/OFF of standalone switch.
- ※ Do not use materials of easily fading, degrading and damage, like plastic container, cloth, wallpaper and so on.
If irradiated by ultraviolet light, these materials would fade, degrade and damage.
- ※ Do not use plant or use unit near plant. The ultraviolet light would make leaves wither and influence the growth of plant.
- ※ Carefully clean the indicator lamp and reflect sheet. The sterilizing effect would sharply decline if they're dirty.
- ※ Please replace the indicator lamp once it reaches the service life. If the indicator lamp lights up beyond its service life, the output of ultraviolet light would weaken and the initial sterilizing effect cannot be achieved. Its service life is 4000hrs (about 6 months) with continuous lighting.
- ※ If use bactericidal lamp to work by necessity, to prevent exposure of skin and eyes, please wear protective mask and gloves for ultraviolet light.

Warning

- ※ When the indicator lamp lighting up, never directly or indirectly use your eyes to look, and never have the reflected light enter your eyes, which may cause eye ache or influence the eyesight.
- ※ Do not have the ultraviolet light radiation (sterilizing ray) directly or indirectly contact your skin, which may cause skin inflammation.
- ※ When install, disassemble or clean the unit, please cut off the power, or else it may cause eye ache, eyesight decline, skin inflammation, etc.

5. Cautions on handling



Warning

1. About handling of flammable or combustible solution



The unit is not explosion proof. Take special care for handling samples on which explosive substances, combustible substances or substances containing them. Flammable or combustible solution will evaporate when left at a room temperature (or at a lower temperature for some types of solutions) and may be ignited and explode from switches, lights and other ignitable sources. Be sure to assure sufficient ventilation when using these materials.

See section “

13. List of dangerous materials” on page 53.

2. Ban on use/countermeasures when an error occurs



If smoke is emerges on the unit or an odd odor is felt, immediately turn the ELB on the main unit off, turn the power supply off and contact your dealer or a Yamato sales office for inspection. Otherwise, a fire or an electrical shock may result. The user shall never attempt to repair the unit to avoid any possible dangers.

3. Ensure adequate ventilation



Do not operate if the vents at the side and back of product are blocked.

If blocked, it may cause the internal temperature rising and performance degradation of product, also result in accident, failure or fire.

4. Do not be splashed by liquid



Do not make liquid enter into the product, especially from the vents at the side and back of product. If liquid enters into product, stop operation. It may cause accident, failure or fire.

5. Do not drop metal object in the chamber



Do not drop metal object in the chamber such as clip, staple and screw.

If metal object drops in the chamber, stop operation. It may cause accident, failure or fire.

6. Do not open...



Do not open the fixed control panel, cover, etc., and do not operate if opened.

It may cause accident, failure or fire.

7. Do not operate without water tray



Do not operate without water tray.

It may cause accident or failure.

8. Do not modify the product



Please do not modify the product as per customer's requirement, it may cause accident, failure or fire.

5. Cautions on handling



Caution

1. Do not step on the unit



Do not step on the unit. Otherwise, the unit may topple over or be damaged resulting a personal injury or a malfunction.

2. Do not put or drop an object on the unit



Do not put or drop an object on the unit. Since the unit contains high precision devices, vibrations or shock may cause a malfunction.

3. When a thunder is heard



When a thunder is heard, turn the ELB on the main unit off then turn the main power off immediately. Otherwise, a lightning strike may result and cause a fire.

4. During night and not to be operated for a long period of time



During the night and when you want to stop the unit for a longer period of time, turn the ELB to "off" and pull out the power cord from the power supply.

5. About recovery from power outage



When the power is applied again after the unit has stopped due to power outage, the unit will automatically return to the status immediately before the power outage and resumes operation. Turn the ELB off if you do not want to resume operation by automatic recovery.

6. Always operate the unit at a correct ambient temperature



Operational temperature range is $RT+5^{\circ}C \sim 70^{\circ}C$.

Never try to operate the unit outside the operating temperature range.

7. When open/close door



- When open/close door, pay attention to the open/close range (space), do not have hand and face approach it, it may hurt hand and face.
- After operation, do not open the door to quickly cool down the samples. Because the heat from chamber may cause deformation of control panel, controller failure and so on.

8. Do not operate the unit with the door open



When the unit is operated with the door open, it cannot control temperature normally, and the heater may overheat result in a possible danger. Be sure to operate the unit with the door closed.

5. Cautions on handling



Caution

9. Set of shelf board and sample



Please correctly set shelf board and sample as per Page 7 Install shelf board. If not, it may cause the product cannot reach its performance and result in accident and failure.

10. Do not operate items not recorded at Instruction Manual



Do not operate items not recorded at Instruction Manual, it may result in unpredictable accident.

6. Maintenance procedures

Daily inspection/maintenance

Be sure to perform daily inspection and maintenance to assure reliable operation of the unit.

Warning

- Be sure to pull out the power cord unless necessary before trying to do inspection and maintenance works.
- Start these works after the device has returned to the normal temperature.
- Never try to disassemble the unit.

Caution

- Wipe off any dirt with a tightly wrung soft cloth. Never try to clean the unit with benzene, thinner or scouring powder, or rub with a scrubbing brush. Deformation, degradation or discoloration may result.


Every month

Please inspect the ELB performance.

Power on to test

Firstly, turn ELB on

Then, use a tip (like ball pen tip) to press the test switch of ELB. If ELB breaks off, it's normal.



Maintenance of inner chamber

Please stop operation, turn off the ELB, and pull out the power cord from distribution board or outlet.

Confirm the temperature in the chamber, and take out the shelf boards and shelf supports.

The inner chamber, shelf boards and shelf rest brackets are made of stainless steel SUS304, and the inner door uses toughened glass. Please use alcohol and cleaning cloth to wipe, and then use dry cloth to wipe.

Pay attention not to use acid or alkali lotion, oil, organic solvent, etc., it may cause corrosion and damage.

There are some sharp protrusions in the chamber, at the shelf boards and shelf supports. Do not get hurt. Especially for novice working, please wear gloves.

6. Maintenance procedures

Daily inspection/maintenance

Drainage of water tray

According to ware moisture and current weather condition, the water collections of water tray are different.

The water accumulation before operation and after operation, take out the tray during operation to drain water.

After draining, put the water tray back.



Maintenance of air suction filter (only for DG450C DG850C)

- Please turn the ELB 「OFF」 .
- Hold the handle of suction filter at the lower right and pull it out. (its temperature might be very high after operation, pay attention not to be scalded)
- Use dust collector to clean the filter screen and revert its installation. (operate the unit with the suction filter being installed)



6. Maintenance procedures

Daily inspection/maintenance

Replacement of germicidal lamp (DG450C DG850C)

※ It's suggested that the germicidal lamp is replaced by YAMATO serviceman

- Please turn the ELB 「OFF」 .

Make sure the power of germicidal lamp is 「OFF」 .

- Use M4 cross screwdriver to screw out the screws at left and right sides of germicidal lamp, and disassemble the protective cover of germicidal lamp (do not lose the screws).



- Hold the lamp and rotate it by 90° .



- Take out the lamp along the gaps of both ends of lamp holder.
- After replacement, revert the installation as per reverse order.

There are sharp bulges in the chamber, at the shelf plate or shelf rest, pay attention not to be injured. It's dangerous to work bare-handed, please wear gloves.



7. When the unit is not to be used for a long time or when disposing

When the unit is not to be used for a long time or when disposing



Caution

When the unit is not going to be used for a long time

- Turn the ELB to off and pull out the power cord.



Warning

When disposing the unit

- Do not leave the unit in the area where children may have access.
- Be sure to remove handles before disposing the unit to prevent the doors from locking.
- In general, dispose the unit as a bulky waste.

Notes about disposition

Always pay attention to the preservation of the global environment.

- We highly recommend taking the unit apart as far as possible for separation or recycling to contribute to the preservation of the global environment. Major components and materials for the unit are as follows:

Names of major components	Major materials
Major mechanism part components	
Enclosure	Steel plate SPCC (powder coating)
Internal bath	Stainless steel
Heat insulator	Rock wool
Door packing	Silicon rubber foam
Nameplates	Polyethylene (PET) resin film
Major electric parts	
Heater	Iron-chrome heater
Boards	Glass fiber and other composite parts
Power cord, wire material and others	Synthetic rubber sheathed and resin sheathed wires

8. Troubleshooting



Safety device and error codes

The unit has the self diagnostic function with a controller and a separate safety device.

Table below shows possible causes and measures when the safety device is triggered.

[Error codes]

When a functional or mechanical abnormality occurs, an error code will be displayed on the control panel. When an abnormality occurs, confirm the error code and immediately stop operation.

Safety device	Symptom	Possible causes and measures
Sensor error	 appears	<ul style="list-style-type: none"> ● Error in the temperature input circuit ● Disconnection or other errors in the temperature sensor. ● Measured temperature is outside the displayable range Contact our service department.
Memory error	 appears	<ul style="list-style-type: none"> ● Memory setting error Contact our service department.
Measured temperature error	— — — — — — — — appears	<ul style="list-style-type: none"> ● When the upper limit alarm of the temperature alarm function is triggered. Contact our service department.

8. Troubleshooting

When a malfunction is suspected

If any of the symptoms below occurs

Symptom	Check
Turning the ELB to on will not activate the unit.	<ul style="list-style-type: none">● If the power cord is connected to the power supply securely.● If power outage is not occurring.● If the standalone overheat prevention device is working.
Temperature does not rise.	<ul style="list-style-type: none">● If the set temperature is below that in the device.● If the power supply voltage has declined.● If the ambient temperature is not low.● If cooling load for inside the bath is not too large.
Temperature fluctuates during operation.	<ul style="list-style-type: none">● If the set temperature is appropriate.● If the power supply voltage has declined.● If ambient temperature fluctuates widely.● If cooling load for inside the bath is not too large.
Displayed temperature differs from the measurement.	<ul style="list-style-type: none">● If the calibration offset setting is not other than "0". Set it to "0." Confirm settings in "Useful functions (calibration offset function)" in page 28.

If power outage occurs

When the power is applied again after the unit has stopped due to power outage, the unit will automatically return to the status immediately before the power outage and resumes operation. Turn the ELB off if you do not want to resume operation by automatic recovery.

- ◆ If the symptom does not match any of the above, immediately turn the ELB on the main unit off, pull out the power cord from the power supply and contact your dealer or one of our sales offices.

9. After sales service and warranty

When requesting a repair

When requesting a repair

If any trouble occurs, immediately stop operation, turn the ELB off, pull out the power plug and contact your dealer or our sales office.

Information necessary for requesting a repair

- | | | |
|---|---|---|
| ◆ Model name of the product | } | Confirm on the warranty card or the nameplate installed on the unit.
See section "Main body" on page 11. |
| ◆ Serial number | | |
| ◆ Date (y/m/d) of purchase | | |
| ◆ Description of trouble (as in detail as possible) | | |

Be sure to indicate the warranty card to our service representative.

Warranty card (attached separately)

- Warranty card is given by your dealer or one of our sales offices and please fill in your dealer, date of purchase and other information and store securely.
- Warranty period is one full year from the date of purchase. Repair service for free is available according to the conditions written on the warranty card.
- For repairs after the warranty period consult your dealer or one of our sales offices. Paid repair service is available on your request when the product's functionality can be maintained by repair.

Minimum holding period of repair parts

The minimum holding period of repair parts for this product is seven years after end of production.

Repair parts here refer to parts necessary for maintaining performance of the product.

10. Specifications

Name	Ware Drying Oven			
Model	DG400C	DG440C	DG800C	DG840C
System	Natural convection	Natural convection Forced circulation air exhaust	Natural convection	Natural convection Forced circulation air exhaust
Temperature control range	RT+5°C~70°C			
Structure	Interior material	Stainless steel SUS304		
	Observation window	250×300mm toughened glass×1	250×700mm toughened glass×1	
	Heater	Iron-chrome heater 1.0KW		Iron-chrome heater 1.34KW
	Air exhaust fan	—	Moisture-proof axial flow fan	— Moisture-proof axial flow fan
	Bactericidal lamp	—	Moisture-proof fluorescent tube 8W	— Moisture-proof fluorescent tube 15W
Control	Control system	PID control of heater output with a micro-computer		
	Setting system	Digital setting using up/down keys		
	Operation mode	Fixed temperature operation, quick auto stop operation Auto stop operation, auto start operation		
	Sensor	K thermocouple		
	Auxiliary functions	Lock function, power outage compensation function, calibration offset function		
Safety device	Controller Self-diagnosis function	Temperature sensor error, display error, temperature input circuit error, measured temperature error		
	Protection device	ELB, hydraulic pressure standalone overheat prevention device		
Spec.	Inner dimension (mm) (w×d×h)	450×450×450		620 × 600 × 1195
	Outer dimension (mm) (w×d×h)	504×562×788	504×562×820	674 × 711 × 1586 674 × 711 × 1618
	Shelf board step, load capacity	10 steps 15kg/pc.		29steps 15kg/pc.
	Shelf support pitch	30mm		
	Internal capacity	92L		445L
	Weight	Approx. 45kg	Approx. 48kg	Approx. 78kg Approx. 83kg
	Power supply (i50/60Hz)	AC115V 9A	AC115V 9A	AC115V 12A AC115V 12A
	Accessories	DG400C/DG440C shelf plate -2pcs, DG800C/840C shelf plate -4pcs, Instruction Manual, Warranty Card		

*DG400C/800C performance: value at power AC115V, DG450C/850C performance: value at power AC220V

*Operating environmental temperature range for this device is 5°C~35°C.

10. Specifications


Name	Ware Drying Oven			
Model	DG410C	DG450C	DG810C	DG850C
System	Natural convection	Natural convection Forced circulation air exhaust	Natural convection	Natural convection Forced circulation air exhaust
Temperature control range	RT+5℃~70℃			
Structure	Interior material	Stainless steel SUS304		
	Observation window	250×300mm toughened glass×1	250×700mm toughened glass×1	
	Heater	Iron-chrome heater 1.0KW		Iron-chrome heater 1.34KW
	Air exhaust fan	—	Moisture-proof axial flow fan	— Moisture-proof axial flow fan
	Bactericidal lamp	—	Moisture-proof fluorescent tube 8W	— Moisture-proof fluorescent tube 15W
Control	Control system	PID control of heater output with a micro-computer		
	Setting system	Digital setting using up/down keys		
	Operation mode	Fixed temperature operation, quick auto stop operation Auto stop operation, auto start operation		
	Sensor	K thermocouple		
	Auxiliary functions	Lock function, power outage compensation function, calibration offset function		
Safety device	Controller Self-diagnosis function	Temperature sensor error, display error, temperature input circuit error, measured temperature error		
	Protection device	ELB, hydraulic pressure standalone overheat prevention device		
Spec.	Inner dimension (mm) (w×d×h)	450×450×450		620 × 600 × 1195
	Outer dimension (mm) (w×d×h)	504×562×788	504×562×820	674 × 711 × 1586 674 × 711 × 1618
	Shelf board step, load capacity	10 steps 15kg/pc.		29steps 15kg/pc.
	Shelf support pitch	30mm		
	Internal capacity	92L		445L
	Weight	Approx. 45kg	Approx. 48kg	Approx. 78kg Approx. 83kg
	Power supply (i50/60Hz)	AC220V 5A	AC220V 5A	AC220V 6.5A AC220V 6.5A
Accessories	DG410C/DG450C shelf plate -2pcs, DG810C/850C shelf plate -4pcs, Instruction Manual, Warranty Card			

*DG410C/450C/810C/850C performance: value at power AC220V.

*Operating environmental temperature range for this device is 5℃~35℃


11. Wiring diagram

DG400C/800C



Symbol	Part name	Symbol	Part name
ELB	Circuit breaker	OH	Thermostat
H	Heater	TH	Temperature sensor (K-thermocouple)
T	Terminal block	CONT	Planar board
SSR	SSR	PIO	Display circuit board
X1, X2	Relay		


DG410C/810C



Symbol	Part name	Symbol	Part name
ELB	Circuit breaker	OH	Thermostat
H	Heater	TH	Temperature sensor (K-thermocouple)
T	Terminal block	CONT	Planar board
SSR	SSR	PIO	Display circuit board
X1, X2	Relay		

11. Wiring diagram

DG450C/850C



Symbol	Part name	Symbol	Part name
ELB	Circuit breaker	X1,X2	Relay
H	Heater	L1	Indicator lamp of air exhaust fan
T	Terminal block	L2	Indicator lamp of sterilization
SSR	SSR	S1	Fan motor switch
OH	Thermostat	FM	Air exhaust fan
TH	Temperature sensor (K-thermocouple)	S2	Door switch
CONT	Planar board	S3	Bactericidal lamp switch
PIO	Display circuit board	SL	Bactericidal lamp
INV	Ballast		

12. List of replacement parts

Common parts

Symbol	Part name	Code No.	Specifications	Manufacturer
TH	Temperature sensor	H010101001	T0304.01-08 Φ 3.2*55*2000	Yamato
CONT	Planar board	B011401002	CN40B-Y	Yamato
PIO	Display circuit board	B011402002	CN40B-Y	Yamato
-	Tough card	B011299001	15P 300mm	Yamato
SSR	SSR	A011006023	KS15/D-38Z25-L	Yamato
OH	Thermostat	B020103002	55. 13225. 070	Yamato

DG400C

Symbol	Part name	Code No.	Specifications	Manufacturer
-	Power cord	A011209001	3*2.0sq 3m	Yamato
ELB	Circuit breaker	A010410004	BV-DN IP+N 16A 30mA	Yamato
X1	Relay	A011002001	HF116F-2/110AL1HSTFW	Yamato
X2	Relay	A011002005	HF13F/A1002Z1D	Yamato
H	Heater	H020601010	DG400A_01_04-14(115V,100W)	Yamato

DG800C

Symbol	Part name	Code No.	Specifications	Manufacturer
-	Power cord	A011209001	3*2.0sq 3m	Yamato
ELB	Circuit breaker	A010410004	BV-DN IP+N 16A 30mA	Yamato
X1	Relay	A011002001	HF116F-2/110AL1HSTFW	Yamato
X2	Relay	A011002005	HF13F/A1002Z1D	Yamato
H	Heater	H020601011	DG800C.01.02-06(115V,1340W)	Yamato

12. List of replacement parts

DG410C

Symbol	Part name	Code No.	Specifications	Manufacturer
-	Power cord	A011209001	3x2.0sq	Yamato
ELB	Circuit breaker	A010410003	BV-DN 1P+N 6A 30mA	Yamato
X1	Relay	A011002002	HF116F-2/220AL1HSTFW	Yamato
X2	Relay	A011002007	HF13F/A2202Z1D	Yamato
H	Heater	H020601004	DG410C_01_04-14(220V 1000W)	Yamato

DG810C

Symbol	Part name	Code No.	Specifications	Manufacturer
-	Power cord	A011209001	3x2.0sq	Yamato
ELB	Circuit breaker	A010410007	BV-DN 1P+N 10A 30mA	Yamato
X1	Relay	A011002002	HF116F-2/220AL1HSTFW	Yamato
X2	Relay	A011002007	HF13F/A2202Z1D	Yamato
H	Heater	H020601005	DG810C.01.02-06(220V,1340W)	Yamato

12. List of replacement parts

DG440C/840C

Symbol	Part name	Code No.	Specifications	Manufacturer
-	Power cord	A011209001	3*2.0sq 3m	Yamato
ELB	Circuit breaker	A010410004	BV-DN IP+N 16A 30mA	Yamato
X1	Relay	A011002001	HF116F-2/110AL1HSTFW	Yamato
X2	Relay	A011002005	HF13F/A1002Z1D	Yamato
S1	Fan motor switch	B011501002	DS-850S-F2-10	Yamato
FM	Air exhaust fan	B080104002	R87F-A1A16H-WR	Yamato
S2	Door switch	A011505004	Z-15GQ22-B	Yamato
S3	Bactericidal lamp switch	B011501002	DS-850S-F2-10	Yamato
L1, L2	Indicator lamp of air exhaust fan	A011102072	AC110V G	Yamato

DG440C

Symbol	Part name	Code No.	Specifications	Manufacturer
H	Heater	H020601010	DG400A_01_04-14(115V,1000W)	Yamato
SL	TUV Lamp	A081702037	TUV 8W	Yamato
INV	UV Ballast (110V)	A010799009	RW12-180-10A	Yamato
SC	T5 Lamp feet	A011399020	T5 Lamp feet	Yamato

DG840C

Symbol	Part name	Code No.	Specifications	Manufacturer
H	Heater	H020601011	DG800C.01.02-06(115V,1340W)	Yamato
SL	TUV Lamp	A081702033	TUV 15W	Yamato
INV	UV Ballast (110V)	A010799004	RW12-425-18A	Yamato
SC	T8 Lamp feet	A011399019	T8 Lamp feet	Yamato

12. List of replacement parts

DG450C/850C

Symbol	Part name	Code No.	Specifications	Manufacturer
-	Power cord	A011209001	3*2.0sq 3m	Yamato
X1	Relay	A011002002	HF116F-2/220AL1HSTFW	Yamato
X2	Relay	A011002007	HF13F/A2202Z1D	Yamato
S2	Door switch	A011505004	Z-15GQ22-B	Yamato
FM	Air exhaust fan	A080104037	KA1238XA2 (IP55) BMT (L)	Yamato
L1, L2	Indicator lamp of air exhaust fan	A011102073	AC220V G	Yamato

DG450C

Symbol	Part name	Code No.	Specifications	Manufacturer
ELB	Circuit breaker	A010410003	BV-DN 1P+N 6A 30mA	Yamato
S1, S3	Switch	A011501003	MR-6-210-C5L-BGAA-220	Yamato
H	Heater	H020601004	DG410C_01_04-14	Yamato
SL	TUV Lamp	A081702037	TUV 8W	Yamato
INV	UV Ballast (220V)	A010799008	RW12-180-10	Yamato
SC	T5 Lamp feet	A011399020	T5 Lamp feet	Yamato

DG850C

Symbol	Part name	Code No.	Specifications	Manufacturer
ELB	Circuit breaker	A010410007	BV-DN IP+N 10A 30mA	Yamato
S1, S3	Switch	B011501002	DS-850S-F2-10	Yamato
H	Heater	H020601005	DG810C.01.02-06	Yamato
SL	TUV Lamp	A081702033	TUV 15W	Yamato
INV	UV Ballast (220V)	A010799003	RW12-425-18	Yamato
SC	T8 Lamp feet	A011399019	T8 Lamp feet	Yamato

13. List of dangerous materials



Never use an explosive substance a flammable substance or a substance containing them for this device.

Explosive substance	Explosive substance	① Nitroglycol, glycerine trinitrate, cellulose nitrate and other explosive nitrate esters
	Explosive substance	② Trinitrobenzen, trinitrotoluenem, picric acid and other explosive nitro compounds
Flammable substances	Explosive substances	③ Acetyl hydroperoxide, methyl ethyl ketone peroxide, benzoyl peroxide and other organic peroxides
	Oxidizing substances	Metal "lithium", metal "potassium", metal "natrium", yellow phosphorus, phosphorus sulfide, red phosphorus, celluloids, calcium carbide (a.k.a, carbide), lime phosphide, magnesium powder, aluminum powder, metal powder other than magnesium and aluminum powder, sodium dithionous acid (a.k.a., hydrosulphite)
		① Potassium chlorate, sodium chlorate, ammonium chlorate, and other chlorates
		② Potassium perchlorate, sodium perchlorate, ammonium perchlorate, and other perchlorates
		③ Potassium perchlorate, sodium perchlorate, ammonium perchlorate, and other inorganic perchlorates
		④ Potassium nitrate, sodium nitrate, ammonium nitrate, and other nitrates
		⑤ Sodium chlorite and other chlorites
		⑥ Calcium hypochlorite and other hypochlorites
	Flammable substances	① Ethyl ether, gasoline, acetaldehyde, propylene chloride, carbon disulfide, and other substances with ignition point at a degree 30 or more degrees below zero.
		② n-hexane, ethylene oxide, acetone, benzene, methyl ethyl ketone and other substances with ignition point between 30 degrees below zero and less than zero.
		③ Methanol, ethanol, xylene, pentyl acetate, (a.k.a.amyl acetate) and other substances with ignition point between zero and less than 30 degrees.
		④ Kerosene, light oil, terebinth oil, isopenthy alcohol (a.k.a. isoamyl alcohol), acetic acid and other substances with ignition point between 30 degrees and less than 65 degrees.
	Combustible gas	Hydrogen, acetylene, ethylene, methane, ethane, propane, butane and other Substance which is a flammable gas at 15 degrees, one air pressure.

14. Standard installation manual

* Install the product according to the following: (Confirm separately for optional items or special specifications)

Model	Serial number	Date	Installation mgr. (company name)	Installation mgr.	Judgment

No	Item	Implementation method	TOC No. Reference page of the operating instruction manual	Judgment
Specifications				
1	Included items	Check for number of staffs against the included item field	10. Specifications field P.44	
2	Installation	· Visual check of environmental conditions Caution: Take care for environment	2. Before operating the unit P.5 · On the installation site	
		· Securing a space		
Operation-related matters				
1	Source voltage	· Measure the user side voltage (outlet) with a tester · Measure voltage during operation (shall meet the standard) Caution: Always use a plug that meets the specification for attaching to the ELB.	2. Before operating the unit P.5 · Be sure to connect the ground wire. P.8 · Power supply is P.44 10.Specifications · Specification-power supply	
2	Operation start	· Starts operation Performs fixed value operation, auto stop operation or auto start operation	2. Before operating the unit P.4~10 · Installation P.17~30 4. Operating procedures	
Description				
1	Operational descriptions	Explain operations of each component according to the operational instructions	4. Operating procedures P.17~30 · Operating procedures P.1~ 1. Safety precautions ~ ~50 13.List of dangerous materials	
2	Error codes	Explain the customer about error codes and procedures for release according to the operational instructions	8. Troubleshooting P.41 ~ 9. After sales service and warranty 43	
3	Maintenance and inspection	Explain operations of each component according to the operational instructions	6. Maintenance procedures P.36 · Daily inspection/maintenance	
4	Completion of installation Entries	· Fill in the installation date and the installation mgr. on the nameplate of the main unit · Fill in necessary information to the warranty card and hand it over to the customer · Explanation of the route for after-sales service	9. After sales service and warranty P.43	

Responsibility

Please follow the instructions in this document when using this unit. Yamato Scientific has no responsibility for the accidents or breakdown of device if it is used with a failure to comply.

Never conduct what this document forbids. Unexpected accidents or breakdown may result in.

Note

- ◆ The contents of this document may be changed in future without notice.
- ◆ Any books with missing pages or disorderly binding may be replaced.

Instruction Manual

Ware Drying Oven

Model DG400C/410C/440C/450C/800C/810C/840C/850C

Second Edition Mar. 22, 2022

Revision Aug. 18, 2023

YAMATO SCIENTIFIC CO., LTD.
Harumi Triton Square Y-36F, 1-8-11 Harumi,
Chuo-ku, Tokyo 104-6136, Japan
Tel : +81-3-5548-7122
Fax : +81-3-5548-0132
<https://www.yamato-scientific.com/>