

## Lab Autoclaves | Automatic Sterilization, Touch Screen Operation

## SJ511C/811C

Operating temp. range 45~138°C

Max. working pressure 0.3MPa

Internal capacity 50L 80L

The steam sterilizer features a 7-inch color touchscreen for simple operation with advanced functions.



## Features

- 7-inch large touchscreen, easy operation with options for Chinese, Japanese and English languages, along with a three-level authorization management system.
- Maximum sterilizing temperature goes up to 138°C, suitable for protein modification.
- By simple setting and operation, these sterilizers manage ordinary sterilization, sterilization of culture mediums and liquids, and dissolution of culture mediums.
- Custom programs can be saved for repeated use.
- Standard equipped with cooling fan to shorten cool down time.
- Equipped with multiple pressure lid opening protection locks and comprehensive safety measures.
- Display power consumption, accumulated power usage, CO<sub>2</sub> emissions, total power-on time, accumulated operating time, and accumulated times of operations.
- Two handles on each side for easy moving.
- Optional automatic water filling function available for purchase.

## Specifications

Model	SJ511C	SJ811C	
System	Automatic high-pressure steam sterilization		
Performance	Operating temp. range	45~138°C: 45~80°C (preheat temp.)/45~60°C (retain temp.)/60~110°C (dissolution)/105~138°C (sterilizing)	
	Temp. unit	0.1°C	
	Max. working pressure	0.3MPa	
Lid structure	Manual upward opening and downward closing with an interlock for safety		
Chamber	Material	3mm thick corrosion-resistant stainless steel, featuring a mirror finish surface	
	Design operational life	20 years	
	Design pressure	0.42MPa	
	Design temp.	151°C	
	Heater	110V, 1500W×2	110V, 2000W×2
	Exhaust valve	Solenoid valves: One for full exhaust and one for slow exhaust	
	Options	Sample sensor interface (1/4), data logger sensor interface (1/4), pressure gauge interface (on solenoid valve piping), automatic water intake (R1/4)	
	Cooling fan blade	Axial fan motor	
Controllers	Temp. control method	PID control by microprocessor	
	Temp. control setting/display method	7-inch color LCD touchscreen for display and settings	
	Timer	1 min~99 h 59 min, resolution: 1 min or 100 h~999 h, resolution: 1 h	
	Operation functions	Instrument sterilization, liquid sterilization, sterilization and retain temp., dissolution and retain temp.	
	Other functions	Power consumption, accumulated power usage, CO <sub>2</sub> emissions, accumulated power-on time, accumulated operating time, and accumulated times of operations, key locking, presetting, saving, preheating, forced cooling, pattern locking, alarm log saving (20 pcs), display of present time, ON-OFF beeping setting, language switch function (Chinese/Japanese/English) and three-level authorization management system, sample sensor (optional)	
Safety device	Sensor failure detection detection, SSR short-circuit detection, broken heater wire detection, prevention of air burning (water level switch, liquid expansion type), alarm against the absence of drain bottle, failure in locking the lid detection, memory error detection, pressure relief valve (0.3 MPa)		
Pressure vessel specification	Class I pressure vessel		
Specification	External dimensions (W×D×H mm)	520×660×881	520×660×1161
	Effective dimensions of chamber (diameter×height mm)	370×465	370×745
	Internal effective capacity of chamber	50L	80L
	Weight	Approx. 95kg	Approx. 115kg
	Power supply rated current	AC220V 14A	AC220 19A
Accessories	Basket	OSQ-90 (diameter 344 × height 200mm), 2 PCS	Basket OSQ-90H (diameter 344 × height 300mm), 2 PCS
		1 steam collection cup, 1 cooling water bottle, 1 heater baffle, 1 set of sterilization test cards (30 cards), 1 filter plug, 1m drainage pipe, 1 hose clamp, long-arm clamp (SJ811C model only)	
Options	Biosafety filtration system, sample temperature sensor, internal chamber temperature measurement sensor, miniprinter, data logger, temperature output terminal (4~20mA), external alarm output terminal, time arrival output terminal, standard basket, extended height basket, automatic water filling system		

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers &amp; Cold Traps

12 Stirrers &amp; Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Disinfection Certificate

Low bed mode

- Convenient retrieval of sterilized items.




Standard Cooling Fan

- Cooling fan operates after sterilization.
- Cooling down to safe temperature.
- Reducing the time before samples can be taken out.
- Free to choose the forced cooling or natural cooling.

Comparison of cooling time with other devices

(Initial water temperature: 25°C, sterilization temperature: 121°C, ambient temperature: 23°C)




Control panel

- 7-inch color LCD touchscreen with operation in Chinese, English and Japanese



Main screen



Sterilization parameter settings (instrument sterilization program)



Operation mode selection (instrument sterilization program)




Operation mode selection (melt and retain temp. program)



Operating status monitoring (instrument sterilization program)

Sterilization starts automatically according to sample temperature sensor



The optional sample temperature sensor controls sterilization time based on the temperature of the sterilized material, ensuring effective sterilization.

Comply with GLP/GMP inspection requirements



Standard equipped with 2 sensor ports on the main unit.



Data logger or built-in micro printer is optional

Safety valve located on the side



Convenient disassembly for safety valve replacement or calibration.

Sterilizers	1
Granulation and Spray Dryers	2
Muffle Furnaces	3
Ovens	4
Incubators	5
Plasma Equipment	6
Water Purifiers	7
Baths	8
Water Circulators	9
Rotary Evaporators	10
Freeze Dryers & Cold Traps	11
Stirrers & Shakers	12
Washers	13
Analysis and Test Devices	14
Options	15

## Lab Autoclaves | Automatic Sterilization, Touch Screen Operation

## SJ530/830

Operating temp. range 45~135°C

Max. working pressure 0.255MPa

Internal capacity 50L 80L

The autoclave features a 7-inch color touchscreen for simple operation with advanced functions.



## Features

- 7-inch large touchscreen, easy operation with options for Chinese, Japanese and English languages, along with a three-level authorization management system.
- Maximum sterilizing temperature goes up to 135°C, suitable for protein modification.
- By simple setting and operation, these sterilizers manage ordinary sterilization, sterilization of culture mediums and liquids, and dissolution of culture mediums.
- Custom programs can be saved for repeated use.
- Standard equipped with cooling fan to shorten cool down time.
- Equipped with multiple pressure lid opening protection locks and comprehensive safety measures.
- Display power consumption, accumulated power usage, CO<sub>2</sub> emissions, total power-on time, accumulated operating time, and accumulated times of operations.
- SJ530/830 has obtained medical device registration and disinfection certificates.

## Specifications

Model	SJ530	SJ830	
System	Automatic high-pressure steam sterilization		
Medical device registration certificate	Yes		
Performance	Operating temp. range	45~135°C: 45~80°C (preheat temp.)/45~60°C (retain temp.)/60~110°C (dissolution)/105~135°C (sterilizing)	
	Temp. unit	0.1°C	
	Max. working pressure	0.255MPa	
Lid structure	Manual upward opening and downward closing with an interlock for safety		
Chamber	Material	3mm thick corrosion-resistant stainless steel, featuring a mirror finish surface	
	Design operational life	20 years	
	Design pressure	0.42MPa	
	Design temp.	151°C	
	Sterilization heater	110V, 1210W×2	110V, 1500W×2
	Exhaust valve	One quick exhaust valve and one slow release valve	
	Options	Sample sensor interface (1/4), data logger sensor interface (1/4), pressure gauge interface (on solenoid valve piping)	
Controllers	Cooling fan blade	Axial fan motor	
	Temp. control method	PID control by microprocessor	
	Temp. control setting/display method	7-inch color LCD touchscreen for display and settings	
	Timer	1 min~99 h 59 min, resolution: 1 min or 100~999 h, resolution: 1 h	
	Operation functions	Instrument sterilization, liquid sterilization, sterilization and retain temp., dissolution and retain temp.	
Other functions	Power consumption, accumulated power usage, CO <sub>2</sub> emissions, accumulated power-on time, accumulated operating time, and accumulated times of operations, key locking, presetting, saving, preheating, forced cooling, sample temperature sensor (optional), pattern locking, error log saving (20 pcs), display of present time, ON-OFF beeping setting, language switch function (Chinese/Japanese/English) and three-level authorization management system		
Safety device	Sensor failure detection detection, SSR short-circuit detection, broken heater wire detection, prevention of air burning (liquid expansion type thermoet), alarm against the absence of drain bottle, failure in locking the lid detection, memory error detection, pressure relief valve (0.255 MPa)		
Pressure vessel specification	Class I pressure vessel		
Specification	External dimensions (W×D×H mm)	520×660×881	520×660×1161
	Effective dimensions of chamber (diameter×height mm)	370×470	370×750
	Internal effective capacity of chamber	50L	80L
	Weight	Approx. 105kg	Approx. 125kg
	Power supply/rated current	AC220V 11.5A	AC220V 14A
Accessories	Basket	Basket OSQ-90 (diameter 344 × height 200mm), 2 PCS	Basket OSQ-90H (diameter 344 × height 300mm), 2 PCS
		1 steam collection cup, 1 cooling water bottle, 1 heater baffle, 1 set of sterilization test cards (30 cards), 1 filter plug, 1m drainage pipe, 1 hose clamp, long-arm clamp (SJ830 model only)	
Options	Biosafety filtration system, sample temperature sensor, internal chamber temperature measurement sensor, miniprinter, data logger, temperature output terminal (4~20mA), external alarm output terminal, time arrival output terminal, standard basket, extended height basket		

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers &amp; Cold Traps

12 Stirrers &amp; Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Medical Device  
Registration  
Certificate

Disinfection  
Certificate

**Low bed mode**


- Convenient retrieval of sterilized items.



**Standard Cooling Fan**

- Cooling fan operates after sterilization.
- Cooling down to safe temperature.
- Reducing the time before samples can be taken out.
- Free to choose the forced cooling or natural cooling.

**Comparison of cooling time with other devices**



**Control panel**


- 7-inch color LCD touchscreen with operation in Chinese, English, and Japanese.



Main screen



Sterilization parameter settings (instrument sterilization program)




Operation mode selection (instrument sterilization program)



Operation mode selection (melt and retain temp. program)



Operating status monitoring (instrument sterilization program)



**Sterilization starts automatically according to sample temperature sensor**

The optional sample temperature sensor controls sterilization time based on the temperature of the sterilized material, ensuring effective sterilization.

**Comply with GLP/GMP inspection requirements**



Standard equipped with 2 sensor ports on the main unit.



Optional data logger and high performance pressure gauge.

**Space-saving top-opening lid design**

Require less setup space compared to other machines.

**Maximum sterilization temperature of 135°C**

Maximum sterilizing temperature goes up to 135°C suitable for protein modification.

**Time-saving**

Standard equipped with automatic start and preheat function for more efficient time using.

Sterilizers	1
Granulation and Spray Dryers	2
Muffle Furnaces	3
Ovens	4
Incubators	5
Plasma Equipment	6
Water Purifiers	7
Baths	8
Water Circulators	9
Rotary Evaporators	10
Freeze Dryers & Cold Traps	11
Stirrers & Shakers	12
Washers	13
Analysis and Test Devices	14
Options	15