

2025-2026

Laboratory Instruments

Corporate Profile

■ Management philosophy

Support R&D and production technology innovation in various fields and contribute to the progress and development of science and technologies.

■ Management principles

As an important member of Yamato Scientific, constantly provide customers the latest and the most valuable products and the service of safety and ease, cooperate with group company, actively expand group business.

To achieve the goal, the company takes laboratory instruments, laboratory facilities and industrial equipment as the products core, engages in design and development, production, sales and service, through the establishment of high-effective team, advanced equipment and perfect system, provide customers with high-quality products and services.

■ Environmental principles

In line with its management philosophy, the company is committed to the proactive promotion of environmental protection activities in all aspects of its business activities, capable of sustaining environmental protection as well as corporate profitability.

■ Environmental policy

1. Compliance with environmental regulations

We will comply with all environment-related regulations and adhere to the provisions of all official agreements and commitments.

2. Improvement of business processes in order to reduce the burden on the environment

We will eliminate all unjustified, unnecessary and inconsistent business processes in order to make more effective use of natural resources and energy and reduce wastes.

3. Promotion of business activities that attempt to reduce the burden on the environment

We will promote activities that will reduce the burden on the environment throughout all of our business processes: design, purchase, production, sales, distribution, use, recycling and disposal.

4. Development of environmentally-conscious products & systems

We will try to improve technologies that promote environmental protection and develop products and systems that remain environmentally friendly throughout their lifecycles.

5. Improvement and development of an environmental management system

We will continue to improve our environmental activities and pollution prevention measures by setting goals and objectives and reviewing them regularly.

6. Notification and publication of our environmental policy

We will educate our employees, and all other people who work for us, to make them familiar with our environmental policy so that they will act in accordance with its principles. We will also make our environmental policy public.

■ Acquisition of ISO certification



Japanese Factory
ISO9001



Japanese Factory
ISO14001



Japanese Factory
ISO13485



Chinese Factory
ISO9001



Chinese Factory
ISO14001



Chinese Factory
OHSAS18001

Design & Development


Japan

The R&D center of Yamato Scientific in Japan undertakes the product development. In a high-efficient development environment, they are committed to new products development contributing to production technology innovation. All the new achievements, not only are provided for enterprise, university and government around the world, but also make a contribution to the world's most cutting-edge research and development (the Nobel Prize) and production technology innovation.



China

Yamato Chongqing has set up a high-efficient R&D team and perfect R&D system, by introducing 3d - CAD and product data management system, the product design is more standard and efficient, and also co-design with Yamato Scientific. Depending on Yamato Scientific's technology accumulation of more than 120 years and its technical support, design and develop all kinds of instruments, equipment and industrial equipment to meet customer needs and ever-changing market needs.



Production & Manufacture

The factory has production lines necessary for metal plate processing, coating, assembling, inspection and packing, and produces various instruments and equipment. Metal plate processing line is 24hr working NCT auto production line; coating line consists of powder coating and paint coating; assembly line consists of multiple production lines. Multi-varieties and small batch production mode is adopted for laboratory instruments and equipment. These production adjustments that we match with the purchase order process enable us to shorten delivery time, plan adequate stock and supply management.

Japanese Factory



Chinese Factory



1. Metal plate processing line

Yamato has NCT automatic processing equipment for punching and cutting raw materials (steel plate, stainless steel, aluminum) and equipment for bending and welding of formed material according to NC programming. The production of processed products from raw materials provides a good preparation for subsequent process.

2. Coating line

In addition to paint coating (30 micrometers thick), we also use a powder coating (60 micrometers thick), which largely improves the production capacity.

3. Assembly line

Flexibly adopt independent assembly production mode and cope with mass production in delivery period according to the orders. Visualize the project process and continue to improve efficiency.

4. Inspection line

Strictly inspect the product's appearance, electrical performance, temp. performance, etc. to ensure and maintain the quality of Yamato products.

Y Laboratory Instruments



Various product lines that precisely capture market demands

Laboratory instruments, as Yamato's main product line, contains sterilizer, spray dryer, plasma cleaner, water purifier, muffle furnace, constant temp. drying oven, constant temp. incubator, constant temp. water bath, constant temp. water circulator, rotary evaporator, freeze dryer, stirrer-shaker, washer, etc. Furthermore, we plan a lineup suitable for customer needs and target at improvement of research environment by providing high-quality products.



Y Industrial Equipment



Contribute to increase productivity

From lithium battery, semiconductor, electric parts to information communication, electronic materials, FPD (LCD, organic EL, PDP and others), precision machine, etc., we provide various equipment and machines to increase the industry productivity.



Laboratory Facilities



Provide suitable environment and new values

We offer equipment needed by various research and development such as fume hood, laboratory tables, clean benches and environmental test rooms. On the basis of standard products, we offer fume hood with built-in water purifier or air conditioners for external atmospheric treatment that meet the requirements of the laboratory.

The laboratory tables are made of pure steel(solid top board), allowing for pharmaceuticals resistance test, as well as special edge processing.

Furthermore, we work constantly with customers to improve the work environment, including ease of movement in the lab and the solution to the storage space and environmental dissatisfaction.

Focus on the ideal next-generation laboratory

Yamato Scientific is able to develop, design and produce laboratory instruments and equipment in its own factory.

With four basic concepts in mind: space saving, efficient work system, storage space and safety measures, freely design to meet the requirements of each customer. Inspection and experimentation devices are chosen for various fields of examination. Provide a free and creative environment based on the theme of people, environment and research and continuously improve to meet the customer's needs.



Support the design, establishment and relocation of laboratory

For the design, establishment and relocation of laboratory, we provide synthetic technical support such as environmental pollution countermeasures and model selection. As per customer's requirements, from basic proposal to enforcement plan, basic design and actual operation, we can satisfactorily solve all the problems about establishment and relocation. Yamato Scientific has professional group and special views about laboratory instruments and equipment, and provides a full range of laboratory services.



Excellent technical strength to help you successfully complete the relocation of laboratories and production lines



The environmental analysis and overall planning of laboratory instruments and equipment when relocation, from laboratory facilities like experiment table and fume hood to experimental data and samples relocation, we provide safe, quick and accurate service. In addition, we plan the layout considering the working efficiency so as to create a high-efficiency and safe laboratory environment.

Laboratory Assistant



Intelligent Robot·Laboratory Assistant

On site of research or quality control, a lot of work need manual operation. Some work does not realize automation, their results depend on human skill and technique. Subdividing these work, tracing the skill and technique, and forming the synthetic system of numerical criteria, these are realized by Intelligent Robot·Laboratory Assistant.

Advantages of Laboratory Assistant

- Different from the former robot system which targeting increase of handling capacity, Yamato robot system can trace the operation procedure confirmed by researcher, realize the accuracy and reproducibility which is hard to achieve by manual operation, and also solve the problem of long-time and continuous operation for manual operation.
- Besides conduct the effective operation layout as per customer's requirement, in order to effectively make use of the laboratory space, the equipment matched with robot also can be customized.
- The robot system, based on human hand operation, saves the training fee of operators and their skill retraining fee.



Characteristics of Laboratory Assistant

- Specially designed arm and auxiliary instrument·equipment
In order to realize the operation demanded by customer, skill imitation is very important. According to motions required by customer, make special arm with additional functions, and it's able to develop, design and produce the synthetic automation of auxiliary instrument·equipment as per operation requirements.
- Combine with analytical and measuring instruments
Proficient in the technique and characteristics of analytical and measuring instruments, choose suitable instrument as needed.
- Pay attention to lead in cost system proposal
Make use of the existing instruments and equipment of customer to reduce the system cost.



Application & Marketing Performance




Cases


Automatic Cell Sample Modulation System (Co-research with Tokyo University)



Application Fields



From Order to Deliver



Sales and Service Network

Domestic Bases and Distributors Network (48 bases / 185 distributors)



**Headquarters of
Yamato Scientific / SUNMEDIX**
Harumi Triton Square
Office Tower Y (36F),
1-8-11 Harumi, Chuo-ku, Tokyo




**Headquarters of
Yamato Material**
24F Kyobashi Edogrand,
2-2-1 Kyobashi, Chuo-ku, Tokyo


Yamato Scientific Group
Gross Sales:
106.50 billion yen
Employees: **1,526**
as of January 2025



**Yamato Scientific Group
Head office**
2-2-2 Nihonbashihoncho,
Chuo-ku, Tokyo



Overseas Bases and Distributors Network (8 bases / 323 distributors)



Laboratory Instruments

| | |
|------------------------------|-------|
| Sterilizers | P.21 |
| Granulation and Spray Dryers | P.37 |
| Muffle Furnaces | P.51 |
| Ovens | P.59 |
| Incubators | P.105 |
| Plasma Equipment | P.127 |
| Water Purifiers | P.137 |
| Baths | P.145 |
| Water Circulators | P.167 |
| Rotary Evaporators | P.181 |
| Freeze Dryers & Cold Traps | P.193 |
| Stirrers & Shakers | P.197 |
| Washers | P.215 |
| Analysis and Test Devices | P.225 |
| Options | P.229 |

A1 | Sterilizers

P.21

SN-SQ



Name: Lab Autoclaves
Model: SN210C/310C/510C, SQ510C/810C
Feature: Standard Type
Page 22

SJ NEW



Name: Lab Autoclaves
Model: SJ511C/811C
Feature: Automatic Sterilization, Touch Screen Operation
Page 24

SJ



Name: Lab Autoclaves
Model: SJ530/830
Feature: Automatic Sterilization, Touch Screen Operation
Page 26

SM



Name: Lab Autoclaves
Model: SM510C/810C/530/830
Feature: Automatic Sterilization, Drying Function
Page 28

SQL




Name: Lab Autoclaves
Model: SQL810C/1010C
Feature: Automatic Sterilization, Large Capacity
Page 30

SJL



Name: Lab Autoclaves
Model: SJL810C/1010C
Feature: Automatic Sterilization, Large Capacity
Page 32

SIT NEW



Name: Dry Heat Sterilizers
Model: SIT410C/610C/810C
Feature: Dry heat sterilization, Standard Type
Page 34

SI



Name: Dry Heat Sterilizers
Model: SI411C/611C
Feature: Dry Heat Sterilization, Gravity Convection
Page 35

SL




Name: Loop Cinerator
Model: SL-21
Feature: Loop Sterilization
Page 36

A2 | Granulation and Spray Dryers

P.37

ADL




Name: Spray Dryer
Model: ADL312SC-A
Feature: Standard Type
Page 38

GB



Name: Spray Dryer
Model: GB211C-A
Feature: High Functionality Type
Page 40

GB



Name: Spray Dryer (For Granulating, Drying, Mixing)
Model: GB211C-B
Feature: Powder Granulating, Drying, Mixing
Page 42

DL



Name: Spray Dryer
Model: DL411C
Feature: Large Capacity
Page 44

GAS

NEW



Name: Organic Solvent Recovery Unit
Model: GAS411C/510C
Feature: N₂ Closed Cycle

Page 46

GWS



Name: Solvent Washing Unit
Model: GWS410
Feature: Corrosive Solvent Wet Washing

Page 48

A3 Muffle Furnaces

P.51

FO



Name: Muffle Furnace
Model: FO111C/211C/311C/411C/511C/611C/711C/811C
Feature: Standard Type

Page 52

FP



Name: Muffle Furnace
Model: FP111C/311C/511C
Feature: Built-in Heating Wires

Page 54

FQ

NEW



Name: Muffle Furnace
Model: FQ310C
Feature: 1700°C

Page 56

A4 Ovens

P.59

DKT

NEW



Name: Forced Convection Ovens
Model: DKT410C/610C/810C
Feature: Programmable, Touch Screen

Page 62

DKN

NEW



Name: Forced Convection Ovens
Model: DKN313C/413C/613C/813C/913C
Feature: Standard Type, Program Operation

Page 64

DKM

NEW



Name: Forced Convection Ovens
Model: DKM311C/411C/611C
Feature: Economical Type, Fixed Temp. Operation

Page 66

DKL

NEW



Name: Forced Convection Ovens
Model: DKL311C/411C/611C
Feature: Basic Function

Page 68

DN-H



Name: Forced Convection Ovens
Model: DN410HC/610HC
Feature: High Temp.

Page 69

DF-DH



Name: Forced Convection Ovens
Model: DF411C/611C, DH411C/611C
Feature: Precision Type

Page 70

DF-DH



Name: Forced Convection Ovens
Model: DF811C/1011C, DH811C/1011C
Feature: Large Capacity, Precision

Page 72

DFS-DHS



Name: Forced Convection Ovens
Model: DFS710C/810C, DHS710C/810C
Feature: Precision Type, Vertical

Page 73

DH



Name: Forced Convection Ovens
Model: DH450C/650C/850C
Feature: Precision Type, 500°C
Page 74

DKG



Name: Forced Convection Ovens
Model: DKG611/611V/811/811V
Feature: Rapid Heating and Cooling
Page 75

DNE




Name: Forced Convection Ovens
Model: DNE411C/611C/811C/911C
Feature: Energy-Saving
Page 76

DNF



Name: Forced Convection Ovens
Model: DNF411C/611C/811C/911C
Feature: Energy-Saving Type, Variable Airflow
Page 78

DVT NEW




Name: Natural Convection Ovens
Model: DVT410C/610C
Feature: Natural Convection, High-Function Type, Touch Screen
Page 80

DVS NEW




Name: Natural Convection Ovens
Model: DVS413C/613C
Feature: Natural Convection, Standard Type, Program Operation
Page 82

DX




Name: Natural Convection Ovens
Model: DX312C/412C/612C
Feature: Natural Convection, Economical Type, Fixed Temp. Operation
Page 84

DY NEW




Name: Natural Convection Ovens
Model: DY311C/411C/611C
Feature: Natural Convection
Page 85

DN-I



Name: Inert Ovens
Model: DN410IC/610IC
Feature: N₂ Replacement, Oxygen-Free Heating
Page 86

DF-I/DH-I



Name: Inert Ovens
Model: DF611IC/DH611IC
Feature: Anaerobic, Oxygen Concentration 500ppm
Page 87

DNN



Name: Inert Ovens
Model: DNN430C/630C/460C/660C
Feature: Anaerobic, Oxygen Concentration 20ppm
Page 88

DF-S·DH-S




Name: Explosion-proof Ovens
Model: DF411SC/611SC, DH411SC/611SC
Feature: Explosion-Proof, Precision Type
Page 89

DF-S·DH-S



Name: Explosion-proof Ovens
Model: DF811SC/1011SC, DH811SC/1011SC
Feature: Explosion-Proof, Precision Type
Page 90

DH-S




Name: Explosion-proof Ovens
Model: DH450SC/650SC/850SC
Feature: Explosion-Proof, 500°C Precision Type
Page 90

DG



Name: Ovens for Labware
Model: DG410C/450C/810C/850C
Feature: Natural Convection, Forced Exhaust
Page 91

DE·DT




Name: Clean Ovens
Model: DE430C/430UC/630C/630UC, DT430C/430UC/630C/630UC
Feature: Dust-Free Environment Heating, Standard Type
Page 92

DES·DTS



Name: Clean Ovens
Model: DES830C/830UC, DTS830C/830UC
Feature: Dust-Free Environment Heating, Large Capacity Type
Page 94

DEC



Name: Clean Ovens
Model: DEC812C/912C
Feature: Dust-Free Environment Heating, Economical
Page 96

DTN




Name: Clean Inert Ovens
Model: DTN430C/630C/450C/650C
Feature: Cleanliness Level 100, Oxygen Concentration 20ppm
Page 97

DIR



Name: Far Infrared Ray Ovens
Model: DIR631C
Feature: Far Infrared
Page 98

DR



Name: High Temp. Ovens
Model: DR210C
Feature: High Temp.
Page 99

DPT NEW




Name: Vacuum Ovens
Model: DPT40C/60C/80C/100C
Feature: High Function Type, Touch Screen
Page 100

DP




Name: Vacuum Ovens
Model: DP23C/33C/43C/63C/83C/103C
Feature: Standard Type
Page 101

ADP



Name: Vacuum Ovens
Model: ADP210C/310C
Feature: Economical Type, Table Type
Page 102

DPF



Name: Vacuum Ovens
Model: DPF43C/63C
Feature: Vertical, High Precision
Page 103

DPH



Name: Vacuum Ovens
Model: DPH43C/63C
Feature: Vertical, 300°C
Page 103

DPHH



Name: Vacuum Ovens
Model: DPHH43C/63C
Feature: Vertical, 400°C
Page 104

DP-HP




Name: Vacuum Ovens
Model: DP43HPC/63HPC
Feature: Vertical, Shelf Heating
Page 104

A5 Incubators

P.105

IS



Name: High-temp. Incubators
Model: IS412C/612C/812C/912C
Feature: Program Operation, Air Jacket
Page 108

IC



Name: High-temp. Incubators
Model: IC412C/612C/812C/912C
Feature: Fixed Temp. Operation, Air Jacket
Page 110

ICF



Name: High-temp. Incubators
Model: ICF410C/610C/810C/910C
Feature: Constant Temp. Operation, Air Jacket Type, High-temp. Uniformity
Page 112

IN NEW



Name: Low-Temp. Incubators
Model: IN613C/613CW/813C/913C
Feature: Program Operations
Page 114

INE




Name: Low-temp. Incubators
Model: INE800
Feature: Inverter Energy-saving, Low-temp.
Page 116

IL



Name: Low-temp. Incubators
Model: IL612C/812C
Feature: Air Jacket Type
Page 117

IJ




Name: Low-temp. Incubators
Model: IJ102/102W/300
Feature: Semiconductor Refrigeration
Page 118

INC NEW




Name: 2-Chamber Incubators
Model: INC822C
Feature: Upper Chamber with High Temp., Lower Chamber with Low Temp.
Page 120

IQ NEW




Name: 2-Chamber Incubators
Model: IQ823C
Name: Upper And Lower Chambers with Low Temp.
Page 121

IPE



Name: CO₂ Incubators
Model: IPE610
Feature: Air Jacket, Dry Heat Sterilization
Page 122

BNA




Name: CO₂ Incubators
Model: BNA610
Feature: Water Jacket Type
Page 123

A6 | Plasma Equipment


P.129

PR




Name: Gas Plasma Reactors
Model: PR500/510/1000
Feature: DP Mode, Barrel
Page 130

PR




Name: Gas Plasma Reactors
Model: PR210C/300/301
Feature: Economic, High-Frequency Plasma
Page 131

PM




Name: Gas Plasma Cleaners
Model: PM110C
Feature: Economical, Low-Frequency Plasma
Page 132

PDC




Name: Gas Plasma Cleaners
Model: PDC200/210/510
Feature: RIE: DP Mode, Parallel Plate Electrode
Page 133

PDC




Name: Gas Plasma Cleaners
Model: PDC610
Feature: RIE: DP Mode, Multi-Stage Parallel Plate Electrode
Page 134

V




Name: Gas Plasma Cleaners
Model: V1000/1000X/1000XS
Feature: RIE: DP Mode, Multi-Stage Parallel Plate Electrode, High-capacity
Page 135

AP



Name: Gas Plasma Cleaners
Model: AP-4000
Feature: Handheld
Page 136

YSP




Name: Gas Plasma Cleaners
Model: YSP62
Feature: Online Type with Automatic Transfer
Page 136

A7 | Water Purifiers


P.137

WGH




Name: Water Purifiers
Model: WGH201
Feature: Ultra-High Performance Type, High-Purity
Page 140

WG




Name: Water Purifiers
Model: WG251/1001
Feature: Standard Type
Page 142

WG




Name: Water Purifiers
Model: WG204
Feature: Economical Type
Page 144

BK·BA




Name: Constant High Temp. Water Baths
Model: BK(A)310C/410C/510C/610C/710C
Feature: Standard Type
Page 148

BF




Name: Immersion Constant Temp. Devices
Model: BF201/401/501/601
Feature: Immersion Type
Page 149

BS



Name: Constant High Temp. Water Baths
Model: BS200/401/601/660
Feature: Liquid Expansion Type
Page 150

BM



Name: Constant High Temp. Water Baths
Model: BM510C/100/401
Feature: Economical Type
Page 151

BOA




Name: Constant High Temp. Oil Baths
Model: BOA311
Feature: Large Capacity
Page 152

BOG·BOS NEW




Name: Constant High Temp. Oil Baths
Model: BOG100/200, BOS100/200
Feature: Combination Type
Page 153

BZ




Name: Constant High Temp. Oil Baths
Model: BZ100/100D/200/300-BF601
Feature: Immersion Thermostat Combination
Page 154

BO



Name: Constant High Temp. Oil Baths
Model: BO500
Feature: Stainless Steel
Page 154

BO



Name: Constant High Temp. Oil Baths
Model: BO510C/601
Feature: Economical Type
Page 155

BN




Name: Oscillating Unit
Model: BN300
Feature: For Constant Temp. Water Baths
Page 155

BW




Name: Oscillating Water Baths
Model: BW101/201/400
Feature: Reciprocating Oscillating
Page 156

BT



Name: Oscillating Water Baths
Model: BT101/311
Feature: Horizontal Reciprocating Oscillating, Pump Stirring
Page 157

BE




Name: Immersion Coolers
Model: BE201/201F/301
Feature: Immersion Type
Page 158

BLG



Name: Constant Low Temp. Baths
Model: BLG100/200
Feature: For Material Synthesis
Page 159

BB



Name: Constant Low Temp. Water Baths
Model: BB311C/411C/611C
Feature: Standard Type
Page 160

BBL



Name: Constant Low Temp. Water Baths
Model: BBL111C/311C
Feature: Desktop Type
Page 162

BL



Name: Constant Low Temp. Water Baths
Model: BL410C/810C
Feature: Large Capacity

Page 163

BV



Name: Constant Low Temp. Water Baths
Model: BV100S
Feature: Semiconductor Cooling

Page 164

HF



Name: Heating Blocks
Model: HF100/200
Feature: Test Tube Heating

Page 165

A9 | Water Circulators

P.167

CF



Name: Cooling Water Circulators
Model: CF313-B/812-B
Feature: External Closed Loop, Standard Type

Page 170

CF



Name: Cooling Water Circulators
Model: CF720C
Feature: External Sealed System, Low Bed Type

Page 171

CFA



Name: Precision Constant Temp. Water Circulators
Model: CFA311C/610C
Feature: External Sealed System, High Precision Temp. Control

Page 172

CLS



Name: Cooling Water Circulators
Model: CLS312C/411C/610C
Feature: External Open System, Standard Type

Page 174

CLH



Name: Precision Constant Temp. Water Circulators
Model: CLH312C/411C/610C
Feature: External Open System, High-Precision Temp. Control

Page 176

CTW-CTA



Name: Precision Constant Temp. Water Circulators
Model: CTW412/812-CTA412/812
Feature: Semiconductor Cooling, Integrated Type

Page 178

CTW-S-CTA-S



Name: Precision Constant Temp. Water Circulators
Model: CTW412S/812S-CTA412S/812S
Feature: Semiconductor Refrigeration

Page 179

A10 | Rotary Evaporators

P.181

RE·REV



Name: Rotary Evaporators
Model: RE212-B/REV212M-B
Feature: Standard Type

Page 184

RE




Name: Rotary Evaporators
Model: RE311
Feature: Automatic Lifting Type

Page 190

A11 | Freeze Dryers & Cold Traps

P.193

DC



Name: Freeze Dryers
Model: DC401/801
Feature: Standard Type
Page 194

CA




Name: Cooling Traps
Model: CA301/801
Feature: Chamber Cooling Type
Page 196

A12 | Stirrers & Shakers

P.197

SFX




Name: Ultrasonic Crushers/Homogenizers
Model: SFX150-Y/250-Y/550-Y
Feature: Cell Destruction, Emulsification, Dispersion, Mixing
Page 198

SA



Name: Shaker
Model: SA300/320/400
Feature: Horizontal-Vertical-Orbital-Dual Vertical Shaking
Page 200

MK



Name: Shaker
Model: MK161
Feature: Reciprocating-Elliptical-Orbital Shaking
Page 202

MK




Name: Shaker
Model: MK201D
Feature: Reciprocating-Orbital Shaking
Page 203

MFD NEW



Name: Magnetic Stirrers
Model: MFD810-B
Feature: Magnetic Stirring
Page 204

MFH NEW



Name: Magnetic Stirrers
Model: MFH810-B
Feature: Magnetic Stirring with Heating Function
Page 205

MA



Name: Magnetic Stirrers
Model: MA300A/300B
Feature: Magnetic Stirring
Page 206

MA



Name: Magnetic Stirrers
Model: MA100/300
Feature: Economical Type
Page 206

MG



Name: Magnetic Stirrers
Model: MG120/600
Feature: Multiple/Simultaneous Stirring
Page 207

MB



Name: Magnetic Stirrers
Model: MB800
Feature: Super Strong Magnetic Stirring
Page 207

MD-MS



Name: Magnetic Stirrers
Model: MD300/500/800-MS500D
Feature: Low Viscosity to High Viscosity Compatible
Page 208

MC-MF



Name: Magnetic Stirrers
Model: MC801/MF820
Feature: Super Strong Magnetic Stirring
Page 209

MH·MG



Name:Magnetic Stirrers
Model:MH301/520/800-MG600H
Feature:With Heating Plate
Page 210

HK·HM



Name:Hot Plates
Model:HK200/300-HM300
Feature:Heating Plate
Page 211

MT



Name:Stirrers For Small Containers
Model:MT-31/51
Feature:Touch-driven Type
Page 211

LM



Name:Laboratory Stirrer
Model:LM100/200
Feature:For Flask Stirring
Page 212

LT



Name:Laboratory Stirrers
Model:LT400A/400B/400C/400D/500A/500B
Feature:Electric Stirring
Page 214

A13 Washers

P.215

M·MH·CPX·CPXH



Name:Ultrasonic Washers
Model:M-Y·MH-Y·CPX-Y·CPXH-Y
Feature:Benchtop
Page 216

AWD



Name:Laboratory Washer for Glassware
Model:AWD510
Feature:Fully-Automatic, Benchtop
Page 218

AW



Name:Laboratory Washer for Glassware
Model:AW83Z
Feature:Fully-Automatic, Large Capacity
Page 220

AW



Name:Laboratory Washer for Glassware
Model:AW62
Feature:Fully-Automatic, Benchtop
Page 221

AW



Name:Laboratory Washer for Glassware
Model:AW47
Feature:Semi-Automatic, Benchtop
Page 222

AW



Name:Ultrasonic Pipet Washer
Model:AW31
Feature:Ultrasonic Cleaning, Benchtop
Page 223

A14 | Analysis and Test Devices

P.225

TE NEW



Name: Thermal Resistance Test System
Model: TE100
Feature: Thermal Resistance Testing for Power Semiconductor Ceramic Substrates

Page [226](#)

A15 | Options

P.229


Fitting for Stacking



Name: Fitting for Stacking
Model: OD·ODN·ODK
Feature: Applicable To Chamber-Type Products

Page [230](#)


Stand



Name: Stand
Model: ON·ONS·OP·OH
Feature: Applicable To Chamber-Type Products

Page [231](#)


Data Logger



Name: Data Logger
Model: SR·AL·EX·KR
Feature: Paper-based, Paperless

Page [232](#)

Vacuum Pump



Name: Vacuum Pump
Model: GLD·NeoDry
Feature: Oil Vane, Dry

Page [233](#)


Air Compressor



Name: Air Compressor
Model: TYW-1/2
Feature: Small, Oil-Free

Page [236](#)


Transformer



Name: Transformer
Model: TR10·TR20
Feature: Standard Type

Page [236](#)

N₂ Generator



Name: N₂ Generator
Model: NiGen LCMS 40-1
Feature: High Capacity

Page [237](#)



1

Sterilizers

Contents

| | |
|----------------------------|----|
| Lab Autoclaves | 22 |
| Dry Heat Sterilizers | 34 |
| Loop Cinerator | 36 |

Lab Autoclaves | Standard

SN210C/310C/510C, SQ510C/810C

Operating temp. range 45~135°C

Max. working pressure 0.255MPa

Internal capacity 20L 32L 47L 50L 80L

High-performance, easy operation and ergonomically designed autoclave, standard equipped with cooling fan to shorten cool down time.




Features

- Maximum sterilizing temperature goes up to 135°C, suitable for protein modification.
- By simple setting and operation, these sterilizers manage ordinary sterilization, sterilization of culture mediums and liquids, and dissolution of culture mediums.
- Custom programs can be saved for repeated use.
- Standard equipped with cooling fan to shorten cool down time.
- The control panel and display are conveniently located on the front of the lid for easy viewing and operation.
- Equipped with multiple pressure lid opening protection locks and comprehensive safety measures.
- The SQ model features a low bed design, making it highly suitable for the retrieval of sterilized items. Compatible for overseas use with a wide power range of AC100-120V and AC200-240V.

Specifications

| Model | | SN210C | SN310C | SN510C | SQ510C | SQ810C |
|-------------------------------|--|--|---|---|---|--|
| System | | Automatic high-pressure steam sterilization | | | | |
| Performance | Operating temp. range | 45~135°C 45~80°C (preheat temp.) /45~60°C (retain temp.) /65~100°C (dissolution)/105~135°C (sterilizing) | | | | |
| | Temp. unit | 0.1°C | | | | |
| | Max. working pressure | 0.255MPa | | | | |
| Lid structure | | Manual upward opening and downward closing with an interlock for safety | | | | |
| Chamber | Material | 3mm thick corrosion-resistant stainless steel, featuring a mirror finish surface | | | | |
| | Design operational life | 20 years | | | | |
| | Design pressure | 0.42MPa | | | | |
| | Design temp. | 151°C | | | | |
| | Heater | 110V, 968W×2 | | 110V, 1150W×2 | 110V, 1210W×2 | 110V, 1500W×2 |
| | Exhaust valve | One quick exhaust valve and one slow release valve | | | | |
| | Options | Sample sensor interface (1/4), data logger sensor interface (1/4), pressure gauge interface (on solenoid valve piping) | | | | |
| Controllers | Cooling fan blade | Axial fan motor | | | | |
| | Temp. control, display/setting method | PID control by microprocessor, digital display/digital setting by ▲/▼ keys | | | | |
| | Timer | Timer: 0 or 1 min~99 h 59 min, timer resolution: 1 min | | | | |
| | Operation functions | Instrument sterilization, liquid sterilization, sterilization and retain temp., dissolution and retain temp., manual mode | | | | |
| Safety device | Other functions | Key locking, presetting, saving, preheating, forced cooling, sample temperature sensor (option), pattern locking, alarm log saving (20 pcs), display of accumulated working time/present time, ON-OFF beeping setting | | | | |
| | Sensor failure detection, SSR short-circuit detection, broken heater wire detection, prevention of air burning (liquid expansion type thermoestat), alarm against the absence of drain bottle, failure in locking the lid detection, memory error detection, pressure relief valve (0.255 MPa) | | | | | |
| Pressure vessel specification | | Class I pressure vessel | | | | |
| Specification | External dimensions (W×D×H mm) | 460×590×848 | 460×590×848 | 460×590×1058 | 520×660×881 | 520×660×1161 |
| | Effective dimensions of chamber (diameter×height mm) | 300×305 | 300×445 | 300×665 | 370×470 | 370×750 |
| | Internal effective capacity of chamber | 20L | 32L | 47L | 50L | 80L |
| | Weight | Approx. 65kg | Approx. 75kg | Approx. 85kg | Approx. 95kg | Approx. 105kg |
| | Power supply (50/60Hz) rated current | AC220V 9A | | AC220V 10.5A | AC220V 11A | AC220V 14A |
| Accessories | | Basket OSM-60(diameter 274 × height 125mm), 2 PCS | Basket OSM-70(diameter 274 × height 200mm), 2 PCS | Basket OSM-70(diameter 274 × height 200mm), 3 PCS | Basket OSQ-90(diameter 344 × height 200mm), 2 PCS | Basket OSQ-90H(diameter 344 × height 300mm), 2 PCS |
| Options | | 1 steam collection cup, 1 drain bottle, 1 heater baffle, 1 set of sterilization test cards (30 cards), 1 filter, 1m drainage pipe 1 hose clamp, long-arm clamp (SQ810C model only), 1 drip tray, user manual, warranty certificate | | | | |
| Options | | Biosafety filtration system, sample temperature sensor, internal chamber temperature measurement sensor, data logger, temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, standard basket, extended height basket | | | | |

CE
Certification

Choose a sterilization program

Instrument sterilize

Sterilization of lab instruments and vessel such as flask, beaker, test tube, scissors.

Fluid sterilize

Sterilization and disinfection of culture media, reagents, or other solutions, and keeps them warm after sterilization.

Sterilize & Retain temp.

Dissolving culture media and maintaining them at a specific temperature to prevent solidification.

Melt & Retain temp.

Sterilization of refining water, purified water and dilution water.

Custom program

Customized temperature and time settings.

Comply with GLP/GMP inspection

Standard equipped with 2 sensor ports on the main unit.



Optional data logger and high performance pressure gauge.

**Multiple safety devices**

- Steam collection cup
- Top lid's interlock mechanism
- Alarm against the absence of drain bottle
- Leakage breaker
- During any malfunction, the autoclave will automatically shutdown
- Alarm buzzer sounding and an error message being displayed
- Alarm recording (20 alarms)

**Easy discharge of sterilization water**



Easier maintenance with larger diameter drain pipe.

Standard cooling fan

- Cooling fan operates after sterilization.
- Cooling down to safe temperature.
- Reducing the time before samples can be taken out.
- Free to choose the forced cooling or natural cooling.

Comparison of cooling time with other devices

(Initial water temperature: 25°C, sterilization temperature: 121°C, ambient temperature: 23°C)

**Sterilization starts automatically according to sample temperature sensor**

The optional sample temperature sensor controls sterilization time based on the temperature of the sterilized material, ensuring effective sterilization.

**Drain bottle located at the front**

The drain bottle is set inside the front door for easy monitoring of the water level.

Space-saving top-opening lid design

Require less setup space compared to other machines.

Maximum set temperature of 135°C

The maximum set temperature of 135°C is suitable for protein modification.

Time-saving

Standard equipped with automatic start and preheat function for more efficient time using.

Universal power supply (selectable power)

All models offer a choice of AC100-120V and AC200-240V (single phase) power supply.

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Lab Autoclaves | Automatic Sterilization, Touch Screen Operation

SJ511C/811C

Operating temp. range 45~138℃

Max. working pressure 0.3MPa

Internal capacity 50L 80L

The steam sterilizer features a 7-inch color touchscreen for simple operation with advanced functions.



Features

- 7-inch large touchscreen, easy operation with options for Chinese, Japanese and English languages, along with a three-level authorization management system.
- Maximum sterilizing temperature goes up to 138℃, suitable for protein modification.
- By simple setting and operation, these sterilizers manage ordinary sterilization, sterilization of culture mediums and liquids, and dissolution of culture mediums.
- Custom programs can be saved for repeated use.
- Standard equipped with cooling fan to shorten cool down time.
- Equipped with multiple pressure lid opening protection locks and comprehensive safety measures.
- Display power consumption, accumulated power usage, CO₂ emissions, total power-on time, accumulated operating time, and accumulated times of operations.
- Two handles on each side for easy moving.
- Optional automatic water filling function available for purchase.

Specifications

| Model | | SJ511C | SJ811C |
|-------------------------------|--|--|---|
| System | | Automatic high-pressure steam sterilization | |
| Performance | Operating temp. range | 45~138℃: 45~80℃ (preheat temp.)/45~60℃ (retain temp.)/60~110℃ (dissolution)/105~138℃ (sterilizing) | |
| | Temp. unit | 0.1℃ | |
| | Max. working pressure | 0.3MPa | |
| Lid structure | | Manual upward opening and downward closing with an interlock for safety | |
| Chamber | Material | 3mm thick corrosion-resistant stainless steel, featuring a mirror finish surface | |
| | Design operational life | 20 years | |
| | Design pressure | 0.42MPa | |
| | Design temp. | 151℃ | |
| | Heater | 110V, 1500W×2 | 110V, 2000W×2 |
| | Exhaust valve | Solenoid valves: One for full exhaust and one for slow exhaust | |
| | Options | Sample sensor interface (1/4), data logger sensor interface (1/4), pressure gauge interface (on solenoid valve piping), automatic water intake (R1/4) | |
| | Cooling fan blade | Axial fan motor | |
| Controllers | Temp. control method | PID control by microprocessor | |
| | Temp. control setting/display method | 7-inch color LCD touchscreen for display and settings | |
| | Timer | 1 min~99 h 59 min, resolution: 1 min or 100 h~999 h, resolution: 1 h | |
| | Operation functions | Instrument sterilization, liquid sterilization, sterilization and retain temp., dissolution and retain temp. | |
| | Other functions | Power consumption, accumulated power usage, CO ₂ emissions, accumulated power-on time, accumulated operating time, and accumulated times of operations, key locking, presetting, saving, preheating, forced cooling, pattern locking, alarm log saving (20 pcs), display of present time, ON-OFF beeping setting, language switch function (Chinese/Japanese/English) and three-level authorization management system, sample sensor (optional) | |
| Safety device | | Sensor failure detection detection, SSR short-circuit detection, broken heater wire detection, prevention of air burning (water level switch, liquid expansion type), alarm against the absence of drain bottle, failure in locking the lid detection, memory error detection, pressure relief valve (0.3 MPa) | |
| Pressure vessel specification | | Class I pressure vessel | |
| Specification | External dimensions (W×D×H mm) | 520×660×881 | 520×660×1161 |
| | Effective dimensions of chamber (diameter×height mm) | 370×465 | 370×745 |
| | Internal effective capacity of chamber | 50L | 80L |
| | Weight | Approx. 95kg | Approx. 115kg |
| | Power supply rated current | AC220V 14A | AC220 19A |
| Accessories | | Basket OSQ-90 (diameter 344 × height 200mm), 2 PCS 1 steam collection cup, 1 cooling water bottle, 1 heater baffle, 1 set of sterilization test cards (30 cards), 1 filter plug, 1m drainage pipe, 1 hose clamp, long-arm clamp (SJ811C model only) | Basket OSQ-90H (diameter 344 × height 300mm), 2 PCS |
| Options | | Biosafety filtration system, sample temperature sensor, internal chamber temperature measurement sensor, miniprinter, data logger, temperature output terminal (4~20mA), external alarm output terminal, time arrival output terminal, standard basket, extended height basket, automatic water filling system | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Low bed mode

- Convenient retrieval of sterilized items.




Standard Cooling Fan

- Cooling fan operates after sterilization.
- Cooling down to safe temperature.
- Reducing the time before samples can be taken out.
- Free to choose the forced cooling or natural cooling.

Comparison of cooling time with other devices

(Initial water temperature: 25°C, sterilization temperature: 121°C, ambient temperature: 23°C)




Control panel

- 7-inch color LCD touchscreen with operation in Chinese, English and Japanese



Main screen

Sterilization parameter settings
(instrument sterilization program)Operation mode selection (instrument
sterilization program)Operation mode selection (melt and
retain temp. program)Operating status monitoring
(instrument sterilization program)Sterilization starts
automatically according
to sample temperature
sensor

The optional sample temperature sensor controls sterilization time based on the temperature of the sterilized material, ensuring effective sterilization.

Comply with GLP/GMP inspection requirements

Standard equipped with 2 sensor ports on
the main unit.Data logger or built-in micro printer is
optional

Safety valve located on the side



Convenient disassembly for safety valve replacement or calibration.

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Lab Autoclaves | Automatic Sterilization, Touch Screen Operation

SJ530/830

Operating temp. range 45~135°C

Max. working pressure 0.255MPa

Internal capacity 50L 80L

The autoclave features a 7-inch color touchscreen for simple operation with advanced functions.



Features

- 7-inch large touchscreen, easy operation with options for Chinese, Japanese and English languages, along with a three-level authorization management system.
- Maximum sterilizing temperature goes up to 135°C, suitable for protein modification.
- By simple setting and operation, these sterilizers manage ordinary sterilization, sterilization of culture mediums and liquids, and dissolution of culture mediums.
- Custom programs can be saved for repeated use.
- Standard equipped with cooling fan to shorten cool down time.
- Equipped with multiple pressure lid opening protection locks and comprehensive safety measures.
- Display power consumption, accumulated power usage, CO₂ emissions, total power-on time, accumulated operating time, and accumulated times of operations.
- SJ530/830 has obtained medical device registration and disinfection certificates.

Specifications

| Model | | SJ530 | SJ830 |
|---|--|--|---------------|
| System | | Automatic high-pressure steam sterilization | |
| Medical device registration certificate | | Yes | |
| Performance | Operating temp. range | 45~135°C: 45~80°C (preheat temp.)/45~60°C (retain temp.)/60~110°C (dissolution)/105~135°C (sterilizing) | |
| | Temp. unit | 0.1°C | |
| | Max. working pressure | 0.255MPa | |
| Lid structure | | Manual upward opening and downward closing with an interlock for safety | |
| Chamber | Material | 3mm thick corrosion-resistant stainless steel, featuring a mirror finish surface | |
| | Design operational life | 20 years | |
| | Design pressure | 0.42MPa | |
| | Design temp. | 151°C | |
| | Sterilization heater | 110V, 1210W×2 | 110V, 1500W×2 |
| | Exhaust valve | One quick exhaust valve and one slow release valve | |
| | Options | Sample sensor interface (1/4), data logger sensor interface (1/4), pressure gauge interface (on solenoid valve piping) | |
| Controllers | Cooling fan blade | Axial fan motor | |
| | Temp. control method | PID control by microprocessor | |
| | Temp. control setting/display method | 7-inch color LCD touchscreen for display and settings | |
| | Timer | 1 min~99 h 59 min, resolution: 1 min or 100~999 h, resolution: 1 h | |
| | Operation functions | Instrument sterilization, liquid sterilization, sterilization and retain temp., dissolution and retain temp. | |
| Specification | Other functions | Power consumption, accumulated power usage, CO ₂ emissions, accumulated power-on time, accumulated operating time, and accumulated times of operations, key locking, presetting, saving, preheating, forced cooling, sample temperature sensor (optional), pattern locking, error log saving (20 pcs), display of present time, ON-OFF beeping setting, language switch function (Chinese/Japanese/English) and three-level authorization management system | |
| | Safety device | Sensor failure detection detection, SSR short-circuit detection, broken heater wire detection, prevention of air burning (liquid expansion type thermoetat), alarm against the absence of drain bottle, failure in locking the lid detection, memory error detection, pressure relief valve (0.255 MPa) | |
| | Pressure vessel specification | Class I pressure vessel | |
| | External dimensions (W×D×H mm) | 520×660×881 | 520×660×1161 |
| | Effective dimensions of chamber (diameter×height mm) | 370×470 | 370×750 |
| Accessories | Internal effective capacity of chamber | 50L | 80L |
| | Weight | Approx. 105kg | Approx. 125kg |
| | Power supply/rated current | AC220V 11.5A | AC220V 14A |
| Options | | Basket OSQ-90 (diameter 344 × height 200mm), 2 PCS Basket OSQ-90H (diameter 344 × height 300mm), 2 PCS 1 steam collection cup, 1 cooling water bottle, 1 heater baffle, 1 set of sterilization test cards (30 cards), 1 filter plug, 1m drainage pipe, 1 hose clamp, long-arm clamp (SJ830 model only) | |
| Options | | Biosafety filtration system, sample temperature sensor, internal chamber temperature measurement sensor, miniprinter, data logger, temperature output terminal (4~20mA), external alarm output terminal, time arrival output terminal, standard basket, extended height basket | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Medical Device
Registration
CertificateDisinfection
Certificate

Low bed mode

- Convenient retrieval of sterilized items.



Control panel

- 7-inch color LCD touchscreen with operation in Chinese, English, and Japanese.



Main screen



Sterilization parameter settings
(instrument sterilization program)Operation mode selection (instrument
sterilization program)Operation mode selection (melt and
retain temp. program)Operating status monitoring
(instrument sterilization program)

Standard Cooling Fan

- Cooling fan operates after sterilization.
- Cooling down to safe temperature.
- Reducing the time before samples can be taken out.
- Free to choose the forced cooling or natural cooling.

Comparison of cooling time with other devices

(Initial water temperature: 25°C, Sterilization temperature: 121°C, Ambient temperature: 23°C)



Sterilization starts automatically according to sample temperature sensor

The optional sample temperature sensor controls sterilization time based on the temperature of the sterilized material, ensuring effective sterilization.

Comply with GLP/GMP inspection requirements



Standard equipped with 2 sensor ports on the main unit.



Optional data logger and high performance pressure gauge.

Space-saving top-opening lid design

Require less setup space compared to other machines.

Maximum sterilization temperature of 135°C

Maximum sterilizing temperature goes up to 135°C suitable for protein modification.

Time-saving

Standard equipped with automatic start and preheat function for more efficient time using.

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Lab Autoclaves | Automatic Sterilization, Drying Function

SM510C/810C/530/830

Operating temp. range Sterilization: 105~135°C Drying: 135~150°C

Max. working pressure 0.255MPa

Internal capacity 50L 80L

High-performance sterilizer with fully automatic operation for sterilization and drying.



Features

- 7-inch large touchscreen, easy operation with options for Chinese, Japanese and English languages, along with a three-level authorization management system.
- Maximum sterilizing temperature goes up to 135°C, suitable for protein modification.
- By simple setting and operation, these sterilizers manage ordinary sterilization, sterilization of culture mediums and liquids, and dissolution of culture mediums.
- Custom programs can be saved for repeated use.
- Standard equipped with cooling fan to shorten cool down time.
- Equipped with multiple pressure lid opening protection locks and comprehensive safety measures.
- Display power consumption, accumulated power usage, CO₂ emissions, total power-on time, accumulated operating time, and accumulated times of operations.
- SM530/830 has obtained medical device registration and disinfection certificates.

Specifications

| Model | | SM510C | SM810C | SM530 | SM830 |
|---|--|---|---|--|---|
| System | | Automatic high-pressure steam sterilization+drying | | | |
| Medical device registration certificate | | None | | Yes | |
| Performance | Operating temp. range | 45~135℃: 45~80℃ (preheat temp.)/45~60℃ (retain temp.)/60~110℃ (dissolution)/105~135℃ (sterilizing)/135~150℃ (drying) | | | |
| | Temp. unit | 0.1℃ | | | |
| | Max. working pressure | 0.255MPa | | | |
| Lid structure | | Manual upward opening and downward closing with an interlock for safety | | | |
| Chamber | Material | 3mm thick corrosion-resistant stainless steel, featuring a mirror finish surface | | | |
| | Design operational life | 20 years | | | |
| | Design pressure | 0.42MPa | | | |
| | Design temp. | 151℃ | | | |
| | Sterilization heater | 110V, 1210W×2 | | | |
| | Dry heat heater | 110V, 295W×2+455W×2 | 110V, 275W×2+625W×2 | 110V, 295W×2+455W×2 | 110V, 275W×2+625W×2 |
| | Drain valve | 1 solenoid valve, 1 hand valve | | | |
| | Exhaust valve | 1 fast exhaust valve, 1 auxiliary exhaust valve, 1 slow exhaust valve | | | |
| | Options | Sample sensor interface (1/4), data logger sensor interface (1/4), pressure gauge interface (on the auxiliary exhaust solenoid valve piping) | | | |
| | Cooling fan blade | Axial fan motor | | | |
| Controllers | Temp. control method | PID control by microprocessor | | | |
| | Temp. control setting/display method | 7-inch color LCD touchscreen for display and settings | | | |
| | Timer | 1 min~99 h 59 min, resolution: 1 min or 100 h~999 h, resolution: 1 h | | | |
| | Operation functions | Instrument sterilization, liquid sterilization, sterilization and retain temp., dissolution and retain temp., sterilization drying, drying | | | |
| | Additional functions | Power consumption, accumulated power usage, CO ₂ emissions, accumulated power-on time, accumulated operating time, and accumulated times of operations | | | |
| | Other functions | Key locking, presetting, saving, preheating, forced cooling, sample temperature sensor (option), pattern locking, alarm log saving (20 pcs), display of present time, ON-OFF beeping setting, language switch function (Chinese/Japanese/English) and three-level authorization management system | | | |
| Safety device | | Sensor failure detection, SSR short-circuit detection, broken heater wire detection, prevention of air burning (liquid expansion type thermoet), alarm against the absence of drain bottle, failure in locking the lid detection, memory error detection, pressure relief valve (0.255 MPa), overpressure alarm | | | |
| Pressure vessel specification | | Class I pressure vessel | | | |
| Specification | External dimensions (W×D×H mm) | 520×660×881 | 520×660×1161 | 520×660×881 | 520×660×1161 |
| | Effective dimensions of chamber (diameter×height mm) | 370×470 | 370×750 | 370×470 | 370×750 |
| | Internal effective capacity of chamber | 50L | 80L | 50L | 80L |
| | Weight | Approx. 113kg | Approx. 137kg | Approx. 113kg | Approx. 137kg |
| | Power supply/rated current | AC220V 13A | | | |
| Accessories | | Basket OSQ-90 (diameter 344 x height 200mm), 2 PCS | Basket OSQ-90H (diameter 344 x height 300mm), 2 PCS | Basket OSQ-90 (diameter 344 x height 200mm), 2 PCS | Basket OSQ-90H (diameter 344 x height 300mm), 2 PCS |
| | | 1 drain bottle, 1 heater baffle, 1 set of sterilization test cards (30 cards), 1 filter plug, 1m drainage pipe, 1 hose clamp, long-arm clamp (SM810C/SM830 model only) | | | |
| Options | | Biosafety filtration system, sample temperature sensor, internal chamber temperature measurement sensor, miniprinter, data logger, temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, standard basket, extended height basket | | | |

Disinfection
CertificateMedical Device
Registration
CertificateCE
Certification

Low bed mode

- Convenient retrieval of sterilized items.



Control panel

- 7-inch color LCD touchscreen with operation in Chinese, English and Japanese



Main screen

Operation mode selection
(Instrument sterilization program)Instrument sterilization program
operationOperation mode selection (Sterilization
& drying program)



Sterilization & drying program operation

Standard Cooling Fan

- Cooling fan operates after sterilization.
- Cooling down to safe temperature.
- Reducing the time before samples can be taken out.
- Free to choose the forced cooling or natural cooling.

Comparison of cooling time with other devices

(Initial water temperature: 25°C, Sterilization temperature: 121°C, Ambient temperature: 23°C)



Sterilization starts automatically according to sample temperature sensor

The optional sample temperature sensor controls sterilization time based on the temperature of the sterilized material, ensuring effective sterilization.

Comply with GLP/GMP inspection requirements

Standard equipped with 2 sensor ports on
the main unit.Optional data logger and high performance
pressure gauge.

Space-saving top-opening lid design

Require less setup space compared to other machines.

Maximum sterilization temperature of 135°C

Maximum sterilizing temperature goes up to 135°C suitable for protein modification.

Time-saving

Standard equipped with automatic start and preheat function for more efficient time using.

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Lab Autoclaves | Automatic Sterilization, Large Capacity

SQL810C/1010C

Operating temp. range 45~135℃

Max. working pressure 0.255MPa

Internal capacity 85L 110L

The autoclave features a 110L capacity for easy operation with advanced functions.



Features

- With 85L and 110L of internal capacity, it can accommodate more items for sterilization in a single cycle.
- Maximum sterilizing temperature goes up to 135℃, suitable for protein modification.
- By simple setting and operation, these sterilizers manage ordinary sterilization, sterilization of culture mediums and liquids, and dissolution of culture mediums.
- Custom programs can be saved for repeated use.
- Standard equipped with cooling fan to shorten cool down time.
- The control panel and display are conveniently located on the front of the lid for easy viewing and operation.
- Equipped with multiple pressure lid opening protection locks and comprehensive safety measures.

Specifications

| Model | | SQL810C | SQL1010C |
|-------------------------------|--|---|---|
| System | | Automatic high-pressure steam sterilization | |
| Performance | Operating temp. range | 45~135℃ 45~80℃ (preheat temp.)/45~60℃ (retain temp.)/65~100℃ (dissolution)/105~135℃ (sterilizing) | |
| | Temp. unit | 0.1℃ | |
| | Max. working pressure | 0.255MPa | |
| Lid structure | | Manual upward opening and downward closing with an interlock for safety | |
| Chamber | Material | 3mm thick corrosion-resistant stainless steel, featuring a mirror finish surface | |
| | Design operational life | 20 years | |
| | Design pressure | 0.42MPa | |
| | Design temp. | 151℃ | |
| | Heater | 220V, 2000W×2 | |
| | Exhaust valve | One quick exhaust valve and one slow release valve | |
| | Options | Sample sensor interface (R1/4), data logger sensor interface (G1), pressure interface (G1/2) | |
| Controllers | Cooling fan blade | Axial fan motor | |
| | Temp. control, display/setting method | PID control by microprocessor, digital display/digital setting by ▲/▼ keys | |
| | Timer | Timer: 0 or 1 min~99 h 59 min, resolution: 1 min | |
| | Operation functions | Instrument sterilization, liquid sterilization, sterilization and retain temp., dissolution and retain temp., manual mode | |
| Safety device | Other functions | Key locking, presetting, saving, preheating, forced cooling, sample temperature sensor (option), pattern locking, alarm log saving (20 pcs), display of accumulated working time/present time, ON-OFF beeping setting | |
| | | Sensor failure detection, SSR short-circuit detection, broken heater wire detection, prevention of air burning (liquid expansion type thermoet), alarm against the absence of drain bottle, failure in locking the lid detection, memory error detection, pressure relief valve (0.255 MPa) | |
| Pressure vessel specification | | Class I pressure vessel | |
| Specification | External dimensions (W×D×H mm) | 680×760×997 | 680×760×1154 |
| | Effective dimensions of chamber (diameter×height mm) | 450×535 | 450×692 |
| | Internal effective capacity of chamber | 85L | 110L |
| | Weight | Approx. 145kg | Approx. 170kg |
| | Power supply (50/60Hz) rated current | AC220V 19A | |
| Accessories | | Basket Diameter 424× Height 200, 2 pcs | Basket Diameter 424× Height 300, 2 pcs |
| | | 1 steam collection cup, 1 cooling water bottle, 1 heater baffle, 1 set of sterilization test cards (30 cards), 1 filter, 1m drainage pipe, 1 hose clamp, long-arm clamp (SQL1010C model only) | |
| Options | | Biosafety filtration system, sample temperature sensor, internal chamber temperature measurement sensor, data logger, temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, standard basket, extended height basket | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators


11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

**Instrument sterilize**

Sterilization of lab instruments and vessel such as flask, beaker, test tube, scissors.

Fluid sterilize

Sterilization and disinfection of culture media, reagents, or other solutions, and keeps them warm after sterilization.

Sterilize & Retain temp.

Dissolving culture media and maintaining them at a specific temperature to prevent solidification.

Melt & Retain temp.

Sterilization of refining water, purified water and dilution water.

Custom program

Customized temperature and time settings.

Comply with GLP/GMP inspection

Standard equipped with 2 sensor ports on the main unit.



Optional data logger and high performance pressure gauge.

**Multiple safety devices**

- Steam collection cup
- Top lid's interlock mechanism
- Alarm against the absence of drain bottle
- Leakage breaker
- During any malfunction, the autoclave will automatically shutdown
- Alarm buzzer sounding and an error message being displayed
- Alarm recording (20 alarms)

**Easy discharge of sterilization water**



Easier maintenance with larger diameter drain pipe.

Standard Cooling Fan

- Cooling fan operates after sterilization.
- Cooling down to safe temperature.
- Reducing the time before samples can be taken out.
- Free to choose the forced cooling or natural cooling.

Comparison of cooling time with other devices

(Initial water temperature: 25°C, sterilization temperature: 121°C, ambient temperature: 23°C)

**Sterilization starts automatically according to sample temperature sensor**

The optional sample temperature sensor controls sterilization time based on the temperature of the sterilized material, ensuring effective sterilization.

**Drain bottle located at the front**

The drain bottle is located in front for easy access and drain water level can be monitored without opening the cabinet door.

Space-saving top-opening lid design

Require less setup space compared to other machines.

Maximum set temperature of 135°C

The maximum set temperature of 135°C is suitable for protein modification.

Time-saving

Standard equipped with automatic start and preheat function for more efficient time using.

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Lab Autoclaves | Automatic Sterilization, Large Capacity

SJL810C/1010C

Operating temp. range 45~135℃

Max. working pressure 0.255MPa

Internal capacity 85L 110L

The steam sterilizer features a 110L large capacity, with touchscreen for easy operation with advanced functions.



Features

- 7-inch large touchscreen, easy operation with options for Chinese, Japanese and English languages, along with a three-level authorization management system.
- Maximum sterilizing temperature goes up to 135℃, suitable for protein modification.
- By simple setting and operation, these sterilizers manage ordinary sterilization, sterilization of culture mediums and liquids, and dissolution of culture mediums.
- Custom programs can be saved for repeated use.
- Standard equipped with cooling fan to shorten cool down time.
- Equipped with multiple pressure lid opening protection locks and comprehensive safety measures.
- Display power consumption, accumulated power usage, CO₂ emissions, accumulated power-on time, accumulated operating time, and accumulated times of operations.

Specifications

| Model | | SJL810C | SJL1010C |
|-------------------------------|--|--|--|
| System | | Automatic high-pressure steam sterilization | |
| Performance | Operating temp. range | 45~135℃: 45~80℃ (preheat temp.)/45~60℃ (retain temp.)/60~110℃ (dissolution)/105~135℃ (sterilizing) | |
| | Temp. unit | 0.1℃ | |
| | Max. working pressure | 0.255MPa | |
| Lid structure | | Manual upward opening and downward closing with an interlock for safety | |
| Chamber | Material | 3mm thick corrosion-resistant stainless steel, featuring a mirror finish surface | |
| | Design operational life | 20 years | |
| | Design pressure | 0.42MPa | |
| | Design temp. | 151℃ | |
| | Sterilization heater | 220V, 2000W×2 | |
| | Exhaust valve | One quick exhaust valve and one slow release valve | |
| | Options | Sample sensor interface (Rc1/4), temperature test interface (G1), pressure test interface (G1/2) | |
| Controllers | Cooling fan blade | Axial fan motor | |
| | Temp. control method | PID control by microprocessor | |
| | Temp. control setting/display method | 7-inch color LCD touchscreen for display and settings | |
| | Timer | 1 min - 99 h 59 min (resolution: 1 min)/100 h - 999 h (resolution: 1 h) | |
| | Operation functions | Instrument sterilization, liquid sterilization, sterilization and retain temp., dissolution and retain temp. | |
| | Additional functions | Power consumption, accumulated power usage, CO ₂ emissions, accumulated power-on time, accumulated operating time, and accumulated times of operations | |
| Safety device | | Sensor failure detection, SSR short-circuit detection, broken heater wire detection, prevention of air burning (liquid expansion type thermoetatt), alarm against the absence of drain bottle, failure in locking the lid detection, memory error detection, pressure relief valve (0.255 MPa) | |
| Pressure vessel specification | | Class I pressure vessel | |
| Specification | External dimensions (W×D×H mm) | 680×760×997 | 680×760×1154 |
| | Effective dimensions of chamber (diameter×height mm) | 450×535 | 450×692 |
| | Internal effective capacity of chamber | 85L | 110L |
| | Weight | Approx. 145kg | Approx. 170kg |
| | Power supply (50/60Hz) rated current | AC220V 18.5A | |
| Accessories | | Basket OSQL-90 (diameter 424 × depth 200mm), 2 PCS 1 steam collection cup, 1 cooling water bottle, 1 heater baffle, 1 set of sterilization test cards (30 cards), 1 filter plug, 1m drainage pipe, 1 hose clamp, long-arm clamp (SJL1010C model only) | Basket OSQL-90H (diameter 424 × depth 300), 2 PCS |
| Options | | Biosafety filtration system, sample temperature sensor, internal chamber temperature measurement sensor, miniprinter, data logger, temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, standard basket, extended height basket | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Low bed mode

- Convenient retrieval of sterilized items.



Control panel

- 7-inch color LCD touchscreen with operation in Chinese, English and Japanese



Main screen



Sterilization parameter



Operation mode selection (instrument sterilization program)



Operation mode selection (melt and retain temp. program)





Sterilization and temp. retain program

Standard Cooling Fan

- Cooling fan operates after sterilization.
- Cooling down to safe temperature.
- Reducing the time before samples can be taken out.
- Free to choose the forced cooling or natural cooling.

Comparison of cooling time with other devices

(Initial water temperature: 25°C, Sterilization temperature: 121°C, Ambient temperature: 23°C)



Sterilization starts automatically according to sample temperature sensor

The optional sample temperature sensor controls sterilization time based on the temperature of the sterilized material, ensuring effective sterilization.

Comply with GLP/GMP inspection requirements



Standard equipped with 2 sensor ports on the main unit.



Optional data logger and high performance pressure gauge.

Space-saving top-opening lid design

Require less setup space compared to other machines.

Maximum sterilization temperature of 135°C

Maximum sterilizing temperature goes up to 135°C suitable for protein modification.

Time-saving

Standard equipped with automatic start and preheat function for more efficient time using.

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Dry Heat Sterilizers | Dry Heat Sterilization, Standard

SIT410C/610C/810C

Operating temp. range RT+10~300°C

Temp. uniformity 5°C (160~200°C) 9°C (201~300°C)

Internal capacity 90L 150L 300L

Forced convection dry heat sterilizer with program function, easy to set.



Features

- Quick, safe and reliable sterilization function.
- Quick and accurate sterilization with HMI interface and preset sterilization programs.
- Temperature monitoring function monitors the temperature during the sterilization process and automatically stops timing and triggers an alarm if the chamber does not reach the sterilization temperature range.
- The sterilizer door is equipped with a safety interlock mechanism.
- Dual use can also be used as a forced convection drying oven.

Safety

- Equipped with self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection switch, three-level authorization management system and other safety functions.

Specifications

| Model | | | SIT410C | SIT610C | SIT810C |
|---------------|--------------------------------------|---|--|---------------|---------------|
| System | | | Forced convection hot air sterilization | | |
| Performance | Operating temp. range | | Room temp. +10~300℃ | | |
| | GB standard | Sterilization temp. range | ±3℃ (160~200℃); ±5℃ (201~300℃) | | |
| | | Sterilization temp. uniformity | 5℃ (160~200℃); 9℃ (201~300℃) | | |
| | Max. temp. reaching time | | Approx. 55 min | | |
| Composition | Interior material | | Stainless steel plate | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | |
| | Insulating material | | Rock wool+centrifugal glass wool | | |
| | Heater | | Stainless steel heating pipe | | |
| | | | 1.6KW | 1.3KW ×2 | 1.8KW×2 |
| | Fan blade/motor | | Centrifugal fan blade, 30W motor×1 | | |
| Exhaust vent | | I.D. 34 mm×2, located at the back | | | |
| Controllers | Temp. control method | | 3-stage PID | | |
| | Temp. setting method | | 7-inch color touch screen setting | | |
| | Temp. display method | | 7-inch color touch screen display | | |
| | Timer | | 1 min~99 h 59 min | | |
| | Operation functions | | Dry heat sterilization operation: Sterilization cycle operation | | |
| | | | Forced air drying operation: fixed temp. operation, program operation, timing operation | | |
| | Additional functions | | Temperature history curve, alarm history records, operation log, data export function, language switching (Chinese, Japanese, English), calibration function, power outage compensation function, screen brightness adjustment, three-level authorization function, micro printer (optional) | | |
| Sensors | | Temperature controller: K-type thermocouple overheating protection: Liquid expansion temperature controller | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating prevention, SSR short-circuit), independent overheating prevention, overcurrent leakage protection switch | | |
| Specification | Internal dimensions (W×D×H mm) | | 450×450×450 | 600×500×500 | 600×500×1000 |
| | External dimensions (W×D×H mm) | | 610×690×960 | 760×740×1010 | 760×740×1780 |
| | Internal capacity | | 90L | 150L | 300L |
| | Shelf load | | 30kg/layer | | |
| | Shelf layers | | 11 layers | 13 layers | 29 layers |
| | Support spacing | | 30mm | | |
| | Power supply (50/60Hz) rated current | | AC220V 7.5A | AC220V 12A | AC220V 16.5A |
| | Weight | | Approx. 80kg | Approx. 100kg | Approx. 140kg |
| Accessories | Shelf | | Stainless steel wire mesh plate | | |
| | | | 2 pcs | | 4 pcs |
| | Supports | 4 pcs | | 8 pcs | |
| Options | Stand | | ONK60C | | — |
| | Stacking fittings | | ODK40C | — | |
| | Others | | Shelf (1 shelf with 2 supports) micro printer, data logger, combination warning light (standby/running/fault), remote monitoring function (network port), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal | | |

Dry Heat Sterilizers | Dry Heat Sterilization, Natural Convection

SI411C/611C

Max. operating temp. 260°C

Temp. distribution accuracy $\pm 5^{\circ}\text{C}$ (at 260°C)

Internal capacity 77L 159L

Gravity convection dry heat sterilizer with program function, easy to set.

Features

- Quick, safe and reliable sterilization function.
- Easy operation, available for fixed temp. operation, program operation, quick auto stop, auto stop and auto start operations.
- Use special function menu keys and up/down keys to realize digital setting. with repeat function of 6 modes in total of 90 segments program controller.
- Use submenu keys to operate overheat protector, deviation correction and key locking.

Safety

- Equipped with self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection switch, key lock and other safety functions.



Specifications

| Model | | | SI411C | SI611C |
|--------------------------|--------------------------------------|--|--|--------------|
| System | | | Natural convection hot air sterilizer | |
| Performance | Operating temp. range | | Room temp. +5~260℃ | |
| | GB standard | Temp. fluctuation | ±0.5℃ (at 260℃) | |
| | | Temp. uniformity | ±3% (at 260℃) | |
| | JTM standard | Temp. adjusting accuracy | ±1.0℃ (at 260℃) | |
| | | Temp. distribution accuracy | ±5.0℃ (at 260℃) | |
| Max. temp. reaching time | | Approx. 70 min | | |
| Composition | Interior material | | Stainless steel plate | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | |
| | Insulating material | | Glass fiber | |
| | Heater | | Nichrome heating wire | |
| | | | 1.2KW | 1.36KW |
| | Fan blade/motor | | — | |
| Exhaust vent | | I.D. 30mm×2 (top) | | |
| Controllers | Temp. control method | | 3-stage PID | |
| | Temp. setting method | | Use special function menu keys and up/down keys to realize digital setting | |
| | Temp. display method | | Measured temp. display: Green 4-digit LED digital display | |
| | | | Setting temp. display: Red 4-digit LED digital display | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | |
| | Operation functions | | Fixed temp. operation, program operation, auto start, quick auto stop, program operation | |
| | Program mode | | Program operation 6 modes in total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | |
| | Additional functions | | Deviation correction, key lock, power outage compensation | |
| Sensors | | K thermocouple (Temp. controller and overheat protector) | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating prevention, SSR short-circuit), overheating prevention, overcurrent leakage protection switch, key lock functions | |
| Specification | Internal dimensions (W×D×H mm) | | 450×430×400 | 600×530×500 |
| | External dimensions (W×D×H mm) | | 550×540×777 | 700×640×877 |
| | Internal capacity | | 77L | 159L |
| | Shelf load | | 15kg/layer | |
| | Shelf layers | | 10 layers | 13 layers |
| | Support spacing | | 30mm | |
| | Power supply (50/60Hz) rated current | | AC220V 6A | AC220V 6.5A |
| | Weight | | Approx. 42kg | Approx. 59kg |
| Accessories | Shelf | | Stainless punching mesh plate | |
| | | | 2 pcs | |
| | Supports | | 4 pcs | |
| Options | Stand | | ONS60C | |
| | Stacking fittings | | ODK82C | ODK84C |
| | Others | | Shelf (1 shelf with 2 supports) micro printer, data logger, combination warning light (standby/running/fault), external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software | |

Loop Cinerator | Loop Sterilization

SL-21

Made in
Japan

Temp. control range 800~850℃

Caliber 15mm

It replaces the traditional sterilization method using gas burners, achieving a safer and more effective sterilization result.



■ Features

The sterilization operation is conducted inside the combustion tube, preventing the spread of bacteria and microorganisms, making it safer and more convenient to use.

- Sterilize platinum inoculation loops and needles inside the combustion tube. Unlike gas flame sterilization, this method uses a heater, allowing for safe operation.
- The heat output is only 1/6 of that of gas flame sterilization, making it suitable for use inside biological safety cabinets or fume hoods.
- Use a highly heat-resistant and evenly heating quartz glass tube. Particularly efficient for high-volume continuous sterilization.

■ Angle adjustment



Whether operating while standing or sitting, it features a convenient structure that allows for the flexible adjustment of the heating port angle based on the operator's position.

■ Stand placement



The platinum inoculation loops, still carrying residual heat after sterilization, it can be placed on the stand. This is especially convenient during high-volume continuous operations, as it prevents contact with cooling objects.

■ Specifications

| Model | SL-21 |
|-----------------------------------|---------------------|
| Internal sterilization temp. | 800~850℃ |
| Sterilization temp. reaching time | 10mins |
| Caliber | 15mm |
| Heater | 130W |
| Power supply | AC100V 50/60Hz 1.5A |
| External dimensions(mm) | W95×D170×H200mm |
| Weight | 1.5kg |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Granulation and Spray Dryers

Contents

| | |
|--------------------------------------|----|
| Spray Dryers | 38 |
| Organic Solvent Recovery Units | 46 |
| Solvent Washing Units | 48 |

Spray Dryers | Standard

ADL312SC-A

Evaporated water capacity 1500ml/h

Temp. adjustment range 40~240℃

Liquid sending pump flow rate range 0~26ml/min

Nozzle options For liquids-air

New intelligent spray dryer with multiple functions and easy operation.




Features

- Adopting the method of instantaneously heating particle samples even if the samples with poor thermal stability are not easily oxidized, and get the uniform fine powder.
- The micro-powder after spraying has low moisture content, will not be oxidized, and no pollution.
- Products are directly dried into fine powder from solution or suspension samples without the need for filtration, separation, grinding, and other pre-processing and post-processing operations associated with traditional drying methods, avoiding contamination.
- It can be connected to the GAS411C organic solvent recovery device to handle samples containing organic solvents.
- The drying chamber and cyclone use a quick assembly and disassembly structure, further improving operability. A standard lift platform is equipped for installing and removing accessories.
- The device is equipped with a power socket for stirrers (200-230V~1A), convenient for stirring the suspension while sampling.
- Use a unique peristaltic sampling pump, nozzle cooling mechanism, pulse nozzle cleaning mechanism, anti-blockage needle, etc., to achieve diversity and stability in spray conditions.

Specifications

| Model | | ADL312SC-A |
|--------------------|------------------------------------|--|
| Applicable samples | | Water-soluble & organic solvent (can be used when connected to an organic solvent recovery device) |
| Performance | Moisture evaporation capacity | Max. 1500ml/h |
| | Temp. regulator setting range | 0~240℃ (inlet temperature), 0~100℃ (outlet temperature) |
| | Temp. adjusting accuracy | Inlet temperature $\pm 1^\circ\text{C}$ |
| | Drying air volume adjustment range | 0.2~0.9 m ³ /min |
| | Spray air flow adjustment range | 0~30L/min |
| | Spray pressure usage range | 0.3~0.6MPa |
| Composition | Nozzle cleaning function | Manual pulse air cleaning from the nozzle front end |
| | External output | Inlet temperature, outlet temperature output (4~20mA) |
| | Temp. regulator | Multi-PID control |
| | Touch screen | temperature adjustment, blower, heater, liquid feed pump, pulse air switch, automatic needle, alarm display, operation curve |
| | Control switch | selectable inlet temperature or outlet temperature control |
| | Temp. sensor | PT100 thermistor |
| | Heater | 3.2KW |
| | Liquid feed pump | Duct type liquid feed pump |
| | Spray air pump | Use spray air compressor (sold separately) or connect to an organic solvent recovery device (sold separately) when using built-in air compressor |
| | Service socket | For stirrers: 200-230V~1A |
| | Suction blower | Tubular blower |
| | Filters | Suction filter, exhaust filter |
| | Solvent recovery | Use a solvent recovery device (sold separately) |
| | Spray nozzle cooling structure | Joints $\times 2$, outer diameter $\Phi 10.5\text{mm}$ (optional cooling water circulation device CF312L-B) |
| | Spray air connection | Joint outer diameter, $\Phi 7\text{mm}$ |
| | Exhaust connection caliber | $\Phi 50\text{mm}$ |
| Safety functions | | Inlet, outlet temperature overheating, liquid feed pump reverse function, overcurrent leakage protection switch, nozzle connection abnormality |
| Specifications | External dimensions | W580×D420×H1355mm |
| | Weight | Approx. 96kg |
| | Power supply | 200-230V~50/60Hz 17-20A |
| Accessories | | Liquid feed hose $\times 2$, exhaust hose (with 1 hose clamp) $\times 1$, exhaust conversion joint, outlet temperature sensor, fuse, anti-static connection wire, intake hose 5m (with 2 hose clamps), stand components, protective cover (COV20), GF301C glass assembly, exhaust joint, air tube A, air tube B, product warranty, user manual |


Control panel




- Brand new 7-inch ultra-large touch screen control panel, available in Chinese/Japanese/English, easy and convenient to operate.
- High-power heaters greatly enhance the time to reach temperature and the temperature setting range is wider to meet more sample experiments.
- It can switch between two hot air circulation systems: hot air inhalation type and ejection type.
- Both two-fluid nozzles and three-fluid nozzles can be used (optional).
- A new automatic nozzle with a cooling mechanism.
- Remote control is possible.
- Experimental data recording and storage can be realized (optional).

System diagram

Hot air inhalation type

GAS411C connection port
(ADL312SC)

Hot air ejection type



Operation




Hot air ejection type connection pipeline



Relationship between blower output power and maximum inlet temperature (reference)

Air volume—maximum inlet temperature curve under 220V power supply conditions



Setting example



Spray dryer ADL312SC+ organic solvent recovery device GAS411C


Spray nozzle



The top part of the spray consists of a liquid nozzle and an air nozzle.

Liquid nozzle (F)

Air nozzle (A)



| Code | Nozzle size | | Nozzle specifications | Applicable models |
|--------|----------------------|--|-----------------------|--|
| 281298 | Liquid cap PF2050-SS | Aperture A=0.508mm Aperture B=1.270mm | 1 | ADL311(S)/GB210-A |
| | Air cap PA64-SS | Aperture C=1.626mm | | |
| 281290 | Liquid cap PF2050-SS | Aperture A=0.508mm Aperture B=1.270mm | 2A | ADL311(S)/GB210-A |
| | Air cap PA70-SS | Aperture C=1.778mm | | |
| 281292 | Liquid cap PF2850-SS | Aperture A=711μm Aperture B=1270μm | 3 | ADL312SC/GB211C/DL411C/DL410 (standard) ADL311(S)/GB210-A |
| | Air cap PA64-SS | Aperture C=1626μm | | |
| 281291 | Liquid cap PF2850-SS | Aperture A=711μm Aperture B=1270μm | 2 | ADL311(S)/ADL312SC/GB210-A/ GB211C/DL410/DL411C |
| | Air cap PA-70-SS | Aperture C=1778μm | | |

Spray Dryers | Multi-functional

GB211C-A

Evaporated water capacity 1500ml/h

Temp. adjustment range 40~240℃

Liquid sending pump flow rate range 0~26ml/min

Nozzle options For liquids-air

New intelligent spray dryer with multiple functions and easy operation.



Features

- Instantaneous heating of granular samples ensures that even thermally stable samples are not oxidized and uniform fine powder can be obtained.
- The micro-powder after spraying has low moisture content, will not be oxidized, and no pollution.
- Since the sample is directly dried from the solvent or suspension into micro-powder, it eliminates the pre-processing and post-processing operations such as filtration, separation, and grinding, and avoids contamination during these operations.
- By connecting the organic solvent recovery device, samples containing organic solvents can be professionally and safely handled.
- After installing the optional granulation component GF200, it can be used as a fluidized bed dryer and granulator.
- Equipped with an electric lift, it is convenient for installing and removing attachments.
- The machine is equipped with power socket (1A) for the stirrer, convenient for stirring the suspension while sampling.
- Use a unique peristaltic sampling pump, nozzle cooling mechanism, pulse nozzle cleaning mechanism, anti-blockage needle, etc., to achieve diversity and stability in spray conditions.
- Brand new 7-inch ultra-large touch screen control panel, available in Chinese/Japanese/English, easy and convenient to operate.
- Both two-fluid nozzles and three-fluid nozzles can be used.
- Experimental data recording and storage can be realized (optional function).
- Remote control is possible (optional function).
- Temperature zoning control is adopted, heating up faster and more stable.
- High-power heaters greatly enhance the time to reach temperature and the temperature setting range is wider, meeting more sample experiments. and the temperature setting range is wider to meet more sample experiments.
- Widely used in the development of food, pharmaceuticals, new materials, and sample coating.


Specifications

| Model | | GB211C-A |
|--------------------------|------------------------------------|--|
| Applicable test material | | Water-soluble & organic solvent (when connected to GAS411C) |
| Performance | Moisture evaporation capacity | Max. 1500ml/h |
| | Temp. regulator setting range | 0~240℃ (inlet temperature), 0~100℃ (outlet temperature) |
| | Temp. adjusting accuracy | Inlet temperature $\pm 1^{\circ}\text{C}$ |
| | Drying air volume adjustment range | 0.2~0.9 m ³ /min |
| | Spray air flow adjustment range | 0~30L/min |
| | Spray pressure usage range | 0.3~0.6Mpa |
| Composition | Nozzle cleaning function | Manual pulse air cleaning from the nozzle front end |
| | External output | Inlet temperature, outlet temperature output (4~20mA) |
| | Temp. regulator | Multi-PID control |
| | Touch screen | Temperature adjustment, blower, heater, liquid feed pump, pulse air switch, automatic needle, alarm display, operation curve |
| | Control switch | Inlet temperature, outlet temperature control switch |
| | Temp. sensor | PT100 thermistor |
| | Heater | 3.2KW |
| | Liquid feed pump | Duct type liquid feed pump |
| | Spray air pump | When using a spray air compressor (sold separately) or connecting to an organic solvent recovery device (sold separately) with a built-in air compressor |
| | Service socket | For stirrer: 200-230V~1A |
| | Suction blower | Tubular blower |
| | Filters | Suction filter, exhaust filter |
| | Solvent recovery | Use a solvent recovery device (sold separately) |
| | Spray nozzle cooling structure | It can connect to CF312L-B: joints $\times 2$, outer diameter $\Phi 10.5\text{mm}$ (sold separately) |
| | Spray air connection | Joint outer diameter, $\Phi 7\text{mm}$ |
| | Exhaust connection caliber | $\Phi 50\text{mm}$ |
| Safety functions | | Inlet, outlet temperature overheating, liquid feed pump reverse function, overcurrent leakage protection switch, nozzle connection abnormality (when connected to GAS411C) |
| Specifications | External dimensions | W760×D420×H1350mm |
| | Weight | Approx. 110kg |
| | Power supply | 200-230V~50/60Hz 17-21A |
| Accessories | | Liquid feed hose $\times 2$, exhaust hose (with 1 hose clamp) $\times 1$, exhaust conversion joint, outlet temperature sensor, fuse, antistatic connection cable, 5m intake hose (with 2 hose ties), nozzle conversion sleeve, stand assembly, protective cover (COV30), GF301C glass components, quick closure, air tube A, air tube B, product warranty, user manual |

Control panel



System diagram




Application



- Food and Pharmaceuticals
Milk powder, egg yolk, soy sauce, coffee, starch, protein, hormones, serum, antibiotics, enzyme flavorings, extracts, etc.
- Organic Chemistry
Paraffin, fuels, detergents, surfactants, pesticides, preservatives, synthetic resins, pigments, etc.
- Inorganic Chemistry
Ferrites, ceramics, toner, magnetic tape materials, photosensitive materials, various industrial chemicals, test sample waste liquids, etc.

Relationship between blower output power and maximum inlet temperature (reference)



Operability



The disassembly or cleaning of the drying chamber, cyclone, and product collection container adopts a quick-insert method, making it very convenient.

Organic Solvent Recovery Device GAS411C




Spray nozzle



The top part of the spray consists of a liquid nozzle and a gas nozzle.

Liquid nozzle (F)

Gas nozzle (A)



| Code | Nozzle size | | Nozzle specifications | Applicable models |
|--------|----------------------|--|-----------------------|--|
| 281298 | Liquid cap PF2050-SS | Aperture A=0.508mm Aperture B=1.270mm | 1 | ADL311(S)/GB210-A |
| | Air cap PA64-SS | Aperture C=1.626mm | | |
| | | | | |
| 281290 | Liquid cap PF2050-SS | Aperture A=0.508mm Aperture B=1.270mm | 2A | ADL311(S)/GB210-A |
| | Air cap PA70-SS | Aperture C=1.778mm | | |
| 281292 | Liquid cap PF2850-SS | Aperture A=711μm Aperture B=1270μm | 3 | ADL312SC/GB211C/DL411C/DL410 (standard) ADL311(S)/GB210-A |
| | Air cap PA64-SS | Aperture C=1626μm | | |
| 281291 | Liquid cap PF2850-SS | Aperture A=711μm Aperture B=1270μm | 2 | ADL311(S)/ADL312SC/GB210-A/ GB211C/DL410/DL411C |
| | Air cap PA-70-SS | Aperture C=1778μm | | |

Repeatability of spray drying tests

| Test No. | Sample name | Sample concentration (%) | Drying conditions | | | | | | | Recovered amount (g) | Recovery rate (%) |
|----------|-----------------|--------------------------|-------------------|-------------------|-------------------------|---------------------------------|------------------------|-----------------------------------|-----------------|----------------------|-------------------|
| | | | Inlet temp. (°C) | Outlet temp. (°C) | Dry air amount (m³/min) | Spray air pressure (Mpa/kg/cm²) | Test sample amount (g) | Sample liquid feed amount (g/min) | Test time (min) | | |
| 1 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147(1.5) | 198 | 6.6 | 30 | 8.1 | 81.8 |
| 2 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147(1.5) | 198.7 | 6.6 | 30 | 8.1 | 81.5 |
| 3 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147(1.5) | 200.6 | 6.7 | 30 | 8 | 79.8 |
| 4 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147(1.5) | 198.1 | 6.6 | 30 | 8.2 | 82.8 |
| 5 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147(1.5) | 199.3 | 6.6 | 30 | 8.4 | 84.3 |

Spray Dryers | Powder Granulating, Drying, Mixing

GB211C-B

Evaporated water capacity 1500ml/h

Temp. adjustment range 40~240℃

Liquid sending pump flow rate range 0~26ml/min

Spray dryer for powder granulation and drying of wet powders, one machine with multi-purpose, save more space.



Features

This device uses a fluidized bed for powder granulation and drying of wet powders. It is a fluidized bed drying and granulating device formed by combining the main unit GB211C with accessory GF200.

The drying chamber is made of ultra-hard glass, allowing observation of the fluidized bed state and spray state. Additionally, it is very convenient for data verification of flow meters, spray pressure gauges, inlet temperature, and outlet temperature. After purchasing the GF301 spray attachment, it can spray as well as fluidized bed granulate, providing strong expandability. It serves multiple purposes, accommodating various experimental objectives, effectively saving laboratory space, and avoiding multiple purchases, thus saving costs.

- 7-inch ultra-large touch screen with options for Chinese/Japanese/English, easy and convenient to operate.
- Automatic lifting function facilitates the installation and removal of the drying chamber.
- Experimental data recording and storage can be realized (optional function).
- Remote control is possible (optional function).
- Temperature zoning control is adopted, heating up faster and more stable.
- High-power heater significantly reduces temperature reaching time and offers the wider temperature setting range, meeting more sample experiments.
- Widely used in the R&D of food, pharmaceuticals, and new materials, as well as sample coating.


Specifications

| Model | | GB211C-B |
|----------------|--------------------------------------|--|
| Performance | Moisture evaporation capacity | Max. 1500ml/h |
| | Temperature regulator setting range | 0~240℃ (inlet temperature), 0~100℃ (outlet temperature) |
| | Temp. adjusting accuracy | Inlet temperature $\pm 1^{\circ}\text{C}$ |
| | Drying air volume adjustment range | 0.2~0.9 m ³ /min |
| | Spray air flow adjustment range | 0~30L/min |
| | Spray pressure usage range | 0.3~0.6Mpa |
| | Nozzle cleaning function | Manual pulse air cleaning from the nozzle front end |
| Composition | External output | Inlet temperature, outlet temperature output (4~20mA) |
| | Temperature regulator | Multi-PID control |
| | Touch screen | Temperature adjustment, blower, heater, liquid feed pump, pulse jet switch, automatic needle insertion, alarm display, operation curve |
| | Control switch | Inlet temperature, outlet temperature control switch |
| | Temperature sensor | PT100 thermistor |
| | Heater | 3.2KW |
| | Liquid feed pump | Duct type liquid feed pump |
| | Spray air pump | Use spray air compressor (sold separately) |
| | Service socket | For stirrer: 200-230V~1A |
| | Suction blower | Tubular blower |
| | Filters | Suction filter, exhaust filter |
| | Spray nozzle cooling structure | It can connect to CF312L-B: joints $\times 2$, outer diameter $\Phi 10.5\text{mm}$ (sold separately) |
| | Spray air connection | Joint outer diameter, $\Phi 7\text{mm}$ |
| | Exhaust connection caliber | $\Phi 50\text{mm}$ |
| | Safety functions | Overheating at inlet and outlet temperatures, liquid feed pump reverse function, overcurrent and leakage protection switch, abnormal nozzle connection (when connecting to GAS411C) |
| Specifications | External dimensions | W760×D420×H1350mm |
| | Weight | Approx. 110kg |
| | Power supply (50/60Hz) rated current | 200-230V~ 50/60Hz 17-21A |
| Accessories | | 2 liquid feed hoses, 1 exhaust hose (with 1 hose tie), exhaust conversion joint, outlet temperature sensor, fuse (250V 2A), antistatic connection cable, 5m intake hose (with 2 hose ties), nozzle conversion sleeve, stand assembly, protective cover (COV30), GF200 glass components |

Control panel



System diagram



Application



- Granulation, drying, mixing of powders
- Pharmaceuticals, food, catalysts, fuels, detergents, ceramics, etc.
- Since it is suitable for about 50~300g samples, it is most suitable for high-value samples or research-level experiments

Operability




The disassembly or cleaning of the drying chamber, cyclone, and product collection container adopts a quick plug-in method, which is very convenient.

Spray nozzle



The top part of the spray consists of a liquid nozzle and an air nozzle.

Liquid nozzle (F) Air nozzle (A)



| Code | Nozzle size | | Nozzle specifications | Applicable models |
|--------|----------------------|---|-----------------------|---|
| 281297 | Liquid cap PF1650-SS | Aperture A = 406μm Aperture B = 1270μm | 1A | ADL311(S)/GB210-A/GB210-B(standard) |
| | Air cap PA64-SS | Aperture C=1626μm | | |
| 281298 | Liquid cap PF2050-SS | Aperture A = 508μm Aperture B = 1270μm | 1 | ADL311(S)/GB210-A |
| | Air cap PA64-SS | Aperture C=1626μm | | |
| 281290 | Liquid cap PF2050-SS | Aperture A = 508μm Aperture B = 1270μm | 2A | ADL311(S)/GB210-A |
| | Air cap PA70-SS | Aperture C=1778μm | | |
| 281292 | Liquid cap PF2850-SS | Aperture A = 711μm Aperture B = 1270μm | 3 | ADL312SC/GB211C/DL411C/DL410 (standard)/ADL311(S)/GB210-A |
| | Air cap PA64-SS | Aperture C=1626μm | | |
| 281291 | Liquid cap PF2850-SS | Aperture A = 711μm Aperture B = 1270μm | 2 | ADL311(S)/ADL312SC/GB210-A/GB211C/DL410/DL411C |
| | Air cap PA-70-SS | Aperture C=1778μm | | |

Implementation case

| No. | Sample | | | Adhesive | | | Measurement conditions | | | | | | Results | |
|-----|--|---------------|--------|------------|-----------------|--------------|------------------------|--------------------------|-----------------------------------|---------------------|-------------|------------------|-----------------------|---------------------|
| | Name | Particle size | Weight | Name | Concentration % | Spray amount | Inlet temp. °C | Liquid feed amount g/min | Spray pressure kg/cm ² | Spray amount g/time | Spray times | Nozzle height cm | Average particle size | 1400~125μm Recovery |
| 1 | Metatitanic acid-lead zirconate titanate | 10 | 280 | PVA | 10 | 42 | 150 | 5 | 0.4 | 6 | 8 | 25 | 180 | 75 |
| 2 | Sodium fluoride | | 200 | PVA + EDTA | 24 | 165 | 150 | 5 | 0.5 | 15 | 11 | 28 | 320 | 97 |
| 3 | Diamond + metal powder | | 538 | PVA | 10 | 108 | 170 | 6 | 0.4 | 5 | 22 | 30 | 260 | 80 |
| 4 | Pharmaceutical | | 200 | Dextrin | 5 | 171 | 150 | 5 | 0.5 | 8~10 | 21 | 27 | 310 | 85 |
| 5 | Flavor | | 300 | Dextrin | 5 | 60 | 150 | 5 | 0.4 | 8 | 12 | 30 | 330 | 98 |
| 6 | Lactose | | 200 | Sorbitol | 10 | 70 | 100 | 14 | 1.0 | 17 | 4 | 25 | 390 | 80 |

Spray Dryers | Large Capacity

DL411C

Evaporated water capacity 3000ml/h

Temp. adjustment range 40~300°C

Liquid sending pump flow rate range 0~70ml/min

Brand-new upgraded large spray dryer that efficiently recovers 100µm micro-particles.



Features

This device can obtain 10~100µm micro-particles, which are considered very difficult to obtain on laboratory spray dryer devices.

Purposes of this device: It can be used in production preparatory experiments, high-value sample experiments, capsule spray drying method application research, and as an alternative method to traditional drying methods in general laboratories.

- It adopts a high-efficiency two-fluid nozzle spray method and a large drying chamber, compared with other products, the sample adhesion amount in the DL411C drying chamber is less, and the drying chamber is designed with a special cold air inlet to effectively prevent spray particles from refluxing and adhering upwards.
- The large hard glass drying chamber is very convenient for observing the spray and particles status.
- The nozzle head can automatically perform pulse cleaning at set intervals.
- The wide range of temperature settings for the hot air inlet can improve processing capacity and experimental efficiency.

Specifications

| Model | | | DL411C | |
|--------------------|------------------------------------|--|---|--|
| Applicable samples | | | Water-soluble & organic solvent (when connected to GAS510C) | |
| Performance | Moisture evaporation capacity | Max. 3000ml/h | | |
| | Spray method | Two-fluid nozzle method (sieve aperture 0.7mm), three-fluid nozzle method (optional) | | |
| | Spray and hot air contact method | Vertical downward spray and flow method | | |
| | Liquid feed pump | Adjustable from 0 to 70 ml/min | | |
| | Temp. regulator setting range | 0~300℃ (inlet temperature), 0~100℃ (outlet temperature) | | |
| | Temp. adjusting accuracy | Inlet temperature ±1℃ | | |
| | Drying air volume adjustment range | 0.12~1.0 m³/min | | |
| | Spray air flow adjustment range | 0~30L/min | | |
| | Spray pressure usage range | 0.3~0.6Mpa | | |
| Composition | Drying chamber | Size | Internal diameter 457×height 975 mm (glass part) | |
| | | Material | The drying chamber, cyclone unit, and product container are made of ultra-hard glass; other piping materials are stainless steel and silicone | |
| | External output | | Inlet temperature, outlet temperature output (4~20mA) | |
| | Temp. regulator | | Multi-PID control | |
| | Touch screen | | Temperature regulation, blower, heater, liquid delivery pump, pulse jet switch, automatic needle cleaning, alarm display, and operational curve | |
| | Control switch | | Inlet temperature, outlet temperature control switch | |
| | Temp. sensor | | PT100 thermistor | |
| | Heater | | 2.0KW×2 | |
| | Spray air pump | | Using a spray air compressor (sold separately) or connecting to an organic solvent recovery device GAS510C (sold separately) employs the built-in air compressor of GAS510C | |
| | Service socket | | For stirrer: 200-230V~1A | |
| | Suction blower | | Brushless blower | |
| | Filters | | Suction filter, exhaust filter | |
| | Solvent recovery | | Using solvent recovery device GAS510C (sold separately) | |
| | Spray nozzle cooling structure | | It can connect CF312L-B: connector ×2, outer diameter Φ10.5mm | |
| | Spray air connection | | Joint outer diameter, Φ7mm | |
| | Exhaust connection caliber | | Φ50mm | |
| | Safety functions | | | Overheating at inlet and outlet temperature, liquid delivery pump reverse function, overcurrent leakage protection switch, independent overheating preventer for the outside of the heater |
| | Specifications | External dimensions | | W1080×D850×H1770mm |
| Weight | | Approx. 180kg | | |
| Power supply | | 200-230V~50/60Hz 23-27A | | |
| Accessories | | | See Accessory list for details | |

The large size of the drying chamber allows the equipment to provide sufficient time for drying fine particles, thus enabling the production of particles ranging 10~100 µm, similar to products manufactured in industries such as ceramics and pharmaceuticals.

- 7-inch ultra-large touch screen with options for Chinese/Japanese/English, easy and convenient to operate.
- The cover of the drying chamber can be raised and rotated, facilitating easy cleaning.
- The system is equipped with a new automatic nozzle featuring a cooling mechanism.
- Both two-fluid and three-fluid nozzles can be used.
- Experimental data recording and storage can be realized.
- Remote control is possible (optional function).
- Temperature zoning control is adopted, heating up faster and more stable.
- A service socket is provided for easy power supply to the stirrer.

Configuration



The cyclone separation unit is equipped with a safety cover to prevent burns and an anti-static brush.

Filtering device



When installed at the air inlet, it effectively filters out particulate impurities from the incoming air, preventing contamination of the spray sample.
When installed at the exhaust port, it can effectively filter uncollected sample matter for secondary recovery purposes.

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps


12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

System diagram



Application cases

① Powder & Particle


Through granulation and spheroidization processing, the flowability of powders is significantly enhanced, allowing for uniform filling during compaction. For this purpose, it is used in aluminum, zirconium, various ceramics, precious metals, super hard alloys, etc.

② Microencapsulation

Mixing and adjusting core substances and membrane materials during spray drying results in encapsulated powder from the primary solution.

[Application Case]

- Ink powder for carbonless copy paper
- Flavor adjustment for pharmaceuticals, adjusting solubility time
- Microencapsulation of fragrances for pharmaceuticals or hygiene-related products
- Microencapsulation of other pigments, fertilizers, oils, adhesives, etc.



0 50 100μm

Powder generated by DL410


③ Spray Cooling Granulation

Wax, oils, fatty acids, etc. which are difficult to dry and powdered.

④ Special Application Cases


Spray concentration, spray reaction, powder granulation, etc.

Drying time



Drying time for complete evaporation of droplets at various hot air temperatures

Airflow and temperature usage curves




Control panel



Spray nozzle

| Code | Nozzle size | Nozzle specifications | Applicable models |
|------------|--|-----------------------|--|
| 281298 | Liquid cap PF2050-SS Aperture A = 0.508mm Aperture B = 1.270mm | 1 | ADL311(S)/GB210-A |
| | Air cap PA64-SS Aperture C = 1.626mm | | |
| 281290 | Liquid cap PF2050-SS Aperture A = 0.508mm Aperture B = 1.270mm | 2A | ADL311(S)/GB210-A |
| | Air cap PA70-SS Aperture C = 1.778mm | | |
| 281292 | Liquid cap PF2850-SS Aperture A = 711μm Aperture B = 1270μm | 3 | ADL312SC/GB211C/DL411C/DL410 (standard) ADL311(S)/GB210-A |
| | Air cap PA64-SS Aperture C = 1626μm | | |
| 281291 | Liquid cap PF2850-SS Aperture A = 711μm Aperture B = 1270μm | 2 | ADL311(S)/ADL312SC/GB210-A/GB211C/ DL410/DL411C |
| | Air cap PA-70-SS Aperture C = 1778μm | | |
| A080999034 | Liquid cap PF40100-SS Aperture A = 1016μm Aperture B = 2540μm | 6 | DL410/DL411C |
| A080999035 | Air cap PA120-SS Aperture C = 3048μm | | |
| A080999049 | Liquid cap PF60100-SS Aperture A = 1524μm Aperture B = 2540μm | 4 | DL410/DL411C |
| A080999035 | Air cap PA120-SS Aperture C = 3048μm | | |
| A080999039 | Liquid cap PF100150-SS Aperture A = 2500μm Aperture B = 3810μm | 5 | DL410/DL411C |
| A080999041 | Air cap PA180-SS Aperture C = 4572μm | | |



Organic Solvent Recovery Units | N₂ Closed Cycle

GAS411C/510C

Circulating airflow 0.12~0.65m³/min

GAS411C recovery capacity Above 1500ml/h

GAS510C recovery capacity Above 3000ml/h

Brand new upgraded organic solvent recovery circulation system.



Features


When using organic solvents in the spray dryer ADL312SC, GB211C-A, DL411C, this recovery unit is used to prevent external discharge.

- Through a nitrogen closed-loop cycle and using a refrigerator and condenser for solvent recovery, it enables the handling of combustible materials.
- It can dry easily oxidizable substances.
- Due to low-temperature drying capabilities, it can also dry materials prone to thermal deformation.
- Under stringent safety measures, spray drying allows for the recovery of products and solvents.

Specifications

| Model | GAS411C | GAS510C |
|-------------------------|--------------------------------------|---|
| Solvent recovery method | Condensation recovery | |
| Performance | Circulating gas | Nitrogen replacement + closed loop |
| | Condensation temp. setting range | -20~30℃ |
| | Circulating airflow | 0.12~0.65 m ³ /min |
| | Compressed air source flow | 16L/min 32L/min |
| Composition | Air compressor (for spraying) | Linear reciprocating air compressor Oil-free piston compressor |
| | Circulation fan | Roots blower |
| | Solvent recovery container | 2L flask (with anti-fall flask holder) |
| | Refrigerator | Air-cooled fully sealed refrigeration unit, 400W R404A Air-cooled fully sealed refrigeration unit, 735W R404A |
| | Solution recovery container | Single condenser cooling structure Double condenser cooling structure |
| | Filters | 100-level high-efficiency filter (corrosion-resistant type) |
| | Control panel | Condensation temperature control and display |
| | | Blower airflow control and display |
| | | Real-time oxygen concentration display |
| | | Switching of spray air sources |
| | Oxygen concentration meter | 0.1~25% Vol |
| | Micro pump | Pump suction type oxygen sensor monitoring |
| Safety functions | | Oxygen concentration upper limit alarm, flammable gas alarm, over-current protection switch, nitrogen forced introduction (when oxygen concentration rises) |
| Specifications | External dimensions | W710×D950×H1450mm |
| | Weight | 210kg 240kg |
| | Power supply (50/60Hz) rated current | 200-230V~50/60Hz 5-6A 200-230V~50/60Hz 5.5-12A |
| Accessories | | Liquid delivery hoses (2 each of silicone and Viton), 2 stainless steel corrugated hoses, PVC exhaust pipe, connecting pipes, recovery flask (2L) |


Control panel



- Brand new 7-inch ultra-large touch screen control panel, available in Chinese/Japanese/English, easy and convenient to operate.
- It can communicate with the control panel of the spray dryer.
- It can automatically select the spray air source.
- Remote control is available.
- All parts have status indicators.
- Circulating airflow display (optional).
- Main functions


- ① Closed system (nitrogen closed-loop type)
- ② Oxygen concentration control function
- ③ Flammable gas detection function
- ④ Inlet temperature overheat detection function
- ⑤ Outlet temperature overheat detection function
- ⑥ Other self-diagnosis functions (temperature sensor disconnection detection/overheating prevention/nozzle)

System diagram



| No. | Part name | No. | Part name |
|-----|---|-----|--------------------------------|
| 1 | Evaporators | A | O-rings |
| 2 | Sensors | B | Seal rings |
| 3 | Ball valves | C | Corrugated tubes |
| 4 | Flask holders | D | Spray nozzle |
| 5 | Recovery flasks | E | Coiled tubes |
| 6 | Filters | F | Aluminum honey-comb rectifiers |
| 7 | Filter chambers | G | Safety helmet |
| 8 | Differential pressure gauges | H | Pressure gages |
| 9 | Flow meter (for N ₂ introduction) | I | Push pin valve |
| 10 | Compressor | J | Three-way valve |
| 11 | Electromagnetic valve (for N ₂ control) | K | Electromagnetic valve |
| 12 | Flow meter (for O ₂ concentration) | L | Seal rings |
| 13 | Activated carbon filter | | |
| 14 | Pump | | |
| 15 | O ₂ Sensors | | |
| 16 | Electromagnetic valve (for exhaust) | | |
| 17 | Blower | | |
| 18 | Electromagnetic valve (for N ₂ introduction) | | |
| 19 | Electromagnetic valve (for air inlet) | | |

Operation



Applicable organic solvents reference

Selection of liquid delivery hoses

Silicone hose: Ethanol, IPA, Methanol, Acetone, Ethyl Acetate

Viton hose: Xylene, Toluene, Benzene, Ethane, Chloroform, Dichloromethane

[Boundary oxygen concentration table]

| Organic Solvent | Boiling Point [°C] | Melting Point [°C] | Boundary Oxygen Concentration [%] |
|-----------------|--------------------|--------------------|-----------------------------------|
| Xylene | (o)-144 | (o)-25 | (o)-10.5 |
| IPA | 82.3 | -88 | 9.0 |
| Benzene | 80.1 | 5.5 | 10.5 |
| Ethanol | 78.4 | -114.3 | 9.9 |
| Ethyl Acetate | 77.1 | -83.6 | 10.0 |
| Ethane | 67.7 | -95.3 | 11.4 |
| Methanol | 64.6 | -97.4 | 9.7 |
| Chloroform | 61.2 | -63.5 | Non-flammable |
| Acetone | 56.2 | -94.6 | 10.4 |
| Dichloromethane | 40.0 | -97.7 | 23.9 |

Usage objects and fields



- Non-oxide ceramics related
- Polymer materials related
- Superconducting materials related
- Pharmaceuticals related
- Food-related

Circulation & airflow monitoring meter (optional)



Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Solvent Washing Units | Corrosive Solvent Wet Washing

GWS410

Maximum flow rate 15L/min

The world's first solvent washing device, mainly used for spray dryers.




GWS410 captures pollutants in harmful gases using water or alkaline solution at atmospheric pressure and room temperature. GWS410 is designed with a washing tank—when solvent vapors enter the washing tank, sprayed water adheres, cleans, and neutralizes solvent particles, then returns to the bottom of the tank.

Features

- Using water or alkaline solution
- Eliminate highly irritating solvents
- Reduce rust and corrosion of equipment
- Easy to operate
- Easy to maintain—just monitor the pH level of water in the water storage tank and the state of the molecular sieve

System diagram




Specifications

| Model | GWS410 |
|-----------------------------|--------------------------------------|
| System | Flow-type spraying washing type |
| Circulating liquid | Water/alkaline solution |
| Circulation pump | Small magnetic circulation pump |
| Maximum flow rate | 15L/min |
| Maximum lift | 8m |
| Harmful gas cleaning method | Pall ring filling + spraying washing |
| Water storage tank capacity | 35L |
| Safety device | Leakage protection circuit breaker |
| Power supply | 200-230V~0.35-1.0A |
| Dimensions (WxDxH) | 800×500×1230 mm |
| Weight | Approx. 120kg |

*Dimensions do not include protrusions.

Connection diagram



Spray Dryer + GWS410 Organic Solvent Cleaning Device + GAS411C Organic Solvent Recovery Device

- (1) Harmful gases exiting the spray dryer ① enter the washing tank unit.
- (2) Through the cleaning tank unit ②, harmful gases come into contact with the filling material ③ and the cleaning liquid ④ that is sprayed from the nozzle, allowing harmful substances to be absorbed by the cleaning liquid.
- (3) Through multiple filling chambers, the gas passes through the smoke collector ⑤ to prevent the cleaning liquid from being discharged.
- (4) The gas, under the influence of the fan, enters the GAS411C as clean air exhaust ⑥.
- (5) Cleaning liquid A enters the washing tank from the water storage tank through circulation pump B, is dispersed onto the filling material ③ via spray nozzle C, and returns to the water storage tank's circulating water tank through piping D.

Reference Materials

Spray drying-spray dryers for organic solvents


| Sieve hole | # 1 | # 2 | # 3 | # 4 |
|----------------------------|-------|-------|-------|-------|
| Above 12 | 5.6 | 0.8 | 1.2 | 1.3 |
| 12~16 | 0.5 | 0.9 | 1 | 1.2 |
| 16~24 | 0.6 | 0.8 | 1.2 | 1.4 |
| 24~32 | 0.7 | 0.8 | 0.9 | 1.1 |
| 32~42 | 1.6 | 1.7 | 1.9 | 1.8 |
| 42~60 | 5.9 | 4.3 | 4.8 | 3.5 |
| 60~80 | 9.6 | 8.5 | 8.5 | 6.6 |
| 80~115 | 13.2 | 15.6 | 13.4 | 12.8 |
| Below 115 | 66.8 | 66.6 | 67 | 70.3 |
| Average particle diameter※ | 135.6 | 135.7 | 138.3 | 136.9 |

Repeatability of granulation test
Due to the numerous operational factors of the mini spray granulation method, the repeatability may vary depending on the operator's skill level. Since the flow state of the particles can greatly affect the test results, adjusting the hot air volume is typically done to achieve consistent flow conditions.

Conditions

| | |
|-----------------------|--|
| Raw material | Sintered Alumina (Average particle size 40) 400g |
| Binder | 5% PVA solution (#500) 25g |
| Inlet temperature | 100℃ |
| Binder delivery speed | 12.4g/min |
| Binder spray times | 6 times |
| Binder spray pressure | 78KPa (0.8kg/cm ²) |
| Nozzle height | Distance from micropore plate 25cm |

Particle size variation




Factors influencing particle size variation include inlet temperature, binder delivery speed, spray pressure, etc., but the greatest influence comes from the amount of binder added relative to the amount of raw material powder. As the binder amount increases, the particle size will naturally increase.

Conditions

| | |
|-----------------------|-------------------------------------|
| Raw material | Lactose (Below 100 sieve hole) 200g |
| Binder | 70% sorbitol solution |
| Inlet temperature | 90℃ |
| Binder delivery speed | 12g/min |
| Binder spray pressure | 98KPa (1.0kg/cm ²) |
| Nozzle height | Distance from micropore plate 25cm |

Repeatability of granulation test



Particles generated by the granulation device typically range in distribution from 0.1 to 1.5μ; they do not achieve a uniform particle size like those produced by extrusion or compression granulation methods, but this size range can be adjusted to some extent based on experimental conditions, etc.

Conditions

| | |
|-----------------------|--------------------------------------|
| Raw material | Lactose (Below 100 sieve hole) 200g |
| Binder | 70% sorbitol solution 7.3g |
| Inlet temperature | 90℃ |
| Binder delivery speed | 12g/min |
| Binder spray times | 7 times |
| Binder spray pressure | 98KPa (1.0kg/cm ²) |
| Nozzle height | Distance from micropore plate 22.5cm |

Implementation cases (spray drying ADL312SC-A)

| Sample | Composition (%) | Inlet temp. (℃) | Outlet temp. (℃) | Dry air amount (m ³ /min) | Spray air pressure KPa(kg/cm ²) | Sample liquid feed amount (g/min) | Sample recovery rate (%) |
|-------------------------------|-----------------|-----------------|------------------|--------------------------------------|---|-----------------------------------|--------------------------|
| Dextrin (solution) | 10 | 150 | 80 | 0.4 | 98 (1.0) | 6.1 | 66 |
| Dextrin (emulsion) | 40 | 150 | 80 | 0.4 | 98 (1.0) | 5.1 | 63 |
| Titanium dioxide (suspension) | 10 | 150 | 85 | 0.42 | 98 (1.0) | 5.3 | 50 |
| Soy sauce | 50 | 130 | 75 | 0.36 | 98 (1.0) | 5.1 | 60 |
| Salt | 10 | 145 | 85 | 0.38 | 98 (1.0) | 5.3 | 52 |

Repeatability of spray drying tests (Spray dryer ADL312SC-A)

| Test No. | Sample | Composition (%) | Drying conditions | | | | | | | Recovery amount (g/min) | Recovery rate (%) |
|----------|-----------------|-----------------|-------------------|------------------|--------------------------------------|---|----------------------------|-----------------------------------|-----------------|-------------------------|-------------------|
| | | | Inlet temp. (℃) | Outlet temp. (℃) | Dry air amount (m ³ /min) | Spray air pressure KPa(kg/cm ²) | Test sample amount (g/min) | Sample liquid feed amount (g/min) | Test time (min) | | |
| 1 | Coffee solution | 5.00 | 150 | 75 | 0.45 | 147 (1.5) | 93.1 | 3.1 | 30 | 4.3 | 92.4 |
| 2 | Coffee solution | 5.00 | 150 | 75 | 0.45 | 147 (1.5) | 93 | 3.1 | 30 | 4 | 86.0 |
| 3 | Coffee solution | 5.00 | 150 | 75 | 0.45 | 147 (1.5) | 91.4 | 2.0 | 30 | 4 | 87.5 |
| 4 | Coffee solution | 5.00 | 150 | 75 | 0.45 | 147 (1.5) | 84.9 | 2.8 | 30 | 3.7 | 87.2 |
| 5 | Coffee solution | 5.00 | 150 | 75 | 0.45 | 147 (1.5) | 83.8 | 2.8 | 30 | 3.7 | 88.3 |

Implementation cases (granulation device GB211C-A)

| Sample | Composition (%) | Inlet temp. (℃) | Outlet temp. (℃) | Dry air amount (m ³ /min) | Spray air pressure KPa(kg/cm ²) | Sample liquid feed amount (g/min) | Sample recovery rate (%) |
|-------------------|-----------------|-----------------|------------------|--------------------------------------|---|-----------------------------------|--------------------------|
| Dextrin | 20% Solution | 140 | 85 | 0.48 | 147 (1.5) | 8.8 | 60 |
| Pharmaceutical | 10% suspension | 145 | 80 | 0.42 | 196 (2.0) | 8.2 | 82 |
| Black tea extract | 20% solution | 155 | 100 | 0.4 | 147 (1.5) | 7.8 | 72 |
| Silica gel | 20% solution | 150 | 75 | 0.48 | 147 (1.5) | 12.6 | 70 |
| Iron oxide | 3% suspension | 175 | 90 | 0.4 | 127 (1.3) | 9.5 | 75 |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Implementation cases (granulation equipment GB211C-B)

| Sample | | Binder | | | Test conditions | | | | | Results | |
|----------------------------------|------------|----------|-------------------|------------------|------------------|------------------------------|---|---------------------|--------------------|-----------------------|--|
| Name | Weight (g) | Name | Concentration (%) | Spray amount (g) | Inlet temp. (°C) | Liquid sending speed (g/min) | Spray pressure KPa(kg/cm ²) | Spray times (times) | Nozzle height (cm) | Average diameter (μm) | Recovery rate for 12~115 sieve holes (%) |
| Silica | 200 | PVA | 5.0 | 77 | 125 | 15 | 59 (0.6) | 4 | 27 | 339 | 58 |
| Iron oxide | 160 | PVA | 2.5 | 50 | 120 | 15 | 98 (1.0) | 4 | 21 | 205 | 62 |
| Ceramics | 200 | PVA | 3.0 | 106 | 120 | 15 | 78 (0.8) | 3 | 22 | 404 | 82 |
| Alumina | 160 | PVA | 3.0 | 60 | 110 | 15 | 59 (0.6) | 4 | 22 | 311 | 88 |
| Silica dioxide | 150 | CMC | 1.0 | 100 | 120 | 15 | 78 (0.8) | 4 | 22 | 306 | 60 |
| Lactose | 200 | Sorbitol | 70.0 | 10 | 100 | 14 | 98 (1.0) | 4 | 25 | 390 | 80 |
| Black tea extract-powdered sugar | 250 | Guar gum | 0.5 | 24 | 85 | 6 | 59 (0.6) | 10 | 28 | 333 | 77 |
| Fat-containing powder | 200 | Glucose | 30.0 | 11 | 85 | 4 | 59 (0.6) | 7 | 22 | 236 | 82 |

Types and features of binders (reference materials)

| Type | Features |
|-------------------------------|---|
| Gelatin | Low concentration solutions have weaker adhesion, so high concentration solutions should be heated before spraying. |
| Dextrin | Although the adhesive strength is weaker, it forms very well during tablet pressing. |
| Potato Starch | Good particle properties and inexpensive, used in medical and food fields. |
| Sodium Aluminum | Due to high viscosity, it is very suitable as a binder, mainly used in the food field. |
| Arabic Gum | Need to be sprayed after heating, requires a large amount of binder. |
| CMC (Carboxymethyl Cellulose) | High viscosity at low temperatures, leading to more powder residue. |
| HPC (Hydroxypropyl Cellulose) | Suitable for hygroscopic materials with good cohesiveness. |
| MC (Methyl Cellulose) | Strong adhesive power, suitable for products that increase particle size. |
| PVA (Polyvinyl Alcohol) | Although it is very good in granulation, it has some difficulty in crushing granulated products. |
| PVP (Polyvinylpyrrolidone) | High molecular weight substances have strong adhesion, suitable for hydrophobic materials. |

Repeatability of spray drying tests (granulation equipment GB211C-A)

| Test No. | Sample | Sample concentration (%) | Drying conditions | | | | | | | Recovery (g/min) | Recovery rate (%) |
|----------|-----------------|--------------------------|-------------------|-------------------|--------------------------------------|---|----------------------------|--------------------------------------|-----------------|------------------|-------------------|
| | | | Inlet temp. (°C) | Outlet temp. (°C) | Dry air amount (m ³ /min) | Spray air pressure KPa(kg/cm ²) | Test sample amount (g/min) | Sample liquid sending amount (g/min) | Test time (min) | | |
| 1 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147 (1.5) | 198.0 | 6.6 | 30 | 8.1 | 81.8 |
| 2 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147 (1.5) | 198.7 | 6.6 | 30 | 8.1 | 81.5 |
| 3 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147 (1.5) | 200.6 | 6.7 | 30 | 8.0 | 79.8 |
| 4 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147 (1.5) | 198.1 | 6.6 | 30 | 8.2 | 82.8 |
| 5 | Coffee solution | 5.00 | 150 | 80 | 0.45 | 147 (1.5) | 199.3 | 6.6 | 30 | 8.4 | 84.3 |

Implementation cases: combination of granulation equipment GB211C-A and organic solvent recovery equipment GAS411C (for organic solvents)

| Sample | Sample concentration (%) | Inlet temp. (°C) | Outlet temp. (°C) | Amount of nitrogen for drying (m ³ /min) | Spraying pressure (kg/cm ²) | Delivery speed (g/min) | Dispersing medium or solvent | Results | | | Others |
|-------------------------------|--------------------------|------------------|-------------------|---|---|------------------------|------------------------------|---------------|-------------------|---------------------------|--------------------------------------|
| | | | | | | | | Powderization | Recovery rate (%) | Solvent recovery rate (%) | |
| Hydroxypropyl Methylcellulose | 10 | 90 | 55 | 0.5 | 1.0 | 9.9 | Notes※ | Excellent | 65.3 | 92.5 | ※ Chloroform 1: Ethanol 1 |
| Cellulose-based Polymers | 5.0 | 70 | 47 | 0.5 | 1.0 | 8.3 | Dichloromethane | Excellent | 72.3 | | |
| Polymer | 2.0 | 100 | 64 | 0.5 | 1.0 | 8.4 | Notes※ | Excellent | 77.8 | 80.7 | ※ Ethanol 95: Water 5 |
| Resins | 23.5 | 80 | 55 | 0.5 | 1.0 | 4.2 | Notes※ | Excellent | 81.9 | 96.7 | ※ Dispersed in (Methanol 4: Water 1) |
| Carbon + Resin | 5.8 | 100 | 70 | 0.5 | 1.0 | 5.3 | IPA | Excellent | 85.1 | 94.1 | |
| Polymer + Inorganic Salt | 10.2 | 140 | 98 | 0.5 | 1.0 | 3.8 | Notes※ | Excellent | 97.6 | 97.4 | ※ Dimethylacetamide |
| Polyvinylpyrrolidone (K30) | 10.0 | 80 | 55 | 0.5 | 1.0 | 7.7 | Ethanol | Excellent | 79.4 | 95.0 | |
| Polyvinylpyrrolidone + Drug | 10.0 | 80 | 70 | 0.5 | 1.0 | 7.7 | Ethanol | Excellent | 75.9 | 95.4 | |
| Plant Extracts | 3.0 | 130 | 71 | 0.5 | 1.0 | 9.1 | Notes※ | Excellent | 96.5 | 91.9 | ※ Ethanol 6: Water 4 |
| Silicon Carbide | 38.5 | 150 | 84 | 0.5 | 1.0 | 12.1 | Ethanol | Excellent | 89.9 | 99.9 | ※ Nozzle used 3S |
| Aluminum Nitride | 13.2 | 150 | 99 | 0.5 | 1.0 | 12.9 | Butyl Acetate | Excellent | 92.2 | 86.7 | ※ Nozzle used 3S |
| Nitride Ceramics | 60.5 | 120 | 83 | 0.5 | 1.0 | 11.3 | MEK | Excellent | 74.7 | 88.7 | |
| Superconducting Materials | 33.3 | 80 | 60 | 0.5 | 1.0 | 15.7 | Acetone | Excellent | 66.6 | 99.6 | |
| Pharmaceutical | 3.61 | 100 | 68 | 0.6 | 1.0 | 10.0 | Notes※ | Good | 73.6 | 87.2 | ※ Ethanol + Dichloromethane |
| Pharmaceutical | 13.2 | 60 | 45 | 0.32 | 1.25 | 6.0 | Notes※ | Good | 87.6 | 94.7 | ※ Dichloromethane + Ethanol |
| W-Cu | 50.0 | 100 | 62 | 0.5 | 0.5 | 20.7 | Ethanol | Good | 60.3 | 91.9 | |
| Modified polystyrene | 48.7 | 140 | 60 | 0.45 | 1.0 | 22.3 | Water | Good | 67.6 | 91.7 | |
| Polymer | 0.5 | 150 | 88 | 0.5 | 1.0 | 8.5 | Notes※ | Good | 83.1 | 97.6 | ※ Ethanol 3 + Water 1 |
| Organic matter | 50.0 | 150 | 88 | 0.4 | 1.0 | 8.3 | Methanol | Good | | | |
| Silica dispersion solution | 10.0 | 100 | 88 | 0.5 | 1.0 | 4.8 | Notes※ | Good | 96.2 | 99.5 | □ Ethanol + Water (A little) |

Muffle Furnaces

Contents

| | |
|-----------------------|----|
| Muffle Furnaces | 52 |
|-----------------------|----|

Muffle Furnaces | Standard

FO111C/211C/311C/411C/511C/611C/711C/811C

Operating temp. range 100~1150°C

Temp. adjusting accuracy $\pm 2^{\circ}\text{C}$

Internal capacity 1.5L 3.75L 7.5L 9L 11.3L 17.5L 23.6L 30L

8 types of internal capacity, equipped with program operation function.

Features

- Support program operation.
- Simple operation with up/down key settings, equipped with a high-operability controller for precise temperature control.
- Temperature, time and indicator values can all be displayed digitally.
- Use a controller with integrated digital settings and overheat protection. Support self-diagnosis functions, overcurrent leakage protection switch and other safety controls in case of anomalies.
- Use R-type thermocouple for long life.
- Equipped with an exhaust port.
- Optional features available, such as exhaust unit, sample trays, nitrogen gas introduction device (with flow meter), temperature output terminal and so on.



Specifications

| Model | | FO111C | FO211C | FO311C | FO411C | FO511C | FO611C | FO711C | FO811C |
|----------------|---|--|--------------|----------------|--------------|----------------|--------------|--------------|--------------|
| Performance | Operating temp. range | 100~1150℃ | | | | | | | |
| | Temp. adjusting accuracy | ±2℃ (at 1150℃) | | | | | | | |
| Composition | Time to reach maximum temp. at room temp. | Approx. 60 min | | Approx. 70 min | | Approx. 80 min | | | |
| | Exterior material | Cold rolled steel plate with chemical proofing coating | | | | | | | |
| | Furnace body | Ceramic fiber | | | | | | | |
| | Sensors | Dual R-type thermocouple | | | | | | | |
| | Heater | FeCrAl alloy heating wire | | | | | | | |
| | | 1KW | 1.5KW | 2KW | 2.2KW | 2.5KW | 3KW | 3.5KW | 4KW |
| Controllers | Exhaust vent | Inner diameter 20mm, upper part | | | | | | | |
| | Cooling fan blade | Axial flow fan | | | | | | | |
| | Temp. control method | Multi-PID control | | | | | | | |
| | Temp. setting/display | Digital setting/digital display | | | | | | | |
| | Overheating protection | Integrated controller | | | | | | | |
| | Overheat protection setting method | Digital setting | | | | | | | |
| | Operation functions | Fixed temp. operation, program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3), auto start and auto stop operation | | | | | | | |
| | Additional functions | Deviation correction function, power failure compensation function, key lock function | | | | | | | |
| | Timer | 1 min 99 h 59 min and 100 h~999 h 50 min | | | | | | | |
| | Safety device | Self-diagnostic circuit (temperature sensor abnormality, automatic overheat prevention), overheat preventer, leakage protection circuit breaker for over-current, key lock function | | | | | | | |
| Specifications | Internal dimensions (W×D×H mm) | 100×150×100 | 100×250×150 | 200×250×150 | 200×300×150 | 300×250×150 | 250×350×200 | 270×350×250 | 300×400×250 |
| | External dimensions (W×D×H mm) | 346×405×516 | 346×505×566 | 446×505×566 | 446×555×566 | 506×505×626 | 506×605×676 | 506×605×726 | 506×655×726 |
| | Internal capacity | 1.5L | 3.75 L | 7.5 L | 9 L | 11.3 L | 17.5 L | 23.6 L | 30 L |
| | Power supply (50/60Hz) rated current | AC220V 5A | AC220V 7.5A | AC220V 9.5A | AC220V 10.5A | AC220V 12A | AC220V 14A | AC220V 16.5A | AC220V 18.5A |
| | Weight | Approx. 24kg | Approx. 30kg | Approx. 37kg | Approx. 38kg | Approx. 44kg | Approx. 52kg | Approx. 58kg | Approx. 62kg |
| Accessories | | Exhaust vent cover 1 pcs | | | | | | | |
| Options | Stand | ON30C | | | ON61C | | | | |
| | Others | Setter plate, exhaust unit (AC220V), sample box, N ₂ introduction device (with flow meter), mini printer, data logger, combination warning light (standby/operation/failure), external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, central monitoring software | | | | | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices


15 Options

Automatic
Overheating
PreventerOverheat
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification

Control panel



Temperature rise characteristic curve (FO311C type)



Exhaust device unit (optional)



Able to quickly exhaust gases generated by rising temperature inside the chamber.
Piping: flexible pipe, aluminum
Length 1.5m
Diameter 50mm
Exhaust unit power:
Single phase AC220V 0.15A

Internal chamber



Improved thermal insulation features and temperature distribution accuracy by adopting a reasonable thermal insulation structure.

External temp. output terminal (optional)





- Time-up output
 - Output at auto stop and program end
 - Measured temperature transmission output
- For recording and monitoring internal chamber temperature.
Temperature output: 4~20mA

Sample tray (optional)



Dimension diagram (mm)



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Muffle Furnaces | Built-in Heating Wires

FP111C/311C/511C

Operating temp. range 100~1150°C

Temp. adjusting accuracy ±1.5°C

Internal capacity 1.5L 7.5L 11.3L

High-function muffle furnace with built-in heater.

Features

- High-quality ceramic furnace body, built-in heater, not contaminated by samples. Program operation is available.
- Simple operation with up/down key settings, equipped with a high-operability controller for precise temperature control.
- Temperature, time and indicated values can be digitally displayed.
- Use a controller with integrated digital settings and overheat protection. Support self-diagnosis functions, overcurrent leakage protection switch, and other safety controls in case of anomalies.
- Use R-type thermocouple for long life.
- Equipped with exhaust vent.
- Optional features available, such as exhaust unit, sample trays, nitrogen gas introduction device (with flow meter), temperature output terminal and so on.



Specifications

| Model | | FP111C | FP311C | FP511C |
|----------------|---|--|--------------|--------------|
| Performance | Operating temp. range | 100~1150°C | | |
| | Temp. adjusting accuracy | ±1.5°C (at 1150°C) | | |
| | Time to reach maximum temp. at room temp. | Approx. 80 min | | |
| Composition | Exterior material | Cold rolled steel plate with chemical proofing coating | | |
| | Furnace body | Ceramic fiber | | |
| | Sensors | Dual R-type thermocouple | | |
| | Heater | FeCrAl alloy heating wire | | |
| | | 1.1KW | 2.4KW | 3.5KW |
| | Exhaust vent | Inner diameter 20mm, upper part | | |
| Controllers | Cooling fan blade | Axial flow fan | | |
| | Temp. control method | Multi-PID control | | |
| | Temp. setting/display | Digital setting/digital display | | |
| | Overheating protection | Integrated controller | | |
| | Overheat protection setting method | Digital setting | | |
| | Operation functions | Fixed temp. operation, program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3), auto start and auto stop operation | | |
| | Additional functions | Deviation correction function, power failure compensation function, key lock function | | |
| | Timer | 1 min 99 h 59 min and 100 h~999 h 50 min | | |
| Safety device | | Self-diagnostic circuit (temperature sensor abnormality, automatic overheat prevention), overheat preventer, leakage protection circuit breaker for overcurrent, key lock function | | |
| Specifications | Internal dimensions (W×D×H mm) | 100×150×100 | 200×250×150 | 300×250×150 |
| | External dimensions (W×D×H mm) | 376×405×516 | 446×505×566 | 507×505×626 |
| | Internal capacity | 1.5L | 7.5L | 11.3L |
| | Power supply (50/60Hz) rated current | AC220V 5.5A | AC220V 11A | AC220V 16A |
| | Weight | Approx. 24kg | Approx. 42kg | Approx. 48kg |
| Accessories | | Exhaust vent cover 1 pcs | | |
| Options | Stand | ON30C | | ON61C |
| | Others | Setter plate, exhaust unit (AC220V), sample box, N ₂ introduction device (with flow meter), mini printer, data logger, combination warning light (standby/operation/failure), external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, central monitoring software | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices


15 Options

Automatic
Overheating
PreventerOverheat
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Control panel



Temperature rise curve (FP311C)



Exhaust unit (optional)



Able to quickly exhaust gases generated by rising temperature inside the chamber.

Piping: flexible pipe, aluminum

Length 1.5m

Diameter 50mm

Exhaust unit power:


Single phase AC220V 0.15A

Internal chamber




Heater not exposed, improved thermal insulation features and temperature distribution accuracy by adopting a reasonable thermal insulation structure.


Dimension diagram (mm)



[FP111C]



[FP311C]



[FP511C]

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Muffle Furnaces | 1700°C

FQ310C

Operating temp. range 700 ~1700°C

Temp. adjusting accuracy ±4°C

Internal capacity 10L

Ultra-high temperature muffle furnace with higher control precision and program operation function.

Features

- 7-inch touch screen for easy operation with three operation modes, including fixed temp. operation, program operation, time start/stop with standard 3-level authorization login.
- 100 programs, each program maximum 30 segments. Standard data storage function, data storage capacity of approx. 6 years. Temperature rise curve tracking function, program operation mode, able to visually display set temperature and actual temperature curve, providing reference for users to set reasonable temperature rise curve.
- The temperature sensor uses a long-life B-type thermocouple.
- Equipped with an exhaust vent.
- Optional features include an exhaust unit, argon gas introduction device (with flow meter) and temperature output terminal.

Safety

- Safety features include self-diagnostic circuit (temperature sensor abnormality, heater disconnection, automatic overheating prevention, SCR short circuit), snap-action temperature switch (outside the furnace body), leakage protection circuit breaker to prevent overcurrent and 3-level authorization login function.



Specifications

| Model | | FQ310C |
|----------------|---|--|
| Performance | Operating temp. range | 700~1700°C |
| | Temp. adjusting accuracy | ±4°C (at 1200°C) |
| | Time to reach maximum temp. at room temp. | Approx. 75 min to 1600°C, approx.95 min to 1700°C |
| Composition | Exterior material | Cold rolled steel plate with chemical proofing coating |
| | Furnace body | Ceramic fiber |
| | Sensors | B-type thermocouple |
| | Heater | Silicon-molybdenum rod 6KW |
| | Exhaust vent | Inner diameter 20mm, upper part |
| | Cooling fan blade | Axial flow fan |
| Controllers | Temp. control method | Multi-PID control |
| | Temp. setting/display | Settings using a 7-inch color touch screen |
| | Overheating protection | Snap-action temperature switch (outside the furnace body) |
| | Operation functions | Fixed temp. operation, program operation, auto start/stop |
| | Additional functions | Temperature history curve, alarm history records, operation log, data export function, language switching (Chinese, Japanese, English), calibration function, power outage compensation function, screen brightness adjustment, three-level authorization function, micro printer (optional) |
| Specifications | Timer | 1 min~99 h 59 min |
| | Safety device | Ceiling plate inspection switch, temperature overheating protection (thermostat), snap-action temperature switch (external of furnace), heater overcurrent protection, heater disconnection alarm, sensor disconnection, SCR or temperature controller communication interruption, SCR core temperature extra overheat |
| | Internal dimensions (W×D×H mm) | 200×250×200 |
| | External dimensions (W×D×H mm) | 650×703×1425 |
| | Internal capacity | 10L |
| Accessories | Power supply (50/60Hz) rated current | AC220V 29A |
| | Weight | Approx. 130kg |
| Options | Others | Exhaust unit (AC220V), argon gas introduction device (with flow meter), mini printer, data logger, combination warning light (standby/operation/failure), external communication function (network port), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal |

*For long-term operation, keep the temperature 100°C below the maximum working temperature.

*A chemical reaction between the heating element and zirconia may cause the zirconia sample to discolor. We can provide process guidance or replace with other models of heating elements. Please consult us for details.

Automatic
Overheating
PreventerOverheat
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Exhaust device unit (optional)



Able to quickly exhaust gases generated by rising temperature inside the chamber.

Piping: flexible pipe, aluminum

Length 1.5m

Diameter 50mm


Exhaust unit power:

Single phase AC220V 0.15A

Control panel




Temperature rise characteristic curve (FQ310C)



Internal chamber



Dimension diagram (mm)



[FQ310C]

Sterilizers 1

Granulation
and Spray
Dryers 2

Muffle
Furnaces 3

Ovens 4

Incubators 5

Plasma
Equipment 6

Water
Purifiers 7

Baths 8

Water
Circulators 9

Rotary
Evaporators 10

Freeze Dryers
& Cold Traps 11

Stirrers &
Shakers 12

Washers 13

Analysis and
Test Devices 14

Options 15

Related products

Crucible

| Product name | Specifications/quantity | Dimensions |
|----------------|-------------------------|----------------------------------|
| Crucible | SSA-H B1/1 pcs | Outer Diameter 46×Height 36 (mm) |
| | SSA-H B2/1 pcs | Outer Diameter 54×Height 43 (mm) |
| | SSA-H B3/1 pcs | Outer Diameter 66×Height 54 (mm) |
| | SSA-H B4/1 pcs | Outer Diameter 73×Height 59 (mm) |
| | SSA-H B5/1 pcs | Outer Diameter 88×Height 72 (mm) |
| Crucible Lid | SSA-H/1 pcs | For B1 |
| | SSA-H/1 pcs | For B2 |
| | SSA-H/1 pcs | For B3 |
| | SSA-H/1 pcs | For B4 |
| | SSA-H/1 pcs | For B5 |
| Crucible Tongs | 1 pcs | 180mm |
| | 1 pcs | 210mm |
| | 1 pcs | 240mm |
| | 1 pcs | 300mm |
| | 1 pcs | 360mm |
| | 1 pcs | 450mm |
| | 1 pcs | 500mm |
| | 1 pcs | 600mm |



Crucible/Crucible Lid



Crucible Tongs



Combustion Boat



Damar Resin Tube

Combustion boat/damar resin tube

| Product name | Specifications/quantity | Dimensions |
|------------------|-------------------------|---|
| Combustion Boat | #01/50 pcs | Width 6×Height 4×Length 30 (mm) |
| | #10/10 pcs | Width 14×Height 12×Length 140 (mm) |
| | #8/50 pcs | Width 15×Height 13×Length 90 (mm) |
| Damar Resin Tube | ACE-T2/1 pcs | Outer Diameter 15×Inner Diameter 12×Length 100 (mm) |
| | ACE-T3/1 pcs | Outer Diameter 21×Inner Diameter 17×Length 100 (mm) |
| | ACE-T4/1 pcs | Outer Diameter 25×Inner Diameter 20×Length 120 (mm) |
| | ACE-T5/1 pcs | Outer Diameter 30×Inner Diameter 24×Length 150 (mm) |
| | ACE-T6/1 pcs | Outer Diameter 40×Inner Diameter 34×Length 150 (mm) |
| | ACE-T7/1 pcs | Outer Diameter 50×Inner Diameter 40×Length 150 (mm) |
| | ACE-T8/1 pcs | Outer Diameter 60×Inner Diameter 50×Length 180 (mm) |

Furnace plate/Heat-resistant gloves

| Product name | Specifications/quantity | Dimensions |
|-----------------------|-------------------------|---|
| Furnace Bed Plate | SSA-S/1 piece | Length 50×Width 50×Thickness 2.5 (mm) |
| | SSA-S/1 piece | Length 100×Width 100×Thickness 2.5 (mm) |
| | SSA-S/1 piece | Length 150×Width 150×Thickness 2.5 (mm) |
| Heat-resistant Gloves | EGM48/1 pair | 2-finger, Length 350mm, Heat-resistant Temperature 1000°C |
| | EGT27/1 pair | 3-finger, Length 350mm, Heat-resistant Temperature 1000°C |
| | EGT47/1 pair | 5-finger, Length 350mm, Heat-resistant Temperature 1000°C |

Setter Plate

| Product name | Specifications/quantity | Dimensions |
|--------------|------------------------------------|--|
| Setter Plate | FP111C_01_02-01(140×90×t8)/1 pcs | Length 140×Width 90×Thickness T8 (mm) |
| | FO210C_01_02-20(240×90×t8)/1 pcs | Length 240×Width 90×Thickness T8 (mm) |
| | FO310C_01_02-17(240×190×t8)/1 pcs | Length 240×Width 190×Thickness T8 (mm) |
| | FO410C_01_02-07(290×190×t8)/1 pcs | Length 290×Width 190×Thickness T8 (mm) |
| | FO510C_01_02-07(290×240×t8)/1 pcs | Length 290×Width 240×Thickness T8 (mm) |
| | FO610C_01_02-06(340×240×t9)/1 pcs | Length 340×Width 240×Thickness T9 (mm) |
| | FO710C_01_02-07(340×260×t9)/1 pcs | Length 340×Width 260×Thickness T9 (mm) |
| | FO810COP01_01-02(360×290×t9)/1 pcs | Length 360×Width 290×Thickness T9 (mm) |



Furnace Plate Corner Plate



Heat-resistant Gloves

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Contents

| | |
|--------------------------------|-----|
| Forced Convection Ovens | 62 |
| Natural Convection Ovens | 80 |
| Inert Ovens | 86 |
| Explosion-proof Ovens | 89 |
| Ovens for Labware | 91 |
| Clean Ovens | 92 |
| Far Infrared Ray Ovens | 98 |
| High Temp. Ovens | 99 |
| Vacuum Ovens | 100 |



Provide safe, environment friendly and energy-saving products.

Multiple choices are available due to various models with different purposes. Meet the demands of temperature range, temperature distribution accuracy, size, price, program operation, special use type and so on.

Function · Safety Device

Forced convection circulation

By means of forced circulation of high-performance centrifugal fan blades, stir the chamber inside to obtain uniform temperature distribution accuracy.

Natural convection

By means of heater or air jacket natural convection in chamber.

Auto overhear prevention function

Built-in overhear prevention function of controller, usually when constant temp. oven or drying oven chamber temp. reaches set temp. + 12°C, and constant temp. bath temp. reaches set temp. + 6°C, the heater cuts off (auto recovery).

Overheat protector

Controller and overheat protector are all-in-one. They share the same power source, but other circuits are standalone. If overheat abnormally, cut off (manual recovery) heater circuit.

Standalone overheat protector

Different circuit from controller, if overheat abnormally, cut off (manual recovery) heater circuit. According to different models, respectively carry digital, liquid-expansion, thermometal type, etc.

Self-diagnosis function

By means of the microcomputer carried in controller, detect the abnormality of circuit, if abnormal occurs, control the device within safety range and the alarm sounds.

Key lock function

Prevent the misoperation during operation.

Backup function

Even in power outage or cut off power source, the set value can be memorized.

Power outage compensation

When the power recovers, it's able to select to stop or continue operation.

Overcurrent breaker

When abnormal current passes, cut off the power to protect the unit.

Overcurrent electric leakage breaker

The breaker has both overcurrent cut-off function and electric leakage cut-off function.

Emergency stop button

If need to emergency stop, press the emergency stop button to cut off.

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Model selection list

| Category | Max. temp. | Temp. fluctuation (°C) | Temp. deviation (°C) or temp. uniformity (%) | Series | Models | Features | Internal capacity (L) | | | | | | | | Page No. |
|--------------|-----------------------|------------------------|--|--------|------------------------|--|---|-------|--------|---------|---------|---------|---------|----------|----------|
| | | | | | | | 0-20 | 20-50 | 50-100 | 100-150 | 150-200 | 200-400 | 400-600 | 600-1000 | |
| General type | 500°C | ±0.5°C | ±6°C | DH | DH450C/650C | Precision type, 500°C | | | 91 | | | 216 | | | 74 |
| | | ±0.5°C | ±8°C | DH | DH850C | Precision type, 500°C, large capacity | | | | | | | 512 | | 74 |
| | 360°C | ±0.2°C | ±4°C | DH | DH411C/611C | Precision type, horizontal | | | 91 | | | 216 | | | 70 |
| | | ±0.2°C | ±5°C | DHS | DHS710C/810 C | Precision type, vertical | | | | | | 418 | 558 | | 73 |
| | | ±0.1°C | ±1.5% | DN-H | DN410HC/610HC | High temperature type | | | 95 | | | 223 | | | 69 |
| | 300°C | ±0.3°C | ±4°C | DH | DH811C/1011C | Precision type, horizontal, large capacity | | | | | | | 512 | 1000 | 72 |
| | 300°C | ±0.5°C | ±1.5% | DKT | DKT410C/610C/810C | High-function type, touch screen | | | 90 | | 150 | 300 | | | 62 |
| | 260°C | ±0.5°C | ±2% | DKN | DKN313C/413C/613C | Standard type, program operation | | 27 | 90 | | 150 | | | | 64 |
| | | ±0.5°C | ±2% | DKM | DKM311C/411C/611C | Economical type, fixed temp. operation | | 27 | 90 | | 150 | | | | 66 |
| | | ±0.5°C | ±2% | DKL | DKL311C/411C/611C | Basic function type | | 27 | 90 | | 150 | | | | 68 |
| | | ±0.5°C | ±1.5% | DNF | DNF411C/611C/811C/911C | Energy-saving type, variable wind speed | | | 90 | | 150 | 300 | 540 | | 78 |
| | | - | - | DKG | DKG611/611V/811/811V | Labor-saving, rapid heating and cooling | | | | | 150 | 300 | | | 75 |
| | | ±0.1°C | ±4°C | DFS | DFS710C/810C | Precision type, vertical | | | | | | 418 | 558 | | 73 |
| | | ±0.1°C | ±3°C | DF | DF411C/611C | Precision type | | | 91 | | | 216 | | | 70 |
| | 210°C | ±0.5°C | ±1.5% | DNE | DNE411C/611C/811C/911C | Energy-saving type | | | 90 | | 150 | 300 | 540 | | 76 |
| | | ±0.5°C | ±2.5% | DKN | DKN813C/913C | Standard type, program operation | | | | | | 300 | 535 | | 64 |
| | 200°C | ±0.2°C | ±2°C | DF | DF811C/1011C | Precision type, horizontal, large capacity | | | | | | | 512 | 1000 | 72 |
| | Natural convection | ±1°C | ±3% | DX | DX312C/412C | Economical type, fixed temp. operation | | 28 | 74 | | | | | | 84 |
| | | ±1°C | ±3% | DY | DY311C/411C | Basic function type | | 28 | 74 | | | | | | 85 |
| | | ±1°C | ±3% | DVT | DVT410C/610C | Natural convection, high-function type, touch screen | | | | 117 | 188 | | | | 80 |
| | | ±1°C | ±3% | DX | DX612C | Natural convection, economical type, fixed temp. operation | | | | 153 | | | | | 84 |
| | | ±1°C | ±3% | DY | DY611C | Basic function type | | | | 153 | | | | | 85 |
| | | ±1°C | ±3% | DVS | DVS413C/613C | Natural convection, standard type | | | 99 | | 162 | | | | 82 |
| | | ±1°C | ±3% | DVT | DVT410C/610C | Natural convection, high-function type, touch screen | | | | 117 | 188 | | | | 80 |
| Special type | Anaerobic | 600°C | ±0.5°C | ±8°C | DNN | DNN460C/660C | N ₂ replacement, oxygen concentration 20ppm | | | 91 | | 261 | | | 88 |
| | | | ±0.2°C | ±1.5% | DN-I | DN410IC/610IC | N ₂ replacement, heating in an oxygen-free environment | | | 95 | | 223 | | | 86 |
| | | 360°C | ±0.3°C | ±5°C | DNN | DNN430C/630C | N ₂ replacement, oxygen concentration 20ppm | | | 91 | | 261 | | | 88 |
| | | | ±0.2°C | ±5°C | DH-I | DH611C | Anaerobic, oxygen concentration 5000ppm | | | | | 216 | | | 87 |
| | | 260°C | ±0.1°C | ±4°C | DF-I | DF611C | Anaerobic, oxygen concentration 5000ppm | | | | | 216 | | | 87 |
| | Explosion-proof | 500°C | ±0.5°C | ±8°C | DH-S | DH450SC/650SC | Explosion-proof, 500°C | | | 91 | | 216 | | | 90 |
| | | | ±0.5°C | ±10°C | DH-S | DH850SC | Explosion-proof, 500°C | | | | | | 512 | | 90 |
| | | 360°C | ±0.3°C | ±5°C | DH-S | DH411SC/611SC | Explosion-proof, precision type | | | 91 | | 216 | | | 89 |
| | | 300°C | ±0.3°C | ±5°C | DH-S | DH811SC/1011SC | Explosion-proof, precision type, large capacity | | | | | | 512 | 1000 | 90 |
| | | 260°C | ±0.2°C | ±4°C | DF-S | DF411SC/611SC | Explosion-proof, precision type | | | 91 | | 216 | | | 89 |
| | For labware | 200°C | ±0.2°C | ±3°C | DF-S | DF811SC/1011SC | Explosion-proof, precision type, large capacity | | | | | | 512 | 1000 | 90 |
| | | 70°C | ±0.5°C | ±5°C | DG | DG410C/450C/810C/850C | Natural convection, forced exhaust | | | 92 | | | 445 | | 91 |
| | Clean | 500°C | ±0.5°C | ±8°C | DTN | DTN450C/650C | Cleanliness 100 class, oxygen concentration 20ppm | | | 91 | | 261 | | | 97 |
| | | | ±0.3°C | ±1.5% | DT | DT430C/430UC/630C/630UC | Heating in a dust-free environment | | | 91 | | 216 | | | 92 |
| | | 360°C | ±0.5°C | ±1.5% | DTS | DTS830C/830UC | Heating in a dust-free environment, large capacity | | | | | 300 | | | 94 |
| | | | ±0.3°C | ±5°C | DTN | DTN430C/630C | Cleanliness 100 class, oxygen concentration 20ppm | | | 91 | | 261 | | | 97 |
| | | 260°C | ±0.3°C | ±1.5% | DE | DE430C/430UC/630C/630UC | Heating in a dust-free environment | | | 91 | | 216 | | | 92 |
| | | | ±0.5°C | ±1.5% | DES | DES830C/830UC | Heating in a dust-free environment, large capacity | | | | | 300 | | | 94 |
| | | 150°C | ±0.5°C | ±3.0% | DEC | DEC812C/912C | Heating in a dust-free environment, economical type | | | | | 236 | 535 | | 96 |
| | Far infra-red heating | 360°C | ±0.2°C | ±5°C | DIR | DIR631C | Far infrared heating | | | | | 216 | | | 98 |
| | High temp. | 700°C | ±2°C | ±5% | DR | DR210C | Temperature 700°C | 13.8 | | | | | | | 99 |
| | Vacuum | 400°C | ±0.5°C | ±1.0°C | DPHH | DPHH43C/63C | Vertical, 400°C | | | 91 | | 216 | | | 104 |
| | | 300°C | ±0.5°C | ±1.0°C | DPH | DPH43C/63C | Vertical, 300°C | | | 91 | | 216 | | | 103 |
| | | 240°C | ±0.5°C | ±1.5°C | DP | DP23C/33C | Standard type | 10 | 27 | | | | | | 101 |
| | | | ±0.5°C | ±1.5°C | ADP | ADP210C/310C | Desktop, economical type | 10 | 27 | | | | | | 102 |
| | | 200°C | ±0.5°C | ±2.5% | DPF | DPF43C/63C | Vertical, high precision | | | 91 | | 216 | | | 103 |
| | | | ±0.5°C | ±1.0°C | DP-HP | DP43HPC/63HPC | Vertical, shelf heating | | | 91 | | 216 | | | 104 |
| | | | ±0.5°C | ±1.0°C | DPT | DPT40C/60C/80C/100C | High-function type, touch screen | | | 91 | | 216 | 512 | 1000 | 100 |
| | | | ±0.5°C | ±1.0°C | DP | DP43C/63C/83C/103C | Standard type | | | 91 | | 216 | 512 | 1000 | 101 |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Forced Convection Ovens | Programmable, Touch Screen

DKT410C/610C/810C

Operating temp. range RT+10~300°C

Temp. uniformity $\pm 4.5^{\circ}\text{C}$ (at 300°C)

Internal capacity 90L 150L 300L

Pioneer in standard high-end constant temperature ovens with multiple functions, including program operation.

Features

- 7-inch touch screen for easy operation with three operation modes, including fixed temp. operation, program operation, timing operation. Standard 3-level authorization login.
- 100 programs, each program maximum 30 segments. Standard data storage function, data storage capacity approximately 10 years. Temperature rise curve tracking function, program operation mode, able to visually display set temperature and achieved temperature, providing reference for users to set reasonable temperature rise curve.
- Automatically stop heating and fan rotation upon opening the door, ensuring the safety of user samples.

Safety

- Self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), independent overheating prevention, overcurrent leakage protection and three-level authorization management system, heating and fan stops operating when door opens.




Specifications

| Model | | DKT410C | DKT610C | DKT810C |
|----------------|--|-------------------|---|--------------------------------|
| System | | Forced convection | | |
| Performance | Operating temp. range | | Room temp. +10~300°C | |
| | GB standard | Temp. fluctuation | $\pm 0.5^{\circ}\text{C}$ (at 300°C) | |
| | | Temp. uniformity | $\pm 1.5\%$ (at 300°C) | |
| | Max. temp. reaching time | | $\leq 120\text{min}$ | |
| Composition | Interior material | | Stainless steel plate | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | |
| | Insulating material | | Rock wool+centrifugal glass wool | |
| | Heater | | Stainless steel heating pipe | |
| | | | 1.6KW | 1.8KW $\times 2$ |
| | Fan blade/motor | | Centrifugal fan blade, 30W motor $\times 1$ | |
| | Cable port | | I.D. 34mm $\times 1$, located on the right | |
| | Exhaust vent | | I.D. 34 mm $\times 2$, located at the back | |
| Controllers | Observation window | | None | |
| | Temp. control method | | 3-stage PID | |
| | Temp. setting method | | 7-inch color touch screen | |
| | Temp. display method | | 7-inch color touch screen | |
| | Timer | | 1 min~99 h 59 min | |
| | Operation functions | | Fixed temp. operation, program operation, timing operation | |
| | Program mode | | 100 programs, each program has a maximum of 30 segments | |
| | Additional functions | | Temperature curve, temperature history curve, alarm records, alarm history records, operation log, data export function, language switching (Chinese, Japanese, English), calibration function, power failure compensation function, screen brightness adjustment, three-level authorization function | |
| Specifications | Sensors | | Temperature controller: K-type thermocouple, overheat protection: Liquid expansion temperature controller | |
| | Safety device | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), independent overheating preventer, overcurrent leakage protection, and heating and fan stops operating when door opens | |
| | Internal dimensions (W \times D \times H mm) | | 450 \times 450 \times 450 | 600 \times 500 \times 500 |
| | External dimensions (W \times D \times H mm) | | 610 \times 690 \times 960 | 760 \times 740 \times 1010 |
| | Internal capacity | | 90L | 150L |
| | Shelf load | | 30kg/layer | 300L |
| | Shelf layers | | 11 layers | 13 layers |
| | Support spacing | | 30mm | 29 layers |
| Accessories | Power supply (50Hz) rated current | | AC220V, 7.5A | AC220V, 12A |
| | Weight | | Approx. 80kg | Approx. 100kg |
| | Shelf | | Stainless steel wire shelf | |
| | Supports | | 2 pcs | 4 pcs |
| | | | 4 pcs | 8 pcs |
| | Stand | | ONK60C | |
| | Stacking fittings | | ODK40C | — |
| | Others | | Shelf (1 shelf with 2 supports) cable port (30/50mm), micro printer, data logger, combination warning light (standby/running/fault), remote monitoring (network port), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal | |

Forced
ConvectionAutomatic
Overheating
PreventerIndependent
Overheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Temperature rise curve



Control panel



Internal chamber



Shelf (Support)



Cable port (standard configuration)




Exhaust vent (standard configuration)




Dimension diagram (mm)



[DKT410C]



[DKT610C]



[DKT810C]

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Forced Convection Ovens | Standard, Program Operation

DKN313C/413C/613C/813C/913C

Operating temp. range RT+10~260°C

Temp. distribution accuracy $\pm 2.8^{\circ}\text{C}$ (at 260°C), $\pm 2.5^{\circ}\text{C}$ (at 210°C)

Internal capacity 27L 90L 150L 300L 535L

Pioneer in standard constant temperature ovens with multiple functions, including program operation.

Features

- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with 90 segments program controller.
- Deviation correction and key locking through submenu keys.

Safety

- Self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), independent overheating preventer, overcurrent leakage protection switch, and key lock safety functions.



Specifications

| Model | | | DKN313C | DKN413C | DKN613C | DKN813C | DKN913C |
|--|--------------------------------------|-----------------------------|--|--------------|--------------|-----------------------------------|---------------------|
| System | | | Forced convection | | | | |
| Performance | Operating temp. range | | Room temp. +10~260℃ | | | Room temp. +10~210℃ | |
| | GB standard | Temp. fluctuation | ±0.5℃ (at 260℃) | | | ±0.5℃ (at 210℃) | |
| | | Temp. uniformity | ±2% (at 260℃) | | | ±2.5% (at 210℃) | |
| | JTM standard | Temp.adjusting accuracy | ±1.0℃ (at 210℃) | | | | |
| | | Temp. distribution accuracy | ±2.5℃ (at 210℃) | | | | |
| Max. temp. reaching time | | | Approx. 60 min | | | | |
| Interior material | | | Stainless steel plate | | | | |
| Exterior material | | | Cold rolled steel plate with chemical proofing coating | | | | |
| Insulating material | | | Glass fiber+rock wool | | | | |
| Heater | | | Stainless steel heating pipe | | | | |
| | | | 0.8KW | 1.2KW | 1.5KW | | 1.8KW |
| Fan blade/motor | | | Centrifugal fan blade, 10W motor×1 | | | Motor: 30W×1 | Motor: 10W×2 |
| Cable port | | | Inner diameter: 30mm×1, located on the right | | | | |
| Exhaust vent | | | I.D. 30mm×2, located on the top | | | I.D. 30 mm×2, located at the back | |
| Observation windows (width×height in mm) | | | 180×180 | 250×280 | | None | |
| Controllers | Temp. control method | | 3-stage PID | | | | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | | | | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display | | | | |
| | | | Setting temp. display: Red 4-digit LED digital display | | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | | |
| | Operation functions | | Fixed temperature operation, program operation, quick auto stop operation, auto stop operation, auto start operation, etc. | | | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | | | |
| | Temp. control sensor | | K-type sensor | | | | |
| Overheat sensor | | | Liquid expansion type (independent circuit for safety) | | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), independent overheating preventer, overcurrent leakage protection, and key lock functions | | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 300×300×300 | 450×450×450 | 600×500×500 | 600×500×1000 | 1070×500×1000 |
| | External dimensions (W×D×H mm) | | 410×451×677 | 560×601×827 | 710×651×877 | 710×651×1598 | 1180×651×1606 |
| | Internal capacity | | 27L | 90L | 150L | 300L | 535L |
| | Shelf load | | 15kg/layer | | | | |
| | Shelf layers | | 6 layers | 11 layers | 13 layers | 29 layers | 29 layers×2 columns |
| | Support spacing | | 30mm | | | | |
| | Power supply (50/60Hz) rated current | | AC220V, 4A | AC220V, 6A | AC220V, 7.5A | AC220V, 14.5A | AC220V, 17A |
| | Weight | | Approx. 35kg | Approx. 50kg | Approx. 65kg | Approx. 110kg | Approx. 190kg |
| Accessories | Shelf | | Stainless punching mesh plate | | | | |
| | | | 2 pieces, 1 of which is fixed to the chamber; do not disassemble | | | 4 pcs | 8 pcs |
| Options | Supports | | 2 pcs | | | 8 pcs | 16 pcs |
| | Stand | | ON30C | ON61C | | — | |
| | Stacking fittings | | — | OD40C | OD60C | — | |
| | Others | | Shelf (1 shelf with 2 supports) , cable port (30/50mm), micro printer, data logger combination, warning light (standby/running/fault), external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software | | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification

Control panel



Internal chamber



DKN613C



DKN913C


Cable port (standard configuration)



Exhaust vent (standard configuration)



Dimension diagram (mm)



Sterilizers 1

Granulation
and Spray
Dryers 2Muffle
Furnaces 3

Ovens 4

Incubators 5

Plasma
Equipment 6Water
Purifiers 7

Baths 8

Water
Circulators 9Rotary
Evaporators 10Freeze Dryers
& Cold Traps 11Stirrers &
Shakers 12

Washers 13

Analysis and
Test Devices 14

Options 15

Forced Convection Ovens | Economical, Fixed Temp. Operation

DKM311C/411C/611C

Operating temp. range RT+10~260°C

Temp. distribution accuracy $\pm 2.5^{\circ}\text{C}$ (at 210°C)

Internal capacity 27L 90L 150L

Economical model for fixed temperature operation, a simple-to-operate forced convection constant temperature oven.

Features

Operability & functions

- Easy to operate, fixed temperature operation, quick auto stop operation, auto stop operation, and auto start operation.
- Digital setting through special operation function menu keys and up/down keys. Able to set overheating protection, deviation correction, and key lock.

Safety

- Equipped with self-diagnostic circuits (temperature sensor anomaly), power failure compensation function, deviation correction function, leakage protection against overcurrent, independent overheating preventer and other safety features.



Specifications

| Model | | | DKM311C | DKM411C | | DKM611C |
|--------------------------|--------------------------------------|---|--|--------------|--------------|---------|
| System | | | Forced convection | | | |
| Performance | Operating temp. range | | Room temp. +10~260℃ | | | |
| | GB standard | Temp. fluctuation | ±0.5℃ (at 260℃) | | | |
| | | Temp. uniformity | ±2%℃ (at 260℃) | | | |
| | JTM standard | Temp. adjusting accuracy | ±1.0℃ (at 210℃) | | | |
| | | Temp. distribution accuracy | ±2.5℃ (at 210℃) | | | |
| Max. temp. reaching time | | | Approx. 60 min (room temp 10~210℃), approximately 90 min (room temp 10~260℃) | | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber+rock wool | | | |
| | Heater | Stainless steel heating pipe | | | | |
| | | 0.8KW | 1.2KW | 1.5KW | | |
| | Fan blade/motor | | Centrifugal fan blade, 10W motor×1 | | | |
| | Cable port | | Inner diameter: 30mm×1, located on the right | | | |
| | Exhaust vent | | I.D. 30mm×2, located on the top | | | |
| Controllers | Temp. control method | | PID control | | | |
| | Temp. setting method | | Digital setting through function menu key and up/down keys | | | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display | | | | |
| | | Setting temp. display: Red 4-digit LED digital display | | | | |
| | Timer | | 1 min.~99 h 59 min. and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, quick auto stop | | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | | |
| | Sensors | | Temperature controller: K-type thermocouple, overheating protection: Liquid expansion temperature controller | | | |
| Safety device | | | Self-diagnostic function (temperature sensor abnormality, SSR short circuit, automatic overheating preventer), key lock, independent overheating preventer, leakage protection circuit breaker for overcurrent | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 300×300×300 | 450×450×450 | 600×500×500 | |
| | External dimensions (W×D×H mm) | | 410×451×677 | 560×601×827 | 710×651×877 | |
| | Internal capacity | | 27L | 90L | 150L | |
| | Shelf load | | 15kg/layer | | | |
| | Shelf layers | | 6 layers | 11 layers | 13 layers | |
| | Support spacing | | 30mm | | | |
| | Power supply (50/60Hz) rated current | | AC220V, 4A | AC220V, 6A | AC220V, 7.5A | |
| | Weight | | Approx. 35kg | Approx. 50kg | Approx. 65kg | |
| Accessories | Shelf | | Stainless punching mesh plate 2 pieces, 1 of which is fixed to the chamber; do not disassemble | | | |
| | Supports | | 2 pcs | | | |
| Options | Stand | | ON30C | ON61C | | |
| | Stacking fittings | | — | OD40C | OD60C | |
| | Others | | Shelf (1 shelf plate with 2 shelf supports), cable port (30/50mm) | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification

Control panel



Cable port (standard configuration)



Exhaust vent (standard configuration)




Internal chamber






DKM611C

Temperature rise curve



Dimension diagram (mm)



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Forced Convection Ovens | Basic Function, Economical

DKL311C/411C/611C

Forced
ConvectionAutomatic
Overheating
PreventerIndependent
Overheating
PreventerSelf-diagnostic
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range RT+10~260°C

Temp. distribution accuracy $\pm 2.5^{\circ}\text{C}$ (at 260°C)

Internal capacity 27L 90L 150L

Ensuring the basic functions of fixed temperature operation models.



Various functional settings allow for simpler operations in forced convection constant temperature ovens with fixed operation.

Features

- Easy to operate, fixed temp. operation, auto stop operation, quick auto stop operation, and auto start operation.

Safety

- Equipped with self-diagnostic circuits (temperature sensor anomaly, temperature overheating prevention), leakage protection circuit breaker for overcurrent for overcurrent, and independent overheating preventer and other safety functions.

Specifications

| Model | | | DKL311C | DKL411C | DKL611C |
|--------------------------|--------------------------------------|-----------------------------|---|--------------|--------------|
| System | | | Forced convection | | |
| Performance | Operating temp. range | | Room temp. +10~260℃ | | |
| | GB standard | Temp. fluctuation | ±0.5℃ (at 260℃) | | |
| | | Temp. uniformity | ±2% (at 260℃) | | |
| | JTM standard | Temp. adjusting accuracy | ±1.0℃ (at 260℃) | | |
| | | Temp. distribution accuracy | ±2.5℃ (at 260℃) | | |
| Max. temp. reaching time | | | Approx. 90 min (room temp + 10~260℃) | | |
| Composition | Interior material | | Stainless steel plate | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | |
| | Insulating material | | Glass fiber | | |
| | Heater | | Stainless steel heating pipe | | |
| | | | 0.8KW | 1.2KW | 1.5KW |
| | Fan blade/motor | | Centrifugal fan blade, 10W motor×1 | | |
| | Cable port | | Inner diameter: 30mm×1, located on the right | | |
| Exhaust vent | | | I.D. 30mm×2, located on the top | | |
| Controllers | Temp. control method | | PID control | | |
| | Temp. setting method | | Digital setting through function menu key and up/down keys | | |
| | Temp. display method | | Achieved temperature display: red 4-digit led, digital display | | |
| | | | Set temperature display: green 4-digit led, digital display | | |
| | Timer | | 1 min to 99 h~50 min | | |
| | Operation functions | | Fixed temperature operation, auto stop operation, quick auto stop operation, auto start operation | | |
| | Additional functions | | Deviation correction function | | |
| Sensors | | | Temperature controller: K-type thermocouple Overheating Protection: Liquid expansion temperature controller | | |
| Safety device | | | Self-diagnostic function (temperature sensor abnormal, temperature overheat prevention), independent overheating prevention, leakage protection circuit breaker for overcurrent | | |
| Specifications | Internal dimensions (W×D×H mm) | | 300×300×300 | 450×450×450 | 600×500×500 |
| | External dimensions (W×D×H mm) | | 410×451×680 | 560×601×830 | 720×651×880 |
| | Internal capacity | | 27L | 90L | 150L |
| | Shelf load | | 15kg/layer | | |
| | Shelf layers | | 6 layers | 11 layers | 13 layers |
| | Support spacing | | 35mm | | |
| | Power supply (50/60Hz) rated current | | AC220V, 4A | AC220V, 6A | AC220V, 7.5A |
| | Weight | | Approx. 35kg | Approx. 50kg | Approx. 65kg |
| Accessories | Shelf | | Stainless punching mesh plate | | |
| | | | 2 pcs | | |
| | Supports | | 4 pcs | | |
| Options | Stand | | ON30C | ON61C | |
| | Stacking fittings | | — | OD40C | OD60C |
| | Others | | Shelf (1 shelf plate with 2 shelf supports), cable port (30/50mm) | | |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Forced Convection Ovens | High Temp.

DN410HC/610HC

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range RT+10~360°C

Temp. distribution accuracy $\pm 3^{\circ}\text{C}$ (at 360°C)

Internal capacity 95L 223L

Reliable constant temperature oven that supports high-temperature thermal tests.

Features

- Support performing high-temperature resistance tests and thermal processing at 360°C.
- Wide operating temperature range, high accuracy in temperature control.
- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with 90 segments program controller.


Safety

- Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), independent overheating preventer, overcurrent leakage protection, key lock and other safety functions.




Temperature rise curve


DN410HC



DN610HC



Dimension diagram (mm)



Specifications

| Model | | | DN410HC | DN610HC |
|--------------------------|--------------------------------------|--|--|----------------|
| System | | | Forced convection | |
| Performance | Operating temp. range | | Room temp. +10~360℃ | |
| | GB standard | Temp. fluctuation | ±0.1℃ (at 360℃) | |
| | | Temp. uniformity | ±1.5% (at 360℃) | |
| | JTM standard | Temp. adjusting accuracy | ±0.2℃ (at 360℃) | |
| | | Temp. distribution accuracy | ±3.0℃ (at 360℃) | |
| Max. temp. reaching time | | Approx. 60 min | | |
| Composition | Interior material | | Stainless steel plate | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | |
| | Insulating material | | Rock wool | |
| | Heater | | Stainless steel heating tube, 3.0KW / Stainless steel heating tube, 4.0KW | |
| | Fan blade/motor | | Centrifugal fan blades, high-temperature self-cooling motor 30W | |
| | Cable port | | Inner diameter: 30mm×1, located on the right | |
| Exhaust vent | | I.D. 30 mm×2, located at the back | | |
| Controllers | Temp. control method | | 3-stage PID | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display Setting temp. display: Red 4-digit LED digital display | |
| | Timer | | 1 min.~99 h 59 min. and 100~999 h 50 min (including timer waiting function) | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | |
| Sensors | | K thermocouple (Temp. controller and overheating protector) | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | |
| Specifications | Internal dimensions (W×D×H mm) | 470×450×450 | | 620×600×600 |
| | External dimensions (W×D×H mm) | 640×695×915 | | 790×846×1065 |
| | Internal capacity | 95L | | 223L |
| | Shelf load | 30kg/layer | | |
| | Shelf layers/shelf support spacing | 12 layers/30mm | | 17 layers/30mm |
| | Power supply (50/60Hz) rated current | AC220V 14A | | AC220V 19A |
| | Weight | Approx. 80kg | | Approx. 120kg |
| Accessories | Shelf | Stainless steel wire mesh plate 2 pcs | | |
| | Supports | 4 pcs | | |
| Options | Stand | OH41C | | OH61C |
| | Others | Shelf (1 shelf includes 2 supports) , cable port (30/50mm), micro printer, data logger, combination warning light (standby/running/fault), viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software | | |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Forced Convection Ovens | Precision, Horizontal

DF411C/611C, DH411C/611C

| | | | | | | | | |
|-----------------------|----------------|----------------|-----------------------------|---------------------|---------------------|-------------------|-----|------|
| Operating temp. range | DF: RT+10~260℃ | DH: RT+10~360℃ | Temp. distribution accuracy | DF: ±1.5℃ (at 260℃) | DH: ±2.5℃ (at 360℃) | Internal capacity | 91L | 216L |
|-----------------------|----------------|----------------|-----------------------------|---------------------|---------------------|-------------------|-----|------|

High-precision constant temperature ovens for thermal processing.



■ Features

- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with 90 segments program controller.
- Overheating protection, deviation correction and key locking through submenu keys.
- Quick exhaust and cooling through exhaust dampers.

■ Safety

- Equipped with self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), overheating preventer, overcurrent leakage protection, key lock and other safety functions.

■ Specifications

| Model | | | DF411C | DF611C | DH411C | DH611C |
|-----------------------------|--------------------------------------|---|--|--------------------------|--|---------------|
| System | | | Forced convection and ventilation | | | |
| Performance | Operating temp. range | | Room temp. +10~260℃ | | Room temp. +10~360℃ | |
| | GB standard | Temp. fluctuation | ±0.1℃ (at 260℃) | | ±0.2℃ (at 360℃) | |
| | | Temp. deviation | ±1.5℃ (at 100℃), ±2℃ (at 200℃), ±3℃(at 260℃) | | ±1.5℃ (at 100℃), ±2℃ (at 200℃), ±3℃ (at 300℃), ±4℃ (at 360℃) | |
| | JTM standard | Temp. adjusting accuracy | ±0.1℃ (at 260℃) | | ±0.2℃ (at 360℃) | |
| | | Temp. distribution accuracy | ±1.5℃ (at 260℃) | | ±2.5℃ (at 360℃) | |
| Max. temp. reaching time | | Approx. 60 min (to 260℃) | | Approx. 80 min (to 360℃) | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber | | Rock wool | |
| | Heater | | Stainless steel heating tubes with fins | | | |
| | | | 2.1KW | 3.0KW | 2.7KW | 3.75KW |
| | Fan blade/motor | | Axial fan blades, motor 20W×1 | | | |
| | Cable port | | I.D. 30 mm, located on the back | | | |
| Other additional mechanisms | | Exhaust dampers (manual) | | | | |
| Controllers | Temp. control method | | 3-stage PID | | | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | | | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display | | | |
| | | | Setting temp. display: Red 4-digit LED digital display | | | |
| | Timer | | 1 min.~99 h 59 min. and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | | | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | | |
| Heater circuit control | | SSR drive | | | | |
| Sensors | | K thermocouple (Temp. controller and overheating protector) | | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating protection, automatic overheating preventer, SSR short-circuit), independent overheating preventer, overcurrent leakage protection, key lock functions | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 450×450×450 | 600×600×600 | 450×450×450 | 600×600×600 |
| | External dimensions (W×D×H mm) | | 1050×630×850 | 1200×780×1000 | 1050×630×850 | 1200×780×1000 |
| | Internal capacity | | 91L | 216L | 91L | 216L |
| | Shelf load | | 30kg/layer | | | |
| | Shelf layers/shelf support spacing | | 9 layers/45mm | 9 layers/60mm | 9 layers/45mm | 9 layers/60mm |
| | Power supply (50/60Hz) rated current | | AC220V 10A | AC220V 14A | AC220V 12.5A | AC220V 17.5A |
| | Weight | | Approx. 112kg | Approx. 156kg | Approx. 112kg | Approx. 156kg |
| Accessories | Shelf | | Stainless steel wire mesh plate | | | |
| | | | 2 pcs | 3 pcs | 2 pcs | 3 pcs |
| | Supports | | 4 pcs | 6 pcs | 4 pcs | 6 pcs |
| Options | Stand | | OP42/45C | OP62/65C | OP42/45C | OP62/65C |
| | Others | | Shelf (1 shelf with 2 supports), cable port (30/50mm), micro printer, data logger, combination warning light (standby/running/fault), automatic door, viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software, touch screen controller | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers


13 Washers

14 Analysis and Test Devices

15 Options

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification


Temperature rise curve



Control panel



Structural diagram



Internal chamber





Options



- ① Exhaust pipe
- ② Data logger
- ③ Observation window

Dimension diagram (mm)



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Forced Convection Ovens | Large Capacity, Precision, Horizontal

DF811C/1011C, DH811C/1011C

Forced
ConvectionAutomatic
Overheating
PreventerIndependent
Overheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range DF: RT+10 ~ 200°C DH: RT+10 ~ 300°C

Temp. distribution accuracy DF: $\pm 3.0^{\circ}\text{C}$ (at 200°C) DH: $\pm 5.0^{\circ}\text{C}$ (at 300°C)

Internal capacity 512L 1000L

Large capacity high-precision constant temperature ovens for thermal processing.



■ Features

- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with 90 segments program controller.
- Overheating protection, deviation correction and key locking through submenu keys.
- Quick exhaust and cooling through exhaust dampers.

■ Safety

- With self-diagnostic circuits (temperature sensor anomaly, heater disconnection, automatic overheating prevention, SSR short circuit), overheat prevention device, leakage protection switch for overcurrent prevention, key lock, and other safety features.

■ Specifications

| Model | | | DF811C | DF1011C | DH811C | DH1011C |
|-----------------------------|--------------------------------------|-----------------------------|--|-------------------------------|---|-------------------------------|
| System | | | Forced convection and ventilation | | | |
| Performance | Operating temp. range | | Room temp. +10~200℃ | | Room temp. +10~300℃ | |
| | GB standard | Temp. fluctuation | ±0.2℃ (at 200℃) | | ±0.3℃ (at 300℃) | |
| | | Temp. deviation | ±1.5℃ (at 100℃), ±2℃ (at 200℃) | | ±1.5℃ (at 100℃), ±2℃ (at 200℃), ±4℃ (at 300℃) | |
| | JTM standard | Temp. adjusting accuracy | ±0.2℃ (at 200℃) | | ±0.3℃ (at 300℃) | |
| | | Temp. distribution accuracy | ±3.0℃ (at 200℃) | | ±5.0℃ (at 300℃) | |
| Max. temp. reaching time | | | Approx. 60 min (to 200℃) | | Approx. 80 min (to 300℃) | |
| Interior material | | | Stainless steel plate | | | |
| Exterior material | | | Cold rolled steel plate with chemical proofing coating | | | |
| Composition | Insulating material | | Glass fiber | | Rock wool | |
| | Heater | | Stainless steel heating tubes with fins | | | |
| | | | 4.5KW | 6.0KW | 6.9KW | 9.0KW |
| | Fan blade/motor | | Axial fan blades, motor 20W×1 | Axial fan blades, motor 20W×2 | Axial fan blades, motor 20W×1 | Axial fan blades, motor 20W×2 |
| | Cable port | | I.D. 30 mm, located on the back | | | |
| Other additional mechanisms | | | Exhaust dampers (manual) | | | |
| Controllers | Temp. control method | | 3-stage PID | | | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | | | |
| | Temp. display method | | Achieved Temp. display: Green 4-digit LED digital display | | | |
| | | | Setting temp. display: Red 4-digit LED digital display | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | | | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments ×1, 15 segments ×2, 10 segments ×3) | | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | | |
| Heater circuit control | | | SSR drive | | | |
| Sensors | | | K thermocouple (Temp. controller and overheating protector) | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 800×800×800 | 1000×1000×1000 | 800×800×800 | 1000×1000×1000 |
| | External dimensions (W×D×H mm) | | 1500×1015×1340 | 1700×1215×1540 | 1500×1015×1340 | 1700×1215×1540 |
| | Internal capacity | | 512L | 1000L | 512L | 1000L |
| | Shelf load | | 30kg/layer | | | |
| | Shelf layers/shelf support spacing | | 15 layers/50mm | 19 layers/50mm | 15 layers/50mm | 19 layers/50mm |
| | Power supply (50/60Hz) rated current | | 3 phase AC380V 7.5A | 3 phase AC380V 10A | 3 phase AC380V 11A | 3 phase AC380V 14.5A |
| | Weight | | Approx. 350kg | Approx. 450kg | Approx. 350kg | Approx. 450kg |
| Accessories | Shelf | | Stainless punching mesh plate | | | |
| | | | 3 pcs | | | |
| | Supports | | 6 pcs | | | |
| Options | | | Shelf (1 shelf with 2 supports), cable port (30/50mm), micro printer, data logger, combination warning light (standby/running/fault), automatic vent, viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software, touch screen controller | | | |

- 1 Sterilizers
- 2 Granulation and Spray Dryers
- 3 Muffle Furnaces
- 4 Ovens
- 5 Incubators
- 6 Plasma Equipment
- 7 Water Purifiers
- 8 Baths
- 9 Water Circulators
- 10 Rotary Evaporators
- 11 Freeze Dryers & Cold Traps
- 12 Stirrers & Shakers
- 13 Washers
- 14 Analysis and Test Devices
- 15 Options

Forced Convection Ovens | Precision, Vertical

DFS710C/810C, DHS710C/810C

Forced
ConvectionAutomatic
Overheating
PreventerIndependent
Overheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range DFS: RT+10~260℃ DHS: RT+10~360℃

Temp. distribution accuracy DFS: $\pm 2.0^{\circ}\text{C}$ DHS: $\pm 3.0^{\circ}\text{C}$

Internal capacity 418L 558L

A precision large-capacity constant temperature oven suitable for various temperature characteristic experiments.

Features

Space-saving and large capacity design allows for high-precision temperature control and adaptation to a wide range of temperature variations, applicable for various thermal processing and temperature experiments.

- The space-saving design with the control panel positioned at the top effectively uses the space in experimental rooms and production lines.
- The control panel's key input is dialog-based, making it simple for anyone to complete the settings. The internal temperature during operation can be visually displayed through LED.
- Equipped with various operational functions, primarily program operation. Current functions can be enhanced through options.
- Despite being a large capacity machine, it also features high-performance temperature regulation accuracy. Additionally, control accuracy and drying characteristics are not affected even during venting operation.



Specifications

| Model | | | DFS710C | DFS810C | DHS710C | DHS810C |
|--------------------------|--------------------------------------|--|--|--------------------|-----------------------------|---------------|
| System | | | Forced convection and ventilation | | | |
| Performance | Operating temp. range | | Room temp. +10~260℃ | | Room temp. +10~360℃ | |
| | GB standard | Temp. fluctuation | ±0.1℃ (at 260℃) | | ±0.2℃ (at 360℃) | |
| | | Temp. deviation | ±4℃ (at 260℃) | | ±5℃ (at 360℃) | |
| | JTM standard | Temp. adjusting accuracy | ±0.2℃ (at 260℃) | | ±0.3℃ (at 360℃) | |
| | | Temp. distribution accuracy | ±2.0℃ (at 260℃) | | ±3.0℃ (at 360℃) | |
| Max. temp. reaching time | | Approx. 50min | | Approx. 60 min | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber | | Glass fiber + ceramic fiber | |
| | Heater | | Stainless steel heating tubes with fins | | | |
| | | | 4.4 KW | | 5.4 KW | |
| | | | | | 5.8 KW | |
| Fan blade/motor | | Centrifugal fan blades (with temperature insulation coating)/motor rated at 120W | | | | |
| Cable port | | Inner diameter: 30mm, located on the right | | | | |
| Controllers | Temp. control method | | PID control | | | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | | | |
| | Temp. display method | | Achieved temperature display: Green 4-digit LED, digital display Setting temperature display: Red 4-digit LED, digital display | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h/1 min or 1 h | | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | | | |
| | Program mode | | Program operation: 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | | | |
| | Additional functions | | Timer function (actual timing within 24 h) | | | |
| | Sensors | | K thermocouple (Temp. controller and overheating protector) | | | |
| Safety device | | | Self-diagnostic circuits (temperature sensor anomaly, heater disconnection, automatic overheating prevention, SSR short circuit), independent overheating prevention device, leakage protection, door switch, control cabinet switch | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 620×750×900 | 620×750×1200 | 620×750×900 | 620×750×1200 |
| | External dimensions (W×D×H mm) | | 770×965×1622 | 770×965×1922 | 770×965×1622 | 770×965×1922 |
| | Internal capacity | | 418L | 558L | 418L | 558L |
| | Shelf load | | 30kg/layer | | | |
| | Shelf layers | | 27 layers | 37 layers | 27 layers | 37 layers |
| | Shelf plate spacing | | 30mm | | | |
| | Power supply (50/60Hz) rated current | | Single phase AC220V 21A | 3-phase AC380V 10A | | |
| | Weight | | Approx. 150kg | Approx. 170kg | Approx. 150kg | Approx. 170kg |
| Accessories | | | Shelf plates (2 stainless steel wire mesh plates), shelf supports (4) | | | |
| Options | | | Shelf plates, viewing windows, micro printers, data loggers, N ₂ introduction devices, abnormal alarm indicator lights, emergency stop switch, power failure manual recovery selection, cable ports (50/100mm), exhaust vents (80mm), exhaust flanges, automatic dampers, temperature output terminals, external alarm output terminals, timer output terminals, central monitoring software, touch screen controller | | | |

Internal chamber



Shelf



Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Forced Convection Ovens | Precision, 500°C

DH450C/650C/850C

Operating temp. range RT+10~500°C

Temp. distribution accuracy DH450C/650C: ±5°C

DH850C: ±8°C

Internal capacity 91L

216L

512L

Fine ovens with a maximum working temperature of 500°C.



Features

- Maximum operation temperature 500°C.
- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with 90 segments program controller.
- Overheating protection, deviation correction and key locking through submenu keys.
- Quick exhaust and cooling through exhaust dampers.

Safety

- Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), overheating preventer, overcurrent leakage protection, key lock and other safety functions.

Specifications

| Model | | | DH450C | DH650C | DH850C |
|-----------------------------|--------------------------------------|---|--|--------------------------|---|
| System | | | Forced convection and ventilation | | |
| Performance | Operating temp. range | | Room temp. +10~500℃ | | |
| | GB standard | Temp. fluctuation | ±0.2℃ (at 300℃), ±0.5℃ (at 500℃) | | |
| | | Temp. deviation | ±1.5℃ (at 100℃), ±2℃ (at 200℃), ±3℃ (at 300℃), ±5℃ (at 400℃), ±6℃ (at 500℃) | | ±1.5℃ (at 100℃), ±2℃ (at 200℃), ±4℃ (at 300℃), ±6℃ (at 400℃), ±8℃ (at 500℃) |
| | JTM standard | Temp. adjusting accuracy | ±0.2℃ (at 300℃), ±0.5℃ (at 500℃) | | |
| | | Temp. distribution accuracy | ±3.0℃ (at 300℃), ±5.0℃ (at 500℃) | | ±5.0℃ (at 300℃), ±8.0℃ (at 500℃) |
| Max. temp. reaching time | | Approx. 60 min (to 500℃) | | Approx. 80 min (to 500℃) | |
| Composition | Interior material | | Stainless steel plate | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | |
| | Insulating material | | Aluminum silicate wool, ceramic fiber blocks, insulation blocks | | |
| | Heater | Alloy heating wire | | | |
| | | 6.0KW | 7.8KW | 9.0KW | |
| | Fan blade/motor | | Axial flow fan blades, self-cooling high-temperature motor 20W | | |
| | Exhaust port/cable port | | Inner diameter 50mm, located at the back side/inner diameter 30mm, located at the back side | | |
| Other additional mechanisms | | Exhaust dampers (manual) | | | |
| Controllers | Temp. control method | | 3-stage PID | | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys to realize | | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display | | | |
| | | Setting temp. display: Red 4-digit LED digital display | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | |
| | Heater circuit control | | SSR drive | | |
| Sensors | | K thermocouple (Temp. controller and overheat protector) | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | | |
| Specifications | Internal dimensions (W×D×H mm) | | 450×450×450 | 600×600×600 | 800×800×800 |
| | External dimensions (W×D×H mm) | | 1200×800×1150 | 1350×950×1300 | 1550×1150×1500 |
| | Internal capacity | | 91L | 216L | 512L |
| | Shelf load | | 30kg/layer | | |
| | Shelf layers/shelf support spacing | | 9 layers/45mm | 9 layers/60mm | 12 layers/60mm |
| | Power supply (50/60Hz) rated current | | Single phase AC220V 30A | 3 phase AC380V 12.5A | 3 phase AC380V 15A |
| | Weight | | Approx. 160kg | Approx. 250kg | Approx. 320kg |
| Accessories | Shelf | | Stainless steel wire mesh plate | | |
| | | | 2 pcs | | |
| | Supports | | 4 pcs | | |
| Options | | | Shelf plate (1 shelf plate with 2 shelf supports), cable ports (30/50mm), micro printer, data logger, automatic dampers, combination warning lights (standby/running/fault), external communication function (RS485), temperature output terminals (4~20mA), external alarm output terminals, time-up output terminals, central monitoring software, touch screen controller | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Forced Convection Ovens | Rapid heating and cooling

DKG611/611V/811/811V

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerMade in
Japan

Operating temp. range RT+30~260°C

Temp. distribution accuracy $\pm 2.0^{\circ}\text{C}$ (260°C)

Internal capacity 150L 300L

The temperature rise time and fall time have been significantly reduced, aiming to improve production efficiency with easy-operation.

Features

- Temperature rise time and fall time reduced by over 50% compared to previous products, greatly enhancing work efficiency.
- By operating the manual adjustment valve at the front of the body, the full exhaust method with linked intake and exhaust is used, highly shortening the temperature fall time.
- With horizontal airflow, high-precision temperature performance can be achieved even under load.
- The door sealing strip is fluororubber, achieving a silicone-free design. (DKG611V/811V)
- Designed for AC200~220V wide voltage power supply.



Specifications

| Model | | DKG611 | DKG611V | DKG811 | DKG811V |
|------------------|--------------------------------------|--|--------------|---|--------------|
| System | | Forced convection and ventilation | | | |
| Performance | Temp. control range | Room temp. +30~260℃ | | | |
| | Temp. adjusting accuracy | JTM K05±0.5℃ (260℃) | | | |
| | Time fluctuation | JIS 1℃ (260℃) | | JIS 2℃ (260℃) | |
| | Temp. distribution accuracy | JTM K05±2℃ (200℃)/±2.5℃ (260℃) | | | |
| | Temp. slope | JIS 10℃ (260℃) | | JIS 15℃ (260℃) | |
| | Max. temp. reaching time | Within 45 min (25→260℃) | | Within 50 min (25→260℃) | |
| | Temp. fall time | Approx. 40 min (260→50℃) | | Approx. 50 min (260→50℃) | |
| Composition | Interior material | Stainless steel plate | | | |
| | Door side seal | Silicone rubber | Fluororubber | Silicone rubber | Fluororubber |
| | Heater | 2.6~3.15KW | | 3.6~4.36KW | |
| | Cable port | Inner diameter 30mm×1, located on the right | | | |
| | Intake vent | Bottom of heating chamber (manual damper for opening/closing) | | | |
| | Exhaust vent | Φ80mm (upper part of main unit) | | | |
| | Exhaust damper | Manual damper for linked intake and exhaust | | | |
| Controllers | Controllers | VS6 (PID control) | | | |
| | Operation functions | Fixed temp. operation/program operation/quick auto stop operation/auto stop operation/auto start operation/program auto start operation | | | |
| | Program mode | 6 modes (PrG1:30 segments, PrG 2~3: 15 segments, PrG4~6:10 segments) waiting function, repeat function, hold function, skip function | | | |
| | Additional functions | Deviation correction, key lock, power failure compensation | | | |
| Safety functions | | Self-diagnostic circuit (automatic overheating prevention, temperature sensor abnormality, heater disconnection, SSR short circuit, main relay fault, data logger error, internal communication error, temperature measurement error) overheat preventer/overcurrent leakage protection, door switch, temperature fuse, external alarm output terminal | | | |
| Specifications | Internal dimensions | 600×500×500mm | | 600×500×1000mm | |
| | External dimensions | 770×696×985mm 770×846×985mm (including exhaust pipe) | | 770×696×1674mm 770×846×1674mm (including exhaust pipe) | |
| | Internal capacity | 150L | | 300L | |
| | Shelf layers/shelf support spacing | 7 layers, 60mm | | 15 layers, 60mm | |
| | Shelf load | 15kg/layer (total weight capacity 35kg) | | 15kg/layer (total weight capacity 75kg) | |
| | Power supply (50/60Hz) rated current | 1ΦAC200~220V 13.5~15.0A | | 1ΦAC200~220V 18.5~20.5A | |
| | Circuit breaker capacity | 20A | | 30A | |
| Accessories | Weight | Approx. 110kg | | Approx. 155kg | |
| | Shelf | Stainless steel plate 2 pcs | | Stainless steel plate 4 pcs | |
| | Supports | 4 pcs | | 8 pcs | |


Control panel



Internal chamber



Heating and cooling time comparison



Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Forced Convection Ovens | Energy-Saving

DNE411C/611C/811C/911C

Operating temp. range RT+10~210°C

Temp. distribution accuracy $\pm 2.0^{\circ}\text{C}$ (at 210°C)

Internal capacity 90L 150L 300L 540L

Power consumption reduced by over 30%, an environmentally friendly constant temperature oven that reduces CO₂ emissions.

Features

- Significantly reduces power consumption and features programmable functions in an environmentally friendly constant temperature oven. Achieves a 30% energy-saving design (compared to previous products) during constant temperature operation through airtight thermal insulation in the chamber.
- Maximum temperature reach time shortened by up to 15 min (compared to previous products). Standby time and recovery time are also reduced, achieving operational efficiency.
- Able to implement various operating modes and is equipped with calibration correction functions, power failure reset mode selection, and user-defined information save and recall functions, among other support features.
- Able to repeat program operation for up to 99 segments or 99 modes. Equipped with accumulated function for powered/operating time.
- Various optional features allow for system upgrades based on user needs.

Safety

- Features self-diagnostic circuits, independent overheating preventers with digital settings, overcurrent leakage protection, key locks, and other safety functions.

Specifications

| Model | | | DNE411C | DNE611C | DNE811C | DNE911C |
|--------------------------|--------------------------------------|--|---|----------------|-----------------------------------|--------------------------|
| System | | | Forced convection | | | |
| Performance | Operating temp. range | | Room temp. +10~210℃ | | | |
| | GB standard | Temp. fluctuation | ±0.5℃ (at 210℃) | | | |
| | | Temp. uniformity | ±1.5% (at 210℃) | | | |
| | JTM standard | Temp. adjusting accuracy | ±0.5℃ (at 210℃) | | | |
| | | Temp. distribution accuracy | ±2.0℃ (at 210℃) | | | |
| Max. temp. reaching time | | | Approx. 60 min | Approx. 70 min | Approx. 45 min | Approx. 60 min |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber | | | |
| | Heater | | Stainless steel heating pipe | | | |
| | | | 1.1KW | 1.2KW | 1.2 KW×2 | 1.5 KW×2 |
| | Forced air motor | | High-temperature self-cooling motor | | | |
| | | | 10W | | 30W | 30W×2 |
| | Cable port | | Inner diameter 33mm (one at upper right side) | | | |
| Intake vent | | | Inner diameter 33mm (one at lower right side) | | | |
| Exhaust vent | | | I.D. 33mm×2, located on the top | | I.D. 33 mm×2, located at the back | |
| Controllers | Temp. control method | | PID control | | | |
| | Temp. setting method | | Use special function menu keys and up/down keys to realize digital setting | | | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display | | | |
| | | | Setting temp. display: Orange 5-digit LED digital display | | | |
| | Timer/timer resolution | | 1 min~99 hours 59 min/1 min | | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | | | |
| | Program mode | | Program operation max 99 segments, repeat operation function | | | |
| Additional functions | | Power/accumulated running time function (65535 hours), time display, deviation correction function, cumulative power consumption, CO ₂ emission, heater output, cumulative power-on time, power failure recovery mode selection, user defined information login | | | | |
| Sensors | | K thermocouple (Temp. controller and overheating protector) | | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 450×450×450 | 600×500×500 | 600×500×1000 | 1090×500×1000 |
| | External dimensions (W×D×H mm) | | 580×646×860 | 730×695×910 | 730×695×1675 | 1220×695×1675 |
| | Internal capacity | | 90L | 150L | 300L | 540L |
| | Shelf load | | 15kg/layer | | | |
| | Shelf layers/shelf support spacing | | 11 layers/30mm | 13 layers/30mm | 29 layers/30mm | 29 layers/30mm×2 columns |
| | Power supply (50/60Hz) rated current | | AC220V 5.5A | AC220V 6A | AC220V 11.5A | AC220V 14.5A |
| | Weight | | Approx. 63kg | Approx. 77kg | Approx. 92kg | Approx. 185kg |
| Accessories | Shelf | | Stainless punching mesh plate | | | |
| | | | 2 pcs | | 4 pcs | 8 pcs |
| Options | Supports | | 4 pcs | | 8 pcs | 16 pcs |
| | Stand | | ON61C | | | |
| | Stacking fittings | | ODN26C | ODN28C | — | |
| | Others | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), data logger, combination warning lights (standby/running/fault), viewing window, external communication function (RS485), temperature output terminals (4~20mA), external alarm output terminals, time-up output terminals | | | |



1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices


15 Options

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Control panel



Structural diagram



DNE structural diagram


Internal chamber



[DNE411C]

[DNE911C]


Shelf (Support)



Cable port (standard configuration)



Dimension diagram (mm)



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Forced Convection Ovens | Energy-saving, Variable Airflow

DNF411C/611C/811C/911C

Operating temp. range RT+10~260°C

Temp. distribution accuracy $\pm 2.5^{\circ}\text{C}$ (at 260°C)

Internal capacity 90L 150L 300L 540L

Variable airflow, power consumption reduced by over 30%, an environmentally friendly constant temperature oven that reduces CO₂ emissions.

Features

- Equipped with variable airflow mechanism: 10-step variable airflow (411C/611C) + airflow [0] (stop), 10-step variable airflow (811C/911C).
- Achieve 30% energy-saving during constant temperature operation through airtight thermal insulation in the chamber (compared to previous products).
- Maximum temperature reach time shortened by up to 15 min (compared to previous products). Standby time and recovery time are also reduced, achieving operational efficiency.
- High airtightness prevents dust and debris from entering the chamber.
- Fixed temp. operation, 99 segments program operation, auto stop operation, auto start operation.
- Temperature and time setting displays, deviation corrections, etc., can be achieved through the VFD fluorescent display.
- Various optional features allow for system upgrades based on user needs.

Safety

- Equipped with self-diagnostic circuits, independent overheating preventers with digital settings, overcurrent leakage protection, key locks, and other safety functions.



Specifications

| Model | | | | DNF411C | | DNF611C | | DNF811C | | DNF911C | |
|----------------|--------------------------------------|-----------------------------|-----------------------------|---|--|----------------|--|-------------------|--|--|--|
| System | | | | Forced convection + natural convection | | | | Forced convection | | | |
| Performance | Forced air (airflow 10) | Operating temp. range | | Room temp. +10~260℃ | | | | | | | |
| | | GB standard | Temp. fluctuation | ±0.5℃ (at 260℃ max airflow setting) | | | | | | | |
| | | | Temp. uniformity | ±1.5% (at 260℃ max airflow setting) | | | | | | | |
| | | JTM standard | Temp. adjusting accuracy | ±0.5℃ (at 260℃ max airflow setting) | | | | | | | |
| | | | Temp. distribution accuracy | ±2.5℃ (at 260℃ max airflow setting) | | | | | | | |
| | | Max. temp. reaching time | | Approx. 105 min | | | | Approx. 60 min | | Approx. 100 min | |
| | Natural convection (airflow 0) | Operating temp. range | | Room temp. +25~120℃ | | | | | | | |
| | | Temp. adjusting accuracy | | ±0.3℃ (airflow 0 at 120℃) | | | | | | | |
| | | Temp. distribution accuracy | | ±3℃ (airflow 0 at 120℃) | | | | | | | |
| | | Max.temp. reaching time | | Approx. 25 min | | | | | | | |
| Composition | Interior material | | | Stainless steel plate | | | | | | | |
| | Exterior material | | | Cold rolled steel plate with chemical proofing coating | | | | | | | |
| | Insulating material | | | Glass fiber | | | | | | | |
| | Heater | | | Stainless steel heating pipe | | | | | | | |
| | | | | 0.6KW×2 | | 0.83KW×2 | | 1.35KW×2 | | 1.65KW×2 | |
| | Forced air motor | | | DC brushless motor (600~1500rpm) variable (10 segments) | | | | | | | |
| | | | | 30W | | | | 30W×2 | | | |
| | Cable port | | | Inner diameter 33mm (upper on the right) | | | | | | | |
| | Intake vent | | | Inner diameter 33mm (lower on the right) | | | | | | One on each side | |
| | Exhaust vent | | | Inner diameter 50mm×1, located on the back | | | | | | Inner diameter 50mm×2, located on the back | |
| Controllers | Temp. control method | | | PID control | | | | | | | |
| | Temp. setting method | | | Digital setting through special function menu keys and up/down keys | | | | | | | |
| | Temp. display method | | | Achieved temp. display: Green 4-digit LED digital display | | | | | | | |
| | | | | Setting temp. display: Orange 5-digit LED digital display | | | | | | | |
| | Timer/timer resolution | | | 1 min~99 hours 59 min/1 min | | | | | | | |
| | Operation functions | | | Fixed temp. operation, auto start, auto stop, program operation | | | | | | | |
| | Program mode | | | Program operation max 99 segments, repeat operation function | | | | | | | |
| | Additional functions | | | Power/accumulated running time function (65535 hours), clock display, deviation correction function, cumulative power consumption, CO ₂ emission, heater output, cumulative power-on time, power failure recovery mode selection, user defined information login, variable airflow function (DNF411C/611C) | | | | | | | |
| | Sensors | | | K thermocouple (Temp. controller and overheating protector) | | | | | | | |
| Safety device | | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | | | | | | | |
| Specifications | Internal dimensions (W×D×H mm) | | | 450×450×450 | | 600×500×500 | | 600×500×1000 | | 1090×500×1000 | |
| | External dimensions (W×D×H mm) | | | 580×646×897 | | 730×695×947 | | 730×695×1685 | | 1220×695×1685 | |
| | Internal capacity | | | 90L | | 150L | | 300L | | 540L | |
| | Shelf load | | | 15kg/layer | | | | | | | |
| | Shelf layers/shelf support spacing | | | 11 layers/30mm | | 13 layers/30mm | | 29 layers/30mm | | 29 layers/30mm×2 columns | |
| | Power supply (50/60Hz) rated current | | | AC220V 6A | | AC220V 8A | | AC220V 13A | | AC220V 16A | |
| | Weight | | | Approx. 61kg | | Approx. 90kg | | Approx. 135kg | | Approx. 210kg | |
| Accessories | Shelf | | | Stainless punching mesh plate | | | | | | | |
| | | | | 2 pcs | | | | 4 pcs | | 8 pcs | |
| Options | Supports | | | 4 pcs | | | | 8 pcs | | 16 pcs | |
| | Stand | | | ON61C | | | | — | | | |
| | Stacking fittings | | | ODN26C | | ODN28C | | — | | | |
| | Others | | | Shelf (1 shelf with 2 supports), cable port (30/50mm), data logger, combination warning light (standby/running/fault), viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal | | | | | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Control panel



Internal chamber



[DNF611C]



[DNF911C]


Damper operation




Exhaust vent (back of main unit)




Dimension diagram (mm)




[DNF411C]



[DNF611C]



[DNF811C]



[DNF911C]

- Sterilizers 1
- Granulation and Spray Dryers 2
- Muffle Furnaces 3
- Ovens 4
- Incubators 5
- Plasma Equipment 6
- Water Purifiers 7
- Baths 8
- Water Circulators 9
- Rotary Evaporators 10
- Freeze Dryers & Cold Traps 11
- Stirrers & Shakers 12
- Washers 13
- Analysis and Test Devices 14
- Options 15

Natural Convection Ovens | High-Function, Touch Screen

DVT410C/610C

Operating temp. range RT+10~300°C

Temp. uniformity ±9°C (at 300°C)

Internal capacity 117L 188L

High-end drying oven with programmable control functions allows for quick program settings.

Features

- 7-inch touch screen for easy operation with three operation modes, including fixed temp. operation, program operation, timing operation. Standard 3-level authorization login.
- 100 programs, each program maximum 30 segments. Standard data storage function, data storage capacity approximately 10 years.
- Temperature rise curve tracking function, program operation mode, able to visually display set temperature and achieved temperature curve, providing reference for users to set reasonable temperature rise curve.
- Automatically stops heating upon door opening to protect user samples.

Safety

- Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), overheating preventer, overcurrent leakage protection, and a three-level authorization management system, automatically stops heating upon door opening.



Specifications

| Model | | DVT410C | DVT610C |
|----------------|-----------------------------------|---|---|
| System | | Natural convection | |
| Performance | Operating temp. range | | Room temp. +10~300°C |
| | GB standard | Temp. fluctuation | ±1°C (at 300°C) |
| | | Temp. uniformity | ±3% (at 300°C) |
| | Max. temp. reaching time | | ≤120 min |
| Composition | Interior material | | Stainless steel plate |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating |
| | Insulating material | | Rock wool+centrifugal glass wool |
| | Heater | | Stainless steel heating pipe |
| | | | 1.6KW |
| | | | 1.3KW × 2 |
| Controllers | Cable port | | Inner diameter: 34mm×1, located on the right |
| | Exhaust vent | | I.D. 34mm×2, located on the top |
| | Temp. control method | | 3-stage PID |
| | Temp. setting method | | 7-inch color touch screen |
| | Temp. display method | | 7-inch color touch screen |
| | Timer | | 1 min~99 h 59 min |
| | Operation functions | | Fixed temp. operation, program operation, timing operation |
| | Program mode | | 100 Programs with 30 segments each |
| | Additional functions | | Temperature history curve, alarm history records, operation log, data export function, language switching (Chinese, Japanese, English), calibration function, power failure compensation function, screen brightness adjustment, three-level authorization function |
| | Sensors | | Temperature controller: K-type thermocouple, overheating protection: Liquid expansion temperature controller |
| Safety device | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), independent overheating preventer, overcurrent leakage protection, and automatically stops heating upon door opening | |
| Specifications | Internal dimensions | | 450×520×500 |
| | External dimensions (W×D×H mm) | | 610×690×960 |
| | Internal capacity | | 117L |
| | Shelf load | | 30kg/layer |
| | Shelf layers | | 12 layers |
| | Support spacing | | 30mm |
| | Power supply (50Hz) rated current | | AC220V 7.5A |
| | Weight | | Approx. 78kg |
| Accessories | Shelf | | Stainless steel wire shelf |
| | Supports | | 2 pcs |
| Options | Stand | | ONK60C |
| | Stacking fittings | | ODK40C |
| | Others | | Shelf (1 shelf with 2 supports) cable port (30/50mm), micro printer, data logger, combination warning light (standby/running/fault), remote monitoring function (network port), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Natural
ConvectionAutomatic
Overheating
PreventerIndependent
Overheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Internal chamber (DVT610C)



- Use high-temperature resistant silicone rubber sealing strips. Enhance sealing, ensuring stable performance.
- The internal stainless-steel construction offers strong corrosion resistance and is easy to clean.
- Shelf plates are made of stainless-steel wire mesh, highly increasing strength.


Control panel



Cable port (standard configuration)



Temperature rise curve




Exhaust vent (standard configuration)



Shelf (support)



Dimension diagram (mm)



Sterilizers 1

Granulation
and Spray
Dryers 2Muffle
Furnaces 3

Ovens 4

Incubators 5

Plasma
Equipment 6Water
Purifiers 7

Baths 8

Water
Circulators 9Rotary
Evaporators 10Freeze Dryers
& Cold Traps 11Stirrers &
Shakers 12

Washers 13

Analysis and
Test Devices 14

Options 15

Natural Convection Ovens | Standard, Program Operation

DVS413C/613C

Operating temp. range RT+5~260℃

Temp. distribution accuracy $\pm 5^{\circ}\text{C}$ (at 260℃)

Internal capacity 99L 162L

The natural convection drying oven with program control function allows for quick program settings.



Features

- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with 90 segments program controller.
- Deviation correction and key locking through submenu keys.

Safety

- Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), independent overheating prevention, overcurrent leakage protection, key lock safety and other safety functions.

Specifications

| Model | | | DVS413C | | DVS613C | |
|--------------------------|--------------------------------------|--|---|--|--------------|--------|
| System | | | Natural convection | | | |
| Performance | Operating temp. range | | Room temp. +5~260℃ | | | |
| | GB standard | Temp. fluctuation | ±1℃ (at 260℃) | | | |
| | | Temp. uniformity | ±3% (at 260℃) | | | |
| | JTM standard | Temp. adjusting accuracy | ±1℃ (at 260℃) | | | |
| | | Temp. distribution accuracy | ±5℃ (at 260℃) | | | |
| Max. temp. reaching time | | Approx. 75 min (room temperature + 5℃~260℃) | | | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber+rock wool | | | |
| | Heater | Stainless steel heating pipe | | | | |
| | | 1.2KW | | | | 1.36KW |
| | Viewing window (width×height mm) | | 250×280 | | | |
| | Cable port | | Inner diameter: 30mm×1, located on the right | | | |
| Exhaust vent | | I.D. 30mm×2, located on the top | | | | |
| Controllers | Temp. control method | | 3-stage PID | | | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | | | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display | | | | |
| | | Setting temp. display: Red 4-digit LED digital display | | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, program operation, auto stop operation, auto start of fixed temp. operation, auto start of program operation, quick auto stop operation | | | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments × 1, 15 segments × 2, 10 segments × 3) | | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | | |
| Sensors | | Temperature controller: K-type thermocouple; Overheating protection: Liquid expansion temperature controller | | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), independent overheating preventer, overcurrent leakage protection, and key lock functions | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 450×490×450 | | 600×540×500 | |
| | External dimensions (W×D×H mm) | | 560×601×827 | | 710×651×877 | |
| | Internal capacity | | 99L | | 162L | |
| | Shelf load | | 15kg/layer | | | |
| | Shelf layers | | 11 layers | | 13 layers | |
| | Support spacing | | 30mm | | | |
| | Power supply (50/60Hz) rated current | | AC220V 5.5A | | AC220V 6.5A | |
| | Weight | | Approx. 48kg | | Approx. 63kg | |
| Accessories | Shelf | | Stainless punching mesh plate 2 pieces, 1 of which is fixed to the chamber; do not disassemble | | | |
| | Supports | | 2 pcs | | | |
| Options | Stand | | ON61C | | | |
| | Stacking fittings | | OD40C | | OD60C | |
| | Others | | Shelf (1 shelf with 2 supports), cable port (30/50mm), micro printer, data logger, combination warning light (standby/running/fault), external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options


Natural
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification

Internal chamber (DVS413C)



- Use high-temperature resistant silicone rubber sealing strips. Enhance sealing, ensuring stable performance.
- The internal stainless-steel construction offers strong corrosion resistance and is easy to clean.
- The shelf is made of perforated metal plate, significantly improving strength.

Temperature rise curve



Control panel



Cable port (standard configuration)




Exhaust vent (standard configuration)



Shelf (support)



Dimension diagram (mm)



Sterilizers 1

Granulation
and Spray
Dryers 2Muffle
Furnaces 3

Ovens 4

Incubators 5

Plasma
Equipment 6Water
Purifiers 7

Baths 8

Water
Circulators 9Rotary
Evaporators 10Freeze Dryers
& Cold Traps 11Stirrers &
Shakers 12

Washers 13

Analysis and
Test Devices 14

Options 15

Natural Convection Ovens | Economical Type, Fixed Temp. Operation

DX312C/412C/612C

Natural
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification

| | | | | | | | | |
|-----------------------|-------------------------|--------------------|-----------------------------|-------|-------------------|-----|-----|------|
| Operating temp. range | DX312C/412C: RT+5~300°C | DX612C: RT+5~280°C | Temp. distribution accuracy | ±10°C | Internal capacity | 28L | 74L | 153L |
|-----------------------|-------------------------|--------------------|-----------------------------|-------|-------------------|-----|-----|------|

Maximum temperature of 300°C, highly practical drying oven.



Features

- Easy operation, available for fixed temp., quick auto stop, auto stop and auto start operations.
- Digital setting through special operation function menu keys and up/down keys. Able to set overheating protection, deviation correction, and key lock.

Safety

- Equipped with self-diagnostic circuits (temperature sensor anomaly), power failure compensation function, deviation correction function, leakage protection against overcurrent, independent overheating preventer, and other safety functions.

DX312C




Specifications

| Model | | DX312C | DX412C | DX612C |
|----------------|--------------------------------------|--|------------------|---------------------|
| System | | Natural convection | | |
| Performance | Operating temp. range | Room temp. +5~300°C | | Room temp. +5~280°C |
| | GB standard | Temp. fluctuation | ±1°C (at 300°C) | |
| | | Temp. uniformity | ±3% (at 300°C) | |
| | JTM standard | Temp. adjusting accuracy | ±1°C (at 300°C) | |
| | | Temp. distribution accuracy | ±10°C (at 300°C) | |
| Composition | Max. temp. reaching time | Approx. 50min | Approx. 100 min | Approx. 110 min |
| | Interior material | Stainless steel plate | | |
| | Exterior material | Cold rolled steel plate with chemical proofing coating | | |
| | Insulating material | Rock wool | | |
| | Heater | Nichrome heating wire | | |
| Controllers | Exhaust vent | I.D. 30mm×2, located on the top | | |
| | Temp. control method | PID control | | |
| | Temp. setting method | Digital setting through function menu key and up/down keys | | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display | | |
| | | Setting temp. display: Red 4-digit LED digital display | | |
| | Timer | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | |
| | Operation functions | Fixed temp. operation, auto start, quick auto stop, auto stop operation | | |
| | Additional functions | Deviation correction, key lock, power failure compensation | | |
| | Heater circuit control | SSR drive | | |
| | Sensors | Temperature controller: K-type thermocouple, overheating protection: Liquid expansion temperature controller | | |
| Safety device | | Self-diagnostic function (temperature sensor anomaly, automatic overheating prevention), key lock function, independent overheating preventer, overcurrent leakage protection switch | | |
| Specifications | Internal dimensions (W×D×H mm) | 300×310×300 | 450×410×400 | 600×510×500 |
| | External dimensions (W×D×H mm) | 400×440×620 | 550×540×720 | 700×640×820 |
| | Internal capacity | 28L | 74L | 153L |
| | Shelf load | 15kg/layer | | |
| | Shelf layers | 6 layers | 9 layers | 12 layers |
| | Support spacing | 35mm | | |
| Accessories | Power supply (50/60Hz) rated current | AC220V 5A | AC220V 7.5A | AC220V 7.5A |
| | Weight | Approx. 23kg | Approx. 38kg | Approx. 56kg |
| Options | Shelf | Stainless punching mesh plate | | |
| | Supports | 2 pcs | | |
| | Stand | ONS30C | ONS60C | |
| | Stacking fittings | ODK80C | ODK82C | ODK84C |
| Others | | Shelf (1 shelf plate with 2 shelf supports), cable port (30/50mm) | | |

Internal chamber (DX412C)



Dimension diagram (mm)



Natural Convection Ovens | Basic Function Type

DY311C/411C/611C

Natural
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range RT+5~280°C

Temp. gradient 30°C (at 280°C)

Internal capacity 28L 74L 153L

Ensures the basic functions of fixed temperature operation models.


Various function settings allow for simpler operation.

Features

- Easy operation, fixed temp. operation, auto stop operation, quick auto stop operation, and auto start operation.

Safety

- Equipped with self-diagnostic circuits (temperature sensor anomaly, temperature overheating prevention), leakage protection against overcurrent, and independent overheating preventer and other safety functions.



Specifications

| Model | | | DY311C | DY411C | DY611C |
|-------------------------|--------------------------------------|--|---|-----------------|--------------|
| System | | | Natural convection | | |
| Performance | Operating temp. range | | Room temp. +5~280℃ | | |
| | GB standard | Temp. fluctuation | ±1℃ (at 280℃) | | |
| | | Temp. uniformity | ±3% (at 280℃) | | |
| | JIS standard | Temp. variation | ±2℃ (at 280℃) | | |
| | | Temp. gradient | 30℃ (at 280℃) | | |
| | Temp. rise speed | | 6.9℃/min | 4.3℃/min | 1.7℃/min |
| Maximum temp. rise time | | Approx. 60 min | Approx. 70 min | Approx. 120 min | |
| Composition | Interior material | | Stainless steel plate | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | |
| | Insulating material | | Rock wool | | |
| | Heater | Nichrome heating wire | | | |
| | | 0.9KW | 1.36KW | | |
| Exhaust vent | | I.D. 34mm×2, located on the top | | | |
| Controllers | Temp. control method | | PID control | | |
| | Temp. setting method | | Digital setting through function menu key and up/down keys | | |
| | Temp. display method | Achieved temperature display: red 4-digit led, digital display | | | |
| | | Setting temperature display: green 4-digit led, digital display | | | |
| | Timer | | 1 min~999 h 50 min | | |
| | Operation functions | | Fixed temp. operation, quick auto stop operation, auto stop operation, auto start operation | | |
| | Additional functions | | Deviation correction function | | |
| | Heater circuit control | | SSR drive | | |
| Sensors | | Temperature controller: K-type thermocouple, overheating protection: Liquid expansion temperature controller | | | |
| Safety device | | | Self-diagnostic function (temperature sensor abnormal, temperature overheating prevention), independent overheating prevention, overcurrent leakage protection switch | | |
| Specifications | Internal dimensions (W×D×H mm) | | 300×310×300 | 450×410×400 | 600×510×500 |
| | External dimensions (W×D×H mm) | | 400×440×630 | 550×540×730 | 700×640×830 |
| | Internal capacity | | 28L | 74L | 153L |
| | Shelf load | | 15kg/layer | | |
| | Shelf layers | | 6 layers | 9 layers | 12 layers |
| | Support spacing | | 35mm | | |
| | Power supply (50/60Hz) rated current | | AC220V 5A | AC220V 7.5A | AC220V 7.5A |
| | Weight | | Approx. 23kg | Approx. 38kg | Approx. 56kg |
| Accessories | Shelf | | Stainless punching mesh plate | | |
| | Supports | | 2 pcs | | |
| Options | Stand | | ONS30C | | |
| | Stacking fittings | | ODK80C | | ODK84C |
| | Others | | Shelf (1 shelf plate with 2 shelf supports), cable port (30/50mm) | | |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Inert Ovens | N₂ Replacement, Oxygen-free Heating

DN410IC/610IC

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification

Operating temp. range RT+15~360°C

Temp. distribution accuracy ±3°C (at 360°C)

Internal capacity 95L 223L

Specific constant temperature ovens for thermal treatment experiments in an oxygen-free environment.



Features

- Support performing high-temperature resistance tests and thermal processing at 360°C.
- Wide operating temperature range, high accuracy in temperature control.
- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with a total of 90 segments program controller.
- Flow setting and introduction through N₂ flow meter.



Safety

- Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), overheating preventer, overcurrent leakage protection, key lock and other safety functions.

Specifications

| Model | | DN410IC | DN610IC |
|----------------|--------------------------------------|--|---|
| System | | Forced convection | |
| Performance | Operating temp. range | Room temp. +10~360°C | |
| | GB standard | Temp. fluctuation ±0.2°C (at 360°C) | Temp. uniformity ±1.5% (at 360°C) |
| | JTM standard | Temp. adjusting accuracy ±0.2°C (at 360°C) | Temp. distribution accuracy ±3.0°C (at 360°C) |
| | Max. temp. reaching time | Approx. 60 min | |
| | Interior material | Stainless steel plate | |
| Composition | Exterior material | Cold rolled steel plate with chemical proofing coating | |
| | Insulating material | Rock wool | |
| | Heater | Stainless steel heating tube, 3.0KW | Stainless steel heating tube, 4.0KW |
| | Fan blade/motor | Centrifugal fan blades, high-temperature self-cooling motor 30W | |
| | Flow meter | Maximum flow: 40L/min | |
| Controllers | N ₂ inlet interface | Φ8 | |
| | Temp. control method | 3-stage PID | |
| | Temp. setting method | Digital setting through special function menu keys and up/down keys | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display Setting temp. display: Red 4-digit LED digital display | |
| | Timer | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | |
| | Operation functions | Fixed temp. operation, auto start, auto stop, program operation | |
| | Program mode | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | |
| | Additional functions | Deviation correction, key lock, power failure compensation | |
| | Sensors | K thermocouple (temp. controller and overheating protector) | |
| | Safety device | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | |
| Specifications | Internal dimensions (W×D×H mm) | 470×450×450 | 620×600×600 |
| | External dimensions (W×D×H mm) | 640×695×915 | 790×846×1065 |
| | Internal capacity | 95L | 223L |
| | Shelf load | 30kg/layer | |
| | Shelf layers/shelf support spacing | 12 layers/30mm | 17 layers/30mm |
| | Power supply (50/60Hz) rated current | AC220V 14A | AC220V 19A |
| | Weight | Approx. 80kg | Approx. 120kg |
| Accessories | Shelf | Stainless steel wire mesh plate | |
| | Supports | 2 pcs | |
| | Stand | OH41C | OH61C |
| Options | Others | Shelf plate, cable port (30/50mm), micro printer, data logger, combined warning light (standby/operation/fault), viewing window, exhaust vent, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, timer output terminal, central monitoring software | |

Temperature rise curve

N₂ replacement performance curveN₂ inlet (outer diameter 9mm)

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Inert Ovens | Anaerobic, Oxygen Concentration 5000ppm

DF611IC/DH611IC

Operating temp. range DF: RT+10~260°C DH: RT+10~360°C

Temp. distribution accuracy DF: $\pm 2.0^{\circ}\text{C}$ (at 260°C) DH: $\pm 3.0^{\circ}\text{C}$ (at 360°C)

High-precision constant temperature ovens for thermal treatment in an oxygen-free environment.

Features

- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with a total of 90 segments program controller.
- Overheat protection, deviation correction and key locking through submenu keys.
- Flow setting and introduction through N_2 flow meter.
- Real-time display of oxygen concentration in the chamber (optional).

Safety

- Equipped with self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection, key lock safety and other functions.



Specifications

| Model | | DF611IC | DH611IC |
|-----------------------------|--------------------------------------|--|--|
| System | | Forced convection | |
| Performance | Operating temp. range | | Room temp. +10~260°C |
| | GB standard | Temp. fluctuation | $\pm 0.1^{\circ}\text{C}$ (at 260°C) |
| | | Temp. deviation | $\pm 2^{\circ}\text{C}$ (at 100°C), $\pm 3^{\circ}\text{C}$ (at 200°C), $\pm 4^{\circ}\text{C}$ (at 260°C) |
| | JTM standard | Temp. adjusting accuracy | $\pm 0.1^{\circ}\text{C}$ (at 260°C) |
| | | Temp. distribution accuracy | $\pm 2^{\circ}\text{C}$ (at 260°C) |
| | Max. temp. reaching time | | Approx. 60 min (to 260°C) |
| Lowest oxygen concentration | | $\leq 5000\text{ppm}$ | |
| Composition | Interior material | | Stainless steel plate |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating |
| | Insulating material | | Glass fiber |
| | Heater | | Stainless steel heating tubes with fins |
| | Fan blade/motor | | Axial fan blades, motor 20W×1 |
| | Cable port | | I.D. 30 mm, located at the back |
| | Exhaust vent | | 8mm tapered joint |
| | Intake port | | Rc3/8 threaded interface |
| Controllers | Temp. control method | | 3-stage PID |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display Setting temp. display: Red 4-digit LED digital display |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) |
| | Additional functions | | Deviation correction, key lock, power failure compensation |
| | Sensors | | K thermocouple (temp. controller and overheating protector) |
| Safety device | | Equipped with self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection, key lock safety and other functions | |
| Specifications | Internal dimensions | | 600×600×600 |
| | External dimensions (W×D×H mm) | | 1200×780×1000 |
| | Internal capacity | | 216L |
| | Shelf load | | 30kg/layer |
| | Shelf layers/shelf support spacing | | 9 layers/60mm |
| | Weight | | Approx. 109kg |
| | Power supply (50/60Hz) rated current | | AC220V 14A |
| Accessories | Shelf | | Stainless steel wire mesh plate |
| | Supports | | 3 pcs |
| | Stand | | 6 pcs |
| Options | Stand | | OP62C |
| | Others | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), micro printer, data logger, combined warning light (standby/operation/fault), viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, timer output terminal, central monitoring software, touchscreen controller |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Inert Ovens | Anaerobic, Oxygen Concentration 20ppm, High Temp.

DNN430C/630C/460C/660C

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionOvercurrent
Leakage Circuit
Breaker

Max. operating temp. DNN430C/630C: 360°C

DNN460C/660C: 600°C

Oxygen concentration 20ppm

Internal capacity 91L 261L

An atmosphere furnace with an oxygen concentration of 20ppm and a maximum operating temperature of 600°C, suitable for high-temperature thermal treatment in an oxygen-free environment.



Features

- Maximum operating temperatures of 360°C and 600°C, with oxygen concentrations below 20ppm, primarily applied in electronic component annealing, LTCC removal, glass substrate annealing, ceramic debinding, etc.
- High gas-tightness pressure chamber construction, short oxygen concentration reaching time, low N₂ consumption.
- Maintaining high gas-tightness through magnetic sealing and a water-cooling mechanism to protect sealing components from thermal effects.
- Quick heating and cooling achievable, with adjustable heating and cooling rates.
- Standard equipped with waste liquid recovery devices to cool and recover gases.
- Liquid crystal touchscreen operation, with both automatic and manual operation modes.

Safety

- Equipped with door detection switch, temperature overheating preventer, oxygen concentration abnormality detection, nitrogen pressure detection, nitrogen flow detection, cooling water flow detection, leakage sensor, overcurrent leakage protection and other devices.

Specifications

| Model | | | DNN430C | DNN630C | DNN460C | DNN660C |
|------------------------------------|--------------------------------------|--|--|-----------------------------|---|--------------------|
| System | | | Nitrogen replacement + Forced convection | | | |
| Performance | Operating temp. range | | Room temp. +50~360℃ | | Room temp. +50~600℃ | |
| | GB standard | Temp. fluctuation | ±0.3℃ (at 360℃) | | ±0.5℃ (at 600℃) | |
| | | Temp. deviation | ±1.5℃ (at 100℃), ±2℃ (at 200℃), ±4℃ (at 300℃), ±5℃ (at 360℃) | | ±1.5℃ (at 100℃), ±2℃ (at 200℃), ±4℃ (at 300℃), ±6℃ (at 400℃), ±8℃ (at 600℃) | |
| | JTM standard | Temp. adjusting accuracy | ±0.3℃ (at 360℃) | | ±0.5℃ (at 600℃) | |
| | | Temp. distribution accuracy | ±4℃ (at 360℃) | | ±6℃ (at 600℃) | |
| | Max. temp. reaching time | | Approx. 80 min | | Approx. 110 min | |
| | Oxygen concentration | | Below 20ppm | | | |
| Oxygen concentration reaching time | | ≤20 min (250L/min nitrogen) | ≤25 min (250L/min nitrogen) | ≤20 min (250L/min nitrogen) | ≤25 min (250L/min nitrogen) | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Aluminosilicate cotton | | | |
| | Heating method | | Stainless steel heating pipe | | Alloy heating wire | |
| | Heating power | Main heater | 6.0KW | 9KW | 10.8KW | 14.2KW |
| | | N ₂ preheating | 1KW | | | |
| Cooling mechanism | | Stainless steel cooling water heat exchanger | | | | |
| Controllers | Temp. control method | | PID control | | | |
| | Temp. controller | | Liquid crystal touchscreen + PLC | | | |
| | Operation functions | | Manual fixed temp. operation, automatic program operation | | | |
| | Program mode | | 9 program operations, each program can set 9 segments (expandable) | | | |
| | Additional functions | | 3-level authorization login | | | |
| | Heater circuit control | | SSR drive | | | |
| | Sensors | | K-type thermocouple (temperature control and overheating prevention) | | | |
| Safety device | | | Door detection switch, overheating preventer, oxygen concentration abnormality, temperature overheating of N ₂ preheating, low water pressure detection, low N ₂ pressure detection, cooling water flow detection, N ₂ flow detection, fan overload detection, overcurrent leakage protection, etc. | | | |
| Specifications | Internal dimensions | | 450×450×450 | 660×600×660 | 450×450×450 | 660×600×660 |
| | External dimensions (W×D×H mm) | | 1090×1170×1570 | 1300×1320×1780 | 1240×1010×1750 | 1300×1320×1780 |
| | Internal capacity | | 91L | 261L | 91L | 261L |
| | Shelf layer/spacing | | 12 layers/30mm | 19 layers/30mm | 12 layers/30mm | 19 layers/30mm |
| | Exhaust vent | | 40KF flange | | | |
| | Nitrogen inlet | | Rc1/2 | | | |
| | Cooling water inlet | | Rc3/4 | | | |
| | Hot water outlet | | Rc3/4×2 | | | |
| | Power supply (50/60Hz) rated current | | 3 phase AC380V 16A | 3 phase AC380V 19A | 3 phase AC380V 20A | 3 phase AC380V 27A |
| Weight | | Approx. 800kg | Approx. 1000kg | Approx. 800kg | Approx. 1000kg | |
| Options | | | Shelf plate, vacuum pump, temperature data logger, multicolor light tower, Ethernet communication, UPS backup power supply, scanner, PC components, temperature measurement components, etc. | | | |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Explosion-proof Ovens | Explosion-Proof, Precision

DF411SC/611SC, DH411SC/611SC

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Max. operating temp. DF: RT+10~260°C

DH: RT+10~360°C

Temp. distribution accuracy

DF: $\pm 2.5^{\circ}\text{C}$ (at 260°C)DH: $\pm 3.0^{\circ}\text{C}$ (at 360°C)

Internal capacity

91L

216L

Precision constant temperature ovens equipped with explosion pressure relief device.

Features

- Featuring a relief safety door (explosion port) and a robust interlocked door knob for secure body and door fastening.
- Contain the following interlock functions:
The device operates under locked door conditions via the door locking mechanism.
An explosion detection switch is installed on the ceiling; once activated, operation will be interrupted.
- Easy to operate, fixed temp. operation, program operation, quick auto stop operation, auto stop operation, and auto start operation.
- Quick exhaust and cooling through exhaust dampers.

Safety

- Equipped with self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection, key lock and other safety functions.



Specifications

| Model | | | DF411SC | DF611SC | DH411SC | DH611SC |
|-----------------------------|--------------------------------------|-----------------------------|--|---------------|--|---------------|
| System | | | Forced convection and ventilation | | | |
| Performance | Operating temp. range | | Room temp. +10~260℃ | | Room temp. +10~360℃ | |
| | GB standard | Temp. fluctuation | ±0.2℃ (at 260℃) | | ±0.3℃ (at 360℃) | |
| | | Temp. deviation | ±2℃ (at 100℃), ±3℃ (at 200℃), ±4℃ (at 260℃) | | ±2℃ (at 100℃), ±3℃ (at 200℃), ±4℃ (at 300℃), ±5℃ (at 360℃) | |
| | JTM standard | Temp. adjusting accuracy | ±0.1℃ (at 260℃) | | ±0.2℃ (at 360℃) | |
| | | Temp. distribution accuracy | ±2.5℃ (at 260℃) | | ±3.0℃ (at 360℃) | |
| Max. temp. reaching time | | | Approx. 60 min (to 260℃) | | Approx. 80 min (to 360℃) | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber | | Rock wool | |
| | Heater | | Stainless steel heating tubes with fins | | | |
| | | | 2.1KW | 3.0KW | 2.7KW | 3.75KW |
| | Fan blade/motor | | Axial fan blades, motor 20W×1 | | | |
| | Cable port | | I.D. 30 mm, located at the back | | | |
| Other additional mechanisms | | | Relief safety door (explosion port), exhaust vent (manual) | | | |
| Controllers | Temp. control method | | 3-stage PID | | | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | | | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display | | | |
| | | | Setting temp. display: Red 4-digit LED digital display | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | | | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | | |
| Heater circuit control | | | SSR drive | | | |
| Sensors | | | K thermocouple (temp. controller and overheating protector) | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 450×450×450 | 600×600×600 | 450×450×450 | 600×600×600 |
| | External dimensions (W×D×H mm) | | 1050×630×1250 | 1200×780×1400 | 1050×630×1250 | 1200×780×1400 |
| | Internal capacity | | 91L | 216L | 91L | 216L |
| | Shelf load | | 30kg/layer | | | |
| | Shelf layers/shelf support spacing | | 9 layers/45mm | 9 layers/60mm | 9 layers/45mm | 9 layers/60mm |
| | Power supply (50/60Hz) rated current | | AC220V 10A | AC220V 14A | AC220V 13A | AC220V 17.5A |
| Weight | | | Approx. 78kg | Approx. 109kg | Approx. 78kg | Approx. 109kg |
| Accessories | Shelf | | Stainless steel wire mesh plate | | | |
| | | | 2 pcs | 3 pcs | 2 pcs | 3 pcs |
| | Supports | | 4 pcs | 6 pcs | 4 pcs | 6 pcs |
| Options | Stand | | OP42C | OP62C | OP42C | OP62C |
| | Others | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), micro printer, data logger, combined warning light (standby/operation/fault), viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, timer output terminal, central monitoring software, touchscreen controller | | | |

Explosion-proof Ovens | Explosion-proof, Precision High-capacity

DF811SC/1011SC, DH811SC/1011SC

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Specifications

| Model | | DF811SC | DF1011SC |
|--------------------------------|-----------------------------|-----------------------------------|-------------------------------|
| System | | Forced convection and ventilation | |
| Operating temp. range | | Room temp. +10~200℃ | |
| GB standard | Temp. fluctuation | ±0.2℃ (at 200℃) | |
| | Temp. deviation | ±2℃ (at 100℃), ±3℃ (at 200℃) | |
| JTM standard | Temp. adjusting accuracy | ±0.2℃ (at 200℃) | |
| | Temp. distribution accuracy | ±3.0℃ (at 200℃) | |
| Max. temp. reaching time | | Approx. 60 min (to 200℃) | |
| Heater | | 4.5KW | 6.0KW |
| Fan blade/motor | | Axial fan blades, motor 20W×1 | Axial fan blades, motor 20W×2 |
| Internal dimensions (W×D×H mm) | | 800×800×800 | 1000×1000×1000 |
| External dimensions (W×D×H mm) | | 1500×1015×1742 | 1700×1215×1932 |
| Internal capacity | | 512L | 1000L |
| Power supply (50/60Hz) | | 3 phase AC380V 7.5A | 3 phase AC380V 10A |
| Weight | | Approx. 210kg | Approx. 280kg |

| Model | | DH811SC | DH1011SC |
|--------------------------------|-----------------------------|---|-------------------------------|
| System | | Forced convection and ventilation | |
| Operating temp. range | | Room Temp. +10~300℃ | |
| GB standard | Temp. fluctuation | ±0.3℃ (at 300℃) | |
| | Temp. deviation | ±2℃ (at 100℃), ±3℃ (at 200℃), ±5℃ (at 300℃) | |
| JTM standard | Temp. adjusting accuracy | ±0.3℃ (at 300℃) | |
| | Temp. distribution accuracy | ±5.0℃ (at 300℃) | |
| Max. temp. reaching time | | Approx. 80 min (to 300℃) | |
| Heater | | 6.9KW | 9.0KW |
| Fan blade/motor | | Axial fan blades, motor 20W×1 | Axial fan blades, motor 20W×2 |
| Internal dimensions (W×D×H mm) | | 800×800×800 | 1000×1000×1000 |
| External dimensions (W×D×H mm) | | 1500×1015×1742 | 1700×1215×1932 |
| Internal capacity | | 512L | 1000L |
| Power supply (50/60Hz) | | 3 phase AC380V 11A | 3 phase AC380V 14.5A |
| Weight | | Approx. 210kg | Approx. 280kg |

Explosion-proof Ovens | Explosion-proof, 500℃, Precision

DH450SC/650SC/850SC

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Specifications

| Model | | DH450SC | DH650SC | DH850SC |
|--------------------------------|-----------------------------|---|--------------------|--------------------|
| System | | Forced convection and ventilation | | |
| Operating temp. range | | Room temp. +10~500℃ | | |
| GB standard | Temp. fluctuation | ±0.5℃ (at 500℃) | | |
| | Temp. deviation | ±2℃ (at 100℃) | | ±2℃ (at 100℃) |
| | | ±3℃ (at 200℃) | | ±3℃ (at 200℃) |
| | | ±4℃ (at 300℃) | | ±5℃ (at 300℃) |
| JTM standard | Temp. deviation | ±6℃ (at 400℃) | | ±8℃ (at 400℃) |
| | | ±8℃ (at 500℃) | | ±10℃ (at 500℃) |
| | Temp. adjusting accuracy | ±0.5℃ (at 500℃) | | |
| | Temp. distribution accuracy | ±5.0℃ (at 500℃) | | ±8.00℃ (at 500℃) |
| Max. temp. reaching time | | Approx. 60 min | | |
| Heater | | 6.0KW | 7.8KW | 9.0KW |
| Fan blade/motor | | Axial flow fan, self-cooling high-temp. motor 20W × 1 | | |
| Internal dimensions | | 450×450×450 | 600×600×600 | 800×800×800 |
| External dimensions (W×D×H mm) | | 1200×800×1550 | 1350×950×1700 | 1550×1150×1900 |
| Internal capacity | | 91L | 216L | 512L |
| Power supply (50/60Hz) | | Single phase AC220V 30A | 3 phase AC380V 12A | 3 phase AC380V 15A |
| Weight | | Approx. 210kg | Approx. 300kg | Approx. 370kg |

Ovens for Labware | Natural Convection, Forced Exhaust

DG410C/450C/810C/850C

Maximum operating temp. RT+5~70°C

Internal capacity 92L 445L

User-friendly and highly praised labware drying oven.

Features

- Equipped with a large observation window for easy observation.
- Capable of drying multiple instruments for storage or use immediately after drying. The interior and exterior are designed with a high thermal insulation double structure.
- The internal heater is a stainless steel heating tube, With a water basin at the bottom, making it safe for instruments with high water droplet levels.
- Interior is made of stainless steel for easy cleaning.
- Models DG450C and DG850C have more multiple functions compared with DG410C and DG810C. An intake filter is installed at the air intake, and a forced exhaust fan is installed on the exhaust side. Additionally, standard sterilization lamps are included for rapid drying and clean drying functions.

Safety

- Equipped with self-diagnostic circuits (temperature sensor anomaly), power failure compensation function, deviation correction function, leakage protection switch against overcurrent, independent overheating preventer, and other safety features.



Specifications

| Model | | DG410C | DG450C | DG810C | DG850C | |
|---------------------------|--------------------------------------|--|---|-------------------------------|---------------------------------|-------------------------------|
| System | | Natural convection | | | | |
| Performance · Composition | Operating temp. range | | Room temp. +5~70℃ | | | |
| | GB standard | Temp. fluctuation | ±0.5℃ (at 70℃) | | | |
| | | Temp. uniformity | ±5.0℃ (at 70℃) | | | |
| | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Heater | Stainless steel heating pipe | | | | |
| | | 1KW | 1.34KW | | | |
| | Viewing window (width×height mm) | | 250×300 | 250×700 | | |
| | Exhaust vent | | Internal diameter: 34mm × 2 | Axial flow fan forced exhaust | Internal diameter: 34mm × 2 | Axial flow fan forced exhaust |
| Air intake | | Internal diameter: 30mm × 2 | Equipped with air intake filter | Internal diameter: 30mm × 2 | Equipped with air intake filter | |
| Sterilization lamp | | None | 1 × 15W sterilization lamp | None | 1 × 15W sterilization lamp | |
| Controllers | Temp. control method | | PID control | | | |
| | Temp. setting method | | Digital setting via function menu key and ▲ ▼ keys | | | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display | | | | |
| | | Setting temp. display: Red 4-digit LED digital display | | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, quick auto stop, program operation | | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | | |
| | Heater circuit control | | SSR drive | | | |
| Sensors | | Temperature controller: K-type thermocouple, overheating protection: Liquid expansion temperature controller | | | | |
| Safety device | | Self-diagnostic function (temperature sensor abnormality, SSR short circuit, automatic overheating preventer), key lock, independent overheating preventer, leakage protection circuit breaker against overcurrent | | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 450×450×450 | | 620×600×1195 | |
| | External dimensions (W×D×H mm) | | 504×562×788 | 504×562×820 | 674×711×1586 | 674×711×1618 |
| | Internal capacity | | 92L | | 445L | |
| | Shelf load | | 15kg/layer | | | |
| | Shelf layers/shelf support spacing | | 10 layers/30mm | | 29 layers/30mm | |
| | Power supply (50/60Hz) rated current | | AC220V 5A | | AC220V 7A | |
| | Weight | | Approx. 45kg | Approx. 48kg | Approx. 78kg | Approx. 83kg |
| Accessories | Shelf | | Stainless punching mesh plate | | | |
| | | | 2 pcs | | 4 pcs | |
| | Supports | | 4 pcs | | 8 pcs | |
| | Water collection tray | | 1 pcs | | | |
| Options | Stand | | ON61C | | — | |
| | Others | | Shelf plate (1 shelf plate with 2 shelf supports) | | | |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Clean Ovens | Dust-free, Standard

DE430C/430UC/630C/630UC·DT430C/430UC/630C/630UC

| | | | | | | | |
|-----------------------|----------------|----------------|-------------|-----------|-------------------|-----|------|
| Operating temp. range | DE: RT+30~260℃ | DT: RT+30~360℃ | Cleanliness | Class 100 | Internal capacity | 91L | 216L |
|-----------------------|----------------|----------------|-------------|-----------|-------------------|-----|------|

Precision clean constant temperature oven with cleanliness level 100.



Features

This product suppresses particle generation and is suitable for constant temperature control in environments that require cleanliness, such as semiconductors, LCDs, electronic products, and precision instruments.

- Featuring a high-temperature resistant HEPA filter, achieving high-performance temperature control through horizontal circulation. The C-type maintains cleanliness level 100 during constant temperature while the UC-type maintains cleanliness level 100 throughout heating, constant temperature, and cooling.
- A sensor is installed on the inner side of the HEPA filter to control the temperature in proximity to the item.
- Designed for easy placement and retrieval of items. Prevent dust from entering the chamber during door opening and closing.
- Flow setting and introduction through N₂ flow meter.
- Equipped with self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), overheating preventer, overcurrent leakage protection, key lock and other safety functions.

Specifications


| Model | | | DE430C/430UC | DE630C/630UC | DT430C/430UC | DT630C/630UC |
|--------------------------------|--------------------------------------|---|--|--------------------------|------------------------|-------------------|
| System | | | Forced convection | | | |
| Performance | Operating temp. range | | Room temp. +30~260℃ | | Room temp. +30~360℃ | |
| | GB standard | Temp. fluctuation | ±0.3℃ (at 260℃) | | ±0.3℃ (at 360℃) | |
| | | Temp. uniformity | ±1.5% (at 260℃) | | ±1.5% (at 360℃) | |
| | JTM standard | Temp. adjusting accuracy | ±0.3℃ (at 260℃) | | ±0.3℃ (at 360℃) | |
| | | Temp. distribution accuracy | ±2.5℃ (at 260℃) | | ±4℃ (at 360℃) | |
| | Max. temp. reaching time | | Approx. 70 min | Approx. 70 min (to 260℃) | | Approx. 80 min |
| Cleanliness | | | C: cleanliness level 100 during constant temperature UC: cleanliness level 100 throughout heating, constant temperature, and cooling | | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber | | Aluminosilicate cotton | |
| | Heater | Stainless steel heating pipe | | | | |
| | | 2.5KW | 3.6KW | 5.2KW | | |
| | Fan blade/motor | | Centrifugal fan, high-temperature self-cooling motor 370W | | | |
| | Differential pressure gauges | | Analog type 0~500Pa | | | |
| | Cable port | | Internal diameter: 30mm×1, located on the right | | | |
| | Additional mechanism | | Exhaust vent (manual) outer diameter: 61mm | | | |
| | HEPA filter | | Dust collection efficiency: 99.97% for 0.3μm particles | | | |
| N ₂ inlet interface | | Outer diameter: 8mm tapered joint | | | | |
| Controllers | Temp. control method | | 3-stage PID | | | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | | | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display | | | | |
| | | Setting temp. display: Red 4-digit LED digital display | | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | | | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | | |
| Sensors | | K thermocouple (temp. controller and overheating protector) | | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 450×450×450 | 600×600×600 | 450×450×450 | 600×600×600 |
| | External dimensions (W×D×H mm) | | 700×1000×1738 | 850×1150×1738 | 700×1000×1738 | 850×1150×1738 |
| | Internal capacity | | 91L | 216L | 91L | 216L |
| | Shelf load | | 30kg/layer | | | |
| | Shelf layers/shelf support spacing | | 12 layers/30mm | 17 layers/30mm | 12 layers/30mm | 17 layers/30mm |
| | Power supply (50/60Hz) rated current | | 3-phase AC380V 5A | 3-phase AC380V 6.5A | | 3-phase AC380V 9A |
| | Weight | | Approx. 220kg | Approx. 270kg | Approx. 220kg | Approx. 270kg |
| Accessories | Shelf | | Stainless steel wire mesh plate | | | |
| | | | 2 pcs | 3 pcs | 2 pcs | 3 pcs |
| | Supports | | 4 pcs | 6 pcs | 4 pcs | 6 pcs |
| Options | | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), micro printer, data logger, combined warning light (standby/operation/fault), viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, timer output terminal, central monitoring software, touchscreen controller | | | |

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Internal chamber



Structural diagram




- ① Manual or automatic (optional) exhaust vent
- ② Cable port ③ Heater ④ Centrifugal fan ⑤ Exhaust vent
- ⑥ Aging judgment via HEPA filter differential pressure meter
- ⑦ Sensor (temperature control at a position close to the sample)


Control panel



Dimension diagram (mm)



[DE430C, DT430C, DE430UC, DT430UC]



[DE630C, DT630C, DE630UC, DT630UC]

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Clean Ovens | Dust-free Environment Heating, Large Capacity

DES830C/830UC·DTS830C/830UC

| | | | | | | |
|-----------------------|-----------------|-----------------|-------------|-----------|-------------------|------|
| Operating temp. range | DES: RT+50~260℃ | DTS: RT+50~360℃ | Cleanliness | Class 100 | Internal capacity | 300L |
|-----------------------|-----------------|-----------------|-------------|-----------|-------------------|------|

Space-saving, large capacity clean constant temperature chamber.



■ Features

With space-saving and large capacity design, it can suppress dust generation during temperature rise and fall. The controller allows for multifunctional, dialog-type input with simple operation, equipped with an easy-to-confirm graphical monitor.

- With large capacity of 300L for handling bulk workpieces, significantly reducing the footprint while effectively utilizing research, development, and production line space.
- Using high-temperature resistant HEPA filters, reaching up to 360℃; type C maintains cleanliness level 100 during constant temperature while type UC maintains cleanliness level 100 throughout heating, constant temperature, and cooling.

■ Specifications

| Model | | DES830C/830UC | DTS830C/830UC |
|----------------|--------------------------------------|---|---------------------|
| System | | Front horizontal air supply circulation | |
| Performance | Operating temp. range | Room temp. +50~260℃ | Room temp. +50~360℃ |
| | Cleanliness | C: cleanliness level 100 during constant temperature UC: cleanliness level 100 throughout heating, constant temperature, and cooling | |
| | GB standard | Temp. fluctuation ±0.5℃ | |
| | | Temp. uniformity ±1.5% (at 260℃) | ±1.5% (at 360℃) |
| | JTM standard | Temp. adjusting accuracy ±0.5℃ | |
| | | Temp. distribution accuracy ±2.0℃ (at 260℃) | ±5.0℃ (at 360℃) |
| Composition | Max. temp. reaching time | Approx. 70 min | Approx. 90 min |
| | Interior material | Stainless steel plate | |
| | Exterior material | Cold rolled steel plate with chemical proofing coating | |
| | Insulating material | Ceramic fiber | |
| | Heater | Stainless steel heating tube with fins | |
| | | 6KW | 9KW |
| Controllers | Fan blade/motor | Centrifugal fan, high-temperature self-cooling motor 200W × 2 | |
| | Differential pressure gauges | Analog type 0~500Pa | |
| | Cable port | Internal diameter: 30mm×1, located on the right | |
| | HEPA filter | Dust collection efficiency: 99.97% for 0.3μm particles | |
| | N ₂ inlet interface | Outer diameter: 8mm tapered joint | |
| | Temp. control method | 3-stage PID | |
| Safety device | Temp. setting method | Digital setting through special function menu keys and up/down keys | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display | |
| | | Setting temp. display: Red 4-digit LED digital display | |
| | Other displays | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | |
| | Timer | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | |
| | Operation functions | Fixed temp. operation, automatic start operation, automatic stop operation, program operation (up to 16 segments, can repeat, be separated, etc.) | |
| Specifications | Additional functions | Timing function (actual time timing within 24 hours), accumulated timing function (can accumulate up to 49,999 hours), external communication function (RS422A) | |
| | Sensors | K thermocouple (temp. controller and overheating protector) | |
| | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | |
| | Internal dimensions | 610×480×1100 | |
| | External dimensions (W×D×H mm) | 830×1078×1955 | |
| | Internal capacity | 300L | |
| Accessories | Shelf load | 30kg/layer | |
| | Shelf layers/shelf support spacing | 35 layers/30mm | |
| | Power supply (50/60Hz) rated current | 3 phase AC380V 10.5A | 3 phase AC380V 15A |
| | Weight | Approx. 350kg | |
| | Shelf | Stainless steel wire mesh plate | |
| | Supports | 3 pcs | |
| Options | | Shelf plate (1 shelf plate with 2 shelf supports), automatic vent, N ₂ introduction device, micro printer, data logger, combined warning light (standby/operation/fault), temperature output terminal (4~20mA), external alarm output terminal, timer output terminal, central monitoring software, touchscreen controller | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers


13 Washers

14 Analysis and Test Devices

15 Options

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Structural diagram







Control panel



Shelf



Cleanliness data during temperature rise and fall



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Clean Ovens | Dust-Free Environment Heating, Economical

DEC812C/912C

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range RT+10~150°C

Cleanliness Class 100

Internal capacity 236L 535L

Economical clean constant temperature oven with a maximum temperature of 150°C and cleanliness level 100.



Features

Suitable for drying or thermal treatment of semiconductors, LCDs, electronic products, precision instruments, etc., in a dust-free environment.

- Using a high-temperature resistant HEPA filter, achieving high-performance temperature control through horizontal circulation, maintaining cleanliness level 100 during constant temperature.
- A sensor is installed on the inner side of the HEPA filter to control the temperature in proximity to the item.
- Equipped with self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), overheating preventer, overcurrent leakage protection, key lock and other safety functions.

Specifications

| Model | | | DEC812C | | DEC912C | |
|----------------|--------------------------------------|---|--|--|--------------------------|--|
| System | | | Forced convection | | | |
| Performance | Operating temp. range | | Room temp. +10~150℃ | | | |
| | GB standard | Temp. fluctuation | ±0.5℃ (at 150℃) | | | |
| | | Temp. uniformity | ±3.0% (at 150℃) | | | |
| | JTM standard | Temp. adjusting accuracy | ±0.5℃ (at 150℃) | | | |
| | | Temp. distribution accuracy | ±3.0℃ (at 150℃) | | | |
| | Max. temp. reaching time | | Approx. 50min | | | |
| Cleanliness | | Cleanliness maintained at level 100 during constant temperature | | | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber | | | |
| | Heater | | Stainless steel heating pipe | | | |
| | | | 2.4KW | | 3.2KW | |
| | Fan blade/motor | | Centrifugal fan, 30W × 2 | | Centrifugal fan, 30W × 4 | |
| | Cable port | | Internal diameter: 30mm×1, located on the right | | | |
| HEPA filter | | Dust collection efficiency: 99.97% for 0.3μm particles | | | | |
| Controllers | Temp. control method | | 3-stage PID | | | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | | | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display | | | |
| | | | Setting temp. display: Red 4-digit LED digital display | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | | | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | | | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | | | |
| Sensors | | K thermocouple (temp. controller and overheating protector) | | | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 500×450×1050 | | 1072×500×1000 | |
| | External dimensions (W×D×H mm) | | 710×720×1600 | | 1170×720×1600 | |
| | Internal capacity | | 236L | | 535L | |
| | Shelf load | | 30kg/layer | | | |
| | Shelf layers/shelf support spacing | | 29 layers/30mm | | 2×29 layers/30mm | |
| | Power supply (50/60Hz) rated current | | AC220V 12A | | AC220V 16A | |
| | Weight | | Approx. 110kg | | Approx. 190kg | |
| Accessories | Shelf | | Stainless punching mesh plate | | | |
| | | | 4 pcs | | 8 pcs | |
| | Supports | | 8 pcs | | 16 pcs | |
| Options | | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), micro printer, data logger, combined warning light (standby/operation/fault), viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, timer output terminal, central monitoring software, touchscreen controller | | | |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Clean Inert Ovens | Class 100, Oxygen Concentration: 20ppm

DTN430C/630C/450C/650C

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionOvercurrent
Leakage Circuit
Breaker

Maximum operating temp. DTN430C/630C: 360°C DTN450C/650C: 500°C

Cleanliness Class 100

Oxygen concentration 20ppm

Internal capacity 91L

Atmosphere furnace with cleanliness level 100 and oxygen concentration below 20ppm, mainly used for curing semiconductor wafers.

Features

- Maximum operating temperatures of 360°C and 500°C, with cleanliness level 100 and oxygen concentration below 20ppm, primarily used for curing semiconductor wafers (photoresist PI, PBO curing), glass substrate baking, high-precision annealing, etc.
- Using a high-performance filter with high-temperature resistant, maintaining cleanliness level 100 inside the chamber, allowing for high-temperature baking in a clean environment.
- High gas-tightness pressure chamber construction, with short oxygen concentration reaching time and extremely low N₂ consumption. Maintaining high gas-tightness through magnetic sealing and a water-cooling mechanism to protect sealing components from thermal effects.
- Quick heating and cooling, with adjustable heating and cooling rates. Standard equipped with waste liquid recovery devices to cool and recover gases.

Safety

- Equipped with door detection switch, overheating preventer, oxygen concentration abnormality, nitrogen pressure detection, nitrogen flow detection, cooling water flow detection, leakage sensor, overcurrent leakage protection and other devices.



Specifications

| Model | | | DTN430C | DTN630C | DTN450C | DTN650C |
|------------------------------------|--------------------------------------|--|--|-----------------------------|---|--------------------|
| System | | | Nitrogen replacement + Forced convection | | | |
| Performance | Operating temp. range | | Room temp. +50~360℃ | | Room temp. +50~500℃ | |
| | GB standard | Temp. fluctuation | ±0.3℃ (at 360℃) | | ±0.5℃ (at 500℃) | |
| | | Average temp. deviation | ±1.5℃ (at 100℃), ±2℃ (at 200℃), ±4℃ (at 300℃), ±5℃ (at 360℃) | | ±1.5℃ (at 100℃), ±2℃ (at 200℃), ±4℃ (at 300℃), ±5℃ (at 400℃), ±8℃ (at 500℃) | |
| | JTM standard | Temp. adjusting accuracy | ±0.3℃ (at 360℃) | | ±0.5℃ (at 500℃) | |
| | | Temp. distribution accuracy | ±4℃ (at 360℃) | | ±6℃ (at 500℃) | |
| | Max. temp. reaching time | | Approx. 80 min | | Approx. 110 min | |
| | Oxygen concentration | | Below 20ppm | | | |
| Oxygen concentration reaching time | | ≤20 min (250L/min nitrogen) | ≤25 min (250L/min nitrogen) | ≤20 min (250L/min nitrogen) | ≤25 min (250L/min nitrogen) | |
| Cleanliness | | Constant temperature maintained at level 100 | | | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Aluminosilicate cotton | | | |
| | Heating method | | Stainless steel heating pipe | | Alloy heating wire | |
| | Heating power | Main heater | 6.0KW | 9KW | 10.8KW | 14.2KW |
| | | N ₂ preheating | 1KW | | | |
| | Cooling mechanism | | Stainless steel cooling water heat exchanger | | | |
| HEPA filter | | Dust collection efficiency: 99.97% for 0.3μm particles | | | | |
| Controllers | Temp. control method | | PID control | | | |
| | Temp. controller | | Liquid crystal touchscreen + PLC | | | |
| | Operation functions | | Manual fixed temp. operation, automatic program operation | | | |
| | Program mode | | 9 program operations, each program can set 9 segments (expandable) | | | |
| | Additional functions | | Operator level certification | | | |
| | Heater circuit control | | SSR drive | | | |
| Sensors | | K-type thermocouple (temperature control and overheating prevention) | | | | |
| Safety device | | | Door detection switch, overheating preventer, oxygen concentration anomaly, overheating of N ₂ preheating, low water pressure detection, low N ₂ pressure detection, cooling water flow detection, N ₂ flow detection, fan overload detection, overcurrent leakage protection, etc. | | | |
| Specifications | Internal dimensions | | 450×450×450 | 660×600×660 | 450×450×450 | 660×600×660 |
| | External dimensions (W×D×H mm) | | 1090×1170×1570 | 1300×1320×1780 | 1090×1170×1570 | 1300×1320×1780 |
| | Internal capacity | | 91L | 261L | 91L | 261L |
| | Shelf layer/spacing | | 12 layers/30mm | 19 layers/30mm | 12 layers/30mm | 19 layers/30mm |
| | Exhaust vent | | KF40 | | | |
| | Nitrogen inlet | | Rc1/2 | | | |
| | Cooling water inlet | | Rc3/4 | | | |
| | Hot water outlet | | Rc3/4, 2 pcs | | | |
| | Power supply (50/60Hz) rated current | | 3 phase AC380V 16A | 3 phase AC380V 19A | 3 phase AC380V 20A | 3 phase AC380V 27A |
| Weight | | Approx. 900kg | Approx. 1100kg | Approx. 900kg | Approx. 1100kg | |
| Options | | | Shelf plate, vacuum pump, temperature data logger, multicolor light tower, ethernet communication function, UPS backup power supply, scanner, PC components, temperature measurement components, etc. | | | |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Far Infrared Ray Ovens | Far Infrared Heating

DIR631C

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range RT+10~360°C

Temp. distribution accuracy $\pm 3.0^{\circ}\text{C}$ (at 360°C)

It uses the properties of far infrared ray heaters (IR heaters) for thermal treatment of polymer materials and other applications.



Features

- High temperature natural convection ovens with maximum operating temperature 360°C.
- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with a total of 90 segments program controller.
- Overheating protection, deviation correction and key locking through submenu keys.
- Quick exhaust and cooling through exhaust dampers.

Safety

- Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), independent overheating preventer, overcurrent leakage protection, key lock and other safety functions.

Control panel



Main body controller

Independent overheating preventer

IR heater controller

Internal chamber



Specifications

| Model | | | DIR631C |
|--------------------------|-----------------------------------|-----------------------------|--|
| System | | | IR radiation + Forced convection and ventilation |
| Performance | Operating temp. range | | Room temp. +10~360℃ |
| | GB standard | Temp. fluctuation | ±0.2℃ (at 360℃, IR heater off) |
| | | Temp. deviation | ±5℃ (at 360℃, IR heater off) |
| | JTM standard | Temp. adjusting accuracy | ±0.2℃ (at 360℃, IR heater off) |
| | | Temp. distribution accuracy | ±3℃ (at 360℃, IR heater off) |
| Max. temp. reaching time | | | 100 min (from room temperature to 360℃, IR heater off) |
| Composition | Main body heater | | Fin-type stainless steel heating tube 3.75KW |
| | IR heater | | 0.2KW×16 pieces×upper and lower surfaces for a total of 6.4KW |
| | Motor · fan | | Capacitor-type motor 20W + axial flow fan |
| | Sensors | | K-type thermocouple double probe×1 (for main body temperature controller and overheating preventer) K-type thermocouple×2 (for IR heater, centrally embedded in IR heater) |
| | Cable port | | I.D. 30 mm, located at the back |
| | Other additional functions | | Exhaust vent (manual type) |
| | Control method | | PID control (main body) |
| Controllers | Operating functions | | Fixed temp. operation, quick automatic stop, automatic stop, automatic start |
| | Additional functions | | Deviation correction function, key lock function, power failure compensation function, door switch detection |
| Safety device | Leakage circuit breaker | | Leakage, short circuit, and overcurrent protection |
| | Overheating preventer | | Automatically cut off the heater circuit when overheating occurs |
| Specifications | Self-diagnostic function | | Temperature sensor abnormality, heater disconnection, SSR short circuit, automatic overheating prevention function |
| | Internal dimensions (W×D×H mm) | | 600×600×600 |
| | Effective dimensions (W×D×H mm) | | 600×600×200 |
| | External dimensions (W×D×H mm) | | 1200×780×1000 |
| | Internal capacity | | 216L |
| | IR heater upper and lower spacing | | 200mm |
| | Weight | | Approx. 230kg |
| Accessories | Power supply | | 3-phase 380V 50/60Hz 18A |
| | Shelf | | Stainless steel wire mesh plate 1 pcs |
| | Supports | | 2 pcs |
| Options | Stand | | OP62C |
| | Others | | Shelves (1 shelf with 2 shelf supports), micro printer, data logger, combined warning light (standby/running/fault), external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

High Temp. Ovens | High Temp., 700°C

DR210C

Natural
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range 300~700°C

Temp. distribution accuracy $\pm 25^{\circ}\text{C}$ (at 700°C)

Internal capacity 13.75L

High-temperature drying oven that combines electric furnace, constant temperature chamber, and drying chamber.

Features

- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with a total of 90 segments program controller.
- Overheating protection, deviation correction and key locking through submenu keys.

Safety

- Equipped with self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit protection), overheating preventer, overcurrent leakage protection, key lock and other safety functions.




Specifications

| Model | | | DR210C | |
|--------------------------|--------------------------------------|---|--|--|
| System | | | Natural convection | |
| Performance | Operating temp. range | | 300~700℃ | |
| | GB standard | Temp. fluctuation | ±2℃ (at 700℃) | |
| | | Temp. uniformity | ±5% (at 700℃) | |
| | JTM standard | Temp. adjusting accuracy | ±2℃ (at 700℃) | |
| | | Temp. distribution accuracy | ±25℃ (at 700℃) | |
| Max. temp. reaching time | | Approx. 70 min (from room temperature + 5 to 700℃) | | |
| Composition | Interior material | | Stainless steel plate | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | |
| | Insulating material | | Ceramic fiber | |
| | Heater | | Iron-chromium alloy heating wire 1.3KW | |
| Controllers | Temp. control method | | 3-stage PID | |
| | Temp. setting method | | Digital setting through special function menu keys and up/down keys | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display | |
| | | | Setting temp. display: Red 4-digit LED digital display | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | |
| | Additional functions | | Deviation correction, key lock, power failure compensation | |
| Sensors | | K thermocouple (temp. controller and overheating protector) | | |
| Safety device | | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | |
| Specifications | Internal dimensions (W×D×H mm) | | 250×250×220 | |
| | External dimensions (W×D×H mm) | | 520×443×612 | |
| | Internal capacity | | 13.75L | |
| | Shelf load | | 15kg/layer | |
| | Shelf layers | | 3 layers | |
| | Support spacing | | 33mm | |
| | Power supply (50/60Hz) rated current | | AC220V 7A | |
| Weight | | | Approx. 36kg | |
| Accessories | Shelf | | Stainless punching mesh plate 2 pcs | |
| | Supports | | Integrated structure with the inner box (33mm spacing) | |
| | Stand | | ON61C | |
| Options | Others | | Shelves, micro printer, data logger, combined warning light (standby/running/fault), external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software | |

Control panel



Dimension diagram (mm)



Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Vacuum Ovens | Multi-function, Touch Screen

DPT40C/60C/80C/100C

Automatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range RT+20~200°C

Operating vacuum range 101~0.1KPa

Internal capacity 91L 216L 512L 1000L

Touch screen and vacuum degree control, various operating modes, easy operation.



■ Features

- 7-inch touch screen for easy operation with fixed temp. operation, program operation, and timer operation. Standard 3-level authorization login.
- Multiple pressure operation mode options combined with constant temperature and program operations, allowing for various programs to be designed.
- Gas exchange can be conducted, effectively preventing oxidation and eliminating impurities in the chamber by replacing the chamber air with N₂ gas.
- The box has space for installing a vacuum pump, maximizing effectiveness in limited spaces.
- Quick-connect flange piping method simplifies connections.
- In emergency situations, manual air intake is available.

■ Safety

- Equipped with self-diagnostic circuit (temperature sensor abnormality, heater disconnection, automatic overheating prevention, SSR short circuit), independent overheating alarm, three-level authority function, and leakage protection switches to prevent overcurrent.

■ Specifications

| Model | | DPT40C | DPT60C | DPT80C | DPT100C |
|----------------|--|--|----------------|----------------|--|
| System | | Pressure reduction, wall heating | | | |
| Performance | Operating temp. range | Room temperature + 20~200°C | | | |
| | Operating pressure range | 101~0.1KPa (760~1Torr) | | | |
| | Temp. fluctuation (GB standard) | ±0.5°C (at 200°C) | | | |
| | Temp. adjustment accuracy (JTM standard) | ±1°C (at 200°C) | | | |
| Composition | Interior material | Stainless steel plate | | | |
| | Exterior material | Cold-rolled steel plate with surface chemical proofing coating | | | |
| | Insulating material | Glass fiber cotton | | | |
| | Heating method | Heating on the external wall | | | |
| Controllers | Heating power | 2.25KW | 4.48KW | 6.5KW | 14.4KW |
| | Viewing window | Tempered glass + safety protection resin panel | | | |
| | Vacuum pump installation space | Yes | | | |
| | Temp. control method | PID control | | | |
| Specifications | Pressure control method | Deviation control | | | |
| | Temp. /pressure setting/display method | 7-inch color LCD touchscreen for display and settings | | | |
| | Timer | 1 min~99 h 59 min | | | |
| | Operation functions | Fixed temp. operation, program operation, timing operation | | | |
| Safety device | Program mode | 100 programs, each with a maximum of 15 segments | | | |
| | Additional functions | Temperature curve, temperature history curve, alarm records, alarm history records, operation log, data export function, Chinese, Japanese, English three-language switching, calibration function, screen brightness adjustment, three-level authority function | | | |
| | Sensors | K-type thermocouple | | | |
| | | Self-diagnostic circuit (temperature sensor abnormality, heater disconnection, automatic overheating, SSR short circuit), independent overheating alarm, leakage protection switches, three-level authority function, safety valve | | | |
| | Internal dimensions (W×D×H mm) | 450×450×450 | 600×600×600 | 800×800×800 | 1000×1000×1000 |
| | External dimensions (diameter×height mm) | 670×670×1524 | 820×820×1675 | 1020×1020×1875 | 1300×1230×2110 |
| | Internal capacity | 91L | 216L | 512L | 1000L |
| | Shelf layers/shelf support spacing | 4 layers/105mm | 4 layers/140mm | 4 layers/190mm | 8 layers/120mm |
| | Exhaust vent | 25KF flange | | | 40KF flange |
| | Air intake | Rc3/8 | | | Rc1/2 |
| | Power supply (50/60Hz) rated current | AC220V 10.5A | AC220V 14.5A | AC380V 18.5A | AC380V 22A |
| | Weight | Approx. 190kg | Approx. 290kg | Approx. 450kg | Approx. 1000kg |
| Accessories | | Stainless steel punched mesh plate, 2 pieces | | | Stainless steel punched mesh plate, 4 pieces |
| Options | | Shelves, vacuum pump, vacuum pump oil, vacuum pump pre-filter, vacuum pump oil-mist filter, vacuum pump exhaust vent, cold trap, N ₂ pressure reduction valve, N ₂ flow meter, air intake filter, air intake silencer, data logger, combined warning light (standby/running/fault), temperature output terminal (4~20mA), pressure output terminal (4~20mA), external alarm output terminal, pirani vacuum gauge | | | |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Vacuum Ovens | Vacuum Standard

DP23C/33C/43C/63C/83C/103C

Automatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionCE
Certification

DP23C/33C Operating temp. range RT+20~240°C

DP43C/63C/83C/103C Operating temp. range RT+20~200°C

Operating vacuum range 101~0.1KPa

Vacuum drying ovens designed for various applications of vacuum drying.



■ Features

- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Digital setting through special function menu keys and up/down keys. 6 modes with a total of 90 segments program controller.
- Overheating protection, deviation correction and key locking through submenu keys.

■ Safety

- Equipped with self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection, key lock and other safety functions.
- For safety, a protective panel made of resin is installed on the viewing window.

■ Specifications

| Model | | DP23C | DP33C | DP43C | DP63C | DP83C | DP103C | |
|------------------------|--|--|-----------------------|------------------------|------------------------|------------------------|------------------------|--|
| System | | Pressure reduction · wall heating | | | | | | |
| Performance | Operating temp. range | Room temp. +20~240℃ | | | Room temp. +20~200℃ | | | |
| | Operating pressure range | 101~0.1KPa (760~1Torr) | | | | | | |
| | Max. temp. reaching time | Approx. 60 min | Approx. 90 min | Approx. 80 min | Approx. 120 min | | | |
| | Temp. fluctuation (GB standard) | ±0.5℃ (at 240℃) | | | ±0.5℃ (at 200℃) | | | |
| | Temp. adjusting accuracy (JTM standard) | ±1.5℃ (at 240℃) | | | ±1.0℃ (at 200℃) | | | |
| Composition | Interior material | Stainless steel plate | | | | | | |
| | Exterior material | Cold rolled steel plate with chemical proofing coating | | | | | | |
| | Insulating material | Glass fiber | | | | | | |
| | Heating method | Heating on the external wall | | | | | | |
| | Heating power | 0.68KW | 1.05KW | 2.25KW | 3.15KW | 3.96KW | 14.4KW | |
| | Viewing window | Tempered glass + safety protection resin panel | | | | | | |
| | Vacuum gauge | Pointer type, -100~0 KPa | | | | | | |
| Controllers | Vacuum pump installation space (W×D×H mm) | - | - | 320×600×540 | 470×750×540 | 635×930×540 | 915×1120×570 | |
| | Temp. control method | 3-stage PID | | | | | | |
| | Temp. setting method | Digital setting through special function menu keys and up/down keys | | | | | | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display | | | | | | |
| | | Setting temp. display: Red 4-digit LED digital display | | | | | | |
| | Timer | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | | | | |
| | Operation functions | Fixed temp. operation, auto start, auto stop, program operation | | | | | | |
| | Program mode | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | | | | | | |
| | Additional functions | Deviation correction, key lock, power failure compensation | | | | | | |
| Heater circuit control | SSR drive | | | | | | | |
| Sensors | K thermocouple (temp. controller and overheating protection) | | | | | | | |
| Safety device | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | | | | | | |
| Specifications | Internal dimensions (W×D×H mm) | 200×250×200 | 300×300×300 | 450×450×450 | 600×600×600 | 800×800×800 | 1000×1000×1000 | |
| | External dimensions (W×D×H mm) | 400×410×682 | 510×460×782 | 670×670×1524 | 820×820×1675 | 1020×1020×1875 | 1300×1228×2110 | |
| | Internal capacity | 10L | 27L | 91L | 216L | 512L | 1000L | |
| | Shelf layer/spacing | 3 layers (fixed)/63mm | 4 layers (fixed)/71mm | 4 layers (fixed)/105mm | 4 layers (fixed)/140mm | 4 layers (fixed)/190mm | 8 layers (fixed)/120mm | |
| | Exhaust vent | External diameter 18mm | | | KF25 | | KF40 | |
| | Air intake | External diameter 18mm | | | Rc3/8 | | Rc1/2 | |
| | Power supply (50/60Hz) rated current | AC220V 3.5A | AC220V 5A | AC220V 10.5A | AC220V 14.5A | AC220V 18.5A | AC380V 22A | |
| | Weight | Approx. 43kg | Approx. 69kg | Approx. 190kg | Approx. 290kg | Approx. 450kg | Approx. 1000kg | |
| Accessories | | Stainless steel punched mesh plate, 2 pieces | | | | | | Stainless steel punched mesh plate, 4 pieces |
| Options | Stand | ONS10C | | | — | | | |
| | Others | Shelves, vacuum pump, N ₂ introduction device, data logger, combined warning light (standby/running/fault), external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software | | | | | | |

Vacuum Ovens | Economical, Table Type

ADP210C/310C

Automatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification

Operating temp. range 40~240℃

Operating vacuum range 101~0.1KPa

Internal capacity 10L 27L

Vacuum drying oven with enhanced and cost-effectiveness features.

Features

- The door sealing strip is made of silicone in one piece, ensuring sealing performance. Easy operation allows for fixed temp. operation, program operation, rapid auto stop, auto stop, and auto start.
- Digital settings can be achieved through a dedicated operation function menu key and ▼▲ keys. Able to set overheating protection, deviation correction, and key lock.

Safety

- Equipped with self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection, key lock safety and other safety functions.
- For safety, a protective panel made of resin is installed on the viewing window.



Specifications

| Model | | ADP210C | ADP310C |
|----------------|---|--|---------------------------------------|
| System | | Pressure reduction · wall heating | |
| Performance | Operating temp. range | 40~240℃ | |
| | Operating pressure range | 101~0.1KPa (760~1Torr) | |
| | Max. temp. reaching time | Approx. 70 min | Approx. 100 min |
| | Temp. fluctuation (GB standard) | ±0.5℃ (at 200℃) | |
| | Temp. adjusting accuracy (JTM standard) | ±1.5℃ (at 200℃) | |
| Composition | Interior material | Stainless steel plate | |
| | Exterior material | Cold rolled steel plate with chemical proofing coating | |
| | Insulating material | Glass fiber | |
| | Heating method | Heating on the external wall | |
| | Heating power | 0.68KW | 1.05KW |
| | Viewing window | Tempered glass + safety protection resin panel | |
| | Vacuum gauge | Pointer type, -100~0 KPa | |
| Controllers | Temp. control method | PID control | |
| | Temp. setting method | Digital setting via function menu key and ▲ ▼ keys | |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display | |
| | | Setting temp. display: Red 4-digit LED digital display | |
| | Timer | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | |
| | Operation functions | Fixed temp. operation, auto start, quick auto stop, program operation | |
| | Additional functions | Deviation correction, key lock, power failure compensation | |
| | Heater circuit control | SSR drive | |
| Sensors | Temperature controller: Pt100 thermistor; overheating protection: liquid expansion temperature controller | | |
| Safety device | | Self-diagnostic circuit (temperature sensor anomaly, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection, key lock functions | |
| Specifications | Internal dimensions | 200×250×200 | 300×300×300 |
| | External dimensions (W×D×H mm) | 400×412×603 | 500×465×705 |
| | Internal capacity | 10L | 27L |
| | Shelf layers | 2 layers (fixed) | 3 layers (fixed) |
| | Shelf spacing | 65mm | 75mm |
| | Exhaust vent | External diameter 18mm | |
| | Air intake | External diameter 18mm | |
| | Power supply (50/60Hz) rated current | AC220V 3.5A | AC220V 5A |
| | Weight | Approx. 30kg | Approx. 55kg |
| Accessories | | Aluminum stamped mesh plate, 2 pieces | Aluminum stamped mesh plate, 3 pieces |
| Options | | Stand (ONS10C), shelves, vacuum pump, micro printer, data logger, centralized monitoring software | |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Vacuum Ovens | Vertical, High Precision

DPF43C/63C

Automatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionCE
Certification

Specifications

| Model | DPF43C | DPF63C |
|--------------------------------------|--|-------------------|
| System | Pressure reduction · internal wall heating | |
| Operating temp. range | 40~200°C | |
| Operating pressure range | 101~0.1KPa (760~1Torr) | |
| Max. temp. reaching time | Approx. 45 min | Approx. 60 min |
| GB standard | Temp. fluctuation | ±0.5°C (at 200°C) |
| | Temp. uniformity | ±2.5% (at 200°C) |
| JTM standard | Temp. adjusting accuracy | ±1.0°C (at 200°C) |
| | Temp. distribution accuracy | ±5°C (at 200°C) |
| Heating method | Heating on the internal wall | |
| Heating power | 4.2KW | 5.4KW |
| Vacuum gauge | Pointer type, -100~0KPa | |
| Internal dimensions (W×D×H mm) | 380×390×380 | 530×540×530 |
| Effective volume | 56L | 151L |
| External dimensions (W×D×H mm) | 670×670×1524 | 820×820×1675 |
| Internal capacity | 91L | 216L |
| Exhaust vent | KF25 | |
| Air intake | Rc3/8 | |
| Power supply (50/60Hz) rated current | AC220V 20A | AC220V 26A |
| Weight | Approx. 290kg | Approx. 450kg |

Vacuum Ovens | Vertical, 300°C

DPH43C/63C

Automatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Specifications

| Model | DPH43C | DPH63C |
|---|--|-----------------|
| System | Pressure reduction · external wall heating | |
| Operating temp. range | 40~300°C | |
| Operating pressure range | 101~0.1KPa (760~1Torr) | |
| Max. temp. reaching time | Approx. 80 min | Approx. 120 min |
| Temp. fluctuation (GB standard) | ±0.5°C (at 300°C) | |
| Temp. adjusting accuracy (JTM standard) | ±1.0°C (at 300°C) | |
| Heating method | Heating on the external wall | |
| Heating power | 3.2KW | 4KW |
| Vacuum gauge | Pointer type, -100~0 KPa | |
| Internal dimensions (W×D×H mm) | 450×450×450 | 600×600×600 |
| External dimensions (W×D×H mm) | 670×670×1524 | 820×820×1675 |
| Internal capacity | 91L | 216L |
| Exhaust vent | KF25 | |
| Air intake | Rc3/8 | |
| Power supply (50/60Hz) rated current | AC220V 15A | AC220V 18.5A |
| Weight | Approx. 190kg | Approx. 290kg |

Vacuum Ovens | Vertical, 400°C

DPHH43C/63C

Automatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Specifications

| Model | DPHH43C | DPHH63C |
|---|--|----------------|
| System | Pressure reduction · internal wall heating | |
| Operating temp. range | 40~400°C | |
| Operating pressure range | 101~0.1KPa (760~1Torr) | |
| Max. temp. reaching time | Approx. 75 min | Approx. 90 min |
| Temp. fluctuation (GB standard) | ±0.5°C (at 400°C) | |
| Temp. adjusting accuracy (JTM standard) | ±1.0°C (at 400°C) | |
| Heating method | Heating on the internal wall | |
| Heating power | 4KW | 5.2KW |
| Vacuum gauge | Pointer type, -100~0KPa | |
| Internal dimensions (W×D×H mm) | 450×450×450 | 600×600×600 |
| External dimensions (W×D×H mm) | 820×820×1675 | 1020×1020×1875 |
| Internal capacity | 91L | 216L |
| Exhaust vent | KF25 | |
| Air intake | Rc3/8 | |
| Power supply (50/60Hz) rated current | AC220V 18.5A | AC220V 24A |
| Weight | Approx. 290kg | Approx. 450kg |

Vacuum Ovens | Vertical, Shelf Heating

DP43HPC/63HPC

Automatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Specifications

| Model | | DP43HPC | DP63HPC |
|---|--------------|--|--------------|
| System | | Pressure reduction · external wall heating · shelf heating | |
| Operating temp. range | | 40~200℃ | |
| Operating pressure range | | 101~0.1KPa (760~1Torr) | |
| Temp. fluctuation (GB standard) | | ±0.5℃ (at 200℃) | |
| Temp. adjusting accuracy (JTM standard) | | ±1.0℃ (at 200℃) | |
| Heating method | | Heating on the external wall and the shelf | |
| Heating power | Wall surface | 2.25KW | 3.15KW |
| | Shelf | 1KW/pcs | 1.2KW/pcs |
| Vacuum gauge | | Pointer type, -100~0 KPa | |
| Control method | | Independent PID control of wall heating and each shelf heating | |
| Internal dimensions | | 450×450×450 | 600×600×600 |
| External dimensions (W×D×H mm) | | 670×669×1500 | 820×819×1650 |
| Internal capacity | | 91L | 216L |
| Exhaust vent | | KF25 | |
| Air intake | | Rc3/8 | |
| Power supply (50/60Hz) | | Single phase AC220V | |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options



Contents

| | |
|-----------------------------------|-----|
| High-temperature Incubators | 108 |
| Low-temperature Incubators | 114 |
| CO ₂ Incubators | 122 |

5 | Incubators



1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Incubator product series to meet various customer needs

The series of constant temperature incubators perform excellently in terms of high-precision temperature control and cultivation stability. The various product series not only control the pH of the culture medium but also meet various requirements of customers for capacity and usage temperature range.

Function · Safety Device

| | | | |
|--|--|---|--|
| Forced convection circulation | By means of forced circulation of high-performance centrifugal fan blades, stir the chamber inside to obtain uniform temperature distribution accuracy. | Self-diagnosis function | By means of the microcomputer carried in controller, detect the abnormality of circuit, if abnormal occurs, control the device within safety range and the alarm sounds. |
| Natural convection | By means of heater or air jacket natural convection in chamber. | Key lock function | Prevent the misoperation during operation. |
| Auto overheat prevention function | Built-in overheat prevention function of controller, usually when constant temp. oven or drying oven chamber temp. reaches set temp. + 12°C, and constant temp. bath temp. reaches set temp. + 6°C, the heater cuts off (auto recovery). | Backup function | Even in power outage or cut off power source, the set value can be memorized. |
| Overheat protector | Controller and overheat protector are all-in-one. They share the same power source, but other circuits are standalone. If overheat abnormally, cut off (manual recovery) heater circuit. | Power outage compensation | When the power recovers, it's able to select to stop or continue operation. |
| Standalone overheat protector | Different circuit from controller, if overheat abnormally, cut off (manual recovery) heater circuit. According to different models, respectively carry digital, liquid-expansion, thermometal type, etc. | Overcurrent breaker | When abnormal current passes, cut off the power to protect the unit. |
| | | Overcurrent electric leakage breaker | The breaker has both overcurrent cut-off function and electric leakage cut-off function. |
| | | Emergency stop button | If need to emergency stop, press the emergency stop button to cut off. |

Model selection list

| Category | Temp. control range | Temp. fluctuation | Temp. uniformity | Series | Models | Features | Capacity (L) | | | | | | | Page No. | | |
|----------------------------|---------------------|-----------------------------------|-------------------------|--------|------------------------|---|---|-------------------------------------|--------|---------|---------|---------|---------|----------|-----|-----|
| | | | | | | | 0-20 | 20-40 | 40-100 | 100-150 | 150-200 | 200-400 | 400-600 | | | |
| High-temp. incubation | Room temp. +5~80°C | ±0.1°C | ±1.5°C | IS | IS412C/612C/812C/912C | Programmable Incubators with Air Jacket | | | 97 | | 159 | 318 | 567 | 108 | | |
| | Room temp. +5~80°C | ±0.1°C | ±1.5°C | IC | IC412C/612C/812C/912C | Constant temperature operation, air Jacket type | | | 97 | | 159 | 318 | 567 | 110 | | |
| | Room temp. +5~80°C | ±0.1°C | ±1.0°C | ICF | ICF410C/610C/810C/910C | High temperature uniformity | | | 89 | | 147 | 294 | 524 | 112 | | |
| Low-temp. incubation | 1-chamber | -10~60°C | ±0.5°C | ≤2.0°C | IN | IN613C/613CW/813C/913C | Program operation | | | | | 143 | 286 | 624 | 114 | |
| | | 0~60°C | - | - | INE | INE800 | Inverter energy-saving, low-temperature type | | | | | | 286 | | 116 | |
| | | 0~60°C | ±0.3°C | ±1.0°C | IL | IL612C/812C | Air Jacket type | | | | | 159 | 300 | | 117 | |
| | | 0~60°C | - | - | IJ | IJ102/102W/300 | Semiconductor cooling | 15.6 | | 43 | | | | | 118 | |
| | 2-chamber | Upper chamber: Room temp. +5~80°C | ±0.1°C | ±1.5°C | INC | INC822C | High-temperature + low-temperature with 2-chamber | | | | 143 | | | | 120 | |
| | | Lower chamber: -10~60°C | ±0.5°C | ≤2.0°C | | | | | | | 150 | | | | | |
| | | | Upper chamber: -10~60°C | ±0.1°C | ±1.0°C | IQ | IQ823C | Low-temperature type with 2-chamber | | | | 143 | | | | 121 |
| | | | Lower chamber: -10~60°C | ±0.1°C | ±1.0°C | | | | | | | 143 | | | | |
| CO ₂ incubation | Room temp. +5~60°C | ±0.1°C | ±0.3°C | IPE | IPE610 | Air Jacket type, dry heat sterilization | | | | | 185 | | | | 122 | |
| | Room temp. +5~50°C | ±0.1°C | ±0.2°C | BNA | BNA610 | Water jacket type | | | | | 167 | | | | 123 | |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

High-temp. Incubators | Air Jacket

IS412C/612C/812C/912C

Operating temp. range RT+5~80°C

Temp. distribution accuracy $\pm 1.0^{\circ}\text{C}$ (at 37°C)

Internal capacity 90L 159L 318L 567L

From small capacity to large capacity, high-function, long-term hot sale constant temperature incubators.

Features

- Due to the air jacket heat transfer method, the temperature distribution is uniform.
- Double door structure reduces heat loss, maintaining a stable constant temperature state.
- Inner door is made of reinforced glass, allowing observation of the interior.
- Easy operation, available for fixed temp., program, quick auto stop, auto stop and auto start operations.
- Use special function menu keys and up/down keys to realize digital setting. With repeat function, 6 modes in total of 90 segments.
- Use submenu keys to operate overheat protector, deviation correction and key locking.

Safety

- Equipped with self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection switch, key lock safety and other safety functions.



Specifications

| Model | | | IS412C | IS612C | IS812C | IS912C |
|----------------|--------------------------------------|--|--|--|---------------|---------------|
| System | | | Airflow natural convection | | | |
| Performance | Operating temp. range | | Room temp. +5~80℃ | | | |
| | GB standard | Temp. fluctuation | ±0.5℃ (at 37℃) | | | |
| | | Temp. uniformity | ±1.5℃ (at 37℃) | | | |
| | JTM standard | Temp. adjusting accuracy | ±0.2℃ (at 37℃) | | | |
| | | Temp. distribution accuracy | ±1.0℃ (at 37℃) | | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber | | | |
| | Heater | | Nichrome heating wire | | | |
| | | | 0.3KW | 0.4KW | 0.7KW | 2.2KW |
| Exhaust vent | | I.D. 30mm×2 (top) | | I.D. 30mm×2 (one at the right and one at the left) | | |
| Controllers | Temp. control method | | PID control | | | |
| | Temp. setting method | | Use special function menu keys and up/down keys to realize digital setting | | | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display | | | |
| | | | Setting temp. display: Red 4-digit LED digital display | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, auto stop, program operation | | | |
| | Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | | | |
| | Additional functions | | Deviation correction, key lock, power outage compensation | | | |
| Sensors | | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple | | | | |
| Safety device | | | Self-diagnostic function (temperature sensor anomaly, SSR short circuit, heater disconnection, automatic overheat prevention), key lock function, overheat prevention device, overcurrent leakage protection switch | | | |
| Specifications | Internal dimensions | | 450×480×450 | 600×530×500 | 600×530×1000 | 1070×530×1000 |
| | External dimensions (W×D×H mm) | | 560×606×820 | 710×656×870 | 710×656×1619 | 1180×655×1619 |
| | Internal capacity | | 97L | 159L | 318L | 567L |
| | Shelf load | | Approx. 15kg/shelf | | | |
| | Shelf layers | | 11 layers | 13 layers | 29 layers | 29 layers×2 |
| | Support spacing | | 30mm | | | |
| | Power supply (50/60Hz) rated current | | AC220V 1.5A | AC220V 2A | AC220V 3.5A | AC220V 6.5A |
| Weight | | Approx. 45kg | Approx. 65kg | Approx. 102kg | Approx. 166kg | |
| Accessories | Shelf | | Stainless punching mesh plate | | | |
| | | | 2 pcs | | 4 pcs | 8 pcs |
| | Supports | | 4 pcs | | 8 pcs | 16 pcs |
| Options | Stand | | ON61C | | — | |
| | Stacking fittings | | OD40C | OD60C | — | |
| | Others | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), micro printer, data logger, combined warning light (standby/operation/fault), viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, timer output terminal, central monitoring software, touchscreen controller | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Natural
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Control panel



Exhaust vent (IS412C/612C)




Internal chamber (IS612C)



Shelf · supports



Dimension diagram (mm)



Sterilizers 1

Granulation
and Spray
Dryers 2Muffle
Furnaces 3

Ovens 4

Incubators 5

Plasma
Equipment 6Water
Purifiers 7

Baths 8

Water
Circulators 9Rotary
Evaporators 10Freeze Dryers
& Cold Traps 11Stirrers &
Shakers 12

Washers 13

Analysis and
Test Devices 14

Options 15

High-temp. Incubators | Fixed Temp. Operation, Air Jacket

IC412C/612C/812C/912C

Operating temp. range RT+5~80°C

Temp. distribution accuracy $\pm 1.0^{\circ}\text{C}$ (at 37°C)

Internal capacity 97L 159L 318L 567L

From small capacity to large capacity, long-term hot sale constant temperature operation type incubators.

Features

- Due to the air jacket heat transfer method, the temperature distribution is uniform.
- Double door structure reduces heat loss, maintaining a stable constant temperature state.
- The inner door is made of reinforced glass, allowing observation of the interior.
- Easy to operate, fixed temp. operation, quick auto stop operation, auto stop operation, and auto start operation are achievable.
- Deviation correction setting and key lock setting are available.

Safety

- Equipped with self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection switch, key lock and other safety functions.



Specifications

| Model | | IC412C | IC612C | IC812C | IC912C | |
|-----------------------------|--------------------------------------|---|---|--|---------------|---------------|
| System | | Airflow natural convection | | | | |
| Performance | Operating temp. range | | Room temp. +5~80℃ | | | |
| | GB standard | Temp. fluctuation | ±0.5℃ (at 37℃) | | | |
| | | Temp. uniformity | ±1.5℃ (at 37℃) | | | |
| | JTM standard | Temp. adjusting accuracy | ±0.2℃ (at 37℃) | | | |
| Temp. distribution accuracy | | ±1.0℃ (at 37℃) | | | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber | | | |
| | Heater | | Nichrome heating wire | | | |
| | | 0.3KW | 0.4KW | 0.7KW | 1.2KW | |
| Exhaust vent | | I.D. 30mm×2 (top) | | I.D. 30mm×2 (one at the right and one at the left) | | |
| Controllers | Temp. control method | | PID control | | | |
| | Temp. setting method | | Use special function menu keys and up/down keys to realize digital setting | | | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display Setting temp. display: Red 4-digit LED digital display | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, quick auto stop | | | |
| | Additional functions | | Deviation correction, key lock, power outage compensation | | | |
| | Sensors | | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple | | | |
| Safety device | | Self-diagnostic function (temperature sensor anomaly, SSR short circuit, heater disconnection, automatic overheat prevention), key lock function, overheat prevention device, overcurrent leakage protection switch | | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 450×480×450 | 600×530×500 | 600×530×1000 | 1070×530×1000 |
| | External dimensions (W×D×H mm) | | 560×606×820 | 710×656×870 | 710×656×1600 | 1180×655×1619 |
| | Internal capacity | | 97L | 159L | 318L | 567L |
| | Shelf load | | Approx. 15kg/shelf | | | |
| | Shelf layers | | 11 layers | 13 layers | 29 layers | 29 layers×2 |
| | Support spacing | | 30mm | | | |
| | Power supply (50/60Hz) rated current | | AC220V 1.5A | AC220V 2A | AC220V 3.5A | AC220V 6.5A |
| Accessories | Weight | | Approx. 45kg | Approx. 65kg | Approx. 102kg | Approx. 166kg |
| | Shelf | | Stainless punching mesh plate | | | |
| | | | 2 pcs | 4 pcs | 8 pcs | |
| Options | Supports | | 4 pcs | 8 pcs | 16 pcs | |
| | Stand | | ON61C | | | |
| | Stacking fittings | | OD40C | OD60C | — | |
| | Others | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), data logger, combination warning lights (standby/running/fault), viewing window, external communication function (RS485), temperature output terminals (4~20mA), external alarm output terminals, time-up output terminals | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Natural
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification

Control panel



Exhaust vent (IC412C/612C)




Interior chamber (IC612C)



Shelf · supports



Dimension diagram (mm)



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

High-temp. Incubators | Fixed Temp. Operation, Air Jacket Type, High Temp. Uniformity

ICF410C/610C/810C/910C

Operating temp. range RT+5~80°C

Temp. distribution accuracy $\pm 0.5^{\circ}\text{C}$ (at 37°C)

Internal capacity 89L 147L 294L 524L

Higher temperature uniformity, long-term hot sale constant temperature operation type incubators.



Features

- Achieved higher precision temperature control and balanced internal chamber distribution through forced circulation by an internal fan.
- Double door structure reduces heat loss, maintaining a stable constant temperature state.
- Inner door is made of reinforced glass, allowing observation of the interior.
- Easy to operate, fixed temp. operation, quick auto stop operation, auto stop operation, and auto start operation are achievable.
- Deviation correction setting and key lock setting are available.

Safety

- Equipped with self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating prevention, SSR short-circuit protection), overheating prevention, overcurrent leakage protection switch, key lock and other safety functions.

Specifications

| Model | | | ICF410C | ICF610C | ICF810C | ICF910C |
|-----------------------------|--------------------------------------|--|--|--|---------------|---------------|
| System | | | Forced convection | | | |
| Performance | Operating temp. range | | Room temp. +5~80℃ | | | |
| | GB standard | Temp. fluctuation | ±0.1℃ (at 37℃) | | | |
| | | Temp. uniformity | ±1℃ (at 37℃) | | | |
| | JTM standard | Temp. adjusting accuracy | ±0.1℃ (at 37℃) | | | |
| Temp. distribution accuracy | | ±0.5℃ (at 37℃) | | | | |
| Composition | Interior material | | Stainless steel plate SUS304 | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Insulating material | | Glass fiber | | | |
| | Heater | | Nichrome heating wire | | | |
| | | | 0.36KW | 0.48KW | 0.88KW | 1.45KW |
| | Fan | | 1 pcs | | 2 pcs | 4 pcs |
| Exhaust vent | | I.D. 30mm×2 (top) | | I.D. 30mm×2 (one at the right and one at the left) | | |
| Controllers | Temp. control method | | PID control | | | |
| | Temp. setting method | | Use special function menu keys and up/down keys to realize digital setting | | | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display Setting temp. display: Red 4-digit LED digital display | | | |
| | Timer | | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) | | | |
| | Operation functions | | Fixed temp. operation, auto start, quick auto stop | | | |
| | Additional functions | | Deviation correction, key lock, power outage compensation | | | |
| Sensors | | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple | | | | |
| Safety device | | | Self-diagnostic function (temperature sensor anomaly, SSR short circuit, heater disconnection, automatic overheat prevention), key lock function, overheat prevention device, overcurrent leakage protection switch | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 450×440×450 | 600×490×500 | 600×490×1000 | 1070×490×1000 |
| | External dimensions (W×D×H mm) | | 560×606×820 | 710×656×870 | 710×656×1619 | 1180×655×1619 |
| | Internal capacity | | 89L | 147L | 294L | 524L |
| | Shelf load | | Approx. 15kg/shelf | | | |
| | Shelf layers | | 9 layers | 12 layers | 29 layers | 29 layers×2 |
| | Support spacing | | 30mm | | | |
| | Power supply (50/60Hz) rated current | | AC220V 2A | AC220V 2.5A | AC220V 4.5A | AC220V 7A |
| Weight | | Approx. 45kg | Approx. 65kg | Approx. 102kg | Approx. 166kg | |
| Accessories | Shelf | | Stainless punching mesh plate | | | |
| | | | 2 pcs | | 4 pcs | 8 pcs |
| | Supports | | 4 pcs | | 8 pcs | 16 pcs |
| Options | Stand | | ON61C | | — | |
| | Stacking fittings | | OD40C | OD60C | — | |
| | Others | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), micro printer, data logger, combined warning light (standby/operation/fault), viewing window, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, timer output terminal, central monitoring software, touchscreen controller | | | |

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerCE
Certification

Control panel



Exhaust vent (ICF410C/610C)




Interior chamber (ICF810C)



Fan switch (built-in)



Dimension diagram (mm)



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Low-temp. Incubators | Program Operations

IN613C/613CW/813C/913C

Operating temp. range -10~60°C

Temp. fluctuation ±0.5°C

Temp. uniformity ≤2.0°C (at 37°C)

Internal capacity 143L 286L 624L

Widely applicable for various constant temperature tests and environmental experiments.

Features

- Widely applicable for various constant temperature tests and environmental tests.
- Achieved high precision temperature control and balanced internal chamber temperature distribution through forced circulation by a fan blade.
- The outer door uses a large hollow glass and thickened tempered glass viewing window structure design, better observation of the interior while improving insulation. (IN613CW type)
- The inner glass door enhances the insulation effect. Observing the sample through the inner glass door minimizes temperature changes.
- Defrosting can be set to cycle defrost based on the continuous operation time of the compressor.
- After power outage recovery, the operation can automatically restart under set temperature conditions through the power outage compensation function.
- Equipped with safety functions such as compressor overload relay, overheat prevention device, overcurrent leakage protection switch, self-diagnostic circuit (temperature sensor anomaly, heater disconnection, SSR short circuit, automatic overheat prevention), and buzzer alarm during anomalies.




Specifications

| Model | | | IN613C | IN613CW | IN813C | IN913C |
|----------------------|--------------------------------------|--|---|--|--|--|
| System | | | Forced convection | | | |
| Performance | Operating temp. range | | -10~60℃ | | | |
| | GB standard | Temp. fluctuation | ±0.5℃ | | | |
| | | Temp. uniformity | ≤2.0℃ (at 37℃), ≤3.0℃ (other working temperatures) | | | |
| | | Heating time | 20~52℃≤20min | | | |
| | | Cooling time | 20~-4℃≤60min | | | |
| Composition | Interior material | | Stainless steel plate | | | |
| | Exterior material | | Cold rolled steel plate with chemical proofing coating | | | |
| | Viewing window | | - | Width 516×height 416 mm | - | |
| | Insulating material | | Foamed polystyrene | | | |
| | Refrigerator | | 158W | | 300W | 519W |
| | Refrigerant | | R134A | | | |
| | Defrost structure | | Cycle operation | | | |
| | Fan blade | | Crossflow fan | | | |
| | Heater | | Nickel-chromium alloy heating wire, 550W | | Nickel-chromium alloy heating wire, 750W | Nickel-chromium alloy heating wire, 750W×2 |
| | Sensors | | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple | | | |
| | Cable port | | Internal diameter 30mm (one on the right side) | Internal diameter 50mm (one on the right side) | Internal diameter 30mm (one on the right side) | |
| | Temp. control method | | PID control | | | |
| | Temp. setting method | | Use special function menu keys and up/down keys to realize digital setting | | | |
| | Temp. display method | | Achieved temp. display: Green 4-digit LED digital display Setting temp. display: Red 4-digit LED digital display | | | |
| | Timer/timer resolution | | 1 min~99 h 59 min or 100~999 h 50 min | | | |
| Operating functions | | Fixed temp. operation, auto start, auto stop, program operation | | | | |
| Program mode | | Program operation 6 modes with a total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) | | | | |
| Additional functions | | Temperature compensation function, key lock function, power outage compensation function | | | | |
| Safety device | | | Compressor overload relay, overheat prevention device, overcurrent leakage protection switch, self-diagnostic circuit (temperature sensor anomaly, heater disconnection, SSR short circuit, automatic overheat prevention), buzzer alarm during anomalies | | | |
| Specifications | Internal dimensions (W×D×H mm) | | 600×477×500 | | 600×477×1000 | 1310×477×1000 |
| | External dimensions (W×D×H mm) | | 710×645×915 | | 710×645×1630 | 1420×645×1630 |
| | Internal capacity | | 143L | | 286L | 624L |
| | Shelf load | | 15kg/layer | | | |
| | Shelf layers/shelf support spacing | | 13 layers/30mm | | 23 layers/30mm | |
| Accessories | Power supply (50/60Hz) rated current | | AC220V 50Hz 3.5A | | AC220V 50Hz 5A | AC220V 10A |
| | Weight | | Approx. 89kg | | Approx. 115kg | Approx. 230kg |
| | Shelf | | Stainless punching mesh plate 3 pcs | | 5 pcs | 8 pcs |
| Options | Supports | | 6 pcs | | 10 pcs | 16 pcs |
| | Stand | | ON61C | | - | - |
| | Stacking fittings | | OD60C | | - | - |
| Others | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), micro printer, data logger, combined warning light (standby/operation/fault), viewing window, external communication function (RS485), temperature output terminal (4-20mA), server socket, central monitoring software, external alarm output terminal, timer output terminal (one-out-of-two), touchscreen controller | | | | |

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker**Interior chamber (IN613C)**

[IN613C]

[IN813C]

Control panel**Temperature fall characteristic curve (IN613C)****Micro printer (optional)****Shaker can be installed (IN613CW)****External output****Shelf · supports****Central monitoring software (optional)**

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Low-temp. Incubators | Inverter Energy-saving, Low-temp. Type

INE800

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerMade in
Japan

Operating temp. range 0~60°C

Temp. distribution accuracy $\pm 0.5^{\circ}\text{C}$ (at 37°C)

Internal capacity 286L

Achieve approximately 44% energy savings through frequency conversion, with a significant reduction in frosting.



Features

Achieve approximately 44% energy savings through frequency conversion (compared with our company's previous products), with a significant reduction in frosting.

- Equipped with frequency conversion function, the refrigerator power can be adjusted to minimize unnecessary power consumption during operation.
- By reducing unnecessary cooling capacity through frequency conversion, it adjusts to the most suitable cooling capacity, thereby highly avoiding the frosting phenomenon.
- The temperature range is 0~60°C, providing the optimal temperature environment for various cultures and preservation.
- To create a more ideal cultivation temperature environment, the temperature distribution accuracy within chamber is improved.
- Standard features include program operation, automatic operation stop, automatic cooling operation start, temperature deviation automatic correction function, and so on. The controller can display power consumption and CO₂ emissions.
- Standard features include temperature analog output (4~20mA), external communication terminal (RS485).

Specifications

| Model | | INE800 |
|----------------|--|---|
| System | | Forced convection |
| Performance | Operating temp. range | 0~60°C |
| | Set temp. range | -5~65°C |
| | Temp. adjusting accuracy | $\pm 0.2^{\circ}\text{C}$ (at 37°C with continuous operation of the refrigerator) $\pm 0.5^{\circ}\text{C}$ (at 37°C with cycle operation of the refrigerator) |
| | Temp. distribution accuracy | $\pm 0.5^{\circ}\text{C}$ (at 37°C with continuous operation) |
| | Max. temp. reaching time | 35min (20~65°C) |
| | Min. temp. reaching time | 50min (20~0°C) |
| | Refrigerator operation modes | Continuous operation, cycle operation, stop operation |
| Composition | Interior material | Stainless steel plate |
| | Exterior material | Electro-galvanized steel plate with chemical resistance coating on the surface |
| | Insulating material | Foamed polyethylene |
| | Refrigerator/refrigerant/action range | 200W rotary compressor/R134A/set below 40°C |
| | Defrost function | Manual defrost (at will), automatic defrost (time/moment) |
| | Air blower fan | DC axial flow fan |
| | Heater | 750W |
| | Sensors | Pt resistance thermometer (temperature controller), K-type thermocouple (overheat prevention) |
| | Cable port | Internal diameter $\Phi 50$ (right side of the body) |
| | Temp. control method | PID control |
| Controllers | Temp. regulator | V-type control |
| | Temp. control method | PID control |
| | Temp. setting method | Digital settings achieved through function menu key and \blacktriangle \blacktriangledown keys |
| | Operation functions | Constant operation, quick automatic stop operation, automatic stop operation, automatic start operation, program operation (30 segments \times 1, 15 segments \times 2, 10 segments \times 3) |
| | Additional functions | Deviation correction, key lock, power outage |
| Safety device | Sensors | Pt resistance thermometer (temperature controller), K-type thermocouple (overheat prevention) |
| | | Self-diagnosis function (temperature sensor abnormality, heating wire disconnection, SSR short circuit, main relay malfunction, automatic overheat prevention function), key lock function |
| Specifications | Internal dimensions (W \times D \times H mm) | 600 \times 477 \times 1000 |
| | External dimensions (W \times D \times H mm) | 710 \times 645 \times 1730 |
| | Internal capacity | 286L |
| | Shelf load | Approx. 15kg/shelf |
| | Shelf layers/shelf support spacing | 23 layers/30mm |
| | Power supply (50/60Hz) rated current | AC100V 10A |
| | Weight | Approx. 135kg |
| | Shelf/supports | 5 stainless steel stamping mesh plates/10 pieces |
| | Others | 2 keys, 1 silicone plug for the cable port |
| | Options | Shelf plate (1 shelf plate with 2 shelf supports), external alarm output terminal, time arrival output terminal |


Control panel



Independent overheating preventer



Heating and cooling curve



1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Low-temp. Incubators | Air Jacket Type

IL612C/812C

Natural
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range 0~60°C

Temp. distribution accuracy ±1.0°C (at 37°C)

Internal capacity 159L 300L

Multi-functional, low-temperature incubator with air jacket.

Features

- Using environmentally friendly refrigerant R134A.
- Achieved high-precision control and uniform temperature distribution accuracy within the chamber through air jacket heat transfer and computer PID control.
- Include a controller with an independent overheat prevention circuit, achieving a variety of functions centered on constant value operation, automatic stop operation, automatic start operation, defrost function, RS485 communication function, temperature output terminal (4~20mA), and alarm output terminal. (IL612C)
- Equipped with an electronic independent overheat prevention device, featuring program operation up to 99 segments, constant value operation, automatic stop operation, automatic start operation, CR5 type controller with RS485 communication function, temperature output terminal (1~5V), and alarm output terminal. (IL812C)
- Remote operation can be achieved by using an optional communication cable to separate the control panel from the main body.

Safety

- Equipped with safety functions such as self-diagnosis function, refrigerator thermal overload protector, refrigerator delay start protection function, overcurrent leakage protection switch, key lock, etc.



Specifications

| Model | | | IL612C | IL812C |
|----------------|--------------------------------------|---|---|--|
| System | | | Conduction, radiation | |
| Performance | Operating temp. range | | 0~60℃ | |
| | GB standard | Temp. fluctuation | ±0.3℃ (at 37℃) | |
| | | Temp. uniformity | ±1.0℃ (at 37℃) | |
| | JTM standard | Temp. adjusting accuracy | ±0.3℃ (at 37℃) | |
| | | Temp. distribution accuracy | ±1.0℃ (at 37℃) | |
| Composition | Inner door | | Tempered glass 5mm | |
| | Insulating material | | Polyurethane foam | |
| | Refrigerator/refrigerant | | 300W/R134A | 350W/R134A |
| | Defrost structure | | Manual defrost/cycle defrost | |
| | Heater | | Nickel-chromium alloy heating wire, 800W | Nickel-chromium alloy heating wire, 850W |
| | Sensors | | Double sensor cable, platinum resistance thermometer PT100 (for temperature control) + K-type thermocouple (for overheat prevention) | |
| | Overheat preventer | | Electronic integrated control type | Electronic independent type |
| | Door lock | | 1 pcs | |
| | Cable port | | Internal diameter 50mm (right side of the body) | |
| | Temp. regulator | | VS6P | CR5 PID control |
| Operation | | Constant value, timed | | |
| Safety device | | Self-diagnosis function, overcurrent leakage circuit breaker, refrigerator overload relay, refrigerator protection delay timer function, overheat prevention device (IL612C), independent overheat prevention device (IL812C) | | |
| Specifications | Internal dimensions (W×D×H mm) | | 600×530×500 | 600×530×1000 |
| | External dimensions (W×D×H mm) | | 710×645×1008 | 710×645×1600 |
| | Internal capacity | | 159L | 300L |
| | Shelf load | | 15kg/layer | |
| | Shelf layers/shelf support spacing | | 12 layers/30mm | 24 layers/30mm |
| | Power supply (50/60Hz) rated current | | AC220V 6.5A | AC220V 7A |
| | Weight | | Approx. 90kg | Approx. 150kg |
| Accessories | Shelf | | Stainless punching mesh plate | |
| | | | 3 pcs | 5 pcs |
| | Supports | | 6 pcs | 10 pcs |
| Door key | | 2 pcs | | |
| Options | Stand | | ON61C | — |
| | Others | | Shelf plate (1 shelf plate with 2 shelf supports), cable port (30/50mm), data logger, combination warning lights (standby/running/fault), viewing window, external communication function (RS485), temperature output terminals (4~20mA), external alarm output terminals, time-up output terminals | |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Low-temp. Incubators | Semiconductor Refrigeration

IJ102/102W/300

Operating temp. range 5~60°C

Temp. distribution accuracy $\pm 1.0^{\circ}\text{C}$ (at 37°C)

Internal capacity 15.6L 43L

Low vibration, CFC-free, low-temperature incubator equipped with semiconductor cooling chips.

Features

- Ideal for constant temperature experiments and sample preservation in a compact design.
- Use semiconductor cooling chips, achieving CFC-free, low vibration.
- Equipped with a programmable temperature controller.
- Dedicated air jacket can be added (optional).
- Double-layered (optional: stacking stand).
- Equipped with rich safety functions such as self-diagnosis function, automatic overheat prevention function, overheat prevention device, overcurrent leakage protection switch, etc.
- The viewing window of IJ102W uses energy-saving double-layer semi-tempered glass, effectively preventing unnecessary heat transfer.



Specifications

| Model | | IJ102 | IJ102W | IJ300 |
|----------------|--------------------------------------|--|--------------|------------------|
| System | | Forced convection | | |
| Performance | Operating temp. range | 5~60°C (condition: room temperature 25°C) | | |
| | Temp. adjusting accuracy | $\pm 0.3^{\circ}\text{C}$ (at 37°C) | | |
| | Temp. distribution accuracy | $\pm 1.0^{\circ}\text{C}$ (at 37°C) | | |
| | Max. temp. reaching time | Approx. 60 min (20~60°C) | | |
| | Min. temp. reaching time | Approx. 100 min (20~5°C) | | |
| Composition | Interior material | Stainless steel plate | | |
| | Exterior material | Electro-galvanized steel plate with chemical resistance coating on the surface | | |
| | Insulating material | Foamed polyurethane | | |
| | Heater | Mica heater 120W | | Tube heater 300W |
| | Cooler | Semiconductor cooling chip/forced convection | | |
| | Viewing window | — | 180×180mm | — |
| | Cable port | — | Φ30mm | — |
| Controllers | Temp. control method | PID control | | |
| | Temp. setting method | Digital settings achieved through function menu key and ▼▲ keys | | |
| | Operation functions | Constant operation, rapid automatic stop operation, automatic stop operation, automatic start operation, program operation (30 segments×1/15 segments×2/10 segments×3) | | |
| | Additional functions | Deviation correction, key lock, power outage compensation | | |
| | Sensors | Pt resistance thermometer (temperature controller), K-type thermocouple (overheat prevention) | | |
| Safety device | | Self-diagnosis function (temperature sensing abnormality, heater disconnection, SSR short circuit, automatic overheat prevention function), overheat prevention device, main relay abnormality, overcurrent leakage protection switch | | |
| Specifications | Internal dimensions (W×D×H mm) | 250×250×250 | | 350×350×350 |
| | External dimensions (W×D×H mm) | 350×399×565 | | 470×496×665 |
| | Internal capacity | Approx. 15.6L | | Approx. 43L |
| | Shelf load | Approx. 15kg/shelf | | |
| | Shelf layers | 7 layers | | 10 layers |
| | Power supply (50/60Hz) rated current | AC100V 4A | | AC100V 9A |
| | Weight | Approx. 20kg | Approx. 22kg | Approx. 37kg |
| Accessories | Shelf | Stainless steel wire mesh plate | | |
| | | 2 pcs | | |
| | Supports | 4 pcs | | |
| | Others | Drain pan 1 | | |
| Options | | Shelf plate (1 shelf plate with 2 shelf supports), air jacket, stacking stand, shaker setting platform, inner door, air jacket inner door, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time arrival output terminal | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionKey Lock
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
BreakerMade in
Japan

Control panel



Drain pan



Internal chamber



Air jacket (optional)
Internal dimensions: W200×D200×H200mm/8L

Stacking stand installation example



| | |
|------------------------------|----------|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

2-Chamber Incubators | Air Jacket Type

INC822C

Natural
ConvectionForced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range Upper section: RT+10~80°C Lower section: -10~60°C

Temp. uniformity $\leq 2.0^{\circ}\text{C}$ (at 37°C)

Internal capacity Upper section: 150L Lower section: 143L

1 product can simultaneously perform high temperature incubation and low temperature incubation.



■ Features

- The upper chamber is a constant temperature chamber with a constant value operation, timer function and overheat prevention device, and the lower chamber is a low temperature incubator with 6 programs setting of 10 segments to 30 segments.
- The lower chamber achieves high-precision temperature control and balanced internal temperature distribution through forced air circulation by fan blades.
- The inner glass door enhances the insulation effect. Observation of samples through the inner glass door minimizes temperature changes.
- Defrosting can be set to cycle defrost based on the continuous operation time of the compressor.
- After power outage, the operation can automatically restart under the set temperature conditions through the power outage compensation function.
- Both upper and lower chambers are equipped with door locks.

■ Safety

- Equipped with safety functions such as refrigerator overload relay, independent overheat prevention device, overcurrent leakage protection switch, self-diagnosis circuit (temperature sensor abnormality, heater disconnection, SSR short circuit, automatic overheat prevention), buzzer alarm in case of abnormality, etc.

■ Specifications

| Model | | INC822C | |
|-----------------|--------------------------------------|--|--|
| Basic structure | | Lower chamber: low temperature incubator | Upper chamber: constant temperature incubator |
| Performance | Operating temp. range | -10~60°C | Room temp. +5~80°C |
| | GB standard | Temp. fluctuation | ±0.1°C (at 37°C) |
| | | Temp. uniformity | ±1.5°C (at 37°C) |
| | Heating time | 20~52°C ≤ 20min | — |
| | Cooling time | 20~-4°C ≤ 60min | — |
| System | | Forced convection | Natural convection |
| Composition | Interior material | Stainless steel (SUS304) | |
| | Heater, insulation material | Nickel-chromium alloy heating wire, foamed polystyrene | Nickel-chromium alloy heating wire, glass wool |
| | Heater power | 550W | 400W |
| | Air blower | Crossflow fan | — |
| | Refrigerator, refrigerant | 158W R134A | — |
| | Defrost mechanism | Cycle operation | — |
| | Cable port | Internal diameter 30 mm right side of the body | Internal diameter 30 mm right side of the body |
| | Exhaust vent | — | 30 mm roof of the body |
| | Controllers | VS6 type program operation temperature controller | VS6 type constant value operation temperature controller |
| | Sensors | Temperature controller: Pt resistance thermometer | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple |
| Specifications | Heater control | SSR control | |
| | Additional functions | Deviation correction, key lock, power outage compensation | |
| | Operation functions | Program operation (30 segments×1, 15 segments×2, 10 segments×3), automatic stop operation, automatic start operation | Automatic stop operation, automatic start operation |
| | Safety device | Refrigerator overload relay, independent overheat prevention device, overcurrent leakage protection switch, self-diagnosis circuit (temperature sensor abnormality, heater disconnection, SSR short circuit, automatic overheat prevention), buzzer alarm in case of abnormality | |
| | Power supply (50/60Hz) rated current | AC220V 6A | |
| | External dimensions (W×D×H mm) | 710×656×1792 | |
| | Internal dimensions (W×D×H mm) | 600×477×500 Internal volume: 143L | 600×530×500 Internal volume: 150L |
| Accessories | Shelf layers/shelf support spacing | 13 layers/30mm | |
| | Shelf load | Approx. 15 kg/piece | |
| | Product weight | Approx. 160 kg | |
| | Shelf | Stainless punching mesh plate | |
| | Supports | 3 pcs | 2 pcs |
| Options | Door key | 6 pcs | 4 pcs |
| | Options | Shelf (1 shelf with 2 supports), micro printer, data logger, combination warning light (standby/running/fault), external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time-up output terminal, centralized monitoring software, and touch screen controller | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

2-Chamber Incubators | Low Temp.

IQ823C

Forced
ConvectionAutomatic
Overheating
PreventerOverheating
PreventerSelf-diagnostic
FunctionPower Failure
Compensation
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range -10~60°C

Temp. uniformity $\leq 2.0^{\circ}\text{C}$ (at 37°C)

Internal capacity 143L×2

A two-chamber low-temperature incubator capable of independent program operation.

Features

- The controller panel is set on the door, improving operability.
- The process of incubation and preservation, etc., can log 6 programs from 10 segments to 30 segments according to experimental conditions.
- Defrosting can be set to cycle defrost based on the continuous operation time of the compressor. Both chambers are equipped with large cable ports with an internal diameter of 50mm.
- Equipped with safety functions such as refrigerator overload relay, independent overheat prevention device, overcurrent leakage protection switch, self-diagnosis circuit (temperature sensor abnormality, heater disconnection, SSR short circuit, automatic overheat prevention), buzzer alarm in case of abnormality, etc.
- Equipped with door lock.

Internal chamber




Cable port



Shelf



Heating and cooling curve



Specifications

| Model | | IQ823C |
|----------------|--|--|
| System | | Forced convection |
| Performance | Operating temp. range | -10~60°C (both upper chamber and lower chamber) |
| | GB standard | Temp. fluctuation $\pm 0.1^{\circ}\text{C}$ (at 37°C) |
| | | Temp. uniformity $\pm 1.0^{\circ}\text{C}$ (at 37°C) |
| | Temp. adjusting accuracy | 20~52°C $\leq 20\text{min}$ |
| | Temp. distribution accuracy | 20~4°C $\leq 60\text{min}$ |
| Composition | Interior material | Stainless steel (SUS304) |
| | Exterior material | Cold rolled steel plate with chemical proofing coating |
| | Insulating material | Foamed polystyrene |
| | Refrigerator | 158W×2 |
| | Refrigerant | R134A |
| | Defrost mechanism | Cycle operation |
| | Heater | Nickel-chromium alloy heating wire 550W×2 |
| Controllers | Cable port | Inner diameter: 50mm (right side) |
| | Temp. control method | PID control |
| | Temp. setting method | Digital settings achieved via function menu key and $\blacktriangle/\blacktriangledown$ keys |
| | Temp. display method | Achieved temp. display: Green 4-digit LED digital display |
| | | Setting temp. display: Red 4-digit LED digital display |
| | Timer | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) |
| | Operation functions | Automatic stop operation, automatic start operation, program operation (30 segments×1, 15 segments×2, 10 segments×3) |
| | Additional functions | Deviation correction, key lock, power outage compensation |
| | Safety device | |
| | Refrigerator overload relay, independent overheat prevention device, overcurrent leakage protection switch, self-diagnosis circuit (temperature sensor abnormality, heater disconnection, SSR short circuit, automatic overheat prevention), buzzer alarm in case of abnormality, etc. | |
| Specifications | Internal dimensions (W×D×H mm) | 600×477×500 each |
| | External dimensions (W×D×H mm) | 710×656×1792 |
| | Internal capacity | 143L×2 |
| | Shelf load | Approx. 15kg/shelf |
| | Shelf layers/shelf support spacing | 13 layers×2 chambers/30mm |
| | Power supply (50/60Hz) rated current | AC220V 7A |
| Accessories | Weight | Approx. 165kg |
| | Shelf/supports | Stainless punching mesh plate: 2 chambers×3 pieces/2 |
| Options | Key | 2 chambers×2 pieces |
| | Shelves, micro printer, data logger, external communication function (RS485), temperature output terminal (4~20mA), external alarm output terminal, time arrival output terminal, centralized monitoring software, touch screen controller | |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

CO₂ Incubators | Air Jacket Type, Dry Heat Sterilization

IPE610

Overheating
PreventerSelf-diagnostic
FunctionOvercurrent
Leakage Circuit
Breaker

Operating temp. range RT+5~60℃

CO₂ concentration adjustment range 0~20%

Temp. uniformity ±0.3℃


Internal capacity Approx.185L

Considered the least contamination risk, air jacket CO₂ incubator with HEPA filter and dry heat sterilization function.



Features

- Dry heat sterilization at 160℃, 2 h according to pharmaceutical standards. Equipped with a program that can perform dry heat sterilization at 160℃ for 2 h. HEPA filter continuously filters the air inside the chamber, reducing the risk of contamination to samples. Inner door is made of reinforced glass, allowing observation of the interior.
- High-precision, long-life CO₂ sensor can perform dry heat sterilization at 160℃ without removal. CO₂ sensor contamination is no longer an issue.



Specifications

| Model | | IPE610 |
|----------------|--|---|
| System | | Air jacket type with dry heat sterilization function |
| Performance | Temp. control range | Room temp. +5~60℃ |
| | Temp. fluctuation | ±0.1℃ (at 37℃) |
| | Temp. uniformity | ±0.3℃ (at 37℃) |
| | Operating humidity range | Above 90%RH |
| | CO ₂ concentration adjustment range | 0~20.0% |
| | CO ₂ concentration control accuracy | ±0.2% (at 37℃, CO ₂ concentration: 5%) |
| Composition | Touch screen | 7-inch color touch screen |
| | Exterior/door | Cold rolled steel plate with chemical proofing coating |
| | Internal chamber/inner door | Stainless steel plate/tempered glass |
| | Door opening method | Standard right opening door (can be changed to left opening door) |
| | Insulating material | Glass fiber |
| | HEPA filter | Dust collection efficiency: 99.97% for 0.3μm particles |
| | Chamber temp. /water bath water temp. sensor | Pt100 thermistor |
| | CO ₂ sensors | Infrared method (new single light, dual wavelength method, no drive part) |
| | Independent overheat prevention sensor | K-type thermocouple |
| | Heater | 900W |
| | Cable port | Φ32mm back side |
| | Others | Adjustable feet, CO ₂ sample port (OUT), temperature output terminal, CO ₂ concentration output terminal, external alarm output terminal, various records display, trend graph, USB data export port |
| Safety device | | Overcurrent leakage protection switch, independent overheat prevention device, temperature upper and lower limit abnormalities, CO ₂ concentration upper and lower limit abnormalities, self-diagnosis function, door open alarm, key lock |
| Specifications | External dimensions | W660×D650×H1020mm |
| | Internal dimensions | W540×D520×H665mm |
| | Internal capacity | Approx. 185L |
| | Shelf dimension/layers/spacing | W471×D540/15 layers/30mm |
| | Shelf material/load capacity | Stainless steel shelves/about 7kg/layer |
| | Power supply (50/60Hz) rated current | AC220V 4.2A |
| Accessories | | Shelf support installation plate 2 pieces, shelf support hanging bars 4 pieces, shelves 4 pieces, shelf supports 8 pieces, humidifying pan 1 piece, silicone stopper for cable port 2 pieces, CO ₂ supply hose 2m, hose clamp 2 pieces, humidifying pan cover 1 piece, sampling tube |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

CO₂ Incubators | Water Jacket Type

BNA610

Overheating
PreventerSelf-diagnostic
FunctionOvercurrent
Leakage Circuit
BreakerMade in
Japan

Operating temp. range RT+5~50°C

CO₂ concentration adjustment range 0~20%


Temp. distribution accuracy ±0.2°C

Internal capacity 167L

Provide a good culture environment and operability in water jacket CO₂ incubator.

Features

- Anti-pollution
Seamless welded design inside the chamber, all corners are arc-shaped, and the four corners on the back are spherical corners. Additionally, the shelf support is integrated with the wall, serving as a dust-proof measure and easy to clean.
- Stable culture operation
Using the water jacket method to make the internal chamber wall temperature uniform, and micro air circulation method of up and down airflow to minimize wall frost.
- Stability and reliability of CO₂ control
Significantly improved durability and stability in continuous operation by using a new single light, dual wavelength NDIR sensor.



[Cross-sectional diagram]



Specifications

| Model | | | BNA610 |
|----------------|--|-----------------------------|---|
| System | | | Water jacket type humidification method: natural evaporation through the humidifying pan; stirring: micro air circulation through the fan blade |
| Performance | Temp. control range | | Room temperature +5 ~ 50℃ (requires room temperature below 32℃ when operating at 37℃) |
| | JIS stand-ard | Temp. adjusting accuracy | ±0.1℃ (at 37℃) |
| | | Temp. distribution accuracy | ±0.2℃ (at 37℃) |
| | Operating humidity range | | Above 95%RH |
| | CO ₂ concentration adjustment range | | 0 ~ 20.0% |
| | CO ₂ concentration control accuracy | | ±0.1% (at 37℃, CO ₂ concentration: 5%) |
| Composition | Touch screen | | 7-inch color touch screen |
| | Exterior/door | | Chemical-resistant coating on the surface of electro-galvanized steel plate |
| | Internal chamber/inner door | | Stainless steel plate/tempered glass |
| | Door opening method | | Standard right opening door (can be changed to left opening door) |
| | Insulating material | | Glass fiber |
| | Chamber temp. /water temp. sensor | | Pt100 thermistor |
| | CO ₂ sensors | | Infrared method (new single light, dual wavelength method, no drive part) |
| | Independent overheat prevention sensor | | K-type thermocouple |
| | Water bath heater | | 130W×2 |
| | Door/door frame heater | | 98W/25W (output amount arbitrarily set) |
| | Cable port | | Φ32mm back side |
| | Others | | Adjustable feet, CO ₂ sample port (OUT), temperature output terminal, CO ₂ concentration output terminal, external alarm output terminal, various records display, trend graph, USB data export port, drain valve (Φ10.5) |
| Safety device | | | Overcurrent leakage protection switch, independent overheat prevention device, temperature upper and lower limit abnormalities, CO ₂ concentration upper and lower limit abnormalities, self-diagnosis function, water bath water level abnormality (electrode type), door open alarm, humidifying water shortage alarm (thermistor), key lock |
| Specifications | External dimensions | | W600×D664 (765)×H880mm |
| | Internal dimensions | | W485×D540×H640mm |
| | Internal capacity | | Approx. 167L |
| | Shelf dimension/layers/spacing | | W471×D451/12 layers/40mm |
| | Shelf material/load capacity | | Shelves (aluminum)/about 5kg/layer |
| | Power supply (50/60Hz) rated current | | AC220V 2A |
| Accessories | | | Shelves 4 pieces, humidifying pan (stainless steel), silicone stopper for cable port 2 pieces, CO ₂ supply hose 2m, hose clamp, water supply and drain hose 2m |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

IPE610 & BNA610 Outline

1. High reliability CO₂ sensor
Long-life sensor with new infrared light source.

2. CO₂ consumption IPE610
Improved airtightness significantly reduces CO₂ consumption.
Stable: about 0.34L/h
1 h after door closes: about 8L/h


3. Internal chamber
Smooth operation of the shelves.
Considering daily sanitation, all corners are R rounded, and the four corners on the back are spherical corners.



4. Power cord hole (Φ32mm)
Convenient for running measurement sensor wiring.
5. Independent overheat prevention device IPE610
Provides dual protection against overheating inside the chamber.
6. Shelves (stainless steel) IPE610
Load capacity about 7kg.
7. Shelves (aluminum) BNA610
Using lightweight, good thermal conductivity aluminum material, the temperature distribution on the shelf is excellent, effectively maintaining the sample temperature.



- Resistive touch screen (operable with gloves)
- Temperature/CO₂ concentration value/setting value
- Temperature/CO₂ concentration chart display
- Calendar function (date, time)
- Accumulated time function
- Timer alarm function (operation continues)
- Various data and record management
- Various data can be exported via USB



Simple screen (can be turned off)

- The advantage of a large screen is not just the large amount of information displayed.
Display mode that is easy to view even from a distance. Three languages: Chinese/English/Japanese.

BNA test data

BNA antibacterial shelf, antibacterial humidifying pan, and clean circulation structure with high-efficiency filter are as options.

Antibacterial shelf antibacterial capability test Result: Antibacterial when it is above 2.0

| Test bacteria | Test piece | Antibacterial activity value |
|-----------------------|---------------------------------------|------------------------------|
| Escherichia coli | SUS304 | 0.9 |
| | Antibacterial shelf (stainless steel) | 3.7 |
| | Antibacterial shelf (aluminum) | 6.3 |
| Staphylococcus aureus | SUS304 | 0.9 |
| | Antibacterial shelf (stainless steel) | 4.9 |
| | Antibacterial shelf (aluminum) | 5.4 |

JIS Z 2801: 2000 "Antibacterial Processed Products - Antibacterial Test Method - Antibacterial Effect"

Test result: measurement result of the number of live bacteria on the test piece (at 35°C, after 24 h)

Unprocessed test piece: Polyethylene sheet

Test piece: SUS304 antibacterial treatment, antibacterial shelf (stainless steel), antibacterial shelf (aluminum)

Escherichia coli: Escherichia coli NBRC 3972

Staphylococcus aureus: Staphylococcus aureus aubsp aureus NBRC 12732

Antibacterial humidifying pan water antibacterial test Result: Antibacterial

After sterilizing and drying each pan, pure water is added, and placed in a 37°C chamber for 2 h, measuring the number of live bacteria in the humidifying water. (Opening the door 5 times a day×30 seconds/time, measuring the number of live bacteria in 1ml of humidifying water diluted 1000 times)

| Specimen | Measurement | Number of live bacteria | |
|---|--|-------------------------|---------------|
| | | When adding water | After 2 weeks |
| Humidifying pan (stainless steel) | Chamber temperature: 37°C Chamber humidity: 95±5%RH | 0 pcs | 21 pcs |
| Antibacterial humidifying pan (antibacterial stainless steel) | | 0 pcs | 1 pcs |



When adding humidifying water



Antibacterial humidifying pan (stainless steel)




Humidifying pan (stainless steel)

Antibacterial humidifying pan Escherichia coli proliferation test Result: Antibacterial

| Test bacteria | Specimen | Measurement | Number of live bacteria |
|------------------|-------------------------------|---------------|-------------------------|
| Escherichia coli | Stainless steel | 6,000,000 pcs | 12,000 pcs |
| | Antibacterial stainless steel | | Below 10 pcs |
| | Plastic plate | | 10,000,000 pcs |

After pre-cleaning and sterilizing the test piece, soaking it in the Escherichia coli test bacteria solution for 4 h, wiping and measuring the number of colonies according to the usual method.

Clean circulation structure cleanliness test (HEPA filter)



Result: FED specified class 100


Measurement condition: After opening the door for 30 seconds, activating the HEPA high-efficiency filter clean circulation structure, eliminating dust and miscellaneous bacteria from the source of bacteria. It only takes about 40 min to reach the FED specified 100 level.

Antibacterial shelf mold resistance test

| Test piece | Mold resistance | | | |
|---------------------------------------|-----------------|---------|---------|---------|
| | 1 week | 2 weeks | 3 weeks | 4 weeks |
| SUS304 | 2 | 2 | 2 | 2 |
| Antibacterial shelf (stainless steel) | 0 | 0 | 0 | 0 |

JISZ2911: 2000 "Mold Resistance Test"
Test bacteria strain
Aspergillus niger IFO 6341
penicillium citrinum IFO 6352
Rhizopus oryzae IFO 31005
Cladosporium cladosporioides IFO 6348
Chaetomium globosum IFO 6347

Recovery performance after door opening for 30 seconds




Selectable door opening direction (left/right) IPE610

Right door opening
(factory setting)

Left door opening

Independent overheat prevention device and terminal location IPE610



Independent overheat prevention device and terminal location BNA610

Terminal block:
Upper section: temperature output, CO₂ output, external alarm output
Lower section (optional): Communication terminal

Standard accessories and optional items selection table IPE610

| Option name | Option model No. | Option name | Option model No. |
|-----------------------------------|------------------|------------------------------------|------------------|
| Standard shelves | IPE610_01_02-10 | Rust inhibitor for humidifying pan | OBN34C |
| Standard humidifying pan | IPE610_01_02-29 | External communication adapter | OBN24C |
| Stacking fittings | OBNE62 | Pressure regulator set | OBN30-1C |
| Heat insulating stacking fittings | OIPE60 | Pressure regulator set | OBN30-2C |
| Low bed stand with casters | OBNE60 | CO ₂ cylinder switcher | OIT28C |
| Stand | OBNE61 | HEPA filter | IPE610_01_02-09 |

Option BNA610

| No. | Product name | Specifications | Product number |
|-----|---|--|----------------|
| ① | Standard shelves | Standard aluminum shelves for BNA, also usable for IPE type. W471×D457mm | 213747 |
| ② | Antibacterial shelf (aluminum) | Antibacterial treated aluminum shelves to inhibit microbial growth. W471×D457mm | 213748 |
| ③ | Antibacterial shelf (stainless steel) | Antibacterial treated stainless steel shelves to inhibit microbial growth, also usable for BNE type. W471×D457mm | 211253 |
| ④ | Standard humidifying pan | Standard humidifying pan for BNE (stainless steel) Approx. 4L, W350×D265×H60mm | 213749 |
| ⑤ | Antibacterial humidifying pan | Stainless steel humidifying pan for IP type to inhibit microbial growth, also usable for BNE type. Approx. 4L, W350×D265×H60mm | 213750 |
| ⑥ | Stacking fittings | Used when stacking two BNE units, placing one on the lower layer, and connecting the upper and lower layers. | 213751 |
| ⑦ | Heat insulating stacking fittings | Special connecting fittings for IP type, allowing simultaneous dry heat sterilization on the lower layer and culture operation on the upper layer. | 211254 |
| ⑧ | Low bed stand with casters | Casters for easy product mobility, usable for two stacked units. W600×D601×H715mm | 213752 |
| ⑨ | Stand | Fixed stand for one product. W600×D601×H715mm | 213753 |
| ⑩ | Rust inhibitor for water jacket | Rust inhibitor to prevent water jacket from rusting. Capacity 50ml (10ml/time ×5 times) | 213758 |
| ⑪ | Antiseptic for humidifying pan | Added to humidifying water to inhibit microbial growth. 25g/unit (approximately 9g/time for 3L humidifying water) | 213759 |
| ⑫ | CO ₂ sample port (IN) | Used to return gas extracted from the sample port (OUT) back into the chamber. | 213760 |
| ⑬ | External communication (RS485) terminal | Allows monitoring running status and remote operation. | 213756 |
| ⑭ | External communication adapter | RS485 to USB converter. | 213754 |
| ⑮ | CO ₂ relief valve | Regulate CO ₂ supply pressure | 213757 |
| ⑯ | CO ₂ cylinder switcher | Automatically switch to the reserve cylinder when the CO ₂ cylinder is empty. | 213753 |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Plasma Equipment

Contents


| | |
|-------------------------------|-----|
| Gas Plasma Reactors | 130 |
| Gas Plasma Dry Cleaners | 132 |

From research and development to industrial production, from surface micro-machining to surface treatment modification, plasma cleaners and reactors have very extensive social application needs.

Plasma cleaners and reactors are increasingly used in dry cleaning applications in semiconductors and electronic materials, such as photoresist stripping on silicon wafers, organic film removal, surface activation, micro polishing, and carbon film removal. YAMATO Scientific's plasma products play an active role in these fields.



YAMATO Scientific's plasma surface treatment equipment is not only used for cleaning applications related to semiconductors but also has broad applications in various other fields.



Plasma treatment has significant effects on enhancing the strength of various bonding, coating, and plating applications.

DP Mode Cylindrical Type PR/PM Series

Use in photoresist stripping from wafers, organic film removal using oxygen or argon gas, surface activation, micro polishing, or carbon film removal in a wide range of fields.




RIE DP mode parallel plate type PDC/V series


RIE/DP 2 modes of plasma processing, it can be used for etching and dry cleaning of silicon wafers, and COB surface activation of sensors and others, improving the stability of wire bonding, and etching of oxides and hydroxides, removal bonding barriers from electrode surfaces.



Cylindrical Structure (DP)



Control panel




Plasma super cleaner

IC packaging technology

- BGA (/Ball Grid Array)
- CSP/(Chip Size Package)
- Imaging equipment assembly


Plasma
Cleaning

- Die Attach
- Wire Bonding
- Molding



Plasma device etching and ashing

Etching




Ashing




Plasma device etching and ashing modes

IC packaging has evolved from ceramic to molded resin, with wire material transitioning from Ni, Al to Cu. Plasma etching and cleaning technology can effectively meet the demands of these product advancements.

Ar plasma for metal oxide removal (PDC/V Series)




Substrate cleaning (V series)



Plasma etching and ashing modes


Organic film removal

- Molding release agents
- Solder residues
- Oil and grease (thin films)
- Organic films (buffing)




Surface modification

- Polyester
- Polypropylene
- Teflon
- Ceramic
- Polyimide




Surface cleaning

- Material surface




Surface etching

- Glass epoxy
- Resin polyimide



Metal oxide removal


- Cu_2O, CuO
- Al_2O_3



Dry etching

- Barrel
- DP


- Parallel plate
- Sputtering



Measurement example

Measurement position

silicon wafer 8 inch



Parallel plate electrode experiment 1

Test material: silicon wafer 8 inch
Photoresist: positive

| Electrode | Parallel plate electrode working surface $\Phi 350$ electrode distance 60mm | | | | |
|-----------------------------|---|--|--|--|--|
| RF power | 500W | | | | |
| Processing time | 5min | | | | |
| Gas flow rate | 500sccm | | | | |
| Pressure | 100Pa | | | | |
| Working surface temperature | 180°C | | | | |

| Measurement position | ① | ② | ③ | ④ | ⑤ | Average value |
|----------------------|-------|-------|-------|-------|-------|---------------|
| ini (Å) | 40290 | 40540 | 40930 | 40620 | 40670 | 40610 |
| after (Å) | 18350 | 19980 | 15930 | 16150 | 16690 | 17420 |
| ini-after (Å) | 21940 | 20560 | 25000 | 24470 | 23980 | 23190 |
| reta (Å/min) | 4388 | 4112 | 5000 | 4894 | 4796 | 4638 |

| Maximum value | Minimum value | Uniformity |
|---------------|---------------|------------|
| 5000 | 4112 | 9.573092 |

Parallel plate electrode experiment 2

Test material: silicon wafer 8 inch
Photoresist: positive

| Electrode | Parallel plate electrode working surface $\Phi 350$ electrode distance 60mm | | | | |
|-----------------------------|---|--|--|--|--|
| RF power | 500W | | | | |
| Processing time | 5min | | | | |
| Gas flow rate | 500sccm | | | | |
| Pressure | 100Pa | | | | |
| Working surface temperature | 160°C | | | | |

| Measurement position | ① | ② | ③ | ④ | ⑤ | Average value |
|----------------------|-------|-------|-------|-------|-------|---------------|
| ini (Å) | 40450 | 40480 | 40820 | 40780 | 40560 | 40618 |
| after (Å) | 25290 | 24780 | 24250 | 24640 | 23890 | 24570 |
| ini-after (Å) | 15160 | 15700 | 16570 | 16140 | 16670 | 16048 |
| reta (Å/min) | 3032 | 3140 | 3314 | 3228 | 3334 | 3209.6 |

| Maximum value | Minimum value | Uniformity |
|---------------|---------------|------------|
| 3334 | 3032 | 4.704636 |

- Sterilizers 1
- Granulation and Spray Dryers 2
- Muffle Furnaces 3
- Ovens 4
- Incubators 5
- Plasma Equipment 6
- Water Purifiers 7
- Baths 8
- Water Circulators 9
- Rotary Evaporators 10
- Freeze Dryers & Cold Traps 11
- Stirrers & Shakers 12
- Washers 13
- Analysis and Test Devices 14
- Options 15

Gas Plasma Reactors | DP Mode, Barrel

PR500/510/1000

Made in
Japan

High frequency output 50~500W/1000W

Reaction chamber $\Phi 215 \times 305 \text{mm} / \Phi 400 \times 500 \text{mm}$

Low-temperature plasma, classic barrel.



PR500



PR510

PR500: manual control

PR510/1000: automatic touch screen control

The classic model of low-temperature plasma is very historic. With automatic coordination function and touch screen operation (PR510/1000), it improves operability and safety.

Usage

- Photoresist removal
- Component cleaning
- Interface activation
- Micro polishing
- Compatible with silicon wafers and glass substrates

Features

- Large diameter ($\Phi 215 \text{mm}$) quartz reaction chamber
- Compact integrated design
- Simple touch screen operation (PR510/1000)

Internal chamber



Silicon wafer ashing state



Specifications

| Model | | PR500 | PR510 | PR1000 |
|--------------------|-----------------------------|--|--------------------------------|-------------------|
| Mode | | DP mode, barrel | | |
| Control section | High frequency output power | Max. 500W | | Max. 1000W |
| | Frequency | 13.56MHz | | |
| | Tuning method | Automatic coordination | | |
| Reaction | Reaction chamber | Quartz Φ215×305mm | | Φ400×500mm |
| | Reaction gas | 2System (O ₂ /CF ₄) | | |
| | Control | Manual | Touch screen automatic control | |
| | Piping material | Stainless steel, teflon | | |
| External dimension | | W438×D520×H760mm | W520×D630×H760mm | W800×D900×H1706mm |
| Weight | | Approx. 60kg | Approx. 70kg | |
| Power supply | | Single phase AC100V 9A | 3-phase AC200V 8A | |

Gas Plasma Reactors | Economic, High-frequency Plasm

PR210C/300/301

Made in Japan

| | | | | | | |
|-----------------------|-------------|-----------------|------------------|---------------------|------------------------------|--------------------|
| High frequency output | PR210C:200W | PR300/301: 300W | Reaction chamber | PR210C: Φ100×L160mm | PR300: Φ64×L160mm×3 chambers | PR301: Φ118×L160mm |
|-----------------------|-------------|-----------------|------------------|---------------------|------------------------------|--------------------|

From sample surface modification to cleaning, a small desktop plasma cleaner.

Usage

- Imparting functional groups to the surface of polymer materials
→Improving adhesion and adhesion. Through oxidation reaction, functional groups such as -OH, >C=O, -COOH are generated on the surface (affected by trace amounts of water and carbon dioxide).
- Similarly, in nitrogen plasma, nitrogen atoms are absorbed on the surface, generating functional groups such as -NH₂.
- Silicon wafer etching.
- Material surface modification (metal, polymer, film, ceramic, etc.).
- Asbestos pre-treatment (film filtration ashing).
- From electronic equipment to the biochemical market.
- Low-temperature ashing (polymer materials, coal, food, etc.).
- Using glass and PDMS plates to bond PDMS chips.
- Surface organic matter removal.

Features

- High-performance output type with 200w rf (PR210C).
- Due to its compact desktop size, it can be installed in very small spaces.
- Synchronized, automatic tuning, easy operation (PR210C).
- Plasma reactors (ashing equipment) (PR210C).



Reaction chamber




Single chamber (PR210C)



Test material rack

Example: asbestos analysis pre-treatment



| Model | | PR210C | PR300 | PR301 |
|---------------------|-----------------------|---|--|---------------------------------|
| Mode | | DP mode, barrel | | |
| Control | High frequency output | PR210C: Max. 200W, PR300/301: Max. 300W | | |
| | Frequency | 13.56MHz | | |
| | Tuning method | Manual 2-axis tuning method | | |
| Reaction | Reaction chamber | Heat-resistant glass Φ100×160mm | Ultra-hard glass Φ64×160mm×3 chambers | Heat-resistant glass Φ118×160mm |
| | Reaction gas | 1 system (oxygen) flow meter | | |
| | Control | Manual shut-off valve | Automatic pressure reduction, automatic shut-off valve | |
| | Piping material | Stainless steel, teflon, copper, copper alloy | | Stainless steel, teflon |
| External dimensions | | W350×D400×H500mm | W438×D520×H556mm | W438×D520×H630mm |
| Weight | | Approx. 10kg | Approx. 36kg | Approx. 34kg |
| Power supply | | AC220V 10A | AC100V 7A | |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Gas Plasma Cleaners | Economic, Low-frequency

PM110C

Low-frequency output Low-frequency high voltage power supply

Reaction chamber Chamber dimensions $\Phi 100 \times L160\text{mm}$

Economic, surface modification and cleaning of target objects by the active action of ionized gas.



Features

- Economic plasma cleaner.
- Simple plasma treatment operation.

Control panel



Viewing window



Promising applications

applications <cleaning>

- Removal of organic films (oil, wax, etc.).
- Cleaning of electronic parts, optical parts, and small components, pre-treatment before electroplating, coating.
- Bacteria disinfection (cell membrane etching by oxygen radicals or disinfection by UV).

<Surface modification>

- Imparting functional groups to the surface of polymer materials → Improving adhesion and adhesion. Through oxidation reaction, functional groups such as -OH, >C=O, -COOH are generated on the surface (affected by trace amounts of water and carbon dioxide).
- Improving the adhesion of various polymers and the biocompatibility of medical and dental materials.
- By treating the culture dish, it can be evenly coated with medium, etc.

<Etching>

- Etching of organic polymer materials. For organic polymers, oxygen can become the etching gas, allowing for resist etching.
- Promoting seed germination experiments.

Specifications

| Model | PM110C |
|--------------------------------|---|
| Plasma generation source | Low-frequency high voltage power supply, output approx. 50w, frequency 60KHz, output voltage 10kV |
| Electrode structure | Capacitive type, 2 segments |
| Gas flow meter | For oxygen, flow rate 30~300ml/min |
| Vacuum adjustment valve | For vacuuming, returning to atmosphere |
| Reaction chamber dimensions | Chamber dimensions $\Phi 100 \times L160\text{mm}$ |
| External dimensions (W×D×H mm) | 310×300×448 |
| Power supply 50/60Hz | AC220V 10A |
| Vacuum pump interface | KF25 (recommended vacuum pump exhaust speed 30L/min) |
| Gas interface | Outer diameter $\Phi 6.35\text{mm}$ |

Plasma experiment examples



No plasma treatment plasma treatment 10 min

Radish seed germination promotion experiment results
Germination status 1 day after plasma treatment (treatment gas: atmosphere)

Gas Plasma Cleaners | RIE-DP Mode, Parallel Plate Electrode

PDC200/210/510

Made in
Japan

High frequency output 300W (PDC200model)/500W (PDC210/510model)

Electrode dimensions 250×170mm (PDC200/210model)/410×170mm (PDC510model)

Small plasma cleaners suitable for research and development.

Using plasma surface treatment equipment.
Focusing on “easy to use” high-performance
desktop models, easy to retrieve data, centered on
electronic materials, suitable for various purposes.

Usage

- CSP, BGA, COB substrate plasma cleaning
- Removal of organic films and metal oxide films
- Dry cleaning of printed circuit boards
- Interface activation treatment
- LED packaging

Features


- DP mode (standard on PDC510, optional on PDC200/210)
- RIE mode
- High uniformity electrode structure
- Simple touch screen controller operation



Internal chamber



Configuration diagram



Specifications

| Model | PDC200 | PDC210 | PDC510 |
|-----------------------------|---------------------------------|--------------------------------------|--------------------------------------|
| Plasma mode | RIE (DP mode optional) | | RIE/DP |
| Electrode structure | Parallel plate (fixed) | | |
| High frequency output | Max300W | Max500W | |
| Frequency | 13.56MHz | | |
| Control/display | LCD touch screen | | |
| Reaction chamber dimensions | W400×D250×H150mm | | W500×D300×H200mm |
| Electrode dimensions | W250×D170mm | | W400×D200mm |
| Reaction chamber material | Aluminum A5052 | | |
| Reaction gas | 2 Systems (Ar, O ₂) | | |
| Filling gas | N ₂ or dry air | | |
| Reaction gas flow control | Flow meter | Polymer controller | |
| Vacuum pump | Optional | Oil pump (Approx. 345L/min) standard | Oil pump (Approx. 500L/min) standard |
| Gas inlet | 2 reaction gas, 1 filling gas | | |
| External dimensions | W540×D600×H600mm | W540×D600×H600mm | W700×D700×H1285mm |
| Weight | Approx. 100kg | Approx. 105kg | Approx. 180kg |
| Power supply | Single phase AC100V 10A | 3-phase AC200V 6A | 3-phase AC200V 8A |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Gas Plasma Cleaners | RIE·DP Mode, Multi-stage Parallel Plate Electrode

PDC610

Made in
Japan

Maximum output 600W

Electrode dimensions $\Phi 215 \times 305 \text{ mm}$

Small, multi-stage plasma cleaners.



Usage

- Improving the adhesion of various materials, surface modification treatment
- Ashing, etching treatment
- Plastic packaging pre-treatment
- Printed circuit board coating pre-treatment
- LED related parts processing
- Cleaning of electronic parts
- Photoresist stripping or surface residue removal
- Cleaning of precision parts such as optical, metal, mechanical parts
- Surface modification of resin products including fluororesin
- FCC compatible

Features

- 600w small high output power
- Electrode can be selected for 1, 2, 3 stages (specify at the time of purchase)
- It can process vertical sample storage box
- It can select RIE/DP mode (specify at the time of purchase)
- Data logger (optional)
- Adaptation point memory function (optional)

Plasma discharge



Internal chamber



Specification

| Model | PDC610 |
|-----------------------------|--|
| Plasma mode | RIE/DP (fixed mode) |
| Electrode structure | 3-stage independent parallel plates (fixed) |
| High frequency output | Max600W |
| Frequency | 13.56MHz |
| Control and display method | PLC · touch screen |
| Reaction chamber dimensions | W350×D270×H300mm |
| Reaction chamber material | A5052 |
| Reaction gas | 2 Systems (Ar, O ₂) |
| Filling gas | N ₂ or dry air |
| Vacuum pump | Oil pump (Approx. 345L/min) - standard equipment |
| External dimensions | W600×D722×H700mm |
| Exterior material | Stainless steel |
| Power supply | 3 Phase AC200 50/60Hz 15A |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Gas Plasma Cleaners | RIE-DP Mode, Multi-stage Parallel Plate Electrode

V1000/1000X/1000XS

Made in Japan

High frequency output 100~1000W (V1000)/100~1000W&1500W (V1000X/1000XS)

Electrode dimensions 280×280mm (V1000) 300×300mm (V1000X) 400×375mm (V1000XS)

Plasma cleaning machine with 13.56MHz high frequency output, featuring RIE-DP dual modes.

Using plasma surface treatment equipment.
Focusing on “easy to use” high-performance desktop models, easy to retrieve data, centered on electronic materials, suitable for various purposes.

Usage


- Ashing and etching of silicon wafers
- Cleaning of substrate surface dirt
- IC/LED packaging, BGA/CSP substrate processing
- Dry cleaning of electronic materials for semiconductor-related parts
- Removal of natural oxide films and organic substances
- Interface active treatment
- Surface treatment of contaminants

Features

- Electrode structure can be changed to enhance plasma effect.
- Equipped with a high-precision adapter and high-performance RF power supply.
- Custom specifications available as needed.

Internal chamber**Specification**

| Model | V1000 | V1000X | V1000XS |
|----------------------------------|---|----------------------------------|----------------------|
| Plasma mode | DP/RIE | | |
| High-frequency power supply | 13.56MHz | | |
| Output power | Max1000W | Max1000 & 1500W | |
| Effective electrode dimensions | W280×D280mm | W300×D300mm | W400×D375mm |
| Electrode | Parallel plates | 2-stage electrodes (independent) | Plate (multipurpose) |
| Reaction gas | Polymer control system 2 systems (filling gas route installed separately) | | |
| Vacuum pump (standard equipment) | Approx. 1000L/min | Approx. 1500L/min | 1000 & 1500L/min |
| Safety device | Door detection switch, interlock mechanism, emergency stop switch, temperature overheat preventer, etc. | | |

Piping connection diagram

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Gas Plasma Cleaners | Handheld

AP-4000

High frequency output 500W/750W



Usage

- Enhance adhesion of various materials and surface modification treatment.
- Clean various substrates.
- Electronic components, surface modification of materials.
- Remove organic matter and natural oxide films.

Features

- Compact design can be embedded in production lines.
- ECO series comes with fixed power (optional), simple operation, and convenient use.
- AP-4000/AP-4000R allows power adjustment via touchscreen, suitable for laboratory or production lines.

Specifications

| Products | | AP-4000ECO | AP-4000RECO | AP-4000 | AP-4000R |
|------------------------------|------------------|--|-------------|--|-------------|
| Basic Specifications | Power supply | AC220V single phase $\pm 10\%$ (10A) 50/60Hz Ground wire must be grounded | | | |
| | Air-gas | CDA (99%), N ₂ , Ar, synthetic gas (H ₂ 5%+N ₂ 95%) | | | |
| | Air-gas pressure | 0.4MPa | | | |
| Air volume: pressure control | | Pressure regulator | | | |
| Pressure setting | | 0.08-0.2MPa | | | |
| | | Adjusted according to nozzle shape | | | |
| Processing width | | ~15mm | ~60mm | ~15mm | ~60mm |
| | | Based on nozzle shape and distance to target | | | |
| External dimensions | | 357×534×186 | | | |
| Weight | | 24kg | 30kg | 25kg | 31kg |
| Operation | | Button switch or external switch (toggle) | | Touch screen switch or external switch (toggle) | |
| Output power | | Fixed 500W or 750W | | 500-1000W (Touch screen adjustable ± 50 KW) | |
| Power consumption | | Under 1KW | Under 1.2KW | Under 1KW | Under 1.2KW |
| Others | | Various I/O available | | Touchscreen control, temperature monitoring function, various I/O, other settings and alarms | |

Gas Plasma Cleaners | Online Type with Automatic Transfer System

YSP62

Made in Japan

High frequency output 600W

Reaction stand 385×362mm

Customizable sample transfer system.



Features

- Quick and convenient operation: Insert sample storage box, automatic processing, remove storage box.
- Flexible substrates can also be transferred.
- Transfer speed can be adjusted between Hi-Low.
- Remote operation available via network (optional).

Specifications

| Model | YSP62 |
|-----------------------------|--|
| Mode | RIE (DP mode optional) |
| High frequency output power | 600W |
| Reaction stage dimensions | W385×D362mm |
| Internal cavity dimensions | W395×D378×H45mm (inner size) |
| Reaction gas | 2 Inlet gas systems (oxygen, argon) |
| Reaction gas flow control | Mass flow controller |
| Vacuum pump | Rotary vacuum pump (Approx. 500L/min) |
| Gas inlet | 2 for reaction gas, 1 for filling gas, 1 for air |
| External dimensions | W1860×D1100×H1800mm |

Water Purifiers

Contents

| | |
|----------------------------------|-----|
| Water Purification Systems | 140 |
|----------------------------------|-----|



Reduces TOC value, high-performance pure water production device meeting JIS K0557 A4 standard.

■ Water quality specifications japan industrial standard JIS K 0557


| Items ※1 | Type and water quality | | | |
|--|-------------------------|---------------------------|-------------------------|-------------------------|
| | A1 | A2 | A3 | A4 |
| Conductivity mS/m (25℃) | Below 0.5 | 0.1※2,3 or below | 0.1※2 or below | 0.1※2 or below |
| Organic carbon (TOC) mgC/L | Below 1 | Below 0.5 | Below 0.2 | Below 0.05 |
| Zinc µgZn/L | Below 0.5 | Below 0.5 | Below 0.1 | Below 0.1 |
| Silicate µgSi ₂ /L | — | Below 50 | Below 5.0 | Below 2.5 |
| Chloride µgCl ⁻ /L | Below 10 | Below 2 | Below 1 | Below 1 |
| Sulfate µgSO ₄ ²⁻ /L | Below 10 | Below 2 | Below 1 | Below 1 |
| Usage examples | Cleaning of Instruments | General chemical analysis | Trace component testing | Trace component testing |

※1. Choose required items based on testing methods.

※2. Conductivity meter actual measured value at pure water device sampling port.

※3. Final process ion exchange device sampling port connected to precision filter, conductivity meter actual measured value below 0.01ms/(25℃).

■ Classified by JIS K 0557



Main products and uses

| Type | Model | Collection method | Pure water collection | | Category (JIS K 0557) | | Usage examples (JIS K 0557) | | | | Page |
|-----------------------------|------------|--|-----------------------|--------------------|-----------------------|--------------------|-----------------------------|-----|-----|-----|------|
| | | | Distilled water | Ion exchange water | Distilled water | Ion exchange water | A1 | A2 | A3 | A4 | |
| Standard type | WG251/1001 | Ion exchange → distillation → filtration | ○ | ○ | A4 | A4 | ○ ○ | ○ ○ | ○ ○ | ○ ○ | 146 |
| Ultra-high-performance type | WGH201 | Ion exchange → distillation → UV lamp → filtration | ○ | ○ | A4 | A4 | ○ ○ | ○ ○ | ○ ○ | ○ ○ | 144 |
| Economical type | WG204 | Ion exchange → distillation | ○ | ○ | A4 | A4 | ○ ○ | ○ ○ | ○ ○ | ○ | 148 |

○ Distilled water ○ Ion exchange water

Model overview

Automatic distillation series

Ultra-Pure Water Production Device
Water Purification Systems

| Model | Ion exchange collection volume (Approx L/min) | Distilled water production volume (Approx L/h) | Collection method |
|-----------------|---|--|-------------------|
| WG251 WGH201 | 0.5~1.0 | 1.5 | |
| WG204 | 1.0 | 1.5 | |
| WG1001 | 0.5~1.0 1~1.5 1~1.5 | 5 5 10 | |

II. Ion exchange water collection
II. Distilled water collection

Model overview

Automatic distillation series

| | ① Water quality monitoring meter | ② Water quality anomaly Alarm | ③ Residual water drainage Initial flow water drainage | ④ Automatic drainage of distilled water | ⑤ Automatic cleaning and drainage of the distiller | ⑥ Cooling water volume, automatic Adjust-ment | ⑦ Dry burn prevention | ⑧ heater overheat detection | ⑨ Leakage detection | ⑩ Water outage detection | TOC digital display | Leakage protection | Circuit breakers | Water supply pump no-load prevention | Ion exchange water collection capacity Setting function | Distilled water collection capacity setting function | Consumable replacement period disp |
|--------|----------------------------------|-------------------------------|---|---|--|---|-----------------------|-----------------------------|---------------------|--------------------------|---------------------|--------------------|------------------|--------------------------------------|---|--|------------------------------------|
| WG251 | ● | ● | ● | ● | | ● | ● | ● | ● | ● | | ● | | ● | ● | ● | ● |
| WG1001 | ● | ● | ● | ● | | ● | ● | ● | ● | ● | | ● | | ● | ● | ● | ● |
| WGH201 | ● | ● | ● | ● | | ● | ● | ● | ● | ● | | ● | | ● | ● | ● | ● |
| WG204 | ● | | ● | | | ● | ● | ● | ● | ● | | ● | | ● | | | ● |

● Digital water quality meter shows measurement values. (WG204 type, analog display)

● with TOC meter

- ① Water quality monitor: Display the water quality of each part, digital display. (WG204 shows in stages)
- ② Water quality abnormal alarm: Displays when water quality decreases, digital display.
- ③ Initial water drainage: Automatically discharges the initial water of volatile impure mixtures periodically.
- ④ Automatic drainage of distilled water: Automatically distills and drains periodically to prevent scale buildup and concentration of distilled water in the distiller.
- ⑤ Automatic cleaning and drainage of the distiller: Ensure swirling inside the distiller for improved cleaning effect after periodic distillation when power is ON, then drains.
- ⑥ Automatic adjustment of cooling water volume: Automatically adjusts the cooling water volume according to the operating status.
- ⑦ Dry-run prevention: Prevent dry running by controlling the liquid level in the distiller and pressure switch.
- ⑧ Overheating detection of the heater: Detects the surface of the heater through an electronic circuit.
- ⑨ Water leakage detection: Stops operation in case of water leakage.
- ⑩ Water outage detection: Stops operation in case of water interruption.

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Water Purifiers | Ultra-High Performance Type

WGH201

Distillation capacity 1.5L/h

Purified water collection Distilled water/ion-exchange water

Water quality level Distilled water: ASTM D1193 Type1; Ion-exchange water: JIS K0557 A4

Simultaneously collects distilled water and ultrapure water with a resistance of over 18MΩ·cm.

Meet ASTM D1193 type1 and JIS K0557 A4 water quality specifications.



Features

- Collect ultrapure water (resistance over 18MΩ·cm).
- Pre-treatment → Ion exchange → Distillation → Ion exchange → Filtration. The water tank is equipped with a UV sterilization lamp.
- Automatically circulates internally periodically to prevent water degradation of the ultrapure water collection path and water tank.
- Equipped with a manual internal circulation function.
- Conform to the water quality standards for elemental analysis specified by the NADCAP certification body for aircraft component manufacturers, ASTM D1193 TYPE 1.
- Meets a wide range of needs from general analysis to various trace analyses. ※ ASTM D1193 Type1 is the highest grade recognized by ASTM International for chemical analysis water.

Specifications

| Model | | | WGH201 |
|----------------------|---|--------------------|--|
| System · Performance | Water quality level | | Ultrapure water: ASTM D1193 Type1 (over 18MΩ·cm); Ion-exchange water: JIS K0557 A4 |
| | Water collection method | Ultrapure water | Pre-treatment resin → Ion exchange → Distillation → Water tank → Ion exchange → Filtration Automatic circulation (manual circulation) |
| | | Ion exchange water | Pre-treatment resin → Ion exchange → Filtration |
| | Water quality | | Ultrapure water: Over 18MΩ·cm; Ion-exchange water: Approx. 18MΩ·cm; Distilled water: Approx. 1.0MΩ·cm |
| | Distilled water production volume | | Approx. 1.5L/h |
| | Purified water collection volume | | Ultrapure water: Approx. 0.9L/min; Ion-exchange water: 1.0L/min |
| | Water collection volume setting range | | 0.1~26L for fixed quantity collection/continuous collection |
| drainage | | | Left and right selectable connection method/water pipe connection |
| Composition | Pre-treatment column | | 0.1μm hollow fiber membrane Activated carbon (PWF-1) (exchange volume approx. 5000L) |
| | Ion-exchange water resin | | CPS-S 4L 1 piece (exchange volume approx. 700L) |
| | Ultrapure water resin | | CPC-S 4L, CPC-H 3L 1 piece each (exchange every 1 year) |
| | Final filtration | | 0.1μm filtration membrane (installed at the water collection port) (exchange volume approx. 500L) |
| | Distilled water storage tank capacity | | PE 30L (full water setting 26L/10L selectable) |
| | Distilled water storage tank sterilization lamp | | UV sterilization lamp included |
| | Water collection stand | | Drawer type, withstands 10kg |
| | Multi-purpose distilled water collection port | | With connector hose (for water tank collection/drainage) |
| Specifications | External dimensions (mm) | | W600×D660×H1650 |
| | Power supply (50/60Hz) Rated current | | AC220V 8.5A |
| Accessories | | | Supply hose 1 piece, hose clamp 1 piece, scale cleaning agent (1kg), ion exchange resin column (CPC-S) 2 pieces, high purity column (CPC-H) 1 piece, pre-treatment column 1 piece, membrane filter 2 pieces, sealing tape 1 piece, drainage hose 2m × 1, water supply and drainage hose 1 piece, water supply hose set 1 piece, water tank drainage hose set 1 piece |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps


12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Water collection flow chart



Water quality analysis example

| Items | Specifications | Specifications | Water collection |
|---------------------------------|----------------|-----------------|-----------------------|
| | JIS K0557 A4 | ASTM D1193Type1 | Ultrapure water |
| Conductivity mS/m at 25 °C | ≦ 0.1 | ≦ 0.00555 | <0.00555 |
| Total organic carbon (TOC) µg/L | ≦ 0.05 | ≦ 0.05 | 0.02 |
| Zinc µg/L | ≦ 0.1 | — | <0.1 |
| Silicate µg/L | ≦ 2.5 | — | <0.1 |
| Total organosilicone µg/L | — | ≦ 3 | 0.21 |
| Chloride µg/L | ≦ 1 | ≦ 1 | <0.1 |
| Sulfate µg/L | ≦ 1 | — | <0.1 |
| Sodium µg/L | — | ≦ 1 | <0.1 |
| Total water quality | | | ASTM Type1 and JIS A4 |

Options/consumables

Options

| Product name | Specifications | Product number |
|-----------------------------------|-----------------|----------------|
| Stand | AS250 | 253174 |
| Water collection pipe set | OWG24 | 253204 |
| Water supply port set | OWH10 | 253686 |
| Raw water pressure reducing valve | OWG42 | 253769 |
| Connection set G | WL100+WG series | 253668 |
| Drainage component | OWI10 | 253211 |
| Water collection stand | OWL40 | 253266 |
| External alarm output | OWG60 | 262780 |

Consumables

| Product name | Specifications | Product number |
|----------------------------------|----------------|----------------|
| Membrane filter (2 pieces) | MFRL727 | 9020010004 |
| Ion exchange resin | CPC-S | 253080 |
| High purity column | CPC-H | CPCNS30011 |
| Pre-treatment column | PWF-1 | 253099 |
| Air vent filter for storage tank | AVF-1 (4210) | 9020020001 |
| UV sterilization lamp | OWG28 | 253773 |
| Scale cleaning agent | 1kg | 8190010001 |
| Scale cleaning agent | | |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Water Purifiers | Standard Type

WG251/1001

| | | | | | | | | |
|-----------------------|----------------|---------------|-------------------|---|---------------------------|------------------------------------|---------------------|--|
| Distillation capacity | 1.5L/h (WG251) | 5L/h (WG1001) | Collection method | Ion exchange method → Distillation method → Filtration method | Purified water collection | Ion-exchange water/Distilled water | Water quality level | A4 (ion-exchange water)/A4 (Distilled water) |
|-----------------------|----------------|---------------|-------------------|---|---------------------------|------------------------------------|---------------------|--|

Reduces TOC value, high-performance pure water production device meeting JIS K0557 A4 standard.



Specifications

| Model | | WG251 | WG1001 |
|----------------------|---|---|--|
| System · Performance | Collection method | Ion exchange → distillation → filtration | |
| | Water supply method | Quick-connect resin water pipe connection hose/Free hose connection | |
| | Drainage method | Left and right selectable connection method/connection hose | |
| | Purified water collection | Ion-exchange water · Distilled water | |
| | Distilled water production volume | Approx. 1.5L/h | Approx. 5L/h |
| | Distilled water collection volume | 0.5~1L/min | |
| | Ion-exchange water collection volume | 0.5~1L/min | |
| | Water collection capacity setting range | 0.1~30L/continuous collection | 0.1~100L/continuous collection |
| | Condenser | Hard glass | |
| | Heater | Ceramic heater 1.4KW | Ceramic heater 1.9KW × 2 pieces |
| Composition | Pre-treatment column | 0.1μm hollow fiber membrane + Activated carbon | |
| | Ion exchange resin column | CPC-S 4L 1 piece (activated carbon high purity exchange column) | CPC-S 4L 2 pieces (activated carbon high purity exchange column) |
| | Distilled water · Ion water final filtration | 0.1μm × 2 (filters) | |
| | Leakage detection | Water leakage detection supply automatic shutoff | |
| | Distilled water storage tank | 30L resin tank | 100L resin tank |
| | Distilled water UV sterilization lamp | Optional | |
| | Water collection stand | Drawer type · Withstands 10kg 5L | Drawer type · Withstands 20kg 10L |
| | Multi-purpose distilled water collection port | WR type water supply connection at right side of main body | |
| | Water level detection | Water level switch 5-stage detection | |
| | Raw water pressure range | 0.5~5×100KPa (0.5~5kgf/cm ²) | |
| Specifications | Power supply (50/60Hz) | Single phase AC100V 15A | Single phase AC200V 20A |
| | External dimensions · weight | Length 600 × Width 660 × Height 780mm Approx. 55kg | Length 600 × Width 660 × Height 1,650mm Approx. 105kg |
| | Water level display | LED display | |
| Display · Settings | Water quality display | Digital (conductivity · resistance switch) | |
| | Other displays | Consumable replacement display (ion exchange resin column · pre-treatment column, sterilization lamp, filter), text display alarm prompt, consumable replacement history display (up to 20 items), alarm occurrence history display, English/Japanese display switch, maintenance display | |
| | Accessories | Water supply and drainage pipe, connection hose accessories, scale cleaning agent, pre-treatment column, ion exchange resin column, filter, hose clamp, sealing tape | |

Features

- High purity water collected by ion exchange → distillation → filtration.
- By using the pre-treatment column (activated carbon + 0.1μm hollow fiber membrane), raw water impurities such as germs, residual chlorine, and organic matter are removed
- Use high-performance ion exchange resin column (CPC-S, 4L), with high purity conductivity and TOC values.
- Water collection port includes a filter (0.1μm hollow fiber membrane) as standard. Distilled water storage tank capacity is 30L (WG251), 100L (WG1001).
- The display panel indicates the replacement period of consumables.
- The drawer-type water collection stand has a drainage function, making it safe to collect water.

Control panel



Water collection stand




- 1 Sterilizers
- 2 Granulation and Spray Dryers
- 3 Muffle Furnaces
- 4 Ovens
- 5 Incubators
- 6 Plasma Equipment
- 7 Water Purifiers
- 8 Baths
- 9 Water Circulators
- 10 Rotary Evaporators
- 11 Freeze Dryers & Cold Traps
- 12 Stirrers & Shakers
- 13 Washers
- 14 Analysis and Test Devices
- 15 Options

Made in
Japan

Water collection flow chart

WG251/1001 type



Water quality analysis example


Japanese Industrial Standards JIS K 0557

| Items | JIS K 0557 Standard A4 | Measurement value | | | | | |
|--|------------------------|--------------------|------------------|-----------|--------------------|------------------|-----------|
| | | WG251 | | | WG1001 | | |
| | | Ion exchange water | Dis-tilled water | Cate-gory | Ion exchange water | Dis-tilled water | Cate-gory |
| Conductivity (mS/m) | Below 0.1 | 0.0055 | 0.081 | A4 | 0.0056 | 0.07 | A4 |
| Organic carbon (TOC) (μgC/L) | Below 50 | 4 | 33 | A4 | 10 | 20 | A4 |
| Zinc (μgZn/L) | Below 0.1 | <0.01 | <0.01 | A4 | <0.01 | <0.01 | A4 |
| Silicate (μgSiO ₂ /L) | Below 2.5 | <1.0 | <1.0 | A4 | <1.0 | <1.0 | A4 |
| Chloride μgCl/L | Below 1 | <0.1 | <0.1 | A4 | <0.1 | <0.1 | A4 |
| Sulfate μgSO ₄ ²⁻ /L | Below 1 | <0.1 | <0.1 | A4 | <0.1 | <0.1 | A4 |
| Total level | | A4 | A4 | | A4 | A4 | |

Use high-quality ion exchange resin column (CPC-S)

Increases resin volume by approximately 25%.

Adds activated carbon to the high-performance ion exchange resin, reducing TOC values.



- ① High-performance pre-treatment column (PWF-1)
Removes trihalomethane, improving water collection quality.
- ② Ion exchange resin column (CPC-S)
- ③ Boiler for distillation
- ④ Storage tank for distillation

Filters



WG251/1001 type comes with filter

Water collection stand



Easy maintenance
Quick-connect pre-treatment column,
ion exchange resin column, convenient
for replacement work

Options/consumables

Options

| Product name | Specifications | Product number |
|-----------------------------------|--------------------|----------------|
| Stand | AS250 | 253174 |
| Water collection pipe set | OWG24 | 253204 |
| Water supply port set | OWH10 | 253686 |
| Raw water pressure reducing valve | OWG42 | 253769 |
| Sterilization lamp | OWG20 (for WG251) | 253202 |
| Sterilization lamp | OWG22 (for WG1001) | 253203 |
| Connection set G | WL100+WG series | 253668 |

Consumables

| Product name | Specifications | Product number |
|----------------------------------|----------------|----------------|
| Pre-treatment column | PWF-1 | 253099 |
| Ion exchange resin column | CPC-S | 253080 |
| Membrane filter (2 pieces) | MFRL727 | 9020020004 |
| Air vent filter for storage tank | AVF-1 (4210) | 9020020001 |
| UV sterilization lamp | OWG28 | 253773 |
| Scale cleaning agent | 1kg | 8190010001 |

- Sterilizers 1
- Granulation and Spray Dryers 2
- Muffle Furnaces 3
- Ovens 4
- Incubators 5
- Plasma Equipment 6
- Water Purifiers 7
- Baths 8
- Water Circulators 9
- Rotary Evaporators 10
- Freeze Dryers & Cold Traps 11
- Stirrers & Shakers 12
- Washers 13
- Analysis and Test Devices 14
- Options 15

Water Purifiers | Economic

WG204

Made in Japan

Distillation capacity 1.5L/h

Collection method Ion exchange method → Distillation method

Purified water collection Ion-exchange water/Distilled water

Water quality level A4 (Ion-exchange water)/A4 (Distilled water)


Economical type, water quality JIS K 0557 A4 standard.



Features

- Collect high purity water by ion exchange method → distillation method.
- Sterilize the raw water through the pretreatment column (activated carbon + 0.1μm hollow fiber membrane), remove methane, residual chlorine, and impurities.
- Use high-performance ion exchange resin column (CPC-S, 4L), with high purity conductivity and TOC values.
- Optional water collection port filter is available.

Water collection flow chart



Water quality analysis example (Japanese Industrial Standards JIS K 0557)

| Items | Ion exchange water | | Distilled water | |
|--|--------------------|----------|-------------------|----------|
| | Measurement value | Category | Measurement value | Category |
| Conductivity (mS/m) | 0.0055 | A4 | 0.081 | A4 |
| Organic carbon (TOC) (μgC/L) | 4 | A4 | 33 | A4 |
| Zinc (μgZn/L) | <0.01 | A4 | <0.01 | A4 |
| Silicate (μgSiO ₂ /L) | <1.0 | A4 | <1.0 | A4 |
| Chloride μgCl ⁻ /L | <0.1 | A4 | <0.1 | A4 |
| Sulfate μgSO ₄ ²⁻ /L | <0.1 | A4 | <0.1 | A4 |
| Total level | A4 | | A4 | |

Specifications

| Model | | WG204 |
|----------------------|---|--|
| System · Performance | Collection method | Ion exchange → distillation |
| | Water supply method | Quick-connect tap water pipe connection/with stop valve |
| | Drainage method | Left and right analytical method selection/water pipe connection |
| | Purified water collection | Ion-exchange water · Distilled water |
| | Distilled water production volume | Approx. 1.5L/h |
| | Distilled water collection volume | Approx. 1L/min |
| | Ion-exchange water collection volume | Approx. 1L/min |
| | Water collection capacity setting range | Continuous water sampling |
| | Condenser | Hard glass |
| | Heater | Ceramic heater 1.4KW |
| Composition | Pre-treatment column | 0.1μm hollow fiber membrane + Activated carbon |
| | Ion exchange resin column | CPC-S 4L 1 piece (high-purity activated carbon filter column) |
| | Distilled water · Ion water final filtration | Optional |
| | Leakage detection | Water leak detection water supply solenoid valve forced shutoff |
| | Distilled water storage tank | 20L resin box |
| | Distilled water UV sterilization lamp | - |
| | Water collection stand | - |
| | Multi-purpose distilled water collection port | Right side of the main body 1 piece |
| | Water level detection | Two-stage water level switch detection |
| | Raw water pressure range | 0.5~5×100KPa (0.5~5kgf/cm ²) |
| Specifications | Power supply | AC100V 15A |
| | External dimensions · weight | Dimensions: 600mm length×575mm width×780mm height, approx. 48kg |
| | Water level display | Connecting tube water level display |
| Display Settings | Water quality display | 5-stage LED display for conductivity |
| | Other displays | Consumable replacement display (ion exchange resin column) |
| | Accessories | Water supply/drainage pipe, connection hose accessories, scale cleaning agent, pretreatment column, ion exchange resin column, hose clamps |

Options/consumables

Options


| Product name | Specifications | Product number |
|-----------------------------------|-----------------|----------------|
| Stand | AS250 | 253174 |
| Water collection pipe set | OWG24 | 253204 |
| Water supply port set | OWH10 | 253686 |
| Raw water pressure reducing valve | OWG42 | 253769 |
| Connection assembly G | WL100+WG series | 253668 |
| Drainage component | OW10 | 253211 |

Consumables

| Product name | Specifications | Product number |
|----------------------------------|----------------|----------------|
| Pre-treatment column | PWF-1 | 253099 |
| Ion exchange resin column | CPC-S | 253080 |
| Filters (2 pcs) | MFRL727 | 9020010004 |
| Air vent filter for storage tank | AVF-1 (4210) | 9020020001 |
| Scale cleaning agent | 1kg | 8190010001 |

Contents

| | |
|--|-----|
| Constant High Temp. Water Baths | 148 |
| Immersion Constant Temperature Devices | 149 |
| Constant High Temp. Oil Baths | 152 |
| Constant Temperature Water Baths with Oscillating Unit | 156 |
| Immersion Coolers | 158 |
| Constant Low Temp. Water Baths | 159 |
| Heating Blocks | 165 |



1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices


15 Options

High operational efficiency and product diversification.

Constant temperature water baths are widely used in fields ranging from biochemistry to semiconductors, are enhanced when combined with other devices rather than used alone, leading to diverse applications. Yamato Scientific's constant temperature water baths aim to create a product lineup while further improving product functionality to meet various needs.

Function · Safety Device

Electrical circuit breakers



Model selection list

| Category | Product name | Temp. control range | Temp. adjustment accuracy | Temp. distribution accuracy | Series | Model | Features | Internal capacity (L) | Page No. |
|--------------------------------|--|--------------------------------|---------------------------|-----------------------------|--------|----------------------------|--|-----------------------|----------|
| High Constant Temp. Water Bath | Constant high temp. water bath | RT+5~80°C | ± 0.02°C | ± 0.3°C | BA | BA310C/410C/510C/610C/710C | Standard type, fixed temp. operation | 27/42/70/109/144 | 148 |
| | | | | | BK | BK310C/410C/510C/610C/710C | Standard type, program operation | 27/42/70/109/144 | 148 |
| | Immersion constant temp. devices | -20~+100°C | ± 0.05°C | — | BF | BF201/401/501 | Immersion type | — | 149 |
| | | 0~200°C | ± 0.05°C | — | BF | BF601 | Immersion type | — | 149 |
| | Constant high temp. water bath | RT+5~boiling water temperature | ± 3°C | ± 5°C | BS | BS200/401/601/660 | Liquid expansion type | 4.7/9/11/14 | 150 |
| | Constant high temp. water bath | RT+5~95°C | ± 1.5°C | — | BM | BM100/401 | Economical type | 4/7/7 | 151 |
| | | RT+5~90°C | ± 2°C | — | BM | BM510C | Economical type | 4 | 151 |
| | Constant high temp. oil bath | RT+10~270°C | ± 0.1°C | ± 0.3°C | BOA | BOA311 | Large dimensions | 37 | 152 |
| | Constant high temp. oil bath | RT+5~240°C | ± 0.3°C | ± 0.3°C | BOG | BOG100 | Combination type | 1 | 153 |
| | Constant high temp. oil bath | RT+5~240°C | ± 0.3°C | ± 0.3°C | BOG | BOG200 | Combination type | 2.2 | 153 |
| | Constant high temp. oil bath | RT+5~250°C | ± 0.3°C | ± 0.3°C | BOS | BOS100 | Combination type | 1 | 153 |
| | Constant high temp. oil bath | RT+5~240°C | ± 0.3°C | ± 0.3°C | BOS | BOS200 | Combination type | 2.2 | 153 |
| | Constant high temp. oil bath | 0~180°C | ± 0.05°C | — | BZ | BZ100/100D/200/300·BF601 | Immersion type constant temperature device combination | 12/13/20/27 | 154 |
| | Constant high temp. oil bath | RT+5~199°C | ± 0.5°C | — | BO | BO500 | Stainless steel plate | 5.2 | 154 |
| | Constant high temp. oil bath | RT+5~180°C | ± 2°C | — | BO | BO510C/601 | Economical type | 4/7 | 155 |
| | Oscillating unit | — | — | — | BN | BN300 | For constant temperature water baths | — | 155 |
| Low Constant Temp. Water Bath | Oscillating water bath | — | — | — | BW | BW101/201/400 | Reciprocating oscillating | 12/20/30 | 156 |
| | Constant high temp. oscillating water bath | RT+5~80°C | ± 0.02°C | ± 0.1°C | BT | BT101/311 | Reciprocating horizontal oscillating, pump stirring | 19/23/34 | 157 |
| | Immersion Coolers | -20~±35°C | ± 0.1°C | — | BE | BE201/201F/301 | Immersion type | — | 158 |
| | Constant low temp. water bath | -80°C~0°C | ± 0.3°C | — | BLG | BLG100 | For raw material synthesis | 0.3 | 159 |
| | Constant low temp. water bath | -40°C~0°C | ± 0.3°C | — | BLG | BLG200 | For raw material synthesis | 1 | 159 |
| | Constant low temp. water bath | -30°C~+80°C | ± 0.1°C | ± 0.3°C | BB | BB311C/411C/611C | Standard type | 6/13/26 | 160 |
| | Constant low temp. water bath | -10°C~+80°C | ± 0.1°C | ± 0.3°C | BBL | BBL111C/311C | Benchtop type | 8/13 | 162 |
| | Constant low temp. water bath | -15°C~+70°C | ± 0.1°C | ± 0.3°C | BL | BL410C/810C | Large capacity | 36/80 | 163 |
| | Constant low temp. water bath | 0~80°C | ± 0.1°C | ± 0.1°C | BV | BV100S | Semiconductor cooling | 6/10/15 | 164 |
| | Heating Blocks | RT+5~200°C | ± 0.2°C | — | HF | HF100/200 | For test tube heating | — | 165 |

| | |
|------------------------------|----------|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Constant High Temp. Water Baths | Standard Type

BK310C/410C/510C/610C/710C, BA310C/410C/510C/610C/710C

Temp. control range RT+5~80°C

Internal capacity 27L 42L 70L 109L 144L

Easy to use, with a observation window to confirm sample status, the constant high temperature water baths range from 27L to 144L in volume.



BK Series: Fixed temp. operation

BA Series: Program Operation

Features

- Feature both BK fixed temp. operation and BA program operation.
- Achieve very high temperature adjustment precision through high-precision temperature regulation methods and efficient circulation pumps stirring inside the bath.
- Temperature and time settings and indicators are all digitally displayed; quick auto-stop, auto-stop, and auto-start functions can be realized.


Safety

- Equipped with self-diagnostic circuit (temperature sensor abnormality, heater disconnection, automatic overheating prevention, SSR short circuit), an independent overheating alarm, and safety functions such as leakage protection switches that prevent overcurrent and key lock.

Specifications

| Model | | BK/BA310C | BK/BA410C | BK/BA510C | BK/BA610C | BK/BA710C |
|----------------------------|---|---|-----------------|-----------------|-----------------|-----------------|
| System | | High efficiency circulation pump stirring | | | | |
| Performance | Operating temp. range | RT+5~80℃ | | | | |
| | Temp. adjusting accuracy | ±0.02~±0.07℃ | | | | |
| | Temp. distribution accuracy | ±0.3℃ | | | | |
| | Max. temp. reaching time | Approx. 120 min | Approx. 110 min | Approx. 165 min | Approx. 160 min | Approx. 200 min |
| Composition | Internal chamber | Stainless steel plate and glass | | | | |
| | Temp. regulator | Computer PID control (BK: fixed temp. Type, BA: fixed temp. type and program type) | | | | |
| | Sensors | W sensor: platinum temperature measuring resistor Pt100ω (temperature regulator) + K thermocouple (overheating preventer) | | | | |
| | Temp. setting method | Digital Setting | | | | |
| | Temp. display method | Digital display | | | | |
| | Overheating preventer and setting method | Computer ON-OFF control Digital Setting | | | | |
| | Heater (SUS316) | 1.3KW | 2.2KW | 2.4KW | 3.5KW | 4.5KW |
| | Stirring pump | 6w | 30w | | 60w | |
| | Timer | 1 min~99 h 59 min~999 h 59 min digital setting auto start, auto stop, quick auto stop | | | | |
| Program function (BA type) | | 1 mode 30 steps, 2 mode 15 steps, 3 mode 10 steps | | | | |
| Safety device | | Self-diagnostic circuit (temperature sensor error detection, heater disconnection protection, automatic overheating preventer, SSR short-circuit), overheating preventer, overcurrent leakage protection switch and key lock functions. | | | | |
| Specifications | External dimensions (W×D×H mm) | 490×360×367 | 590×410×367 | 690×460×417 | 738×560×467 | 830×560×517 |
| | Water bath dimensions (W×D×H mm) | 300×300×300 | 400×350×300 | 500×400×350 | 548×500×400 | 640×500×450 |
| | Observation window dimensions (Width×height mm) | 240×215 | 350×225 | 440×265 | 350×225 | 440×265 |
| | Internal bath capacity | Approx. 27L | Approx. 42L | Approx. 70L | Approx. 109L | Approx. 144L |
| | Shelf sections | 5 | | 6 | 9 | 10 |
| | Drainage pipe | Φ15×20mm | | | | |
| | Power supply (50/60Hz) | AC220V | | | | |
| | Rated current | 7A | 11A | 12A | 17A | 22A |
| | Weight | Approx. 19kg | Approx. 25kg | Approx. 30kg | Approx. 36kg | Approx. 46kg |
| Accessories | | 1 shelf, 1 set of clamp supports, 1 drain pipe, 1 hose plug | | | | |
| Options | | Container installation bracket, viscosity meter bracket, external circulation pump, cooling pipe, water bath lid, external communication function (RS485) | | | | |

Data (Temperature rise curve)



Control panel



Water bath lid



1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Immersion Constant Temp. Devices | Immersion Type

BF201/401/501/601

Made in
Japan

Temp. control range RT+5~80°C (201/401/501 Type); RT+5~180°C (601 Type)

A multifunctional immersion type constant temperature device with diverse usage methods.



Specifications

| Safety device | | Basic Type | Multifunctional type | | Oil/water dual use | |
|---------------------------|---|--|--|--|----------------------------|--|
| Model | | BF201 | BF401 | BF501 | BF601 | |
| Performance | Temp. setting range | -20.0~90℃ ^{*1} | | | 0~200℃ | |
| | Operating temp. range | RT+5~80℃ ^{*2} | | | RT+5~180℃ ^{*4} | |
| | Temp. adjusting accuracy | ±0.05/0.1℃ ^{*3} | ±0.02/0.05℃ ^{*3} | | ±0.05~±0.2℃ ^{*5} | |
| Composition | Stirring structure | Propeller stirring (fixed) | Jet stirring (10 levels adjustable) | | Propeller stirring (fixed) | |
| | Heater | Stainless steel heating pipe | | | | |
| | | 1.0KW | | | 1.2KW | |
| Control | Temp. control method | PID control | | | | |
| | Temp. setting method | Digital setting | | | | |
| | Temp. display method | Green Led digital display | | | | |
| | | Minimum Display: 0.1℃ | | | Minimum Display: 1℃ | |
| | | Setting temp. and measured temp. conversion | Displayed on main display (auxiliary display shows set temperature, etc.) | | | |
| | Timer | — | | | | |
| | Timer resolution | — | | | | |
| | Operation functions | Fixed temp. operation | Fixed temp. operation: continuous operation based on set temperature, quick auto-stop Program operation: any mode 1~3, maximum 10 parts (can be cyclic, linked operation) | | | |
| | | Additional functions | Temperature pre-setting (recording, adjust 10 temperatures) | Temperature presetting (can record 10 temperatures), cumulative time function (~49999 h), lock function, temperature calibration, power failure compensation mode Variable pump discharge outlet angle function | | |
| | | | Temperature output ^{*6} | | | |
| External alarm output | | | | | | |
| Time limit reached output | | | | | | |
| Heater circuit control | External communication function | | | | | |
| | Bidirectional thyristor | | | | | |
| | Platinum temperature measuring resistor | | | | | |
| Safety device | | Self-diagnostic function (automatic overheating prevention, temperature sensor abnormality, heater disconnection, bidirectional thyristor short circuit, main relay failure), circuit protector, float type boil-dry preventer, buzzer alarm for abnormality | | | | |
| Specifications | External dimensions | W140×D138×H312mm | | | | |
| | Weight | Approx. 4kg | | | | |
| | Clampable thickness | Within 35mm | | | | |
| | Power supply | Single phase 100V 11A | | | Single-phase 100V 13A | |
| Accessories | | Test chamber BY100 Type (Polypropylene) | | | | |

*1. When using at room temperature +5°C or below, please pair with constant low-temperature water baths or immersion type coolers like BE.

*2. When using test chamber BY100.

*3. When setting temperature 37°C/80°C and using test chamber BY100.

*4. When using test chamber BZ100D, viscosity 50cSt silicone oil, and temperature exceeds 81°C, do not use water.

*5. When setting temperature 40°C/180°C, test chamber BZ100D, using silicone oil with viscosity 50cSt.


*6. 1 terminal mixed output.

The performance, operability, and safety of immersion type constant temperature devices are fully equipped, with basic type, high-performance type equipped with multiple functions and oil-related type.

Features

- The beautifully designed BF type without protrusions not only allows for simple assembly but can also be combined with shaker circuits and water baths of different materials, serving as a constant temperature water bath, oscillating incubator, or oil bath.
- From the basic type with only fixed temp. operation to multifunctional, high-performance types with program functions, external output, and external communication functions, there are 4 machine models.
- Jet stirring strength can be varied in 10 levels, and the jet direction can vary approximately 10 degrees horizontally. Selection is based on experimental purposes and corresponding conditions. (Only for 401/501 Type)
- Built-in temperature preset function allows for convenient retrieval and recording of frequently used temperatures.
- In case of abnormal occurrences, safety features such as circuit protectors, boil-dry preventers, and self-diagnostic functions are included.

Jet Stirring Structure Diagram



Constant High Temp. Water Baths | Liquid Expansion Type


BS200/401/601/660

Made in
Japan

Temp. control range RT+5~Boiling RT(Water)

Temp. distribution accuracy ± 3.0 (at 70°C)

Standard configuration includes overheating preventer, overcurrent leakage circuit breaker, ensuring high safety.



Features

BS200 type

- 4 types of sample baskets are included for efficient use according to container size. By removing the top cover, a 5L beaker can also be used.

BS401 type

- 4 types of sample baskets are included for efficient use according to container size. By removing the top cover, a 5L beaker can also be used.

BS601 type

- Equipped with rings and lids for use according to container size, allowing 6 samplers to be used simultaneously.

BS660 type

- Include chains and lids that can be used separately according to the container, allowing simultaneous use of 6 samplers.



Use according to container size.



The top lid is removable.



Simple water level adjustment can be done by moving the overflow up and down



The clamps are die-cast to provide good support for the instruments.

Control panel



Specifications

| Model | | BS200 | BS401 | BS601 | BS660 |
|----------------|--------------------------------------|---|---|---|---|
| Performance | Operating temp. range | RT+5°C to boiling temperature (water) | | | |
| | Temp. adjusting accuracy | $\pm 3.0^\circ\text{C}$ (at 70°C) | | | |
| | Temp. distribution accuracy | $\pm 5.0^\circ\text{C}$ (at 70°C) | | | |
| | Max. temp. reaching time | Approx. 30 min | Approx. 40 min | Approx. 55 min | Approx. 60 min |
| Composition | Internal chamber | Stainless steel plate | | | |
| | Heater | Copper tube heater (nickel plated) | | | |
| | | 900W | 1.4KW | | 1.3KW |
| | Safety device | Overheating preventer (liquid expansion type, self-holding relay, alarm light), overcurrent leakage circuit breaker | | | |
| Specifications | Internal dimensions | $\Phi 200 \times H149\text{mm}$ | $W328 \times D302 \times H99\text{mm}$ | $W790 \times D150 \times H102\text{mm}$ | $W502 \times D302 \times H99\text{mm}$ |
| | Effective internal dimensions | $\Phi 186 \times H115\text{mm}$ | $W298 \times D260 \times H64\text{mm}$ | $W776 \times D140 \times H73\text{mm}$ | $W468 \times D260 \times H66\text{mm}$ |
| | External dimensions | $W300 \times D364 \times H216\text{mm}$ | $W390 \times D438 \times H214\text{mm}$ | $W864 \times D288 \times H215\text{mm}$ | $W565 \times D437 \times H214\text{mm}$ |
| | Internal bath capacity | Approx. 4.7L | Approx. 9L | Approx. 12L | Approx. 14L |
| | Power supply (50/60Hz) rated current | AC100V 9A | AC100V 13A | | |
| | Opening dimensions · number | — | $\Phi 115\text{mm}$ · 4 pieces | $\Phi 115\text{mm}$ · 6 pieces | |
| | Weight | Approx. 5KG | Approx. 9.5KG | Approx. 16KG | Approx. 12KG |
| Accessories | | Chain (excluding BS200), lid, container mounting fixture, stand, horizontal rod for mounting clamp (excluding BS200), connector (for stand), thermometer (0~+100°C alcohol type), base plate, sample basket (BS200) | | | |
| Options | | Cooling water distributor (only for BS401/601) | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices



15 Options

Constant High Temp. Water Baths | Economy Type

BM510C/100/401

Operating temp. range RT+5~90°C (BM510C); RT+5~95°C (BM100/401)

Compact design with better operability for constant temperature water bath.



Features

- The water bath is with protection to prevent burns during contact. (BM100/401)
- Temperature setting uses up and down keys for numeric setting, making it easy to operate. (BM510C/401)
- The water bath has no heater, facilitating water changes and tank removal for cleaning. (BM510C)
- Equipped with an easy-drain drainage valve. (BM401)

Specifications

| Model | | BM510C | BM100 | BM401 |
|----------------|------------------------------|---|---|--|
| Performance | Operating temp. range | RT+5°C~90°C | RT+5°C~95°C | |
| | Temp. adjusting accuracy | ±1.5 (during stirring) | ±2°C (at 60°C) | ±0.3°C (at 60°C) |
| Composition | Temp. control method | PID control | ON/OFF control | PID control |
| | Temp. setting/display method | Digital setting · display | Analog setting (glass thermometer indication) | Digital setting · display |
| | Attached functions | Key lock function, RE signal transmission and reception function, power failure protection operation, deviation correction function | — | Temperature preset (1-point memory · retrieval) |
| | Heater | 1KW (at 200V) 1.44KW (at 240V) | SUS316 heating tube 500W | SUS316 heating tube 1KW |
| | Sensors | K-type thermocouple | Liquid expansion type | Pt100Ω |
| Specifications | Safety device | Water bath protection cover (ABS heat-resistant resin), self-diagnostic function (automatic overheating prevention, sensor anomaly), circuit breaker, temperature fuse, water bath detection switch | Water bath protection cover, overcurrent protection (fuse: 7A), water bath protection cover, temperature fuse | Water bath protection cover, self-diagnostic function (automatic overheating prevention, sensor anomaly, bidirectional thyristor short circuit, heater disconnection, main relay failure), circuit breaker, temperature fuse |
| | Internal bath capacity | Approx. 4L | Approx. 4L | Approx. 7L |
| | Water bath dimensions | Maximum internal diameter Φ240×bottom diameter 165×height 122mm | Internal diameter Φ200×height 120mm | Inner diameter Φ250×height 150mm |
| | External dimensions | W340×D349×H231mm | W240×D300×H150mm | W310×D360×H230mm |
| | Weight | Approx. 5.5kg | Approx. 3.5kg | Approx. 7kg |
| | Power supply | AC200~240V 5.5~6.5A | AC100V 6A | AC100V 13A |
| Accessories | | — | Thermometer (10~110°C) | — |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Constant High Temp. Oil Baths | Large

BOA311

Made in
Japan

Temp. control range RT+10~270°C

Temp. distribution accuracy $\pm 0.2^\circ\text{C}$

Large capacity constant temperature oil bath operating at 270°C.



Features

- Use jet stirring method to achieve high-temperature distribution accuracy.
- Various configuration, standard equipped with external alarm output, temperature output terminal (4~20mA and 1~5V interchangeable), RS485 communication function, key lock function, deviation correction function.


Safety

- Equipped with 3-layer overheating prevention function (when set temperature $+6^\circ\text{C}$, the heater automatically disconnects, overheating protector, independent overheating protector).
- Standard equipped with emergency stop switch.
- Equipped with a float switch for dry-burning prevention.
- Equipped with a float switch to prevent liquid overflow.
- In case of liquid leakage from the operation panel, the glass panel can also serve a protective role.
- A large display light alerts attention to high-temperature operating status.
- Standard features include self-diagnostic function, overheating protector, leakage protection, key lock function, power failure compensation function.

Operation panel



Temperature rise curve



Exhaust device



Specifications

| Model | | BOA311 |
|-----------------|-------------------------------|--|
| Performance | Operating temp. range | RT+10°C~270°C |
| | Temp. adjusting accuracy | $\pm 0.1^\circ\text{C}$ |
| | Temp. distribution accuracy | $\pm 0.2^\circ\text{C}$ |
| | Max. temp. reaching time | Approx. 70 min |
| | Temp. control | Digital controller PID control |
| Composition | Temp. sensor | PT100Ω (temperature controller), K-type thermocouple (overheating protector) |
| | Function | Fixed temp. operation |
| | Stirring method | Jet stirring |
| | Heater | SUS316 heating tube 4.5KW |
| Safety device | | Overcurrent leakage circuit breaker, self-diagnostic function, temperature sensor anomaly, heater anomaly, automatic overheating prevention function, overheating protector, independent overheating protector, oil level overflow detection float switch, detection float switch, emergency stop switch |
| Other functions | | Drain valve, large operating indicator light, external alarm output terminal, temperature output terminal (1~5V, 4~20mA alternating), RS485 communication function, key lock function, power failure compensation function, deviation correction function |
| Specifications | Effective internal dimensions | W296×D340mm×H270mm |
| | Capacity | 37L |
| | External dimensions | W531×D520×H578 (oil bath height 397) mm |
| | Power supply (50/60Hz) | AC200V 23A |
| | Weight | Approx. 37KG |
| Accessories | | Shelf 1 pc, cover 1 pc |

Features of silicone oil

Since silicone oil is specifically designed for high-temperature heat transfer, please select heat-resistant open-type silicone oil (heat-resistant dimethyl silicone oil with a viscosity below 100mm²/s [cSt]).

Recommended silicone oil:

Operating temperature below 200°C: Toshiba TSF458-50

Operating temperature 200~270°C: Toshiba TSF458-100

Features of silicone oil

| Manufacturer | Toshiba Silicone Oil | |
|---------------------------------------|----------------------------|-----------------------------|
| Recommended temperature | Below 200°C | 200°C~270°C |
| Product name | TSF458-50 | TSF458-100 |
| Appearance | Light yellow transparent | |
| Specific gravity (25°C) | 0.961 | 0.963 |
| Viscosity (25°C) | 50mm ² /s (cSt) | 100mm ² /s (cSt) |
| Volatility (150°C 24h) | 0.3% | |
| Viscosity temperature coefficient | 0.59 | |
| Flash point | 325°C | 342°C |
| Pour point | Below -50°C | |
| Viscosity increase rate (300°C, 168h) | 40% | 35% |

Constant High Temp. Oil Baths | Combination Type

BOG100/200, BOS100/200

Made in
Japan

Temp. control range RT+5~240°C (BOG); RT+5~250°C (BOS)

Temp. adjusting accuracy $\pm 0.3^{\circ}\text{C}$ (at 200°C, during stirring)

A combination that maximizes temperature performance and convenient functions.

Features

- There is a transparent glass model available for selection, which offers excellent visibility.
- The controller is separated from the oil bath. It can operate safely at a distance of up to 2m (cable length).
- The bottom of the bath pot is designed with grooves to match various shapes of magnetic stirrers. The top of the stirrer is embedded in the grooves at the bottom of the pot, preventing the risk of misalignment or detachment when using hot magnetic stirrers with the oil bath.
- The recommended oil volume is marked on the bath, helping to reduce the risk of oil overflow while avoiding oil waste.
- Features a safety design with independent overload prevention device and heater protection cover.
- Insert the optional external temperature sensor into the sample container to measure the temperature of the reaction solution and control the overall system to maintain the desired temperature.



Specifications

| Model | | | BOG100 | BOG200 | BOS100 | BOS200 |
|----------------|---|--|--|--|--|--|
| Composition | Liquid bath | | Hard glass | | Stainless steel plate | |
| | Exterior material | | Stainless steel/chrome-free galvanized steel plate sintered coating | | | |
| | Heater | | Stainless steel tube heater | | | |
| | Heater capacity | | 310W at 100V | 425W at 100V | 310W at 100V | 425W at 100V |
| Performance | Set temp. range | | 0~260℃ | | | |
| | Temp. control range | | RT+5~240℃ | | RT+5~250℃ | |
| | Temp. adjusting accuracy | | ±0.3℃ (at 200℃, during stirring) | | | |
| Controllers | Temp. control method | | PID control, temperature difference control | | | |
| | Temp. setting · display method | | ▲ ▼ Digital setting based on keys · display method (display in units of 1℃) | | | |
| | Operation functions | | Fixed temp. operation | | | |
| | Additional functions | | Calibration offset function, power failure compensation selection function, LED brightness switching function | | | |
| | Heater circuit control | | Bidirectional cross control | | | |
| | Sensors | | Pt100 sensor (Class A) | | | |
| Safety device | | | Self-diagnostic function (temperature sensor anomaly detection, automatic overheating prevention function, communication anomaly, outside set temperature range anomaly), fuse (6.3A short circuit, overcurrent protection), independent overheating protector | | | |
| Specifications | External dimensions (diameter × width × height) | | Φ150*205*140 (including heating unit) mm | Φ180*235*180 (including heating unit) mm | Φ155*210*140 (including heating unit) mm | Φ185*240*180 (including heating unit) mm |
| | Internal bath dimensions (diameter × height) | | Φ140×100 mm | Φ170×140 mm | Φ140×100 mm | Φ170×140 mm |
| | Internal bath capacity | No-load recommendation | Approx. 0.8L | Approx. 1.7L | Approx. 0.8L | Approx. 1.7L |
| | | No-load maximum | Approx. 1.0L | Approx. 2.2L | Approx. 1.0L | Approx. 2.2L |
| | Controller unit (width×depth×height) | | 150×90×45mm (excluding protruding parts) | | | |
| | Power supply (frequency) | | Single-phase AC100V (50/60Hz) | | | |
| | Rated current | | 4.0A | 5.0A | 4.0A | 5.0A |
| | Weight | Liquid bath, heating unit (including cables) | Approx. 1.8kg | Approx. 2.3kg | Approx. 2.5kg | Approx. 2.9kg |
| | | Controller unit | | Approx. 0.5kg | | |
| Total weight | | Approx. 2.3kg | Approx. 2.8kg | Approx. 3.0kg | Approx. 3.4kg | |
| Accessories | | | Heating protection cover (1 piece), operation handle (outer diameter Φ12x90mm: 3 pieces), stirring bar (about Φ8x50mm: 3 pieces), plug cap (3 pieces), identification label (1 piece), branded bundling strap (3 pieces), instruction manual, warranty. | | | |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Constant High Temp. Oil Baths | Immersion Thermostat Combination

BZ100/100D/200/300·BF601

Made in Japan

Operating temp. range RT+5~180°C

Immersion thermostat BF601 combined with four types of stainless steel experimental baths.

Through the immersion thermostat BF601 that shares oil/water and stainless steel experimental bath, it can be used as a water bath or oil bath.

Features

- Programs 1-3 can be easily selected through easy-to-set dialog keys.
- It features a jet stirring function.
- Operating temperature range from room temperature +5°C~180°C.
- 12L~27L, 4 options available.



Specifications

| Model | Capacity | Dimensions |
|--------|----------|---------------|
| BZ100 | 12L | 230×390×150mm |
| BZ100D | 13L | 240×300×200mm |
| BZ200 | 20L | 300×500×150mm |
| BZ300 | 27L | 300×500×200mm |

Constant High Temp. Oil Baths | Stainless Steel

BO500

Operating temp. range RT+5~199°C

Stainless steel high-temperature constant temperature oil bath.

Features

- To prevent burns from contact, a protective film is set at the bottom of the water bath.
- This unit has a dedicated magnetic stirrer (MB800 model)

Specifications

| | Model | BO500 |
|----------------|--------------------------|-------------------|
| Performance | Operating temp. range | RT+5°C~199°C |
| | Temp. adjusting accuracy | ±0.5°C |
| Composition | Internal bath | Stainless steel |
| | Temp. regulator | PID control |
| | Temp. sensor | Pt100 |
| | Temp. setting method | Digital setting |
| | Humidity display method | Glass thermometer |
| Specifications | Heater | Heating tube 700W |
| | Internal bath capacity | Approx. 5.2L |
| | Water bath dimensions | Φ240×height 130mm |
| | External dimensions | W250×D290×H130mm |
| | Weight | Approx. 1.4kg |
| | Power supply | AC100V 8A |



1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices



15 Options

Constant High Temp. Oil Baths | Economy Type

BO510C/601

Temp. control range RT+10~180°C (BO510C); RT+5~180°C (BO601)

Compact design with better operability for constant temperature oil bath.



Features

- Temperature setting uses up and down keys for numeric setting, making it easy to operate.
- The oil bath is with protection to prevent burns during contact. (BO601)
- The oil bath has no heater, facilitating oil changes and bath removal for cleaning. (BO510C)
- Equipped with an easy-drain drainage plug. (BO601)

Specifications

| Model | | BO510C | BO601 |
|----------------|------------------------------|---|--|
| Performance | Operating temp. range | RT+5°C~180°C | |
| | Temp. adjusting accuracy | ±2 (during stirring) | ±2°C (at 100°C) |
| Composition | Temp. control method | PID control | |
| | Temp. setting/display method | Digital setting · display | |
| | Attached functions | Key lock function, RE signal transmission and reception function, power failure protection operation, deviation correction function | Temperature preset (1-point memory · retrieval) |
| | Heater | 1KW (at 200V); 1.44KW (at 240V) | SUS316 heating tube 1KW |
| | Sensors | K-type thermocouple | Pt100Ω |
| Specifications | Safety device | Water bath protection cover (ABS heat-resistant resin), self-diagnostic function (automatic overheating prevention, sensor anomaly), circuit breaker, temperature fuse, water bath detection switch | Water bath protection cover, self-diagnostic function (automatic overheating prevention, sensor anomaly, bidirectional thyristor short circuit, heater disconnection, main relay failure), circuit breaker, temperature fuse |
| | Internal bath capacity | Approx. 4L | Approx. 7L |
| | Water bath dimensions | Maximum internal diameter Φ240×bottom diameter 165×height 122mm | Inner diameter Φ250×height 150mm |
| | External dimensions | W340×D349×H231mm | W310×D360×H230mm |
| | Weight | Approx. 5.5kg | Approx. 87kg |
| | Power supply | AC200~240V 5.5~6.5A | AC100V 12A |

Oscillating Unit | For Constant Temp. Water Baths

BN300

Made in Japan

Oscillating amplitude 10~30mm variable

Oscillating frequency 20~160 times/min

Convenient to load and unload oscillating units.



Combination Usage Examples



BN300 model + BK410C model



BN300 model + BZ300 model



BN300 model + BV300 model

Specifications

| Model | | BN300 |
|----------------|--|---|
| Composition | Oscillating frequency | 20~160 times/min |
| | Oscillating amplitude | 10~30mm variable |
| | Oscillating frequency setting and display method | Digital setting · digital display |
| | Shakers | Capacitive motor 25w |
| | Safety device | Motor overload, power fuse |
| Specifications | Container rack spec. | 80 test tubes (diameter 15mm); 6 beakers (300ml); 3 beakers (500ml) |
| | Spring spacing | 15mm spacing (horizontal and vertical the same) |
| | External dimensions | W340×D644×H250mm |
| | Power supply (50/60Hz) rated current | AC100V 1A |
| | Weight | Approx. 10kg |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Oscillating Water Baths | Reciprocating Oscillating

BW101/201/400

Made in
Japan

Oscillating amplitude 10~40mm variable

Oscillating frequency 20~160 times/min

High-precision oscillating constant temperature water bath combined with an immersion thermostat.

The entire series adopt a mechanism with variable oscillating amplitude, ranging from small to large product series, featuring a superior structure of oscillating constant temperature water baths.

Features

- The oscillating amplitude is variable from 10~40mm, oscillating frequency is 20~160 times/min, with stepless speed regulation.
- The entire series is equipped with a digital display function to confirm the oscillating frequency.
- The oscillating rack employs a dual-spring method to ensure the protection of test tubes, beakers, and other vessels.
- The BW201/400 model is equipped with a drainage pump.
- Include a power socket for convenient combination use with immersion thermostat devices.

**Control panel****Power socket****Specifications**

| Model | | BW101 | BW201 | BW400 |
|-----------------------|--|--|---------------|--------------|
| Performance | Oscillating method | Reciprocating oscillating | | |
| | Oscillating amplitude | 10~40mm variable | | |
| | Oscillating frequency | 20~160 times/min (stepless speed regulation) | | |
| | Operating temp. range | Room temperature | | |
| Functions/Composition | Material | Water bath: stainless steel plate | | |
| | Oscillating frequency control method | Feedback phase control | | |
| | Oscillating frequency display method | Digital display | | |
| | Attached functions | Dedicated service socket AC100V 12A | | |
| | Liquid used | Water only | | |
| | Drainage method | Natural drainage | Drainage pump | |
| Safety device | | Leakage circuit breaker, motor overheat prevention (automatic reset thermal protector) | | |
| Specifications | Internal dimensions of the water bath (W×D×H mm) | 230×390×150 | 300×500×150 | 380×535×150 |
| | External dimensions (W×D×H mm) | 295×445×285 | 370×560×285 | 445×585×295 |
| | Water bath capacity | Approx. 12L | Approx. 20L | Approx. 30L |
| | Weight | Approx. 28kg | Approx. 35kg | Approx. 42kg |
| | Power supply (50/60Hz) | AC100V 12A (15A when in BF model combination state) | | |
| Accessories | | Oscillating stand | | |
| Options | | Immersion constant temperature devices | | |

Oscillating rack**Container rack spec.**

| | For BW101 | For BW201 | For BW400 |
|----------------------------|-----------|-----------|-----------|
| Test tubes (diameter 15mm) | 60 pcs | 126 pcs | 190 pcs |
| Triangular beakers (50ml) | 10 pcs | 21 pcs | 30 pcs |



1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Oscillating Water Baths | Horizontal Reciprocating Oscillating, Pump Stirring

BT101/311

Made in
Japan

Temp. control range RT+5°C~80°C

Temp. distribution accuracy $\pm 0.1^\circ\text{C}$

Oscillating amplitude 10~40mm variable

Oscillating frequency 20~160 times/min

Variable oscillating amplitude enhances usability.



Rich in operability and safety features, the oscillating structure is integrated into the bath, providing a superior structure for oscillating constant temperature water baths.

Features

- The oscillating amplitude is variable from 10~40mm, oscillating frequency is 20~160 times/min with stepless speed regulation and digital display.
- Temperature is set using up and down keys for digital setting, with high precision temperature control. Coupled with the stirrer in the bath, high precision distribution accuracy is achieved.
- Standard features include automatic overheat prevention, temperature sensor failure, and overcurrent leakage circuit breaker safety features.

Specifications

| Model | | BT101 | BT311 |
|----------------|-----------------------------------|---|---|
| System | | Horizontal reciprocating vibration, pump stirring | |
| Performance | Operating temp. range | RT+5~80°C | |
| | Temp. adjusting accuracy | $\pm 0.02 \sim \pm 0.08^\circ\text{C}$ | |
| | Temp. distribution accuracy | $\pm 0.1^\circ\text{C}$ | |
| | Max. temp. reaching time | Approx. 95 min | Approx. 120 min |
| | Oscillating frequency | 20~160 times/min, with variable knob and digital display | |
| | Oscillating amplitude | 10~40mm variable type | |
| Composition | Internal bath | Stainless steel plate | |
| | Temp. regulator | Microcomputer PID control | |
| | Sensors | Pt100 ω (temperature regulator) + K-type thermocouple (overheat preventer) | |
| | Temp. setting/display method | Digital setting · digital display | |
| | Overheat preventer | Microcomputer on/off control | |
| | Overheat preventer setting method | Digital Setting | |
| | Heater | Copper tube heater (nickel plated) | |
| | | 1.2KW | 2.5KW |
| | Stirrers | Magnetic pump | |
| | | 6W | 10W |
| Safety device | Shaker | Gear motor | |
| | | 25W | 40W |
| Specifications | Timer | 1 min~99 h 59 min~999 h 50 min | |
| | | Self-diagnostic function (heater disconnection, sensor disconnection, SSR short circuit, automatic overheat prevention function), overcurrent leakage circuit breaker, dry burning prevention, overheat preventer | |
| | Container rack spec. | Test tubes | Diameter 16.5mm, suitable for lengths 150~200mm (JIS) |
| | | Beakers | Round, triangular 100/300/500ml |
| | Internal dimensions | 130 pcs | 260 pcs |
| | | 12/5/3 pieces | 24/10/6 pieces |
| Accessories | Internal dimensions | | 302×350×250 |
| | External dimensions (W×D×H mm) | | 579×414×325 |
| | Internal bath capacity | | 19L (When full of water) |
| | Power supply (50/60Hz) | | AC100V 14A |
| | Weight | | Approx. 27kg |
| | Oscillating rack | | 1 sets |
| Accessories | Spring spacing | | 20mm (Same for vertical and horizontal) |

Control panel



Water bath cover (optional)



Oscillating stand

Container rack spec.

| | For BT101 Model | For BT311 Model |
|------------------------------|-----------------|-----------------|
| Test tubes (Diameter 16.5mm) | 130 pcs | 260 pcs |



Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Immersion Coolers | Immersion Type

BE201/201F/301

Made in
Japan


Operating temp. range -20~+35°C

Versatile immersion coolers.

When combined with a constant temperature water bath with temperature regulation function, it can serve as constant low temperature water bath below room temperature.

Features

- Use environmentally friendly refrigerants instead of freon.
- Just immerse the cooling coil in the liquid to quickly cool the existing water bath.
- Compact and lightweight, but powerful cooling ability, quickly cool liquids.
- The cooling coil of BE201F/301 is made of stainless steel flexible tubing, which not only has excellent corrosion resistance but can also be bent according to the size and shape of the water bath.



Specifications

| Model | | BE201 | BE201F | BE301 |
|-----------------------------|----------------------------------|-----------------------------------|------------------------------|------------------------------------|
| Performance | Operating temp. range | -20~35°C | | |
| | Cooling capacity | 190W at 0°C/Nybrine 60% 5L RT20°C | | 350W at 0°C/Nybrine 60% 10L RT20°C |
| Composition | Refrigerator | Reciprocating 160W | | Reciprocating 400W |
| | Refrigerant | R134a · 110g | R134a · 120g | R404a · 170g |
| | Cooling coil (diameter × length) | SUS304 spiral tube Φ35×170mm | SUS304 spiral tube Φ15×500mm | SUS304 spiral tube Φ15×1000mm |
| External dimensions (L×W×H) | | 410×320×303mm | | |
| Power capacity | | AC100V 3.0A | | AC100V 5.5A |
| Weight | | 24kg | | 25kg |

Cooling coil

BE201
Diameter 30×length 170mmBE201F
Diameter 15×length 500mmBE301
Diameter 15×length 1000mm

Constant Low Temp. Baths | For Material Synthesis

BLG100/200

Made in Japan

Temp. control range -80°C/-40°C

Water bath capacity 300ml/1000ml

Eliminates the trouble of monitoring and replenishing dry ice and liquid nitrogen. Since it is a unit with stirring function, just place a suitably sized beaker inside.



- After starting operation, it reaches the minimum temperature in about 1 h (within 30 min under no load).
- Equipped with a powerful neodymium magnet, capable of stirring multiple reactions simultaneously.
- Includes a timer function (convenient for next-day experiments).
- -80°C Model (BLG100), suitable for small experiments (~100ml container).
- -40°C Model (BLG200), suitable for experiments involving material synthesis (~500ml container).
- Liquid refrigerants (ethanol, fluorinated ethers) or aluminum beads can be used in the water bath.


Specifications

| Model | | BLG100 | BLG200 |
|------------------------|----------------------------------|--|--|
| System | | Stirling cooler 80W helium refrigerant (Freon-free) | |
| Performance | Temp. control range | -80°C~0°C (Surrounding temp. 23°C±5°C) | -40°C~0°C (Surrounding temp. 23°C±5°C) |
| | Temp. adjusting accuracy | ±0.3°C (Cooling bath temperature at stability) | |
| | Time fluctuation | ±0.3°C | |
| | Time to reach the lowest temp. | Approx. 55 min (20°C→-80°C) | Approx. 65 min (20°C→-40°C) |
| Composition | Water bath capacity | Aluminum/Approx. 300ml | Aluminum/Approx. 1000ml |
| | Temp. control method | PID control + frequency control | |
| | Temp. display and setting method | Digital | |
| | Maximum suitable container | 100ml three-necked beaker | 500ml three-necked beaker |
| | Stirring function | Neodymium magnet, 100-1200rpm, rotational number digital display | |
| | Timer | Max. 99 h 50 min | |
| External dimensions | | W210×D425×H295mm | W210×D450×H295mm |
| Weight | | Approx. 15kg | |
| Power supply (50/60Hz) | | AC100V 2A | |
| Accessories | | Insulation cover, stainless steel columns×2 | |

Control panel




Temperature decrease curve



Bath dimension comparison



Dimension diagram



Options

Crank clamp & bracket combination



| Product name | Product number |
|-----------------------------------|----------------|
| Crank clamp & bracket combination | 221774 |

Teflon stirring particles FO type



| Overall length | Diameter | Product number |
|----------------|----------|----------------|
| 10mm | 5mm | 221775 |
| 20mm | 10mm | 221776 |
| 30mm | 15mm | 221777 |

With side pipe pear-shaped beaker (side pipe size Φ8×30)



| Capacity | Main pipe | Product number |
|----------|-----------|----------------|
| 5ml | 15/25 | 221763 |
| 10ml | 15/25 | 221764 |
| 20ml | 15/25 | 221765 |
| 30ml | 15/25 | 221766 |
| 50ml | 15/25 | 221767 |

Three-necked round-bottom flask (side pipe 15/25)



| Capacity | Main pipe | Product number |
|----------|-----------|----------------|
| 50ml | 15/25 | 221768 |
| 100ml | 15/25 | 221769 |
| 100ml | 29/42 | 221770 |
| 200ml | 29/42 | 221771 |
| 300ml | 29/42 | 221772 |
| 500ml | 29/42 | 221773 |

※ BLG100 not applicable.

Constant Low Temp. Water Baths | Standard

BB311C/411C/611C

Temp. control range -30~+80°C

Temp. distribution accuracy $\pm 0.3^{\circ}\text{C}$

Internal capacity 6L 13L 26L

Constant low temperature water baths that can adjust from -30~+80°C with an accuracy of $\pm 0.1^{\circ}\text{C}$.

Features

- Capable of achieving temperature adjustment precision of $\pm 0.1^{\circ}\text{C}$.
- Extend the use of the water bath itself while ensuring effective space inside the bath.
- Equipped with external circulation function.
- Switching between internal and external circulation can be easily achieved with a simple operation of the water valve.
- Shelf height can be adjusted in 3 segments, allowing even small containers to be used.
- Primarily operates with automatic stop and automatic start functions, with RS485 communication function, temperature output terminal (4~20mA) as standard features.

Features

- Self-diagnostic function, leakage circuit breaker to prevent overcurrent, pressure detection for the refrigeration unit, anti-dry heating float switch, overload protection relay for the refrigeration unit, delayed function for refrigeration unit protection.



Specifications

| Model | | | BB311C | BB411C | BB611C |
|------------------------|--|--|--|---|--------------|
| Performance | Operating temp. range | | -30~80℃ | | |
| | Temp. adjusting accuracy | | ±0.1℃ | | |
| | Temp. distribution accuracy | | ±0.3℃ | | |
| | Temp. unit | | 0.1℃ | | |
| Cooling capacity | | Approx. 420W (361Kcal/h), at liquid temp. 15℃ | Approx. 510W (439Kcal/h), at liquid temp. 15℃ | Approx. 730W (628Kcal/h), at liquid temp. 15℃ | |
| Temp. control | | PID control | | | |
| Operation functions | | Fixed operation, automatic stop operation, automatic start operation | | | |
| Temp. sensor | | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple | | | |
| Temp. setting, display | | Digital setting, digital display | | | |
| Composition | Refrigerator, refrigerant | | Air cooling | | |
| | | | 300W, R404A | 350W, R404A | 600W, R404A |
| | Device circulation capacity (50/60Hz) | Max. flow rate | 2.8/3.2l/min (jet stirring) | | |
| | | Max. lift head | 1.1/1.4m | | |
| | Heater | | Stainless steel heating pipe | | |
| | | | 850W | 1.2KW | |
| | Cooling pipe | | Nickel-plated copper | | |
| | External circulation port | | Both the outlet and return ports have an outer diameter of 11mm | | |
| | Operating ambient temp. range | | 5~30℃ | | |
| | Safety device | | Self-diagnosis function, overload protection relay for the refrigeration unit, pressure detection for the refrigeration unit, anti-dry heating float switch, delay function for refrigeration unit protection, key lock function, automatic overheating protection function, leakage protection switches for over-current, overheating protector | | |
| Other functions | | Water valve, condenser filter, RS485 communication function, temperature output terminal (4~20mA), refrigeration unit pressure indicator, deviation correction function, power failure compensation function | | | |
| Specification | Water bath dimensions (W×D×H mm) | | 150×300×170 | 250×315×190 | 330×435×200 |
| | Effective water bath dimensions (W×D×H mm) | | 120×140×140 | 220×150×160 | 300×285×170 |
| | Internal bath capacity | | 6L | 13L | 26L |
| | External dimensions (W×D×H mm) | | 395×520×880 | 420×550×880 | 440×650×880 |
| | Power supply (50/60Hz) | | AC220V 7A | | AC220V 10A |
| | Weight | | Approx. 46kg | Approx. 53kg | Approx. 70kg |
| Accessories | | Top cover, drainage pipe, overflow pipe, shelf | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers


13 Washers

14 Analysis and Test Devices


15 Options

Heating/Cooling · Cooling Capacity Curve · Flow/Head Curve

Heating Cooling Curve




Cooling Capacity Curve



Cooling Capacity Curve

Flow/Head Curve



Control panel



Back side



- Due to the setting of a shut-off valve on the discharge side, external circulation can be switched quickly with a quick-connection. The nozzle diameter is $\Phi 10.5\text{mm}$.
- RS485 external communication
- Temperature output terminal

Usage cases



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Constant Low Temp. Water Baths | Benchtop Type

BBL111C/311C

Temp. control range -10~+80°C

Temp. distribution accuracy $\pm 0.3^\circ\text{C}$

Internal capacity 8L 13L

Benchtop constant low temperature water baths that ensure a wide temperature range and ample water bath space.

Features

- Achieves precise temperature control with an adjustment accuracy of $\pm 0.1^\circ\text{C}$.
- Expands the usability of the water bath itself and fully guarantees effective space within the bath.
- Equipped with an external circulation function.
- Switching between internal and external circulation can be easily achieved by operating the water valve.
- Shelf height can be adjusted in 2 segments, allowing for the use of small containers.
- Features automatic stop operation, automatic start operation, RS485 communication function, and temperature output terminal (4~20mA) as standard.

Safety

- Self-diagnosis function, leakage circuit breaker to prevent overcurrent, pressure detection for the refrigeration unit, anti-dry heating float switch, overload protection relay for the refrigeration unit, delay function for refrigeration unit protection.




Specification


| Model | | BBL111C | BBL311C | |
|-------------------------------|---|--|---|--|
| System | | Pump-circulated cooling/External closed-loop circulation | | |
| Performance | Operating temp. range | -10~80℃ | | |
| | Temp. adjusting accuracy | ±0.1℃ | | |
| | Temp. distribution accuracy | ±0.3℃ | | |
| | Temp. unit | 0.1℃ | | |
| | Cooling capacity | Approx. 300W (260Kcal/h), at liquid temp. 15℃ | Approx. 410W (350Kcal/h), at liquid temp. 15℃ | |
| Composition | Temp. control | | PID control | |
| | Operation functions | | | Fixed operation, automatic stop operation, automatic start operation |
| | Temp. sensor | | | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple |
| | Temp. setting, display | | | Digital setting, digital display |
| | Refrigerator, refrigerant | | | Air cooling |
| | | | 160W, R134A | 350W, R404A |
| | Stirring method | | | Magnetic pump |
| | Device circulation capacity (50/60Hz) | Max. flow rate | 3.7/4.7L/min | |
| | | Pump capacity | 11/12L/min | |
| | | Max. lift head | 1/1.5m | |
| | | Pump capacity | 1.5/2.1m | |
| | Heater | | 700W | 900W |
| | Cooling pipe | | Nickel-plated copper | |
| | External circulation port | | Both the outlet and return ports have an outer diameter of 11mm | |
| Operating ambient temp. range | | 5~30℃ | | |
| Safety device | | Self-diagnosis function, overload protection relay for the refrigeration unit, pressure detection for the refrigeration unit, anti-dry heating float switch, delay function for refrigeration unit protection, key lock function, automatic overheating protection function, leakage protection switch for over-current, overheating protector | | |
| Other functions | | Water valve, condenser filter, RS485 communication function, temperature output terminal (4~20mA), refrigeration unit pressure indicator, deviation correction function, power failure compensation function | | |
| Specifications | Water bath dimensions (W×D×H mm) | 300×150×177.5 | 300×240×177.5 | |
| | Effective water bath dimensions (W×D×H mm) | 238×100×85 | 238×190×85 | |
| | Internal bath capacity (effective capacity) | 8L (6.75L) | 13L (11.5L) | |
| | External dimensions (W×D×H mm) | 500×530×500 | 500×600×500 | |
| | Power supply (50/60Hz) | AC220V 5A | AC220V 6A | |
| | Weight | Approx. 35kg | Approx. 55kg | |
| Accessories | | Shelves, top cover, drainage pipe, overflow pipe | | |

Heating/cooling · cooling capacity curve

Heating Cooling Curve




Cooling Capacity Curve



Flow/head curve

Flow/Head Curve



Constant Low Temp. Water Baths | Large Capacity

BL410C/810C

Temp. control range -15~+70°C

Temp. distribution accuracy $\pm 0.3^\circ\text{C}$

Internal capacity 36L 80L

Large capacity constant low temperature water baths with observation window.



Features

- Wide operating temperature range: $-15\sim+70^\circ\text{C}$.
- Automatic stop operation and automatic start operation are both achievable. RS485 communication function and temperature output terminal (4~20mA) are also standard equipment.
- Set a large observation window.
BL410C: width 250×height 135 mm
BL810C: width 450×height 135 mm
- Shelf height can be adjusted in segments.

Safety


- Self-diagnosis function, leakage circuit breaker to prevent overcurrent, pressure detection for the refrigeration unit, anti-dry heating float switch, overload protection relay for the refrigeration unit, delay function for refrigeration unit protection.

Control panel




BL410C/BL810C

Heating/cooling · cooling capacity curve



Cooling Capacity Curve



Specifications

| Model | | BL410C | BL810C |
|-----------------|----------------------------------|--|--|
| System | | Pump circulated cooling | |
| Performance | Operating temp. range | $-15\sim+70^\circ\text{C}$ | |
| | Temp. adjusting accuracy | $\pm 0.1^\circ\text{C}$ | |
| | Temp. distribution accuracy | $\pm 0.3^\circ\text{C}$ | |
| | Temp. unit | 0.1°C | |
| | Cooling capacity | Approx. 540W (464Kcal/h), at liquid temp. 15°C | Approx. 920W (791Kcal/h), at liquid temp. 15°C |
| Composition | Temp. control | PID control | |
| | Operation functions | Fixed operation, automatic stop operation, automatic start operation | |
| | Temp. sensor | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple | |
| | Temp. setting, display | Digital setting, digital display | |
| | Refrigerator, refrigerant | Air cooling | |
| | | 350W, R404A | 600W, R404A |
| | Stirring method | Jet stirring is realized by magnetic pump | |
| | Heater | 1.3KW | 2KW |
| | Operating ambient temp. range | $5\sim30^\circ\text{C}$ | |
| Safety device | | Self-diagnosis function, overload protection relay for the refrigeration unit, pressure detection for the refrigeration unit, anti-dry heating float switch, delay function for refrigeration unit protection, key lock function, automatic overheating protection function, leakage protection switches for over-current, overheating protector | |
| Other functions | | Water valve, condenser filter, RS485 communication function, temperature output terminal (4~20mA), refrigeration unit pressure indicator, deviation correction function, power failure compensation function | |
| Specifications | Water bath dimensions | 400×300×300 | 600×400×350 |
| | Internal bath capacity | 36L | 80L |
| | Shelf layer quantity and spacing | 2 layers, adjustable 30mm spacing | |
| | External dimensions (W×D×H mm) | 680×390×805 | 880×490×855 |
| | Power supply (50/60Hz) | AC220V 8.5A | AC220V 14A |
| | Weight | Approx. 57kg | Approx. 85kg |
| Accessories | | Bottom shelf, variable mid-shelf | |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Constant Low Temp. Water Baths | Semiconductor Cooling

BV100S

Made in
Japan

Temp. control range 0~+80°C

Temp. distribution accuracy $\pm 0.1^\circ\text{C}$

Internal capacity 6L

Semiconductor cooling, constant low temperature water bath with low oscillation.

Rich in operability and safety features, the oscillating structure is integrated into the bath, providing a superior structure for oscillating constant temperature water baths.

Features

- The water bath is equipped with a magnetic pump with stirring function, achieving high temperature distribution accuracy for precise control.
- The water bath cover features a flat design, significantly improving the overall sealing after installation.
- Easy operation, available for fixed temp. operation, program operation, quick auto stop, auto stop and auto start operations.

Safety

- Safety functions include self-diagnosis function, leakage protection switch for over-current, anti-dry heating float switch, overheating protector, and key lock function.



Control panel



BN300 combination example



Specifications

| Model | | BV100S |
|---------------|---|---|
| System | | Semiconductor cooling + heater + magnetic pump stirring |
| Performance | Operating temp. range | 0~80°C (Ambient temperature 23°C) |
| | Temp. adjusting accuracy | $\pm 0.1^\circ\text{C}$ (Water temperature 20°C, Room temperature 20°C) |
| | Temp. distribution accuracy | $\pm 0.1^\circ\text{C}$ (Water temperature 20°C, Room temperature 20°C) |
| | Max. temp. reaching time | Approx. 90 min |
| | Time to reach the lowest temp. | Approx. 210 min |
| Composition | Temp. control | PID control |
| | Operation functions | Fixed temp. operation, auto start, auto stop, program operation |
| | Program function | Program operation 6 modes in total of 90 segments (30 segments×1, 15 segments×2, 10 segments×3) |
| | Temp. setting method | Use special function menu keys and up/down keys to realize digital setting |
| | Temp. display method | Measured temp. display: Green 4-digit LED digital display |
| | | Setting temp. display: Red 4-digit LED digital display |
| | Timer | 1 min~99 h 59 min and 100~999 h 50 min (including timer waiting function) |
| | Temp. sensor | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple |
| | Heater | 500W |
| | Cooler | Cooling plates ×4 |
| Safety device | Cooling fan | DC fan |
| | Stirrers | Magnetic pump |
| | | 3W |
| | Self-diagnosis function (temperature sensor anomaly, heater disconnection, automatic overheating prevention function, SSR path), leakage protection switches for over-current, anti-dry heating float switch, overheating protector, key lock | |
| | Water bath dimensions | 200×220×150 |
| | Internal bath capacity | 6L |
| | Drainage pipe diameter | Inner diameter 15mm×outer diameter 20mm |
| | External dimensions (W×D×H mm) | 340×538×415 |
| | Power supply (50/60Hz) | AC220V 5A |
| | Weight | Approx. 35kg |
| Accessories | | Drainage pump |

Heating Blocks | Test Tube Heating

HF100/200

Made in Japan

Temp. control range RT+50~+200°C

Detection End Pt100Ω

Rapid heating for test tubes.



No oil or water required, use aluminum blocks for rapid heating, concentration, decomposition, and reaction, achieving higher precision in temperature adjustment.

Features

- Feature ten different types of aluminum blocks (optional) to accommodate various experimental containers.
- Achieve temperature adjustment accuracy of $\pm 0.2^\circ\text{C}$ with deviation correction function for higher precision experiments.
- The controller features convenient automatic stop and automatic start timing functions.
- Operating temperature range from room temperature $+5^\circ\text{C}$ to 200°C allows for wide applications.
- Number of aluminum blocks: HF100 1 piece; HF200 2 pieces.

Control panel





Specifications

| Model | | HF100 | HF200 |
|-----------------|--------------------------------|--|-------------------------|
| Performance | Operating temp. range | RT+5~200°C | |
| | Temp. adjusting accuracy | $\pm 0.2^\circ\text{C}$ (at 200°C , using $\Phi 16.5$ test tube block) | |
| | Max. temp. reaching time | Approx. 30 min (RT23°C) | Approx. 35 min (RT23°C) |
| Composition | Number of aluminum blocks | 1 aluminum block | 2 aluminum blocks |
| | Temp. control | Digital controller PID control | |
| | Sensors | Pt100Ω | |
| | Temp. setting and display | Digital setting/display/minimum setting digit 0.1°C | |
| | Function | Fixed temp. operation, quick auto stop, auto stop, auto start | |
| | Heater | Mica heater 370w | Mica heater 600w |
| Safety device | | Over-current leakage breaker, automatic overheating prevention function (main relay disconnecting when set temperature $+12^\circ\text{C}$, manual reset), independent overheating protector (manual reset type, action temperature approximately 230°C) | |
| Other functions | | Key lock function, deviation correction function, power failure compensation function | |
| Specifications | Bath dimensions (W×D×H mm) | 112×112×70 | 222×112×70 |
| | External dimensions (W×D×H mm) | 230×310×139 | 340×310×139 |
| | Power supply (50/60Hz) | AC100V 3.8A | AC100V 6.1A |
| | Weight | Approx. 5kg | Approx. 6.5kg |
| Accessories | | Aluminum block handling handle | |

Aluminum block (optional)

| Product number | Model | Suitable test tubes | Quantity of placement | Depth for 1 piece |
|----------------|---------|---------------------|-----------------------|-------------------|
| ① 213173 | ① OBH80 | 0.5ml micro tube | 48 pcs | 29mm |
| ② 213174 | ② OBH82 | 1.5ml micro tube | 36 pcs | 37mm |
| ③ 213181 | ③ OBH84 | 12mm test tube | 36 pcs | 50mm |
| ④ 213182 | ④ OBH86 | 15mm test tube | 25 pcs | 50mm |
| ⑤ 213183 | ⑤ OBH88 | 16.5mm test tube | 25 pcs | 65mm |
| ⑥ 213184 | ⑥ OBH90 | 18mm test tube | 20 pcs | 65mm |
| ⑦ 213175 | ⑦ OBH92 | 21mm test tube | 12 pcs | 65mm |
| ⑧ 213176 | ⑧ OBH94 | 24mm test tube | 12 pcs | 65mm |
| ⑨ 213185 | ⑨ OBH96 | 25mm test tube | 12 pcs | 65mm |
| ⑩ 213177 | ⑩ OBH98 | 30mm test tube | 8 pcs | 65mm |

Temperature rise curve



(Note) ⑥ Thermometer installation example (Thermometer sold separately)
⑦ Aluminum block handle installation example



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Water Circulators

Contents

| | |
|--------------------------------|-----|
| Closed Water Circulators | 170 |
| Opened Water Circulators | 174 |

6 | Water Circulators



1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers


14 Analysis and Test Devices

15 Options

Constant temperature water circulators with a powerful product lineup

Widely used in the fields from biochemistry to semiconductors, with models for energy-saving, high power, low bed, benchtop, space-saving, with solvent recovery, and with external linkage, the product variety is extensive.

Function · Safety Device



Model selection list

| Category | Product name | Temp. control range | Max. flow rate (L/min) | Max. lift head | Series | Model | Features | Internal capacity (L) | Page No. |
|----------------------|---|---------------------|------------------------|----------------|--------|---------|--|-----------------------|----------|
| External closed loop | Cooling water circulator | -20°C~30°C | 6.3/7.2 | 4.0/5.6 | CF | CF313-B | Standard type | 4 | 170 |
| | | | 10.0/11.8 | 9.5/12.5 | CF | CF812-B | Standard type | 16 | 170 |
| | Cooling water circulator | -20~RT | 16.9/19.6 | 9.7/13.2 | CF | CF720C | Low bed mode | 16 | 171 |
| | Precision constant temp. water circulator | -10°C~+80°C | 8.9/10.3 | 6.6/9.0 | CFA | CFA311C | Precision type | 13 | 172 |
| External opened loop | Cooling water circulator | -10~RT | 5.4/6.2 | 3.5/5.0 | CLS | CLS312C | Standard type | 1.5 | 174 |
| | | | 5.4/6.3 | 3.7/5.3 | CLS | CLS411C | Standard type | 3 | 174 |
| | | | 6.7/7.8 | 6.2/8.7 | CLS | CLS610C | Standard type | 3 | 174 |
| | Precision constant temp. water circulator | -10°C~+80°C | 5.4/6.2 | 3.5/5.0 | CLH | CLH312C | Precision type | 1.5 | 176 |
| | | | 5.4/6.3 | 3.7/5.3 | CLH | CLH411C | Precision type | 3 | 176 |
| | | | 6.7/7.8 | 6.2/8.7 | CLH | CLH610C | Precision type | 3 | 176 |
| | Precision constant temp. water circulator | -10°C~+70°C | 8 | — | CTW | CTW412 | Semiconductor cooling, integrated type | — | 178 |
| | | | 11 | — | CTW | CTW812 | Semiconductor cooling, integrated type | — | 178 |
| | | 0°C~+70°C | 8 | — | CTA | CTA412 | Semiconductor cooling, integrated type | — | 178 |
| | | | 11 | — | CTA | CTA812 | Semiconductor cooling, integrated type | — | 178 |
| | Precision constant temp. water circulator | -10°C~+70°C | 8 | — | CTW-S | CTW412S | Semiconductor cooling, separated type | — | 179 |
| | | | 11 | — | CTW-S | CTW812S | Semiconductor cooling, separated type | — | 179 |
| | | 0°C~+70°C | 8 | — | CTA-S | CTA412S | Semiconductor cooling, separated type | — | 179 |
| | | | 11 | — | CTA-S | CTA812S | Semiconductor cooling, separated type | — | 179 |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Cooling Water Circulators | External Closed Loop, Standard Type

CF313-B·CF812-B

Operating temp. range -20°C~30°C

Cooling capacity CF313-B: Approx. 450W (at 10°C)

CF812C-B: Approx. 1230W (at 10°C)

Strong cooling capacity, water-saving external closed cooling water circulator.

Cooling water is steadily supplied to the closed cooling parts of research instruments, analytical and measuring instruments, industrial equipment, etc., through refrigeration unit and circulation pump.

Features

- Use environmentally friendly refrigerants R452A/R407C.
- The operating part is waterproof, featuring a simple and easy-to-understand digital setting display.
- CF812-B can connect to 2 sets of rotary evaporators (basic model).
- Equipped with safety designs such as overcurrent leakage protection switch, refrigeration unit overload relay, pump overheating protector, and time delay timer for refrigeration unit protection.

New Control Unit

Use high-brightness white LED to enhance visibility and operability. It can switch display between measuring temperature and set temperature.

Specifications

| Model | | CF313-B | CF812-B |
|----------------|--|--|--|
| System | | External closed circulation by circulation pump | |
| Performance | Operating temp. range | -20°C~30°C | |
| | Temp. adjusting accuracy | ±1.0°C (≥0°C), ±1.5°C (<0°C) | |
| | Cooling capacity | Approx. 450W at liquid temp. 10°C | Approx. 1230W at liquid temp. 10°C |
| | | Approx. 330W at liquid temp. -10°C | Approx. 560W at liquid temp. -10°C |
| Composition | Temp. control | Refrigeration Unit ON-OFF Control | |
| | Temp. sensor | PT100Ω | |
| | Temp. setting, display | White LED digital display, key input | |
| | Refrigerator, refrigerant | Air cooling | |
| | | 450W, R452A | 650W, R407C |
| | Circulation pump | Magnetic pump | |
| | | 15W | 156W |
| | Device Circulation Capacity | Max. flow rate | 9L/min |
| | | Pump capacity | 9L/min |
| | | Max. lift head | Approx. 4.3m |
| | | Pump capacity | Approx. 4.3m |
| Specifications | Cooling pipe | Stainless Steel SUS304 | Stainless Steel SUS304 |
| | External circulation port | Snap-on connector, O.D. 10mm rigid pipe connection (with attached connectors, can connect to I.D. 9mm soft hose) | |
| | Operating ambient temp. range | 5~35°C | |
| | Safety device | Overcurrent leakage protection switch, refrigeration unit overload relay, pump overheating protector, refrigeration unit protection time-delay timer, temperature sensor abnormality, high/low temperature limit warning, fan protection | |
| | Other functions | Circulation port condensation water drainage pipe, suction port dust filter, cooling operation switch, circulation pump switch, deviation correction, power outage reset mode selection function | |
| Specifications | Water bath dimensions | Length 198×width 142×height 195mm | Inner diameter 285×height 260mm |
| | Water bath material | Stainless steel | |
| | Water bath capacity | Approx. 3.9L (liquid volume 3.5L) | Approx. 15.5L (liquid volume 14L) |
| | Power supply | AC220V 3.8A | AC220V 5.5A |
| | External dimensions (W×D×H mm) (including protrusions) | 205×396×535 225×434×564 (Protruding Dimensions) | 340×370×838 340×408×920 (Protruding Dimensions) |
| | Weight | Approx. 30kg | Approx. 43kg |
| Accessories | | 1 condensation water drainage hose (1m), 2 hose clamps for fixing hoses, 2 hose nozzles (for soft hose connections), and 1 instruction manual. | |

Easy Connection of Circulation Pipes

Connection is completed simply by inserting rigid pipes (O.D.10mm). The connector is rotational, allowing free direction changes without burdening the hose. Also it's able to connect soft hoses with I.D. 9mm through the attached hose connector.

Condensation water drainage

A drain port designed for condensation is near the connectors on the top, connecting the attached drainage hose to drain condensate water.

Easy-to-clean filter

The filter mounting plate on the front is fixed with magnets, making it easy to remove for cleaning. The filter is sheet-like, making it easy to replace and clean.

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Cooling Water Circulators | External Closed Loop, Low Bed Type

CF720C

Operating temp. range -20°C~RT

Cooling capacity Approx. 880W (at 10°C)

Low bed type external closed cooling water circulator, can be installed under the fume hood or laboratory bench.


Cooling water is steadily supplied to the closed cooling parts of research instruments, analytical and measuring instruments, industrial equipment, etc., through refrigeration unit and circulation pump.

Equipped with two circulation systems, 880W cooling capacity, capable of simultaneously cooling 2 sets of rotary evaporators (basic model).


Features

- Use environmentally friendly refrigerant R404A.
- The operating part is digital setting and display, with a timing function.
- The circulation water flow adjustment valve, drainage valve, overflow port, and power switch are all installed on the front of the body, allowing all operations are in front of the unit body.
- Equipped with safety designs such as overcurrent leakage protection switch, refrigeration unit overload relay, pump overheating protector, and a time delay timer for refrigeration unit protection.


Cooling curve



Flow/head curve



System configuration example.



Specifications

| Model | | CF720C |
|-----------------|---------------------------------------|--|
| System | | External closed circulation by circulation pump |
| Performance | Operating temp. range | -20°C~room temperature |
| | Temp. adjusting accuracy | ±2°C |
| | Cooling capacity | Approx. 880W (750Kcal/h) at liquid temp. 10°C Approx. 680W (750Kcal/h) at liquid temp. 0°C Approx. 470W (750Kcal/h) at liquid temp. -10°C |
| Composition | Temp. control | Refrigeration Unit ON-OFF Control |
| | Temp. sensor | T-type thermocouple |
| | Temp. setting, display | Digital setting, digital display |
| | Refrigerator, refrigerant | Air-cooled sealed rotary type 600W R404A 680g |
| | Circulation pump | Magnetic pump 127W |
| | Pump capacity (50/60Hz) | Max. flow rate 30L/min |
| | | Max. lift head 10.3/14.3/m |
| | Device circulation capacity (50/60Hz) | Max. flow rate 16.9/19.6L/min |
| | | Max. lift head 9.7/13.2/m |
| | Cooling pipe | Nickel-plated copper |
| Safety device | External circulation port | Outer diameter 10.5mm pipe connector (2 each for inlet and outlet, with water valves) |
| | Operating ambient temp. range | 5~35°C |
| Safety device | | Overcurrent leakage protection switch, refrigeration unit overload relay, pump overheating protector, refrigeration unit protection time delay timer, bypass for pump protection |
| Other functions | | drainage, overflow, condenser filter, key lock function, deviation correction function, power outage compensation function, temperature output terminal |
| Specifications | Water bath material | Stainless steel |
| | Water bath capacity | Approx. 16L (liquid volume 14L) |
| | Power supply (50/60Hz) | AC220V 5A |
| | External dimensions (W×D×H mm) | 465×435×600 |
| | Weight | Approx. 60kg |
| Accessories | | Top cover, 4×1.5m circulation hoses, 4 clamps, 1 drainage hose, 1 overflow hose |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Precision Constant Temp. Water Circulators | External Closed Loop, High Precision Temp. Control

CFA311C/610C

Operating temp. range -10°C~+80°C

Temp. adjusting accuracy ±0.1°C

Precision type, external closed water circulator with a temperature adjusting accuracy of ±0.1°C.

Features

- Provide high-precision circulating water with an operating temperature range of -10~+80°C, and a temperature adjustment precision of ±0.1°C.
- As a cooling device, it has a strong cooling capacity.
- Monitors low water level, refrigeration unit pressure anomalies, refrigeration unit operation, and circulation pump operation, allowing real-time monitoring of the device's operating status.
- Equipped with various functions such as automatic stop operation, automatic start operation, temperature output terminal, deviation correction function, and RS485 external communication function.

Safety

- Feature various protection devices such as overcurrent leakage protection switch, self-diagnosis function, key lock function, refrigeration unit overload protector, refrigeration unit protection time delay timer, refrigeration unit pressure detection, float switch for preventing pump dry running, bypass for circulation pump protection, automatic overheating protection function, and overheating protector.



Specifications

| Model | | CFA311C | CFA610C |
|----------------|---------------------------------------|--|--|
| System | | External closed circulation by circulation pump | |
| Performance | Operating temp. range | -10~80°C | |
| | Temp. adjusting accuracy | ±0.1°C | |
| | Temp. unit | 0.1°C | |
| | Cooling capacity | Approx. 330W (284Kcal/h), at liquid temperature 10°C | Approx. 900W (770Kcal/h), at liquid temperature 10°C |
| Composition | Temp. control | PID control | |
| | Temp. sensor | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple | |
| | Temp. setting, display | Digital setting, digital display | |
| | Refrigerator, refrigerant | Air cooling | |
| | | 300W, R404A | 675W, R404A |
| | Circulation pump | Magnetic pump | |
| | | 45W | 65W |
| | Device circulation capacity (50/60Hz) | Maximum flow rate | 8.9/10.3L/min |
| | | Pump capacity | 15/17L/min |
| | | Maximum lift | 6.6/9.0m |
| Specifications | | Pump capacity | 8.0/11.0m |
| | Heater | Stainless steel heating tube 1.0KW | |
| | Cooling pipe | Nickel-plated copper | |
| | External circulation port | Rc3/8 (CFA301 standard equipment pagoda joint) | |
| | Operating ambient temp. range | 5~35°C | |
| | Safety device | Self-diagnosis function, overcurrent leakage protection switch, refrigeration unit overload protector, refrigeration unit pressure detection, refrigeration unit protection time delay timer, float switch for preventing pump dry running, bypass for circulation pump protection, automatic overheating protection function, overheating protector | |
| | Other functions | Drain valve, condenser filter, RS485 external communication function, refrigeration unit pressure indicator, key lock function, deviation correction function, temperature output terminal (4-20mA) | |
| | Water bath dimensions | Dimensions: width 250×depth 315×height 180mm | Inner diameter 300×height 236mm |
| | Water bath material | Stainless steel | |
| | Water bath capacity | Approx. 13L | Approx. 16L |
| Accessories | Power supply (50/60Hz) | AC220V 8A | AC220V 17A |
| | External dimensions (W×D×H mm) | 380×565×720 | 420×565×1050 |
| | Weight | Approx. 60kg | Approx. 77kg |
| Accessories | | Drainage pipe, overflow pipe | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Internal water bath (CFA610C)



Control panel



Water level monitoring, refrigeration unit monitoring, refrigeration unit action indicator

Back (CFA311C)




- RS485 external communication terminal
- Temperature output terminal
- circulation port

Back (CFA610C)




- RS485 external communication terminal
- Temperature output terminal
- circulation port


Flow/Head Curve




Heating Cooling Curve



Dimension diagram (mm)



Cooling Capacity Curve



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Cooling Water Circulators | External Opened Loop, Standard

CLS312C/411C/610C

Operating temp. range CLS312C: -10°C~RT CLS411C/610C: -15°C~RT

Temp. adjusting accuracy ±1.5~2°C

Strong cooling capacity, external opened cooling water circulator with ON-OFF control of the refrigeration unit.

Features

- External opened system cooling water circulator with temperature ranges of -10°C and -15°C is used.
- Equipped with a flow adjustment valve to freely adjust the circulation pump flow rate according to the size of water bath.
- Use a flow sensor to monitor pipeline flow, with an automatic stop protection device in case of pipeline blockage, indicated by a flashing red light.
- Include quick automatic stop, automatic stop, automatic start, deviation correction, power outage compensation, temperature output terminal, and other functions.

Safety

- Feature self-diagnosis function, overcurrent leakage protection switch, key lock, power outage compensation function, flow monitoring, refrigeration unit monitoring, refrigeration unit protection time delay timer, flow anomalies, high-pressure anomalies, and other safety protection devices.



Specifications

| Model | | | CLS312C | CLS411C | CLS610C |
|-------------------------------|---------------------------------------|----------------|---|---|---|
| System | | | External opened loop circulation by circulation pump | | |
| Performance | Operating temp. range | | -10℃~RT | -15℃~RT | |
| | Temp. adjusting accuracy | | ±1.5~±2℃ | | |
| | Cooling capacity | | Approx. 450W (387Kcal/h), at liquid temp. 15℃ | Approx. 570W (490Kcal/h), at liquid temp. 15℃ | Approx. 820W (705Kcal/h), at liquid temp. 15℃ |
| Composition | Temp. control | | Refrigeration Unit ON-OFF Control | | |
| | Temp. sensor | | T-type thermocouple | | |
| | Temp. setting, display | | Digital setting, digital display | | |
| | Refrigerator, refrigerant | | Air cooling | | |
| | | | 200W, R404A | 350W, R404A | 600W, R404A |
| | Circulation pump | | Magnetic pump | | |
| | Device circulation capacity (50/60Hz) | Max. flow rate | 5.4/6.2L/min | 5.4/6.3L/min | 6.7/7.8L/min |
| | | Pump capacity | 10.0/11.0L/min | 10.0/11.0L/min | 15.0/17.0L/min |
| | | Max. lift head | 3.5/5.0m | 3.7/5.3m | 6.2/8.7L/min |
| | | Pump capacity | 4.9/6.9m | | 8.0/11.0L/min |
| | Cooling pipe | | Nickel-plated copper | | |
| | External circulation port | | Inlet and outlet ports are outer diameter 13mm | | |
| Operating ambient temp. range | | 5~30℃ | | | |
| Safety device | | | Self-diagnosis function, overcurrent leakage protection switch, key lock, power outage compensation function, flow monitoring, refrigeration unit monitoring, refrigeration unit protection time delay timer, flow anomalies, high-pressure anomalies | | |
| Other functions | | | Flow adjustment valve, drainage valve, condenser filter, temperature output terminal, deviation correction function, key lock function, power outage compensation function | | |
| Specifications | Water bath dimensions | | Inner diameter 120× depth 200mm | Inner diameter 150× depth 200mm | |
| | Water bath material | | Stainless steel | | |
| | Water bath capacity | | 1.5L | 3L | |
| | Power supply (50/60Hz) | | AC220V 3A | AC220V 3.5A | AC220V 5.5A |
| | External dimensions (W×D×H mm) | | 410×460×550 | 380×513×725 | 380×615×725 |
| | Weight | | Approx. 40kg | Approx. 45kg | Approx. 60kg |
| Accessories | | | 1 thermal insulation pipe for circulation, 1 pump circulation pipe, 4 clamps, 1 drainage pipe | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Usage cases



The test bath is an optional product

Control panel







Water level monitoring,
refrigeration unit monitoring,
refrigeration unit action indicator

Back




- Due to the installation of a shut-off valve on the discharge side, external circulation can be switched using a quick connector.
- RS485 external communication
- Temperature output terminal




Cooling curve



Flow/head curve



Dimension diagram (mm)



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Precision Constant Temp. Water Circulators | External Opened Loop, High-Precision Temp. Control

CLH312C/411C/610C

Operating temp. range CLH312C: -10°C~80°C CLH411C/610C: -15°C~80°C

Temp. adjusting accuracy ±0.1°C

Precision type, external opened water circulator with a temperature adjusting accuracy of ±0.1°C.

Features

- Operating temperature range of -10~+80°C and -15~+80°C for external opened cooling water circulator.
- Equipped with a flow adjustment valve to freely adjust the circulation pump flow rate according to the size of water bath.
- Use a flow sensor to monitor pipeline flow, with an automatic stop protection device in case of pipeline blockage, indicated by a flashing red light.
- Include quick automatic stop, automatic stop, automatic start, deviation correction, power outage compensation, temperature output terminal, and other functions.

Safety

- Feature self-diagnosis function, overcurrent leakage protection switch, key lock, power failure compensation function, as well as flow monitoring, refrigerating machine monitoring, delayed timer for refrigerating machine protection, abnormal flow, abnormal high pressure, and other safety protection devices.



Specifications

| Model | | | CLH312C | CLH411C | CLH610C | |
|-------------------------------|---------------------------------------|--|---|---|---|--|
| System | | | External opened loop circulation by circulation pump | | | |
| Performance | Operating temp. range | | -10~+80℃ | -15~+80℃ | | |
| | Temp. adjusting accuracy | | ±0.1℃ | | | |
| | Cooling capacity | | Approx. 450W (387Kcal/h), at liquid temp. 15℃ | Approx. 570W (490Kcal/h), at liquid temp. 15℃ | Approx. 820W (705Kcal/h), at liquid temp. 15℃ | |
| Composition | Temp. control | | PID control | | | |
| | Temp. sensor | | Temperature regulator: Pt thermal resistor, for overheat protection: K-type thermocouple | | | |
| | Temp. setting, display | | Digital setting, digital display | | | |
| | Heater | | Stainless steel heating pipe | | | |
| | | | 900W | 900W | 1.5KW | |
| | Refrigerator, refrigerant | | Air cooling | | | |
| | | | 200W, R404A | 350W, R404A | 600W, R404A | |
| | Circulation pump | | Magnetic pump | | | |
| | Device circulation capacity (50/60Hz) | Max. flow rate | 5.4/6.2L/min | 5.4/6.3L/min | 6.7/7.8L/min | |
| | | Pump capacity | 10.0/11.0L/min | 10.0/11.0L/min | 15.0/17.0L/min | |
| | | Max. lift head | 3.5/5.0m | 3.7/5.3m | 6.2/8.7L/min | |
| | | Pump capacity | 4.9/6.9m | 8.0/11.0L/min | | |
| | Cooling pipe | | Nickel-plated copper | | | |
| External circulation port | | Inlet and outlet ports are outer diameter 13mm | | | | |
| Operating ambient temp. range | | 5~30℃ | | | | |
| Safety device | | | Self-diagnosis function, overcurrent leakage protection switch, key lock, power failure compensation function, as well as flow monitoring, refrigerating unit monitoring, delayed timer for refrigerating unit protection, abnormal flow, abnormal high pressure, overheating protector | | | |
| Other functions | | | Flow adjustment valve, drainage valve, condenser filter, temperature output terminal, deviation correction function, key lock function, power outage compensation function | | | |
| Specification | Water bath dimensions | | Inner diameter 120× depth 200mm | Inner diameter 150× depth 200mm | | |
| | Water bath material | | Stainless steel | | | |
| | Water bath capacity | | 1.5L | 3L | | |
| | Power supply (50/60Hz) | | AC220V 6A | AC220V 7.5A | AC220V 11A | |
| | External dimensions (W×D×H mm) | | 410×460×550 | 380×513×725 | 380×615×725 | |
| | Weight | | Approx. 40kg | Approx. 45kg | Approx. 60kg | |
| Accessories | | | 1 thermal insulation pipe for circulation, 1 pump circulation pipe, 4 clamps, 1 drainage pipe | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Usage cases



The test bath is an optional product

Control panel





Water level monitoring,
refrigeration unit monitoring,
refrigeration unit action indicator

Back




- Due to the installation of a shut-off valve on the discharge side, external circulation can be switched using a quick connector.
- RS485 external communication
- Temperature output terminal


Cooling curve




Flow/head curve




Dimension diagram (mm)



[CLH312C]



[CLH411C]



[CLH610C]

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Precision Constant Temp. Water Circulators | Semiconductor Cooling, Integrated Type

CTW412/812·CTA412/812

Made in Japan

Cooling method Water-cooled (CTW type) Air-cooled (CTA type)

Operating temp. range -10~+70°C (CTW type) 0~+70°C (CTA type)

Temp. adjusting accuracy ±1°C

Compact integrated precision constant temperature water circulator, can be placed on a desk or embedded in other devices.



The cooling and heating parts and power control parts are integrated, allowing high-precision control of constant temperature water near room temperature, using a Freon-free cooling method.

Features

- As a circulation type, it can be used as a constant temperature source for cooling, heating and heat preservation of hand-held bath.
- Standard equipped with RS485 communication function, enabling centralized management through a computer.
- There are water-cooled (CTW type) and air-cooled (CTA type) that does not require cooling water.
- Different types with different cooling capacities and circulation pump capacities can be selected based on the usage purpose.
- It has fixed value operation, program operation, automatic stop operation, quick automatic stop operation, automatic start operation, external sensor control, external communication function, alarm output, and other functions.

Specifications

| Model | | CTW412 | CTW812 | CTA412 | CTA812 | |
|----------------------------|---|--|---|---|----------------|------------------|
| System | | Water-cooled semiconductor cooling, external opened circulation | | Air-cooled semiconductor cooling, external opened circulation | | |
| Performance | Operating temp. range | | -10~+70℃ (CTW type) 0~+70℃ (CTA type) | | | |
| | Temp. adjusting accuracy | | ±0.1℃ | | | |
| | Cooling capacity | | 126Kcal/h (147W) | 250kcal/h (291W) | 83kcal/h (97W) | 163kcal/h (189W) |
| | Circulation pump | Maximum flow rate | 8L/min | 11L/min | 8L/min | 11L/min |
| Maximum discharge pressure | | 29.4KPa | 78.4KPa | 29.4KPa | 78.4KPa | |
| Composition | Temp. adjustment | | Microcomputer PID control, dialog-type digital setting display | | | |
| | Control sensor | | Thermoresistor (Pt100Ω) | | | |
| | Timing function | | 1 min~99 h 50 min, 100~9999 h (with hour•moment switching function) | | | |
| | Operation functions | | Fixed value operation (with pause function), fixed value automatic stop operation, quick automatic stop operation, fixed value automatic start operation, program operation (99 segments, with segment function, repeated, gradient operation, etc.), program automatic start operation | | | |
| | Other functions | | Self-diagnosis function, RS485 external communication function, alarm output terminal, deviation correction function, power failure compensation function, electricity consumption and electricity cost record display function, external temperature sensor switching, operating instructions function, etc. | | | |
| | Cooling and heating liquid contact material | | Stainless steel (SUS304) | | | |
| Safety device | | Self-diagnostic function, overcurrent leakage circuit breaker, key lock function | | | | |
| Specifications | External dimensions (W×D×H mm) | | 291×380×360 | 331×480×380 | 291×380×360 | 331×480×380 |
| | Circulation port pipe joint | | Outer diameter Φ12.7 | | | |
| | Number of temp. control modules | | 4 pcs | 8 pcs | 4 pcs | 8 pcs |
| | Power supply (50/60Hz) | | AC200V 1.7A | AC200V 3A | AC200V 1.7A | AC200V 3.3A |
| | Weight | | Approx. 19.5kg | Approx. 29.5kg | Approx. 24kg | Approx. 36kg |
| Accessories | | Insulated hose for circulating water I.D. Φ11.5×1m (2 pieces), insulated hose for discharging hot water I.D. Φ12×length 3m (1 piece, only for CTW) | | | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Precision Constant Temp. Water Circulators | Semiconductor Cooling, Split Type

CTW412S/812S·CTA412S/812S

Made in Japan

Cooling method Water-cooled (CTW type) Air-cooled (CTA type)

Operating temp. range -10~+70°C (CTW type) 0~+70°C (CTA type)

Temp. adjusting accuracy ±1°C

Easy to assemble, split-type precision constant temperature water circulator that can be placed on a desk or embedded in other devices.

The cooling and heating parts and power control parts are split-type, allowing high precision control of constant temperature water near room temperature, using a Freon-free cooling method.

Features

- As a circulation type, it can be used as a constant temperature source for cooling, heating and heat preservation of hand-held bath.
- Standard equipped with RS485 communication function, enabling centralized management through a computer.
- There are water-cooled (CTW type) and air-cooled (CTA type) that does not require cooling water.
- Different types with different cooling capacities and circulation pump capacities can be selected based on the usage purpose.
- It has fixed value operation, program operation, automatic stop operation, quick automatic stop operation, automatic start operation, external sensor control, external communication function, alarm output, and other functions.



Specifications

| Model | | CTW412S | CTW812S | CTA412S | CTA812S | |
|----------------------------|---|--|---|--|----------------|------------------|
| System | | Water-cooled semiconductor cooling, external open system circulation | | Air-cooled semiconductor cooling, external open system circulation | | |
| Performance | Operating temp. range | | -10~+70℃ (CTW type) 0~+70℃ (CTA type) | | | |
| | Temp. adjusting accuracy | | ±0.1℃ | | | |
| | Cooling capacity | | 126Kcal/h (147W) | 250kcal/h (291W) | 83kcal/h (97W) | 163kcal/h (189W) |
| | Circulation pump | Maximum flow rate | 8L/min | 11L/min | 8L/min | 11L/min |
| Maximum discharge pressure | | 29.4KPa | 78.4KPa | 29.4KPa | 78.4KPa | |
| Composition | Temp. adjustment | | Microcomputer PID control, dialog-type digital setting•display | | | |
| | Control sensor | | Thermoresistor (Pt100Ω) | | | |
| | Timing function | | 1 min~99 h 50 min, 100~9999 h (with hour•moment switching function) | | | |
| | Operation functions | | Fixed temp. operation (with pause function), fixed temp. automatic stop operation, quick automatic stop operation, fixed temp. automatic start operation, program operation (99 segments, with segmentation function, repeated, gradient operation, etc.), program automatic start operation | | | |
| | Other functions | | Self-diagnosis function, RS485 external communication function, alarm output terminal, deviation correction function, power failure compensation function, electricity consumption and electricity cost record display function, external temperature sensor switching, operating instructions function, etc. | | | |
| | Cooling and heating liquid contact material | | Stainless steel (SUS304) | | | |
| Safety device | | Self-diagnostic function, overcurrent leakage circuit breaker, key lock function | | | | |
| Specifications | External dimensions (W×D×H mm) | | 291×380×360 | 331×480×380 | 291×380×360 | 331×480×380 |
| | Circulation port pipe joint | | Outer diameter Φ12.7 | | | |
| | Number of temp. control modules | | 4 pcs | 8 pcs | 4 pcs | 8 pcs |
| | Power supply (50/60Hz) | | AC200V 1.7A | AC200V 3A | AC200V 1.7A | AC200V 3.3A |
| | Weight | Heat Exchange Department | Approx. 12.5kg | Approx. 20kg | Approx. 16kg | Approx. 30kg |
| Power Control Department | | Approx. 10kg | | | | |
| Accessories | | Insulated hose for circulating water I.D. Φ11.5×1m (2 pieces), insulated hose for discharging hot water I.D. Φ12×length 3m (1 piece, only for CTW) | | | | |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

10

Rotary Evaporators


Contents

| | |
|--------------------------|-----|
| Rotary Evaporators | 182 |
|--------------------------|-----|



The popular RE series rotary evaporators are integrated with a vacuum controller!

The integration of the rotary evaporator with the vacuum controller improves the overall operability and layout flexibility of the system. All operations can be controlled through the vacuum controller at the center of the main unit, eliminating the need to operate individual devices sequentially. This means that the vacuum pump and organic solvent recovery device can be set up in the space behind the rotary evaporator, increasing the layout flexibility of the operational area.



The VR102S-B vacuum controller is an optional accessory for the RE series, requiring a bracket and wiring for installation.


The REV series maintains the same space-saving design as before.

Automatic distillation, no complicated solvent data required!

The automatic mode of the vacuum controller allows operation to start simply by pressing a button after the sample is placed, without the need for previous solvent data. Distillation can be performed stably with both single and mixed solvents. This mode is particularly convenient for scenarios requiring multiple solvents or different distillation conditions each time.

Auto Mode
Automatic gradient descent •
Set target value for operation

For solvents with unknown pressure set values, the pressure change during pressure reduction is detected, and the pressure reduction gradient is controlled to prevent bumping. The starting point of gradient operation (pressure, time) is automatically calculated according to the set pressure value.





Compact design while maintaining convenience!

All functions and conveniences of the RE series are retained, with an all-around space-saving design in length, width, and height. When placed at an angle, 3 devices^{*1} can be stored in 1 fume hood for simultaneous use. Despite the space-saving design, the lift height of the lifter (200mm travel) is greater than that of other brands, ensuring convenient height for installation and flask replacement. Although the lifter operates manually, it moves very easily with the aid of a spring, making it more convenient than electric lifts when using small distillation flasks.

^{*1} When using vertical glass set



Wireless connection to the vacuum pump via Wi-Fi!

The RUN/STOP and pumping speed adjustment of the vacuum pump is wirelessly controlled by the vacuum controller. The vacuum pump operates at full power only during the initial stage of distillation, then reduces speed for quieter and more energy-efficient operation. The wireless connection eliminates the need for wiring, reducing the amount of wiring in the crossover system.

^{*} When using the specified vacuum pump N820G and vacuum control unit G

^{*} Wi-Fi IEEE802.11b/g/n frequency 2.4GHz



| | |
|------------------------------|-----------|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

The integration of the rotary evaporator and vacuum controller in the REV series

Creates new value through control functionality

Features

- The vacuum controller features a color LCD screen, allowing operation status to be identified by color. There are three operating modes for different applications to reduce the risk of bumping.
- The rotary evaporator can be operated in a one-touch manner with the vacuum pump. Wireless connection allows control of vacuum pump N820G or other specified vacuum pumps without additional wiring, enhancing layout flexibility.

REV Series



REV212MAW-B

REV212MBW-B



Rotary evaporator operation panel



Vacuum controller operation panel

Rotary evaporator Horizontal New release

REV212MA-B

[Main specifications] ◉ main unit + glass set A
◉ Dimensions: W719×D324×H534mm
(Excluding protruding parts) ◉ weight: 10.5kg

REV212MAW-B

[Main specifications] ◉ main unit + glass set A + water bath
◉ Dimensions: W744×D365×H534mm
(excluding protruding parts) ◉ Weight: 15kg

REV212MAO-B

[Main specifications] ◉ main unit + glass set A + oil bath
◉ Dimensions: W744×D365×H534mm
(excluding protruding parts) ◉ Weight: 15kg

REV212MB-B

[Main specifications] ◉ main unit + glass set B
◉ Dimensions: W529×D324×H745mm
(Excluding protruding parts) ◉ weight: 10.5kg

REV212MBW-B

[Main specifications] ◉ main unit + glass set B + water bath
◉ Dimensions: W554×D365×H745mm
(excluding protruding parts) ◉ Weight: 15kg

REV212MBO-B

[Main specifications] ◉ main unit + glass set B + oil bath
◉ Dimensions: W554×D365×H745mm
(excluding protruding parts) ◉ Weight: 15kg

Rotary evaporator Vertical New release

REV212ME-B

[Main specifications] ◉ main unit + glass set E
◉ Dimensions: W529×D324×H753mm
(Excluding protruding parts) ◉ weight: 10.5kg

REV212MEW-B

[Main specifications] ◉ main unit + glass set E + water bath
◉ Dimensions: W554×D365×H753mm
(excluding protruding parts) ◉ Weight: 15kg

REV212MEO-B

[Main specifications] ◉ main unit + glass set E + oil bath
◉ Dimensions: W554×D365×H753mm
(excluding protruding parts) ◉ Weight: 15kg

User-friendly controller

1 Operating mode

Displays the current operating mode

2 Status bar

Displays different colors corresponding to different operating status

3 Measured pressure

Displays current pressure

4 Set vacuum degree

Display the currently set vacuum degree

5 Running time

Displays the running time from start to present



Leak key

- Open the leak valve to restore normal pressure. (Convenient in emergencies, such as bumping)
- Long press for more than 3 seconds to manually activate the cleaning function.

Hold key

Temporarily maintains the pressure during operation

Power key


Operating knob

Basic settings can be completed using just the knob!

All settings can be set by rotating or pressing the knob, no instruction manual required, simple design that can be operated by feeling


Run/stop key

Three operating modes can be selected based on purpose:




-Manual mode

Continuously operate according to the set vacuum level and set hysteresis range.



-Gradient Mode

To prevent bumping, gradually reduce the pressure to the set value during pressure reduction. The starting point of gradient operation (pressure, time) is automatically calculated according to the set pressure value.



-Auto Mode

For solvents with unknown pressure set values, the pressure change during pressure reduction is detected, and the pressure reduction gradient is controlled to prevent bumping. The starting point of gradient operation (pressure, time) is automatically calculated according to the set pressure value.

The easy-to-use, intuitive RE series

Rotary evaporator can be equipped with vacuum control functionality!



RE Series

RE212AW-B

RE212BW-B

Plenty of convenience, widely praised!

- Easily adjust the height by operating the lift.
- The glass set and water/oil bath can be installed and used freely on the left and right.
- AC power adapter and wire are stored in the body for a clean and tidy appearance.
- Easy maintenance for water bath/oil bath.



Control panel

Rotary evaporator Horizontal

RE212A-B

[Main specifications] ○ body + glass set A
○ Dimensions: W719×D324×H534mm
(excluding protruding parts) ○ weight: 9kg

RE212AW-B

[Main specifications] ○ main unit + glass set A + water bath
○ Dimensions: W744×D365×H534mm
(excluding protruding parts) ○ Weight: 13.5kg

RE212AO-B

[Main specifications] ○ main unit + glass set A + oil bath
○ Dimensions: W744×D365×H534mm
(excluding protruding parts) ○ Weight: 13.5kg

Rotary evaporator Vertical

RE212B-B

[Main specifications] ○ body + glass set B
○ Dimensions: W529×D324×H745mm
(excluding protruding parts) ○ weight: 9kg

RE212BW-B

[Main specifications] ○ main unit + glass set B + water bath
○ Dimensions: W554×D365×H745mm
(excluding protruding parts) ○ Weight: 13.5kg

RE212BO-B

[Main specifications] ○ main unit + glass set B + oil bath
○ Dimensions: W554×D365×H745mm
(excluding protruding parts) ○ Weight: 13.5kg

RE212E-B

[Main specifications] ○ body + glass set E
○ Dimensions: W529×D324×H753mm
(excluding protruding parts) ○ weight: 10.5kg

RE212EW-B

[Main specifications] ○ main unit + glass set E + water bath
○ Dimensions: W554×D365×H753mm
(excluding protruding parts) ○ Weight: 13.5kg

RE212EO-B

[Main specifications] ○ main unit + glass set E + oil bath
○ Dimensions: W554×D365×H753mm
(excluding protruding parts) ○ Weight: 13.5kg

Intuitive operability, wide range of rotation speed settings.

Control the operation of the device through the central panel. Slowly rotate the knob to change the speed with each rotation; when rotating quickly, it changes with every 10 rotations. The speed can be set within a wide range of 5~315rpm. Once the speed is determined, press the Run/Stop key, and the distillation flask begins to rotate. The rotation mode can be divided into forward and reverse, with the option to select forward-reverse cycles or automatic reverse within the set time.


※ The current measured speed and set speed will be displayed in the speed display part. Speed settings can also be made during operation. The brightness of the LED can be adjusted in 8 levels.



Prevent liquid retention structure, improve the lifespan of the vacuum seal.

To prevent condensed liquid retained in the condenser from soaking into the vacuum seal part along the inner wall, a liquid retention prevention structure has been designed in front of the vacuum seal part inside the condenser. In addition, Glass set B considers the condensation inside the feed hose glass sleeve and has also designed a liquid retention prevention structure. By preventing the flow of liquid towards the vacuum part, the lifespan of the vacuum seal ring is improved.

※ When using ketones or ether solvents, the standard vacuum seal ring will corrode and swell, please use a PTFE vacuum seal ring (optional).



Prevention of backflow of condensate in the feed hose (Glass set B)

To prevent the condensate inside the condenser pipe from flowing back into the distillation flask through the feed hose, a feed hose glass sleeve is designed inside the condenser. The glass sleeve wraps the feed hose from the outside, successfully creating a complex structure within the glass pipe.



Balancing cooling capacity and space saving

Condensers A and B both ensure a condensation surface area of approximately 0.143m², guaranteeing optimal cooling capacity for condensation. The dimensions of Glass set B/E, including the water bath/oil bath, are more compact overall than previous products, allowing multiple units to be placed in a space-limited fume hood.

● RE212B-B picture



In the fume hood, set up 3 vertical glass sets.

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Model selection overview

| Composition | Model | REV212MA-B | REV212MAW-B | REV212MAO-B | REV212MB-B | REV212MBW-B | REV212MBO-B | REV212ME-B | REV212MEW-B | REV212MEO-B |
|--------------------|-----------|------------|-------------|-------------|------------|-------------|-------------|------------|-------------|-------------|
| Body | REV212M-B | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Glass set A | RG202A-B | ● | ● | ● | — | — | — | — | — | — |
| Glass set B | RG202B-B | — | — | — | ● | ● | ● | — | — | — |
| Glass set E | RG202E-B | — | — | — | — | — | — | ● | ● | ● |
| Water bath | BM312-B | — | ● | — | — | ● | — | — | ● | — |
| Oil bath | BO312-B | — | — | ● | — | — | ● | — | — | ● |

| Composition | Model | RE212A-B | RE212AW-B | RE212AO-B | RE212B-B | RE212BW-B | RE212BO-B | RE212E-B | RE212EW-B | RE212EO-B |
|--------------------|----------|----------|-----------|-----------|----------|-----------|-----------|----------|-----------|-----------|
| Body | RE212-B | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Glass set A | RG202A-B | ● | ● | ● | — | — | — | — | — | — |
| Glass set B | RG202B-B | — | — | — | ● | ● | ● | — | — | — |
| Glass set E | RG202E-B | — | — | — | — | — | — | ● | ● | ● |
| Water bath | BM312-B | — | ● | — | — | ● | — | — | ● | — |
| Oil bath | BO312-B | — | — | ● | — | — | ● | — | — | ● |

Body + glass set (A/B/E) specifications

| Series | | REV212M-B | | | RE212-B | | |
|-------------------|---|--|---|--|--|---|--|
| Vacuum controller | | ●*1 | | | - | | |
| Model | | REV212M-B + RG202A-B | REV212M-B + RG202B-B | REV212M-B + RG202E-B | RE212-B + RG202A-B | RE212-B + RG202B-B | RE212-B + RG202E-B |
| Composition | Condenser pipe | Dual serpentine tubes (Cooling surface area: approx. 0.143m ² , with insulation material) Vacuum sealing ring (FKM, surface PTFE treated) Condensate retention prevention | Longitudinal dual serpentine tubes (Cooling surface area: approx. 0.143m ² , with insulation material) Vacuum sealing ring (FKM, surface PTFE treated) Condensate retention prevention, condensate backflow prevention | Longitudinal dual serpentine tubes (Cooling surface area: approx. 0.149m ² , with insulation material) Vacuum sealing ring (FKM, surface PTFE treated) | Dual serpentine tubes (Cooling surface area: approx. 0.143m ² , with insulation material) Vacuum sealing ring (FKM, surface PTFE treated) Condensate retention prevention | Longitudinal dual serpentine tubes (Cooling surface area: approx. 0.143m ² , with insulation material) Vacuum sealing ring (FKM, surface PTFE treated) Condensate retention prevention, condensate backflow prevention | Longitudinal dual serpentine tubes (Cooling surface area: approx. 0.149m ² , with insulation material) Vacuum sealing ring (FKM, surface PTFE treated) |
| | Rotary joint | TS29/38, inner diameter Φ18×L284mm | TS29/38, inner diameter Φ18×L208mm | | TS29/38, inner diameter 18×L284mm | TS29/38, inner diameter Φ18×L208mm | |
| | Flask | Distillation flask (TS29/38, 1L), flask clip (POM), recovery flask (S35/20, 1L), flask clip (SUS) | | | | | |
| | Motor | DC servo motor | | | | | |
| | Interface size | Outer diameter Φ10mm (cooling water · vacuum) | | | | | |
| Performance | Rotational speed range | 5~315rpm White LED digital setting and display | | | | | |
| | Lifter stroke | Manual balance type (~200mm stepless adjustment, 90° rotation lock/release) | | | | | |
| | Distillation flask evaporation capacity | ~23mL/min (water evaporation rate) | | | | | |
| | Vacuum degree | Below 1000Pa | | | | | |
| Safety functions | | Motor overload stop function, overcurrent protection, overvoltage protection | | | | | |
| Other functions | | Forward/reverse/timer automatic forward-reverse switch, power outage recovery (mode selection) | | | | | |
| Specifications | Distillation flask | 50~2000mL | | | | | |
| | Recovery flask | 100~2000mL | | | | | |
| | Power supply | AC100~230V | | | | | |
| | Current | 1A | | | | | |
| | External dimensions | W719× D324× H534mm (Lifter highest position: H734mm) | W529× D324× H745mm (Lifter highest position: H945mm) | W529× D324× H753mm (Lifter highest position: H953mm) | W719× D324× H534mm (Lifter highest position: H734mm) | W529× D324× H745mm (Lifter highest position: H945mm) | W529× D324× H753mm (Lifter highest position: H953mm) |
| | Weight | Approx. 10.5kg (main body: approx. 8.5kg + glass set: approx. 2kg) | | | Approx. 9kg (main body: approx. 7kg + glass set: approx. 2kg) | | |
| Accessories | Body | Manual, warranty card, AC power adapter, power cord (approx. 2m), guide groove, rear protective plate, hook and loop fastener (single-sided, double-sided), pipe holder | | | | | |
| | Glass set | Warranty card, distillation flask, distillation flask clip, recovery flask, recovery flask clip, vacuum sealing ring, flask propeller (for ST29), cooling pipe insulation material manual, sample introduction tap, connecting rotary shaft, vacuum nozzle (gray), cooling nozzle (black), cooling pipe nut, ring spring, VR102S quick guide | | | Warranty card, distillation flask, distillation flask clip, recovery flask, recovery flask clip, vacuum sealing ring, flask propeller (for ST29), cooling pipe insulation material manual, sample introduction tap, connecting rotary shaft, vacuum nozzle (gray), cooling nozzle (black), cooling pipe nut, ring spring | | |

*1 For vacuum controller specifications, refer to VR102S controller specifications

Water bath/oil bath specifications

| Components | | Water bath | Oil bath |
|------------------|--|--|--------------------------|
| Model | | BM312-B | BO312-B |
| Performance | Temp. control range | Room temp. +10~90℃ | Room temp. +10~180℃ |
| | Temp. adjusting accuracy | ±1℃ | ±1.5℃ (water), ±2℃ (oil) |
| Composition | Exterior material | PBT (with glass fiber) | |
| | Internal bath material | Stainless steel | |
| Safety functions | | Automatic overheating prevention, abnormal temperature upper limit, independent overheating prevention, overcurrent fuse, etc. | |
| Other functions | | Deviation correction, overshoot notification, power outage compensation function (mode selection), 2A service socket | |
| Specifications | Power supply | AC200~230V 5~6A | |
| | Inner dimensions and liquid volume of water bath | Φ240×H119mm about 5L | |
| | External dimensions | Φ262 (maximum depth 286mm)×H240mm | |
| | Weight | Approx. 4.5kg | |
| Accessories | | Manual, warranty card, power cord (approx. 3m), fuse tube | |

Vacuum controller specifications

| Model | | VR102S |
|------------------|---|--|
| Performance | Vacuum setting range | 0~1013hPa |
| | Vacuum measurement range | 0~1100hPa |
| Controllers | Display | Color LCD (2.3 inches) |
| | Display items | Measured vacuum degree, set vacuum degree, operating time, various status displays |
| | Function | Operating modes |
| | | Holding function |
| | | Leak |
| | | Pressure unit switching |
| | | Procedures after operation |
| | | Wireless specifications |
| | | Other functions |
| Components | Pressure sensor | Semiconductor pressure sensor (SUS316L diaphragm) |
| | Leak solenoid valve | For atmospheric introduction |
| | Leak port | Outer diameter Φ10mm interface PP |
| | Vacuum connection port | Hard pipe connection (outer diameter Φ4mm) (included) |
| | Vacuum control solenoid valve connection port | Mini DIN (3pin) |
| | RE212 interface | Mini DIN (8pin) |
| | Power connection port | DC jack |
| Safety functions | | Communication error, pressure sensor error, memory error, Leak error, high-pressure error, Leak function during abnormal stop |
| Specifications | External dimensions | 86mm×113mm×83mm (excluding protrusions) |
| | Power ^{*1} | DC24V (AC100~240V less than 1A) |
| | Weight | 0.5kg |
| Accessories | | Vacuum tube branch joint (vacuum pipeline side: outer diameter Φ10 port, vacuum controller connection side: outer diameter Φ4mm×Φ2mm × 700 mm, with PTFE hard pipe), manual, warranty card, VR102S quick guide |

^{*1} When connected to RE212, power is taken from the optional connection cable. Otherwise, a separately purchased AC adapter and power cord are required for use.

Vacuum pump specifications

| Model | N820G |
|--|---|
| Exhaust speed ^{*1} | 9~20L/min (with manual adjustment knob) |
| Vacuum pump material | PTFE/FFPM |
| External input | Remote operation via external signal ^{*2} |
| Achievable vacuum degree | Lowest speed 600Pa (air ballast valve closed) Highest speed 800Pa (air ballast valve closed) |
| Vacuum hose connection port outer diameter | Φ10.6/Φ9 (2-stage) |
| Power supply | AC100~240V 50/60Hz 1A |
| External dimensions | W163×D259 × H220mm (excluding protrusions) |
| Allowable ambient temp. | +5~+40℃ |
| Allowable solvent temp. | +5~+40℃ |
| Safety functions | Overcurrent protection, temperature protection (motor), lock protection (motor) |
| Weight | Approx. 8.8kg |

^{*1} Measurement conditions: atmospheric pressure, temperature 20℃

^{*2} Separate purchase of VR102S vacuum controller and control unit G required

| Model | N816 |
|-----------------------------|--|
| Exhaust speed | 16L/min |
| Vacuum pump material | PTFE/FFPM |
| External input | None |
| Achievable vacuum degree | 2000 Pa |
| Vacuum hose connection port | ID 6 |
| Power supply | AC100~240V 50/60Hz 1A |
| External dimensions | W90×D361×H141mm |
| Allowable ambient temp. | +5~+40℃ |
| Allowable solvent temp. | +5~+40℃ |
| Safety functions | Overheat protection switch, power fuse |
| Weight | Approx. 3.95kg |

Cooling water circulation device specifications

| Model | | CF313-B | CF812-B |
|-----------------|--------------------------------------|--|--|
| System | | External sealed loop circulation by circulation pump | |
| Performance | Operating temp. range ^{*1} | -20~30℃ | |
| | Temp. control accuracy ^{*2} | ±1.0℃(≧0℃) ±1.5℃(<0℃) | ±1.0℃ |
| | Cooling capacity ^{*2} | Approx. 450W at liquid temperature 10℃ Approx. 330W at liquid temperature -10℃ | Approx. 1230W at liquid temperature 10℃ Approx. 560W at liquid temperature -10℃ |
| | Maximum flow rate ^{*3} | Approx. 9L/min | Approx. 12L/min |
| | Maximum lift ^{*3} | Approx. 4.3m | Approx. 10.3m |
| Composition | Temp. control | Refrigerator ON-OFF control | |
| | Temp. sensor | Pt100Ω | |
| | Refrigerator · refrigerant | Air-cooled 450W R452A | Air-cooled 650W R407C |
| | Circulation pump | 15W | 156W |
| Safety device | | Overcurrent leakage circuit breaker, temperature sensor error, temperature upper/lower limit warning, temperature upper/lower limit abnormality, refrigerator high-pressure switch, fan motor protection, circulation pump protection, refrigerator protection delay timer | |
| Other functions | | Suction port dust filter, cooling operation switch, circulation pump switch, deviation correction, power outage reset mode selection function | |
| | Water bath material | Stainless steel | |
| | Water bath capacity | Approx. 3.9L (liquid volume 3.5L) | Approx. 16L (liquid volume 14L) |
| | External dimensions | W205×D396×H535mm (excluding protrusions) | W340×D370×H838mm (excluding protrusions) |
| | Power supply | AC220V 3.8A | AC220V 5.5A |
| | Weight | Approx. 30kg | Approx. 43kg |
| Accessories | | Soft hose nozzles (for soft hose connection) 2 pcs, 1m condensate drain hose (inner diameter Φ12mm) 1 pc, condensation water hose clamp 2 pcs, manual 1 pc | |

^{*1} No heating function. Depending on the ambient temperature and connection conditions, it may not reach -20℃

^{*2} Performance values at ambient temperature 20℃, power supply voltage AC220V. Temperature control accuracy is the specification value calculated using JTMK05, temperature variation is calculated using JIS

^{*3} Capacity value at tap water, liquid temperature 20℃

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |



Water bath (oil bath)

Removing and cleaning the water bath, changing the water (oil) is very convenient. Digital temperature setting and display make operation very simple. According to the usage temperature range, you can choose either the BM type water bath or the BO type oil bath.



Space-saving

It can also be stored in a fume hood, and the depth of the vertical condenser (glass set B/C) is only 400mm, saving width space.



Realizes wide voltage

The main body can handle 100~240V, and it comes with power cords for different countries. The water bath and oil bath can each handle 100~120V and 200~240V.

Vacuum sealing ring


Made of fluororubber. The sealing ring part adopts a multi-segment shape, which can maintain a very stable decompression state. As an option, we also have a full teflon version. When using ketone and ethers, please use the full teflon version.

Thoughtful safety structure

It can detect abnormalities in the drive part (rotation part, lifting part) and stop safely. In addition, the water bath and oil bath temperature abnormalities can also be automatically detected and stopped.

Stable rotation from low to high speed

Adopt feedback control, combined with digital setting and display, to achieve stable rotation accuracy and improve the reproducibility of experiments.



Standard equipped with electric lifting

All models are standard equipped with an electric lifting device. The device can be moved up and down by a switch.



Digital setting and display

Rotation speed, vacuum degree, and steam temperature are all digitally set and displayed. Experimental conditions can be repeated countless times. In addition to Japanese, the display can also be switched to English. (The digital display of the RE311 model is limited to the rotation speed, vacuum regulator, and steam temperature display can be optionally installed)



Movable rotary joint

The fixed position of the rotating connection part can be changed by 80mm, which is very convenient when using distillation flasks of different capacities or changing the contact position of the vacuum sealing ring. (Japan Patent No.: 3220033)



Distillation flask extraction structure

Distillation flask can be easily removed.



Glass part

The vertical condenser (glass set B/C) is designed to prevent liquid retention. To effectively exhibit the capabilities of the B/C type condenser and improve condensation efficiency, the condensation area has been expanded by approximately 20% (compared to our previous products).



Hose joint

For easy and simple connection of cooling water and vacuum hoses to the glass parts, resin-made hose joints are used.

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Rotary Evaporators | Automatic Lifting Type

RE311

Distillation/recovery flask 1L (standard)

Rotational speed control range 20~250rpm

Water bath/oil bath temp. range RT+5~100°C/RT+5~180°C

The basic rotary evaporator equipped with an electric lifter to improve operability.



RE311A-W

Features

- Equipped with an electric lifter, the device can be moved up and down by switch control.
- The rotational speed can be digitally set and displayed.
- By adopting a feedback control motor, stable rotation from low to high speed can be achieved.
- This product uses wide voltage; please connect it to a power supply of 100~240V.
- Space-saving design, easily stored in a fume hood.
- Equipped with an electric lifter rise function in the event of a power outage.



RE311B-W



RE311C-W

Specifications

| Model | | | RE311 |
|----------------------|-------------------------|-------------------------------|--|
| Composition | Motor | | DC brushless motor (for rotation) DC motor (for lifter) |
| | | | |
| Performance-Function | Rotation speed | | 20~250r/min (rpm) |
| | Electric lifter stroke | | 130mm |
| | Rotation speed setting | | Potentiometer type (Digital display) |
| | Safety functions (Body) | | Manual lower limit range setting, motor overload |
| | Lifter function | | Electric lifter |
| Specifications | External dimensions | When installing glass set A | Width 828×Depth 400×Height 586 (716: when lifter is raised) mm |
| | | When installing glass set B/C | Width 643×Depth 400×Height 727 (857: when lifter is raised) mm |
| | Weight | | Approx. 14.4kg (excluding glass parts, water/oil bath) |
| | Power supply | Body (50/60Hz) rated current | AC100~240V 1.5A (excluding water/oil bath) (3A) |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Movable rotary joint



Due to the variable fixed position of the rotating connection part can be changed into 80mm, it is very convenient when using distillation flasks of different capacities or when changing the contact position of the vacuum seal ring.

Distillation flask extraction structure



Simply by turning the distillation flask extraction nut, the distillation flask can be easily removed.

Hose joint



For easy and simple connection of cooling water and vacuum hoses to the glass parts, resin-made hose joints are used.

Options



Vacuum regulator VR300




Steam temperature indicator TA300



Hose connection fittings



Water bath BM510 (oil bath BO510C)



Clamp KS35B

Clamp KC29

Control panel



Glass parts (anti-liquid accumulation structure)



Anti-liquid accumulation design of the vertical condenser (glass set B/C). To effectively exhibit the capabilities of the B/C type condenser and improve condensation efficiency, the condensation area has been expanded by approximately 20% (compared to our previous products).

| Product number | Product name | Model |
|----------------|--|-------|
| 222121 | Water bath AC200~240V | BM510 |
| 222115 | Oil bath AC200~240V | BO410 |
| 255277 | Vacuum regulator (for RE301) | VR300 |
| 255368 | Teflon vacuum seal ring | ORE70 |
| 255341 | Teflon solenoid valve | ORE80 |
| 255342 | SUS pressure sensor | ORE90 |
| 255317 | Steam temp. display | TA300 |
| 255285 | Water trap for RE311/601/801 (with nozzle) | ORE40 |
| 255284 | Hose connection fittings for RE311/601/801 | ORE30 |
| 255296 | Cooling hose Φ9mm, 2m×1 | |
| 255297 | Vacuum hose Φ6mm, 5m×1 | |
| RE20000100 | Distillation flask 100ml, frosted | ORG20 |
| RE20000200 | Distillation flask 200ml, frosted | ORG22 |
| RE20000300 | Distillation flask 300ml, frosted | ORG24 |
| RE20000500 | Distillation flask 500ml, frosted | ORG26 |
| RG00A30040 | Distillation flask 1l, frosted | ORG16 |
| RE20002000 | Distillation flask 2l, frosted | ORG28 |
| RE47002 | Recovery flask 300ml, frosted | ORG34 |
| RE47001 | Recovery flask 500ml, frosted | ORG36 |
| 2551730413 | Recovery flask 1l, frosted | ORG18 |
| RE47003 | Recovery flask 2l, frosted | ORG38 |
| 255505 | Surface-coated distillation flask 1L | ORG58 |
| 255511 | Surface-coated recovery flask 1L | ORG56 |
| F0410005 | Clamp | KS35B |
| F0410001 | Clamp | KC29 |
| 255281 | Solvent recovery device | RT200 |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Freeze Dryers & Cold Traps

Contents

Freeze Dryers194

Cooling Traps.....196

Freeze Dryers | Standard

DC401/801

Cooling temp. -45°C (401 type)/-85°C (801 type)

Internal capacity 4L

High-operability standard freeze dryer.



- Safety structure automatically discharges when power is OFF or in case of a power outage, preventing vacuum oil backflow.
- Using a hot gas bypass method, it can smoothly defrost and remove ice inside the tank.



- Equipped with a Pirani vacuum gauge for higher precision freeze-drying.
- Equipped with safety valves and power sockets for vacuum pump linkage.
- Use environmentally friendly refrigerants that do not destroy the ozone layer.

Specifications

| Model | | DC401 | DC801 |
|----------------|--|---|---|
| Performance | Cold trap temp. | -45°C | -85°C |
| | Dehumidification amount | 0.6L | 1.0L |
| Composition | Temp. sensor | — | Pt100Ω |
| | Temp. display | — | 7-segment LCD (display temperature is the surface temperature of the cooling pipes) |
| | Refrigerator | Air-cooled 400W | Air-cooled 350W |
| | Refrigerant | R404A | R600a other, mixed refrigerant |
| | Linking gauge | — | High pressure/low pressure indicator |
| | Water bath shape/material | Cylindrical shape, stainless steel plate | |
| | Drain pipe | Vacuum hose with plugged | |
| | Vacuum gauge | Pirani vacuum gauge | |
| | Defrosting structure | Hot gas bypass method | |
| | Exhaust vent (vacuum pump connection port) | Φ17 | |
| Specifications | Power socket | Vacuum pump Max.6A/heating Max.2A | |
| | Usage environment temp. range | 5~30°C | |
| | Function | Cold trap defrosting: hot gas defrost | |
| | Safety device | Overcurrent leakage circuit breaker, overload protection for refrigeration, backflow prevention valve | |
| | Cold trap dimensions | Φ153×height 235mm | |
| Specifications | External dimensions | Width 300×Depth 450×Height 720mm | Width 405×Depth 500×Height 846mm |
| | Internal capacity | Approx. 4L | |
| | Power supply (50/60Hz) rated current | AC100V 6A | AC100V 7A |
| | Applicable vacuum pump exhaust amount | 50/60L/min | |
| Accessories | Weight | Approx. 43kg | Approx. 65kg |
| | | Vacuum silicone lubricant, vacuum hose, fuse 10A | |
| | | Vacuum silicone lubricant, vacuum hose, fuse 10A, 4 casters fixed brackets | |

Beaker chamber, multi-branch pipe, beaker lid, and glassware used in DC401/801 are optional items.

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers


13 Washers

14 Analysis and Test Devices

15 Options

Made in Japan


Cooling curve




Control panel




Sublimation data



When placing 1 pre-frozen item with 0.15L of water into a 0.5L sample beaker (pear-shaped), then installing it onto the cold trap, the sublimation amount is calculated from the residual amount in the beaker after a certain period.



Options



| Product No. | 212560 | 212561 | 212562 | 212563 | 212564 |
|----------------|-----------------|---------------------|-------------------------------|-------------------------------------|-----------------------|
| Product name | Beaker chamber | Multi-branch pipe A | Multi-branch pipe B | Drying chamber | Closed drying chamber |
| Specifications | Shelf sections | | | Shelf layers ×1 | |
| | Mounting rack | | | Inner diameter 18.5mm | |
| | Rack spacing | | | 60mm Petri dish hanging quantity ×7 | |
| Port quantity | 12 | 8 | 16 | Temperature regulation 30°C±2°C | |
| Dimensions | Φ195×height 303 | | Width 304×depth 60×height 263 | Φ252×height 240 | Φ252×height 425 |
| | | | Width 624×depth 60×height 263 | | |

| Product name | Specifications | Product number |
|--|----------------|----------------|
| Valve | | 212565 |
| Dry beaker | 120ml×5 pcs | 212820 |
| | 250ml×5 pcs | 212821 |
| Dry beaker lid (with glass joint) | Five straight | 212570 |
| | Five 45° bend | 212571 |
| Serum bottles (with rubber stopper) | 50ml×10 pcs | 212814 |
| | 30ml×10 pcs | 212815 |
| | 10ml×10 pcs | 212816 |
| Ampoule connector | Single 5 pcs | 212572 |
| | Double 5 pcs | 212573 |
| | Triple 5 pcs | 212574 |
| Test tube connector (with glass interface) | Straight | 212590 |
| | 45° Bend | 212591 |
| Eggplant shaped flask (male) | 200ml | 212594 |
| | 300ml | 212595 |
| | 500ml | 212596 |
| Eggplant shaped flask (male) connector | 29 | 212597 |
| Eggplant shaped flask (female) | 200ml | 212566 |
| | 300ml | 212567 |
| | 500ml | 212568 |
| Eggplant shaped flask (female) connector | 29 | 212569 |
| Microtube rack | 1.5ml×24 pcs | 212580 |
| Glass joint | Straight | 212598 |
| | 45° Bend | 212599 |

Cooling Traps | Chamber Cooling Type

CA301/801

Made in
Japan

Minimum temp. achieved -45°C (Model 301); -85°C (Model 801)

Maximum dehumidification capacity 0.9kg (water) (Model 301); 1.0kg (water) (Model 801)

With two models achieving minimum temperatures of -45°C and -85°C, it is space-saving and comes with casters.



Effectively collects water vapor and harmful vapors discharged from vacuum drying ovens, rotary evaporators, etc., in the vacuum system. Space-saving cold trap compatible with both acid and organic solvents.

Features

- Built-in refrigeration unit with easy movement on casters.
- Significantly reduces the amount of steam inhaled by the vacuum pump.
- The pre-freezing tank can be used as a low-temperature tank.
- Use environmentally friendly refrigerants.
- Specially equipped with glass condensers for acids and organic solvents.

Specifications

| Model | CA301 | CA801 |
|-------------------------------|--|---|
| System | Direct cold trap or glass cold trap | |
| Performance | Dehumidification capacity Maximum 0.9kg (water-based) | Maximum 1.0kg (water-based) |
| Minimum temp. achieved | -45°C | -85°C |
| Refrigerator | Air-cooled 400W | Air-cooled 350W |
| Refrigerant | R404A | R600a other, mixed refrigerant |
| Lid | SUS304 with nozzle (outer diameter 17.6mm) | |
| Water chamber/material | Cylindrical-SUS304 | |
| Usage environment temp. range | 5~35°C | 5~30°C |
| Temp. display | 7-segment LCD ※ Display temperature is the surface temperature of the cooling pipe | |
| Temp. sensor | Pt100Ω | |
| Defrosting structure | None | Hot gas bypass method |
| Safety device | Overcurrent leakage circuit breaker, refrigeration unit overload protection | Overcurrent leakage circuit breaker, refrigeration unit delay timer, refrigeration unit overload protection, refrigeration unit high-pressure abnormality, sensor disconnection abnormality |
| Other functions | Drain valve (with plug vacuum hose) | Drain valve (with plug vacuum hose), combination meter |
| Tank size | Inner diameter 153×height 235mm | |
| Internal capacity | Approx. 4L (liquid volume 3L) | |
| Power supply | AC100V 5A | AC100V 7A |
| External dimensions (mm) | Width 345×depth 475×height 726 | Width 405×depth 500×height 850 |
| Weight | Approx. 47kg | Approx. 65kg |
| Accessories | — | Caster fixing brackets 4 pcs |

Options

| Product name | Specifications | Product number |
|---------------------|----------------|----------------|
| Condenser | OCA10 | 221487 |
| Rubber tube reducer | Brass Φ30×Φ18 | 242185 |
| | Brass Φ30×Φ12 | 242186 |
| | SUS Φ30×Φ18 | 221496 |
| | SUS Φ30×Φ12 | 241497 |
| Stainless steel lid | | 281440 |



Example of combination with glass condenser




Glass condenser




Stainless steel lid

Cooling curve




Combination usage examples

Example of combination with vacuum drying oven



Example of combination with rotary evaporator



Contents

| | |
|--|-----|
| Ultrasonic Crushers/Homogenizers | 198 |
| Shakers | 200 |
| Magnetic Stirrers | 204 |
| Hot Plates | 211 |
| Laboratory Stirrers | 212 |

Ultrasonic Crushers/Homogenizers | Cell Destruction, Emulsification, Dispersion, Mixing

SFX150-Y·SFX250-Y·SFX550-Y

Ultrasonic power 150W (Model 150), 250W (Model 250), 550W (Model 550)

Ultrasonic mode Continuous/Pulse

Widely used for liquid processing such as breaking nano-level biological cells, homogenization, emulsification, accelerating reactions, dispersion, mixing, degassing, etc. No need for materials like broken glass to promote reactions, reducing operating costs.

Features

- 40 kHz ultrasonic frequency, low noise, providing a comfortable operating environment. (Model 150)
- Processing volume from 0.2 ml to 150 ml, suitable for most small sample volumes. (Model 150)
- Processing volume from 0.2 ml to 500 ml (Model 250), and from 0.2 ml to 1000 ml (Model 550).
- 20 kHz ultrasonic frequency, capable of processing the vast majority of samples, with powerful performance and excellent reliability. (Model 250/550)
- Selectable continuous or pulse ultrasonic modes.
- Amplitude setting range from 10 to 100%.
- Specialized microprobe amplitude limiting mode to extend probe lifespan.
- Multifunctional button-type LCD display, easy to operate, and intuitive display.
- High-efficiency transducer with high conversion rate, does not generate heat, no need for cooling, can work continuously for a long time, providing more stable and convenient time, temperature, and energy ultrasonic modes, selectable based on sample features.
- Real-time monitoring of experimental data during the processing, including power, output energy, temperature values, and experiment progress.
- Real-time temperature monitoring system accurately responds to specific user temperature control needs. (Temperature probe required)
- It can store up to 20 sets of programs. Programmable parameters include continuous or pulse mode; time, energy, or temperature control mode; pulse ultrasonic settings for time or energy; amplitude settings; temperature settings, etc.
- Real temperature control achieved by automatically adjusting the pulse time ratio to keep the sample temperature within the user-set range (requires temperature probe).




Handheld transducer



Replaceable microprobes suitable for various sample sizes

- LED indicator lights to alert the user when the device is active.
- The round top ring provides a comfortable grip.
- Embedded start button to prevent accidental operation.
- Ergonomic handle for easy and comfortable control.



- Continuous or pulse ultrasonic
- Time, energy, and temperature control modes
- Experiment progress prompt
- Real-time process monitoring
- Microprobe protection mode
- Temperature mode setting
- Editable time or energy mode
- Store 20 sets of data

Specifications

| Model | SFX150-Y | SFX250-Y | SFX550-Y |
|------------------------------|------------------------|--------------------------------|---|
| Maximum output power | 150W | 250W | 550W |
| Maximum energy setting | 999999J | 999999J | |
| Maximum time setting | 99 h 59 min 59 seconds | 99 h 59 min 59 seconds | |
| Minimum pulse energy setting | 1J | 1J | |
| Minimum pulse time setting | 0.01 s | 0.01 s | |
| Standard frequency | 40KHz | 20KHz | 20KHz |
| Standard probe | 1/8" probe | 1/2" (13mm) 7-4 titanium alloy | 1/2" (13mm) or 3/4" (19mm) 7-4 titanium alloy |
| Maximum current | 1A | 1.5A | 6A |
| Rated voltage | AC 200-240V at 50/60Hz | AC 200-240V at 50/60Hz | |
| Main unit weight | 6.5kg | 6.5kg | |
| Main unit dimensions | 348L*203W*242H (mm) | 348L*203W*242H (mm) | |

Made in
Mexico

■ Sonifier® SFX150-Y Series:



Cup probe: Limited to 3.0" diameter




Microprobes: Dimensions range from 1/8" to 1/4"



40 kHz conventional transducer

■ Sonifier® SFX250-Y/SFX550-Y Series: Optional Accessories



Lotus-shaped cooling cup

transducer support bracket

soundproof box

■ Accessories ordering information

| Catalog No. | Note |
|--------------|---|
| 200-060-022R | Temperature probe |
| 101-063-275 | Soundproof box |
| 101-147-043 | (19mm) non-replaceable tip probe |
| 101-147-047 | 2" cup probe (requires 100-008-249 bolt) |
| 101-147-048 | 3" cup probe |
| 101-148-013 | Replaceable tip |
| 101-147-049 | 1/2" extension probe with replaceable tip 101-148-013 |

| Catalog No. | Note |
|--------------|--|
| 101-148-062 | Tapered 1/8" diameter microprobe for 1/2" threaded probe |
| 101-148-070 | Tapered 1/4" diameter microprobe for 1/2" threaded probe |
| 101-063-212 | Dual-stage 1/8" diameter microprobe with connector, directly connectable to transducer |
| 101-147-050 | Dual-stage microprobe connector part |
| 101-148-063 | Dual-stage microprobe lower end probe part |
| 101-063-1110 | Transducer fixing ring for 1/2" laboratory stand |
| 101-135-066R | 20 KHz 102-C transducer |
| 101-118-039 | Crescent wrench (2 required) |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Shakers | Horizontal • Vertical • Orbital • Dual Vertical Shaking

SA300/320/400

Shaking rate 20~300rpm (SA300/400)

20~210rpm (SA320)

Amplitude 40mm



SA320

SA300

Two-dimensional shaking of horizontal and vertical (SA300), two-dimensional shaking of horizontal orbital and vertical orbital (SA320), dual vertical shaking type (SA400), applicable to large-volume sample extraction, cultivation, mixing, stirring, etc.

■ Features

- Use DC brush motor for stable torque from low to high speed rotation.
- Compact and capable of carrying ample load for shaking.
- With timer, it can set any time.

SA300/320

- Body can lie flat for horizontal shaking, enabling two-dimensional shaking of vertical and horizontal, widely used.
- SA320 orbital shaking model provides good effect as an incubator in horizontal position. Additionally, high shaking effect even at low orbital speed during vertical shaking.
- Equipped with various optional items for different experiments.

SA400

- Simultaneously shakes 6×1l erlenmeyer flasks or 4×2l erlenmeyer flasks. Rotation speed setting and easy-to-confirm scale setting and digital rotation counter.



SA320

SA400

■ Specification

| Model | SA300 | SA320 | SA400 |
|-------------------------|--|--|---|
| Shaking method | Horizontal · vertical reciprocating shaking | | Vertical reciprocating dual shaking |
| Maximum number of tubes | Erlenmeyer flask 1,000 ml×3, 2,000 ml×2 | | Erlenmeyer flask 1,000 ml×6, 2,000 ml×4 |
| Oscillating frequency | Horizontal shaking 20~300 rpm vertical shaking 20~300 rpm | Horizontal orbital shaking 20~210 rpm vertical orbital shaking 20~210 rpm | Vertical shaking 20~300 rpm |
| Speed setting display | Variable scale on dial | | Dial setting, digital display |
| Timer | Dial setting 0~60 min (minimum scale 5 min) with continuous operation switch | | |
| Motor | DC motor 90w | | |
| External dimensions | Width 460×Depth 460×Height 423mm | | W520×D460×H483 mm |
| Body weight | Approx. 55.5kg | | Approx. 39kg |
| Power supply | AC220V 1A | | |
| Accessories | One fuse, one set of carbon brushes | | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Made in
Japan

Horizontal shaking



Orbital shaking



Vertical shaking



Options

Centrifuge tube shaking rack



Horizontal-vertical shaking for SA300/320/400

| Product number | Specifications | Number of stands |
|----------------|------------------------------|------------------|
| 232087 | Φ16~35mm Length 110~130mm | 18 pcs |

Test tube shaking rack



Horizontal shaking for SA300/320

| Product number | Specifications | Number of stands |
|----------------|--------------------------------|------------------|
| 232086 | Φ16.5~18mm Length 160~190mm | 18 pcs |

2L separating funnel shaking rack



Vertical shaking for SA300/320/400

| Product number | Specifications |
|----------------|--------------------------------|
| 232089 | 50ml, 100ml~1,000ml 2,000ml |

Separating funnel shaking rack



Vertical shaking for SA300/320/400

| Product number | Specifications |
|----------------|----------------|
| 232096 | 100ml~1,000ml |

Test tube shaking rack base



Horizontal shaking for SA300/320

| Product number | Specifications |
|----------------|--|
| 232088 | 2 layers, maximum test tube rack width 238×depth 121×height 105 mm |

● Installation Platform To Be Purchased Separately.

Non-slip mat



Horizontal shaking for SA300/320

| Product number | Specifications |
|----------------|---|
| 232071 | Thickness 1 mm Width 450×Depth 396mm |

● Installation Platform To Be Purchased Separately.

Erlenmeyer flask shaking rack



Horizontal shaking for SA300/320

| Product number | Specifications |
|----------------|----------------|
| 232097 | 100ml~1,000ml |

Erlenmeyer flask clamp



Horizontal shaking for SA300/320

| Product number | Specifications | Number of stands |
|----------------|----------------|------------------|
| 232062 | For 100 ml | 10 sets |
| 232063 | For 200 ml | 9 sets |
| 232064 | For 300 ml | 5 sets |
| 232065 | For 500 ml | 4 sets |
| 232066 | For 1,000 ml | 2 sets |

● Installation Platform To be purchased separately

Installation platform



Horizontal shaking for SA300/320

| Product number | Specifications |
|----------------|---|
| 232095 | Number of erlenmeyer flask clamps: 28 pcs for 100 ml, 19 pcs for 200 ml, 14 pcs for 500 ml, 9 pcs for 1,000 ml |

Shaking rack



Horizontal shaking for SA300/320

| Tilted erlenmeyer flask clamp shaking rack | | |
|--|----------------|------------------|
| Product number | Specifications | |
| 232067 | For 100 ml | |
| 232068 | For 200 ml | |
| 232069 | For 300 ml | |
| Tilted test tube shaking rack | | |
| Product number | Specifications | Number of stands |
| 232080 | For Φ12mm | 50 pcs |
| 232081 | For Φ16.5mm | 20 pcs |
| 232082 | For Φ18mm | 20 pcs |

2-layer spring shaking rack



Horizontal shaking for SA300/320

| Product number | Specifications |
|----------------|---|
| 232079 | Φ16 mm test tubes 320 pcs (Grid spacing 20 mm) |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Shakers | Reciprocating • Elliptical • Orbital Shaking

MK161

Made in
Japan

One-touch switching between reciprocating • elliptical • orbital in 3 stages.

Features

- Compact design, combined with incubator as shaking incubator.
- Use DC brushless motor, although small in size, Achieve stable shaking speed and high torque shaking power (with additional weight of 8 kg).
- Shaking speed and time are set by scale, digital display.
- With pause function for temporary stop, switching between timed operation and continuous operation is one-touch.
- Various mixing, extraction, stirring modes can be selected depending on combination of shaking table and special rack (optional).

Specifications

| Model | MK161 |
|--------------------------|--|
| Shaking method | Continuous switching (manual) between orbital • elliptical • reciprocating |
| Amplitude | Orbital movement: 30 mm reciprocating movement: 30 mm |
| Oscillating frequency | 20~200rpm |
| Shaking speed adjuster | Scale setting/digital display |
| Timer | Scale setting/digital display digital 0.1 min (6 sec)~99.9 min |
| Shaking table dimensions | Body: w300×d254 mm mounted: W290×D250 mm |
| External dimensions | Width 350×Depth 300×Height 150mm |
| Weight | Approx. 15kg |
| Power supply | AC100V 0.6A |



Options


Installation stand



Both flask clamp and shaking rack installations require the installation stand.

| Product number | Specifications |
|----------------|----------------|
| 232061 | W290×D250mm |

triangular flask clamp



Example of flask clamp and installation stand setup

| Product number | Capacity | Number of stands |
|----------------|--------------|------------------|
| 232062 | For 100 ml | 10 sets |
| 232063 | For 200 ml | 9 sets |
| 232064 | For 300 ml | 5 sets |
| 232065 | For 500 ml | 4 sets |
| 232066 | For 1,000 ml | 2 sets |

shaking rack



example of shaking rack and installation stand setup

Shaking rack for inclined triangular flask

| Product number | Capacity |
|----------------|------------|
| 232067 | For 100 ml |
| 232068 | For 200 ml |
| 232069 | For 300 ml |

- Able to shake two sets simultaneously.

Shaking rack for inclined tube

| Product number | Diameter | Number of stands |
|----------------|----------|------------------|
| 232080 | Φ12mm | 50 pcs |
| 232081 | Φ16.5mm | 20 pcs |
| 232082 | Φ18mm | 20 pcs |

- Able to shake two sets simultaneously.

Shaking rack



Inclined multi-tube shaking rack

| Product number | Capacity | Number of stands |
|----------------|----------|------------------|
| 232070 | 15ml | 12 |

Inclined 50ml centrifuge tube shaking rack

| Product number | Diameter |
|----------------|----------|
| 232083 | Φ29mm |

Non-slip mat



| Product number | Specifications |
|----------------|-----------------|
| 232084 | W290×D250×H30mm |

Spring shaking rack



| Product number | Specifications |
|--|-----------------|
| 232050 | W290×D250×H66mm |
| Number of racks (test tubes): Φ16mm×64 tubes (when placed at 45 degrees) | |
| Number of racks (triangular flasks): 50ml×20, 100ml×10, 200ml×9, 300ml×5, 500ml×4, 1,000ml×2 | |

2-layer shaking rack



| Product number | Specifications |
|--|------------------|
| 232056 | W290×D250×H110mm |
| Number of racks (test tubes): Φ16mm×64 tubes (when placed at 45 degrees) | |
| Number of racks (triangular flasks): 50ml×20, 100ml×10, 200ml×9, 300ml×5, 500ml×4, 1,000ml×2 | |

Shakers | Reciprocating · Orbital Shaking

MK201D

Made in
Japan**Compact type for both reciprocating and orbital shaking.**

■ Features

- It can mix, extract, and stir with a reciprocating/orbital switching shaker.
- Timer can be set for shaking time.
- Shaking speed is displayed digitally.
- By choosing different shaking tables (sold separately), various containers can be shaken.
- Shaking table is easy to attach and detach, and safe without displacement during vibration.

■ Specifications

| Model | MK201D |
|----------------------------|--|
| Shaking method | Reciprocating orbital (manual switch) |
| Amplitude | Orbital motion: 30mm Reciprocating motion: 30mm |
| Oscillating frequency | 20~200rpm infinitely variable speed |
| Shaking speed display mode | Digital display |
| Timer | 0.5 seconds~100 h (unit and scale conversion type) |
| External dimensions | Width 442×Depth 415×Height 130mm |
| Weight | 19kg |
| Power supply | AC100V 2A |

■ Options

Shaking rack for fixed triangular flask



| Product number | Capacity | Number of stands |
|----------------|--------------|------------------|
| 232170 | For 100 ml | 20 sets |
| 232171 | For 200 ml | 10 sets |
| 232172 | For 300 ml | 6 sets |
| 232173 | For 500 ml | 4 sets |
| 232174 | For 1,000 ml | 2 sets |

Centrifuge tube shaking rack



Installation example of 2 sets (inclined installation stand sold separately)

| Product number | Specifications | Number of stands |
|----------------|--------------------------------------|------------------|
| 232179 | 50ml $\Phi 30 \times 110$ mm or less | 8 pcs |

- Able to shake two sets simultaneously.
- Inclined installation stand sold separately.

Non-slip mat



| Product number | Specifications |
|----------------|---------------------------------------|
| 232182 | W400×D330mm material: silicone rubber |

Shaking rack for inclined triangular flask (angle adjustable)



Installation example of 2 sets (inclined installation stand sold separately)

| Product number | Capacity | Number of stands |
|----------------|------------|------------------|
| 232175 | For 100 ml | 8 sets |
| 232176 | For 200 ml | 4 sets |
| 232177 | For 500 ml | 2 sets |

- Able to shake two sets simultaneously.
- Inclined installation stand sold separately.

Universal spring shaking rack



| Product number | Specifications | Number of stands |
|----------------|----------------|------------------|
| 232180 | 250ml | 8 pcs |

2-layer spring shaking rack for triangular flasks and test tubes



| Product number | Specifications |
|----------------|--|
| 231398 | Test tubes $\Phi 16$ mm (45-degree incline) 105 tubes Triangular flasks 100ml 16 sets triangular flasks 200ml 10 sets triangular flasks 500ml 4 sets triangular flasks 1,000ml 2 sets |

Multi-tube shaking rack



Installation example of 2 sets (inclined installation stand sold separately)

| Product number | Specifications | Number of stands |
|----------------|--|------------------|
| 232178 | 15ml $\Phi 16.5 \times 110$ mm or less | 12 pcs |

- Able to shake two sets simultaneously.
- Inclined installation stand sold separately.

Spring shaking rack for triangular flask and test tube



| Product number | Specifications |
|----------------|---|
| 232181 | 17×15=255 grids (inner dimension 14×14mm square grid) |

Inclined installation platform



| Product number | Specifications |
|----------------|--|
| 232178 | Used in combination with dedicated shaking racks |

Magnetic Stirrers | Magnetic Stirring

MFD810-B

Made in
Japan

Rotation rate 50~1600rpm (set in 10rpm increments)

Heat Resistant Temp. 310°C

A combination that maximizes temperature performance and convenient functions.







Features

- Suitable for synthesis experiments requiring high temperature accuracy and stirring efficiency.
- Due to strong magnetic force, even when moving the flask, stirring continues. The magnetic stirrer continues stirring even when separated from the stirring table.
- Feature a decoupling detection mode (stir bar detachment detection), which allows it to resume if the stir bar detaches from the container upon detection.
- When used with the dedicated oil bath BOG/BOS series, the recessed bottom of the oil bath can be embedded in the top of the MFD810-B to prevent the risk of displacement and falling due to vibration.
- Rich and convenient rotation functions.
- The stirring table (hot plate) is designed to be heat-resistant and chemical-resistant.

Specifications

| Model | | MFD810-B |
|---|-------------------------------|---|
| Performance | Speed Range | 50~1600rpm (set in 10rpm increments) |
| | Motor | DC brushless motor (31W) |
| Composition | Magnet used | Neodymium magnet |
| | Display | White led digital display |
| | Exterior material | Aluminum die-casting |
| | Stirring table | Aluminum die-casting (ceramic coating) |
| Stirring motor functions (rotation modes) | | Constant speed, decoupling detection, intermittent, reverse, slow ramp |
| Safety functions | | Overcurrent fuse |
| Other functions | | Service socket |
| Specifications | External dimensions (l*d*h) | 165×275×90mm stirring table diameter Φ135mm |
| | Load capacity | Less than 30kg |
| | Power supply | 220V±5% 0.15A 50Hz |
| | Overcurrent fuse capacity | For internal circuit: 0.5A for service outlet: 5A |
| | Power supply line | 2m insert type FG plug |
| | Weight | Approx. 2.8kg |
| | Usage environment temp. range | 4~40°C |
| Accessories | | Heat cut plate 1pc, protection cover 1pc, screw caps 2pcs, power cable 1pc, spare fuse (for internal circuit, service outlet) 1pc each, instruction manual 1 copy, warranty card 1 copy |

Related optional items

Laboratory lift
OA147Container overturn
prevention frame
OA146Dual-hand
aluminum block
handle
OA145Single-hand
aluminum block
handle
OA144Strong magnetic
stirring bar
OA149/150/151

Magnetic Stirrers | Magnetic Stirring with Heating Function

MFH810-B

Made in
Japan

Rotation rate 50~1600rpm (set in 10rpm increments)

Heat Resistant Temp. 310°C

Magnetic stirring with heating function.



Features




- Due to strong magnetic force, even when moving the flask, stirring continues. The magnetic stirrer continues stirring even when separated from the stirring table.
- Obtained a patent for aluminum block thermostatic bath.
- Multiple aluminum block systems available.
- Feature a decoupling detection mode (stir bar detachment detection), which allows it to resume if the stir bar detaches from the container upon detection.
- When used with the dedicated oil bath BOG/BOS series, the recessed bottom of the oil bath can be embedded in the top of the MFH810-B to prevent the risk of displacement and falling due to vibration.
- Rich and convenient rotation functions.
- The stirring table (hot plate) is designed to be heat-resistant and chemical-resistant.

Specifications

| Model | | MFH810-B |
|---|--------------------------------|--|
| Hot Plates | | Yes |
| Performance *1 | Speed Range | 50~1600rpm (set in 10rpm increments) |
| | Temp. control range: | Room Temperature +25°C ² ~310°C (set in 1°C increments) |
| | Temp. adjusting accuracy | ±1.5°C at 100°C (internal temperature sensor); ±1.0°C at 50°C (external temperature sensor) |
| | Display | White led digital display |
| Composition | Motor | DC brushless motor (31W) |
| | Magnet used | Neodymium magnet |
| | Temp. control method | PID control |
| | Temp. sensor | PT100 |
| | Heater | 600W Mica heater |
| | Exterior material | Aluminum die-casting |
| | Stirring table | Aluminum die-casting (ceramic coating) |
| Stirring motor functions (rotation modes) | | Constant speed, decoupling detection, intermittent, reverse, slow ramp |
| Safety functions | | Overcurrent fuse |
| Other functions | | Service socket |
| Specifications | External dimensions (L*D*H) *3 | 165×275×90mm stirring table diameter Φ135mm |
| | Load capacity | Less than 30kg |
| | Power supply | 220V±5% 3A 50Hz |
| | Overcurrent fuse capacity | For internal circuit: 7A for service outlet: 5A |
| | Power supply line | 2m insert type FG plug |
| | Weight | Approx. 3.0kg |
| | Usage environment temp. range | 4~40°C |
| Included accessories *4 | | Heat cut plate 1pc, protection cover 1pc, screw caps 2pcs, power cable 1pc, spare fuse (for internal circuit, service outlet) 1pc each, instruction manual 1 copy, warranty card 1 copy, external temperature sensor 1pc |

*1 Power Supply AC 220V±5%, Room Temperature 23°C±5°C, Humidity 65% RH±20%, when using internal sensor without load, when using external temperature sensor value is inserted into the optional aluminum block. *2 At room temperature 23°C, speed 1600rpm (maximum condition) self-heating value. *3 Excluding protruding parts. *4 This product does not include a stir bar.

Related optional items

Protection cover
(for MFH)
OA155Laboratory lift
OA147Container overturn
prevention frame
OA146Dual-hand aluminum
block handle
OA145Single-hand aluminum
block handle
OA144Strong magnetic
stirring bar
OA149/150/151

Magnetic Stirrers | Magnetic Stirring

MA300A/300B

Made in
Japan

Rotation rate 100~1,200rpm

Stirring capacity 50~3,000ml (Model 300A) 50~1,000ml (Model 300B)

Compact digital stirrer



Features

MA300A type

- Highly versatile basic magnetic stirrer.
- Simple and highly functional design makes operation easier.
- The rotation speed is displayed digitally, making it easy to set and reproduce stirring conditions.

MA300B type

- Digital display ultrathin high-quality stirrer.
- Excellent for operation and setting, ultrathin for easy storage and management.

Specifications

| Model | MA300A | MA300B |
|--------------------------------------|-----------------------------------|-----------------------------------|
| Stirring table material | Stainless steel plate | |
| Stirring table dimensions | Width 172×depth 156mm | |
| Stirring capacity | 50~3,000ml | 50~1,000ml |
| Rotation times | 100~1,200rpm | |
| Motor | DC Motor | |
| Power supply (50/60Hz) rated current | AC100V 0.5A | |
| External dimensions | Width 172×depth 156×height 53.7mm | Width 172×depth 156×height 37.7mm |
| Accessories | Stirring rod 30mm 1 unit | |

Magnetic Stirrers | Economic Type

MA100/300

Made in
Japan

Rotation rate 100~1,500rpm (Model 100) 200~1,500rpm (Model 300)

Stirring capacity 50~1,000ml (Model 100) 100~3,000ml (Model 300)

Durability is excellent, capable of stable rotation from low to high speed.

Features

MA100 type

- Small and low-cost type.
- Use DC motor electronic control.

MA300 type

- Low-price, high-durability magnetic stirrer.
- Maximum stirring volume is 3L.



Specifications

| Model | MA100 | MA300 |
|--------------------------------------|----------------------------------|---------------------------------|
| Stirring table material | Stainless steel plate | |
| Stirring table dimensions | Width 77×Depth 135mm | Width 178×Depth 165mm |
| Stirring capacity | 50~1,000mL | 100~3,000mL |
| Rotation times | Approx. 100~1,500rpm | Approx. 400~1,550rpm |
| Motor | DC motor electronic control 0.3W | AC motor shaded pole type 0.7W |
| Power supply (50/60Hz) rated current | AC100V 0.2A | AC100V 0.5A |
| External dimensions | Width 84×Depth 150×Height 60mm | Width 181×Depth 197×Height 83mm |
| Weight | Approx. 1.1kg | Approx. 1.6kg |
| Accessories | Stirring rod 30mm 1 unit | |

Magnetic Stirrers | Simultaneous Stirring

MG120/600

Made in Japan

Rotation rate 100~1,500rpm (Model 120) 200~1,500rpm (Model 600)

Stirring capacity 50~1,00ml (Model 120) 50~2,000ml (Model 600)

Features

Model MG120 (12 simultaneous stirrings)

- Capable of high-load simultaneous stirring of up to 12 samples.
- Small, ultra-thin stirrer with excellent operability.
- The surface is made of chemically resistant ceramic steel plate, with a leak prevention structure.

Model MG600 (6 individual stirrings)

- The surface is made of chemically resistant ceramic steel plate, with good chemical resistance.
- 6 stirring chambers can individually adjust rotation speed.
- Capable of stable rotation from low to high speed.



Specifications

| Model | MG120 | MG600 |
|--------------------------------------|---|---|
| Stirrers | 12-linked type (with OFF knob) | 6 units with individual rotation adjustment (with OFF knob) |
| Stirring table (hot plate) material | Ceramic steel plate (stainless steel plate ceramic coating) | |
| Stirring table dimensions | Width 295×Depth 198mm | Width 449×Depth 300mm |
| Stirring capacity | 50~100ml×12 (maximum) | 50~2,000ml×6 units |
| Rotation times | Approx. 100~1,500rpm | Approx. 200~1,500rpm |
| Motor | DC motor×12 linked electronic control | AC Motor 0.7W×6 units |
| Power supply (50/60Hz) rated current | AC100V 0.15A | AC100V 0.5A |
| External dimensions | Width 303×Depth 234×Height 55mm | Width 456×Depth 335×Height 90mm |
| Weight | Approx. 4kg | Approx. 7.5kg |
| Accessories | Stirring rod 25mm×12 | Stirring rod 30mm×6 |

Magnetic Stirrers | Super Strong Magnetic Stirring

MB800

Made in Japan

Rotation rate 70~1,200rpm

Stirring capacity 100~10,000ml

Detects the temperature at the bottom of the water bath and cuts off the power.

Features

MB800 type

- Ideal for the BO500 type oil bath magnetic stirrer, with a chemically resistant aluminum oxide membrane treated stirring plate.
- The power to the oil bath is automatically cut off when the bottom temperature of the oil bath reaches the set value.

Specifications

| Model | MB800 |
|--------------------------------------|--|
| Stirring table material | Aluminum |
| Stirring table dimensions | Width 250×Depth 220mm |
| Stirring capacity | 100~10,000ml |
| Rotation times | 70~1,200rpm |
| Motor | AC motor electronic control |
| Overheat preventer | 70~200°C |
| Sensors | Thermistor |
| Service socket (for oil bath) | AC100V 10A |
| Safety device | Ground fault circuit interrupter, overheating protection for oil bath power cutoff |
| Power supply (50/60Hz) rated current | AC100V 11A (maximum when connected to oil bath) |
| External dimensions | Dimensions: width 250×depth 270×height 150mm |
| Weight | Approx. 4.2kg |
| Accessories | Stirring rod 40mm×1 |



Combination example: MB800+BO500

12

Stirrers & Shakers

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Magnetic Stirrers | Low Viscosity to High Viscosity Compatible

MD300/500/800·MS500D

Made in
Japan

| | | | | |
|------------------------|------------------|------------------|--------------------|-------------------|
| Rotation rate (rpm) | MD300 (80~1,500) | MD500 (70~1,300) | MS500D (10~1,400) | MD800 (50~1,400) |
| Stirring capacity (ml) | MD300 (50~3,000) | MD500 (50~5,000) | MS500D (50~10,000) | MD800 (50~10,000) |



Excellent chemical resistance, leak prevention.

Magnetic stirrer with excellent chemical resistance ceramic steel plate, prevents controller failure through leak prevention structure.

■ Features

MD300 type

- DC motor electronic control, capable of stable stirring.

MD500 type

- Capable of stable rotation from low to high speed.

Model MS500D (Slow stirring)

- Capable of stable rotation even at low speed of 10rpm.
- Adopt feedback control, stable rotation under various loads from low to high speed.
- Equipped with a digital rotation counter for easy setting and confirmation of rotation speed.
- Maximum stirring volume is 10L.

MD800 type

- Equipped with a 2500G strong magnet, suitable for stirring 10L large capacity and high viscosity liquids.

■ Specifications

| Model | MD300 | MD500 | MS500D | MD800 |
|---|---|---------------------------------|---------------------------------|----------------------------------|
| Stirring table material | Ceramic steel plate (stainless steel plate ceramic coating) | | | |
| Stirring table dimensions | Width 167×Depth 147mm | Width 191×Depth 177mm | Width 190×Depth 180mm | Width 217×Depth 214mm |
| Stirring capacity | 50~3,000ml | 50~5,000ml | 50~10,000ml | |
| Rotation times | Approx. 100~1,500rpm | Approx. 70~1,300rpm | Approx. 10~1,400rpm | Approx. 50~1,400rpm |
| Motor | DC motor electronic control | AC motor 2W | Optical pulse feedback system | AC motor electronic control 3W |
| Power supply (50/60Hz) rated current | AC100V 0.2A | AC100V 0.5A | | |
| External dimensions | Width 174×Depth 180×Height 59mm | Width 199×Depth 225×Height 89mm | Width 196×Depth 200×Height 81mm | Width 224×Depth 263×Height 142mm |
| Weight | Approx. 1.4kg | Approx. 2.2kg | | Approx. 3.2kg |
| Accessories | Stirring rod 30mm×1 | | | Stirring rod 40mm×1 |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Magnetic Stirrers | Super Strong Magnetic Stirring

MC801

Made in
Japan

Excellent chemical resistance, super strong magnetic stirring type.

Features

- Multi-purpose magnetic stirrer, capable of right, left, or reverse rotation.
- Powerful stirring and rotational variations for effective mixing.
- Reversal rotation time can be set to 120 seconds.
- Equipped with a high-precision controlled 30W DC brushless motor for stable stirring.
- Equipped with a large strong ferrite magnet.



Specifications

| Model | MC801 |
|--------------------------------------|---|
| Stirring table material | Ceramic steel plate (stainless steel plate ceramic coating) |
| Stirring table dimensions | Width 270×Depth 243mm |
| Stirring capacity | 100~10,000ml |
| Rotation times | Approx. 100~1,500rpm |
| Motor | DC motor electronic control |
| Timer | Reversal time adjustable between 10~120 seconds |
| Power supply (50/60Hz) rated current | AC100V 0.3A |
| External dimensions | Width 278×Depth 290×Height 100mm |
| Weight | Approx. 4.2kg |
| Accessories | Stirring rod 30mm×1, 40mm×1 |

Magnetic Stirrers | Super Strong Magnetic Stirring

MF820

Features

- Using rare earth magnets, super strong magnetic force of 6000G (Gauss).
- Capable of easy stirring even with high viscosity solutions.
- When a magnetic line cover facing upwards is installed, capable of simultaneous stirring of multiple small containers (stirring table optional).
- Ultra-thin with excellent operability.
- Equipped with a 7A service socket.



Specifications

| Model | MF820 |
|--------------------------------------|--|
| Stirring table material | Stainless steel plate |
| Stirring table dimensions | Width 270×Depth 260mm |
| Stirring capacity | 100~20,000ml |
| Rotation times | Approx. 80~1,500rpm |
| Motor | DC brushless motor 30w high precision feedback control |
| Timer | — |
| Power supply (50/60Hz) rated current | AC100V 8A (including 7A service socket) |
| External dimensions | Width 276×Depth 295×Height 80mm |
| Weight | 3.8kg |
| Accessories | Stirring rod 30mm×1, 40mm×1 |



Combination Example: MF820+Stirring Plate (optional)

12

Stirrers & Shakers

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Magnetic Stirrers | With Hot Plate

MH301/520/800·MG600H

Made in
Japan

| | | | | |
|------------------------|-------------------|-------------------|--------------------|----------------------|
| Rotation rate (rpm) | MH301 (400~1,500) | MH520 (150~1,150) | MH800 (100~1,400) | MG600H (100~2,000) |
| Stirring capacity (ml) | MH301 (100~3,000) | MH520 (50~5,000) | MH800 (200~10,000) | MG600H (100~2,000×6) |



Features

MH301/800 type

- The stirring plate is made of aluminum ceramic coating material with excellent heat conduction, and the ceramic coating also has excellent chemical resistance.
- Equipped with a strong magnet, capable of stirring large capacity and high viscosity samples.
- Temperature adjustment of the hot plate is knob type.

MH520 type

- Stable rotation for better effects.
- Round hot plate with excellent thermal efficiency, high temperature rises speed and stability.
- Ceramic coated stirring plate with excellent chemical resistance.
- Equipped with circuit protection device.
- Maximum stirring capacity of up to 5L with strong power.
- Use stirring rods of 20~60mm.

Specifications

| Model | MH301 | MH520 | MH800 |
|--------------------------------------|-------------------------------|--|--------------------------------|
| Stirring table material | Aluminum ceramic coating | | |
| Stirring table dimensions | Width 176×Depth 151mm | Φ168mm | Width 299×Depth 285mm |
| Stirring capacity | 100~3,000ml | 50~5,000ml | 200~10,000ml |
| Rotation times | Approx. 400~1,500rpm | Approx. 150~1,150rpm | Approx. 100~1,400rpm |
| Heater | 400W | 470W | 1000W |
| Temp. adjustment | Bidirectional thyristor | | |
| Hot plate temp. | Max.300°C (set with OFF knob) | Max.325°C | Max.250°C (set with OFF knob) |
| Motor | AC motor electronic control | Asynchronous motor, phase control/electromagnetic control dual-use | AC motor electronic control 3W |
| Power supply (50/60Hz) rated current | AC100V 3A | AC100V 6A | AC100V 11A |
| External dimensions | 184×202×114 (W×D×H mm) | 190×223×123 (W×D×H mm) | 309×315×151 (W×D×H mm) |
| Weight | Approx. 2.6kg | Approx. 3.1kg | Approx. 6.7kg |
| Accessories | Stirring rod 30mm×1 | Stirring rod 30mm×1 | Stirring rod 40mm×1 |



Combination example: MH520+container tip-over protection cover

6 units individually adjustable (individual stirring and heating)

Features

MG600H type

- 6 units with individually adjustable speed and heating.
- The hot plate has a chemically resistant ceramic coating.



Specifications

| Model | MG600H |
|--------------------------------------|--|
| Stirrers | 6 units with individual rotation adjustment |
| Stirring plate (hot plate) material | Aluminum ceramic coating |
| Stirring table dimensions | Φ126mm/6 pcs |
| Stirring capacity | 100~2,000mL×6 units |
| Rotation times | Approx. 300~1,500rpm |
| Hot Plates | Width 230mm×6 units with individual temperature adjustment |
| Heater | 230W×6 |
| Temp. adjustment | Bidirectional thyristor |
| Hot plate temp. | Max. 250°C |
| Motor | AC shaded pole motor |
| Power supply (50/60Hz) rated current | AC100V 15A |
| External dimensions | Width 606×Depth 420×Height 122mm |
| Weight | Approx. 14kg |
| Accessories | Stirring rod 30mm×6 |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Hot Plates

HK200/300·HM300

Made in
Japan

Features

Model HK200/300 (sand chamber type)

- Temperature accuracy remains stable even at high temperatures. (Considering the arrangement of heater and temperature sensor)
- Carefully designed to reduce heat dissipation to the casing.
- The hot plate is sturdy and level, can be used as a sand chamber, or to place sample containers for processing and experiments.

Model HM300 (medium temperature type)

- Stainless steel heating surface with excellent chemical resistance.
- Wider heating surface for more efficient use.
- Specifically for medium temperatures, temperature adjustment accuracy can be maintained within $\pm 0.5^{\circ}\text{C}$ at 40°C .
- Low and stable type, lightweight, and easy to move.



Specifications

| Model | HK200 | HK300 |
|--------------------------------------|---|-------------|
| Operating temp. range | 50~250°C | |
| Temp. adjusting accuracy | $\pm 5^{\circ}\text{C}$ (center of the hot plate at 250°C, no load) | |
| Max. temp. reaching time | Approx. 40 min (center of the hot plate, no load) | |
| Temp. adjustment method | Liquid expansion type | |
| Hot plate material | Aluminum (hard acid-resistant aluminum processing) | |
| Power supply (50/60Hz) rated current | AC100V 7.5A | AC100V 9A |
| Heater capacity | 750W | 900W |
| Hot plate dimensions | 338×238×25 | 388×288×25 |
| External dimensions (W×D×H mm) | 392×330×160 | 442×380×160 |
| Weight | 6.2kg | 7.8kg |
| Safety device | Circuit protection device | |

Specifications

| Model | HM300 |
|--------------------------------------|-------------------------------------|
| Operating temp. range | Room temp. +5~80°C |
| Temp. adjusting accuracy | $\pm 0.5^{\circ}\text{C}$ (at 40°C) |
| Temp. adjustment method | PID control |
| Hot plate material | Stainless steel plate |
| Power supply (50/60Hz) rated current | AC100V 2A |
| Heater capacity | 160W (80W×2 units) |
| Hot plate dimensions | Width 450×Depth 300mm |
| External dimensions | Width 450×Depth 305×Height 90mm |
| Weight | 5.0kg |

Stirrers for Small Containers | Touch-driven Type

MT-31/51

Made in
Japan

Rotation rate 2,800rpm (Model 31) 600~3,000rpm (Model 51)

Vibration surface diameter 70mm

Features

- MT-31 is a fixed speed type, MT-51 is a variable speed type.
- A compact stirrer that can operate in both continuous and touch modes.
- Use an independent touch switch, driving the vibration surface by simply pressing the sample container and stopping when released.
- MT-51 model comes with a rotational speed adjustment knob.

Specifications

| Model | MT-31 | MT-51 |
|--------------------------------------|---|----------------------------------|
| Rotation rate | 2,800rpm | 600~3,000rpm |
| Motor | AC motor 2W | DC motor 10W |
| External dimensions | Width 104×Depth 155×Height 123mm | Width 128×Depth 165×Height 125mm |
| Weight | 1.6kg | 1.8kg |
| Switch | ON1 touch switch, ON2 continuous switch | |
| Applicable containers | Test tubes, centrifuge tubes, cuvettes, erlenmeyer flasks (under 100ml) | |
| Vibration surface diameter | 70mm | |
| Vibration surface material | Polyurethane | |
| Power supply (50/60Hz) rated current | AC100V 1A | |



12

Stirrers & Shakers

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Laboratory Stirrers | For Flask Stirring

LM100/200

Rotation rate 50~1000rpm



LM100

LM200

※ Stirring shaft and wings sold separately

It can be directly mounted on flasks.

■ Features

- Integrated stirring seal and drive unit, can be directly inserted into the flask, allowing stirring under vacuum or sealed conditions without time-consuming shaft adjustments.
- Capable of high-speed rotation from 50rpm~1000rpm.
- Use a 30W DC brushless motor for simple, high-torque operation.
- It can be vacuumed to 399.9Pa without loss of rotational efficiency.
- Compact and lightweight design allows easy operation even with droppers and reflux condensers mounted on 3-neck flasks.
- The LM200 model features digital display for rotational speed, making verification and reproduction very convenient.
- The LM100 model is an economical analog display type, allowing control and adjustment of rotational speed while monitoring the liquid state.
- Equipped with various stirring shafts and wings for small to large capacity.
- Shaft seal uses an oil seal, standard material is fluoro rubber, but can also use more corrosion-resistant PTFE (optional).
- The O-rings for the stirring seal are made of fluoro rubber as standard, but more corrosion-resistant perfluoro rubber is available (optional).
- Stirring seal gaskets are available in TS29/42 and TS24/40 sizes (optional).
- AC100~240V can be used.

■ Anaerobic reaction system example with nitrogen gas




Water bath system example

| Product name | Model | Product code | Quantity |
|---|-----------------|--------------|----------|
| • Flask stirrer | LM200 | 231528 | 1 |
| • Oil bath | BO601 | 222118 | 1 |
| • PTFE stirring shaft (with stirring wings) | OLM14 | 231617 | 1 |
| • Fixing clamp ×2 | OLM46 | 231633 | 2 |
| • Flask holder | OLM50 | 231635 | 1 |
| • Cooling water circulation device | CF301 | 221512 | 1 |
| • Diaphragm pump | PG301 | 255158 | 1 |
| • 3-neck flask 500mL | Sold separately | | |
| • Glass burette | Sold separately | | |
| • Thermometer, thermometer stand | Sold separately | | |
| • Piping | Sold separately | | |

- 1 Sterilizers
- 2 Granulation and Spray Dryers
- 3 Muffle Furnaces
- 4 Ovens
- 5 Incubators
- 6 Plasma Equipment
- 7 Water Purifiers
- 8 Baths
- 9 Water Circulators
- 10 Rotary Evaporators
- 11 Freeze Dryers & Cold Traps
- 12 Stirrers & Shakers
- 13 Washers
- 14 Analysis and Test Devices
- 15 Options


Made in
Japan

Options



| No. | Product name | Model | Product code | Specification |
|-----|--|-------|--------------|--|
| ① | FKM oil seal | OLM60 | 231640 | FKM Black (standard installed) 2 pieces |
| ② | PTFE oil seal | OLM62 | 231641 | 2 pieces |
| ③ | PTFE stirring shaft (with stirring wings) | OLM10 | 231615 | L350×Φ30mm 100~300mL |
| | | OLM12 | 231616 | L450×Φ30mm 100~500mL |
| | | OLM14 | 231617 | L450×Φ50mm 200~1000mL |
| | | OLM16 | 231618 | L450×Φ80mm 300~1000mL |
| | | OLM18 | 231619 | L600×Φ80mm 1000~1500mL |
| | | OLM20 | 231620 | L600×Φ100mm 1000~5000mL |
| | | OLM22 | 231621 | L600×Φ120mm 3000~5000mL |
| ④ | Glass stirring shaft | OLM24 | 231622 | L350mm |
| | | OLM26 | 231623 | L400mm |
| | | OLM28 | 231624 | L530mm |
| ⑤ | PTFE half-moon stirring wing | OLM30 | 231625 | Φ40×16×t3mm 50~300mL |
| | | OLM32 | 231626 | Φ50×17×t3mm 100~500mL |
| | | OLM34 | 231627 | Φ60×17×t4mm 500~5000mL |
| | | OLM36 | 231628 | Φ100×17×t4mm 1000~5000mL |
| | | OLM38 | 231629 | Φ125×30×t5mm for beakers |
| ⑥ | Glass half-moon stirring wing | OLM40 | 231630 | Φ50×17×t3.3mm 100~5000mL |
| | | OLM42 | 231631 | Φ80×17×t3.8mm 500~5000mL |
| ⑦ | Fixing clamp | OLM44 | 231632 | For Φ5~Φ13mm |
| ⑧ | Fixing clamp | OLM46 | 231633 | For Φ6~Φ17mm |
| ⑨ | Fixing clamp | OLM48 | 231634 | For Φ9.5~Φ29mm |
| ⑩ | Flask clamp | OLM50 | 231635 | Adjustment range 3~55mm, shaft diameter 10mm Suitable flasks: 50~3000mL 50mL~3000mL |
| ⑪ | Flask clamp | OLM52 | 231636 | Adjustment range 3~80mm, shaft diameter 12mm Suitable flasks: 50~5000mL |
| ⑫ | Y-type stand | - | 231086 | H725×Φ25 Foot size W400 (inner 310) ×420mm |
| ⑬ | FFKM O-ring | OLM64 | 231642 | For TS29/42 |
| | | OLM66 | 231643 | For TS24/40 |
| | FKM O-ring | OLM64 | LT00038897 | Standard installation for TS29/42 |
| | | OLM66 | LT00038898 | Optional for TS24/40 |
| ⑭ | 24/40 stirring seal set TS29/42 stirring seal set | OLM58 | 231639 | assembly kit including oil seal, o-ring, etc. |
| | | OLM68 | 231644 | |
| ⑮ | Flask clamp | OLM54 | 231637 | For TS29/42, 2 pieces for TS24/40, 2 pieces |
| | | OLM56 | 231638 | |

Related products



Specifications

| Model | LM100 | LM200 |
|--|---|--|
| Product code | 231527 | 231528 |
| Usage environment temperature range | 5°C~35°C | |
| Performance | Rotational speed adjustment range | 50~1000rpm |
| | Maximum torque | Approximately 0.1N·m |
| | Maximum attainable vacuum | Below 399.9Pa |
| Housing | Exterior material | PBT/ADC12 (surface coated) |
| | Motor | DC brushless motor Output power 30W |
| | Rotation adjustment | Power ON/OFF knob |
| | Rotation speed display | - 3-digit×10 rpm digital display |
| | Operation indicator light | LED light (green) Rotation speed display |
| Safety functions | | Overload stop function, slow start function, anti-entanglement protection cover |
| Specifications | Target containers | 50~5000mL (flasks) |
| | Suitable ground joints | TS29/42 (TS24/40 optional) |
| | Compatible stirring shaft | Dedicated stirring shaft |
| | Compatible stirring shaft outer diameter | Φ8mm (Φ7.9~7.95) |
| | Wet parts material | PTFE (glass joint) FKM (oil seal, O-ring) |
| | External dimensions (mm) | W69×D108×H222 |
| | Power supply voltage | Single-phase AC100~240V/1A 50/60Hz |
| Accessories | | Instruction manual, warranty card, shaft, M5 hex bolt, hex wrench, AC adapter, power cord, protective cover, flask clamp |

O-ring and oil seal corrosion resistance comparison

A: Good B: Usable under conditions C/D: Not usable

| Chemical name | FKM O-ring (standard) FKM oil seal (stan- dard) | FFKM O-ring (op- tional) PTFE oil seal (optional) |
|------------------------------|---|---|
| Acetone | D | A |
| Acetone 60°C | D | A |
| Hydrochloric acid (10%) | A | A |
| Muriatic acid (10%) 70°C | A | A |
| Hydrochloric acid (20%) | A | A |
| Muriatic acid (20%) 80°C | A | A |
| Hydrochloric acid (36%) | A | A |
| Muriatic acid (36%) 70°C | A | A |
| Xylene | A | A |
| Chloroform | B | A |
| Acetic acid (10%) | D | A |
| Acetic acid (100%) | D | A |
| Acetic acid (25%) | D | A |
| Acetic acid (50%) | D | A |
| Acetic acid (50%) 70°C | D | A |
| Acetic acid (anhy- drous) | D | A |
| Tetrahydrofuran | D | A |
| Toluene | C | A |
| Pyridine | D | A |
| Hexane | A | A |
| Benzene | C | A |
| Benzene 70°C | C | A |

* Corrosion resistance test by immersion.

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Laboratory Stirrers

LT400A/400B/400C/400D/500A/500B

Made in
Japan

Highly reliable stirrers covering a wide range of viscosities from high to low. Use DC brushless motor and feedback control for stable high-torque rotation. Additionally, it supports a wide voltage range of AC100~125V.

■ Features

- Choose the most suitable model for samples from 6 types covering a wide range of rotational speeds.
- Available in standard torque LT400 and high torque LT500 models.
- Use a DC brushless motor, so there is no need to replace carbon brushes. Safe design that does not produce sparks.
- Equipped with a digital tachometer for easy setting and verification of rotational speed.
- Come with a torque indicator to monitor the load on the stirring shaft.

■ Safety

- Equipped with overcurrent protection circuit, motor overheat protection, and safety cover.



Operation section

■ Specifications

| Model | LT400A | LT400B | LT400C | LT400D | LT500A | LT500B |
|--------------------|---|----------------------------------|--------------------------|------------------------|------------------------|-----------------------|
| Applicable samples | High viscosity use | Medium viscosity use | Medium-low viscosity use | Low viscosity use | High viscosity use | Medium viscosity use |
| Rotation range | 10~300rpm | 15~600rpm | 25~1,200rpm | 60~3,000rpm | 15~600rpm | 25~1,200rpm |
| Rated torque | 0.9N·m (9.0kgf·cm) | 0.5N·m (5.0kgf·cm) | 0.3N·m (3.0kgf·cm) | 0.1N·m (1.0kgf·cm) | 1.0N·m (10.0kgf·cm) | 0.6N·m (6.0kgf·cm) |
| Motor | DC brushless motor 30W | | | DC brushless motor 70w | | |
| Rotation control | Feedback control | | | | | |
| Display panel | Digital display of rotation speed, overload display, torque indicator (1 scale 20%) | | | | | |
| Chuck | Drill chuck for 8mm diameter | | | | | |
| Safety device | Overcurrent protection circuit, motor overheat protection, safety cover | | | | | |
| Specifications | External dimensions | Width 146×Depth 154×Height 165mm | | | | |
| | Power supply | AC100V~AC125V | | | | |
| | Power supply line | With 2m power cord with plug | | | | |
| | Weight | 2.4kg | | | | |
| | Accessories | Clamp, safety cover, wrench | | | | |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Washers

Contents

| | |
|--|-----|
| Ultrasonic Washers | 216 |
| Laboratory Washers for Glassware | 218 |
| Ultrasonic Pipet Washers | 223 |

Ultrasonic Washers | Benchtop Type

M-Y·MH-Y·CPX-Y·CPXH-Y

Washing method | Ultrasound/timing/heating/defoaming

Ultrasonic power (W) | 70/110/160/280

Tank capacity (L) | 1.9/2.8/5.7/9.5/20.8

Spherical waves, powerful and uniform ultrasonic vibrations ensure no blind spots in cleaning. In addition to degreasing and removing adhesive substances, it can also be used for dispersing, dissolving, and extracting applications.



■ Features

Quality Assurance

- High equipment reliability with a failure rate of only 0.6%.

Cleaning Effect

- Feedback load variation ensures consistently stable and uniform operation.
- Different applications can switch between HI and LO power. (CPX-Y, CPXH-Y) .
- Uniformly generated ultrasound can be used for dispersing and dissolving applications.

Alarm Function

- Equipped with an alarm function to ensure equipment performance. (CPX-Y, CPXH-Y)

Waterproof Effect

- Control panel located at the top, vibration board set vertically to improve waterproof effect.

Sweep Function

- The frequency fluctuates periodically within a range of $\pm 1,000\text{Hz}$ for uniform cleaning.

Degassing Function

- New cleaning fluid contains gases (mainly air) that affect the cleaning effect; removing them improves the ultrasonic cleaning effect. CPX-Y/CPXH-Y models come standard with degassing function; for M-Y/MH-Y models, operate the ultrasound for 5~10 min before cleaning when replacing cleaning fluid to achieve the same effect.

Spherical wave technology

- Using a separate vibrator to generate spherical waves in the liquid, it achieves powerful and uniform cleaning throughout the chamber.

■ Common specifications

| | |
|-------------------|--|
| Baths | Stainless steel SUS304 |
| Power supply | AC220V 50/60Hz |
| Power supply line | 1.5m with plug |
| Accessories | Stainless steel tray Water bath cover (1800/2800/3800/5800: resin; 8800: stainless steel) |

■ Specifications

| Model | Function | Dimensions (mm) (L×W×H) | Water tank dimensions (mm) (L×W×H) | Capacity (L) | Weight (kg) | Ultrasonic (W) | Heater (W) | Current (A) |
|------------|---|-------------------------|------------------------------------|--------------|-------------|----------------|------------|-------------|
| M1800-Y | Mechanical timer (5-min increments 0~60 min, continuous mode) | 251×302×303 | 150×137×100 | 1.9 | 3.5 | 70 | None | 0.9 |
| M1800H-Y | Mechanical timer (5-min increments 0~60 min, continuous mode), heating switch | | | | | | 45 | 1.35 |
| CPX1800-Y | Digital timer (1~99 min, continuous mode), degassing mode, power HI/LO switching | | | | | | None | 0.9 |
| CPX1800H-Y | Digital timer (1~99 min, continuous mode), heating adjustment (room temperature (69°C), degassing mode, power HI/LO switching, AUTO mode) | 336×305×304 | 238×138×100 | 2.8 | 4.4 | 110 | 45 | 1.4 |
| M2800-Y | Mechanical timer (5-min increments 0~60 min, continuous mode) | | | | | | None | 1.3 |
| M2800H-Y | Mechanical timer (5-min increments 0~60 min, continuous mode), heating switch | | | | | | 65 | 2.05 |
| CPX2800-Y | Digital timer (1~99 min, continuous mode), degassing mode, power HI/LO switching | 336×305×304 | 238×138×100 | 2.8 | 4.4 | 110 | None | 1.3 |
| CPX2800H-Y | Digital timer (1~99 min, continuous mode), heating adjustment (room temperature (69°C), degassing mode, power HI/LO switching, AUTO mode) | | | | | | 65 | 2.1 |
| M3800-Y | Mechanical timer (5-min increments 0~60 min, continuous mode) | 397×318×381 | 297×152×150 | 5.7 | 5.9 | 110 | None | 1.3 |
| M3800H-Y | Mechanical timer (5-min increments 0~60 min, continuous mode), heating switch | | | | | | 130 | 2.75 |
| CPX3800-Y | Digital timer (1~99 min, continuous mode), degassing mode, power HI/LO switching | | | | | | None | 1.3 |
| CPX3800H-Y | Digital timer (1~99 min, continuous mode), heating adjustment (room temperature (69°C), degassing mode, power HI/LO switching, AUTO mode) | 398×398×381 | 301×240×150 | 9.5 | 6.9 | 160 | 130 | 2.8 |
| M5800-Y | Mechanical timer (5-min increments 0~60 min, continuous mode) | | | | | | None | 1.8 |
| M5800H-Y | Mechanical timer (5-min increments 0~60 min, continuous mode), heating switch | | | | | | 200 | 4.05 |
| CPX5800-Y | Digital timer (1~99 min, continuous mode), degassing mode, power HI/LO switching | 596×466×391 | 495×292×150 | 20.8 | 11.9 | 280 | None | 1.8 |
| CPX5800H-Y | Digital timer (1~99 min, continuous mode), heating adjustment (room temperature (69°C), degassing mode, power HI/LO switching, AUTO mode) | | | | | | 200 | 4.1 |
| M8800-Y | Mechanical timer (5-min increments 0~60 min, continuous mode) | | | | | | None | 3.2 |
| M8800H-Y | Mechanical timer (5-min increments 0~60 min, continuous mode), heating switch | 596×466×391 | 495×292×150 | 20.8 | 11.9 | 280 | 400 | 7.55 |
| CPX8800-Y | Digital timer (1~99 min, continuous mode), degassing mode, power HI/LO switching | | | | | | None | 3.2 |
| CPX8800H-Y | Digital timer (1~99 min, continuous mode), heating adjustment (room temperature (69°C), degassing mode, power HI/LO switching, AUTO mode) | | | | | | 400 | 7.6 |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Made in
Mexico

Operation section




Drain port



With drain valve
3800/5800/8800 comes with 90 cm hose.

Ultrasound



Optional accessories

| | |
|--|-----------------------------|
| | Water bath lid (standard) |
| | Inner frame (without holes) |
| | Inner frame (with holes) |
| | Mesh frame |
| | Cleaning rack (standard) |
| | Beaker frame |
| | Beakers |



Two effective cleaning methods


Direct cleaning
Clean using a stainless steel tray, mesh frame, and inner frame (with holes). Suitable for general cleaning with high efficiency and simplicity.

※ Do not use the stainless steel tray for heavy objects.



Indirect cleaning
It can use different cleaning liquids simultaneously, allowing the use of both cleaning fluid and cleaning agent. Dirt does not remain in the bath, and less cleaning agent is used.

※ Strong acids or strong alkaline cleaning liquids can be used.



| Type | Product name | Product number |
|----------------|------------------------------------|----------------|
| For model 1800 | Mesh frame A12-M | 291815 |
| | Inner frame (without holes) A12-2 | 291812 |
| | Inner frame (with holes) A12-3 | 291813 |
| | Beaker frame (600ml/400ml×1) | 291839 |
| For model 2800 | Mesh frame A22-A | 291820 |
| | Inner frame (without holes) A22-2 | 291817 |
| | Inner frame (with holes) A22-3 | 291818 |
| | Beaker frame (250ml×2) | 291840 |
| | Beaker frame (600ml×2) | 291841 |
| For model 3800 | Mesh frame A32-M | 291825 |
| | Inner frame (without holes) A32-2 | 291822 |
| | Inner frame (with holes) A32-3 | 291823 |
| | Beaker frame (600ml/400ml×2) | 291843 |
| | Beaker frame (250ml×3) | 291842 |
| For model 5800 | Inner frame (without holes) A52-2 | 291827 |
| | Inner frame (with holes) A52-3 | 291828 |
| | Beaker frame (600ml×4) | 291844 |
| For model 8800 | Inner frame (without holes) A82-2 | 291832 |
| | Inner frame (with holes) A82-3 | 291833 |
| | Beaker frame A82-4 (600ml/400ml×6) | 291834 |

Special beaker types

| Capacity | Material | Product number |
|----------|-----------------|----------------|
| 250mL | Glass | 291801 |
| 400mL | Polypropylene | 291802 |
| 600mL | Glass | 291803 |
| 600mL | Stainless steel | 291804 |



Special detergent

Detergent matched to the cleaning object can achieve better cleaning results.

| Detergent | Liquid type | Main cleaning object | Capacity | Product number |
|-----------|-------------|---|----------|----------------|
| BC-III | Neutral | Abrasive powder, lenses, glass, aluminum products, fingerprints | 4L×4 | 001900113 |
| | | | 4L×1 | 001900116 |
| OR-IV | Acidic | Rust, metal oxides | 4L×4 | 001900119 |
| | | | 4L×1 | 001900118 |
| IS-III | Alkaline | Cutting oil for metal and plastic, mineral oil for metal processing, carbon | 4L×4 | 001900111 |
| | | | 4L×1 | 001900115 |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Laboratory Washers for Glassware | Fully-Automatic, Benchtop

AWD510

Washing method Upper, middle, and lower rotating spray nozzles

Hot water sterilization ~93°C arbitrary setting

Conforms to ISO15883 standards, cleaning process can be identified by color change in bath lighting.



Features

Reliable cleaning function

- Able to safely clean and sterilize.
- ~93°C arbitrary rinse temperature setting.
- Sliding two-layer washing, which can meet the needs of large-scale washing at one time.
- Thorough cleaning with rotating nozzles on the upper, middle, and lower levels.
- Automatic dosing and dispensing of cleaning agents (liquid detergent).
- By using jet washing (optional), it is possible to clean the inside of containers such as beakers that are difficult to clean.
- Drying (optional) uses hot air filtered by a high-efficiency filter.
- Since hot water is cooled to below 60°C before discharge, the drainage pipeline does not need to consider heat-resistant circuits.

User-friendly operation

- The Chinese color LCD touch screen is easy to operate.
- With LED lights on the front panel and illumination colors in the washing bath, the progress of the cleaning process can be confirmed even from a distance.
- Using double-layered hollow glass windows, the surface will not heat up during hot water disinfection.
- An electronic lock is used to avoid accidentally opening the door during cleaning.
- Sample vials can also be easily cleaned with a cleaning frame specifically for sample vials (optional).
- According to the cleaning utensils, various cleaning racks (optional) and special cleaning agents are equipped.

Specifications

| Model | | AWD510 | AWD510DRY |
|-----------------------------|----------------------------------|---|---|
| System | Washing method | Upper, middle, and lower rotating nozzles spray | |
| | Cleaning process | Pre-wash, wash, rinse, pure water rinse | Pre-wash, wash, rinse, pure water rinse, dry |
| | | ※ Pre-wash and wash/rinse temperature setting range OFF~80°C/93°C | |
| | Hot water supply method | 2KW heating | |
| | Cleaning agent supply | Automatic dosing and dispensing (liquid detergent) | |
| | Drainage method | Forced drainage after hot water cooling | |
| | Drying method (optional) | — | HEPA hot air drying (about 60°C) |
| Operating | Usage environment temp. range | | 5~40°C (indoor use only) |
| | Water supply connection | Tap water (required) | 5~25°C, water pressure 0.1~0.5MPa, main body interface G3/4, water supply side interface G1/2 |
| | | Warm water (optional) | 5~60°C, water pressure 0.1~0.5MPa, main body interface G3/4, water supply side interface G1/2 |
| | Cleaning agent | Special cleaning agent (purchased separately) | |
| Structure | Pure water connection (optional) | | 0.02~0.1MPa (when using pure water rinse), hose connection outer diameter Φ10.5mm |
| | Exterior material | Stainless steel plate (SUS304) | |
| | Internal chamber material | Stainless steel plate (SUS316) | |
| | Feet | Adjustable leveling feet | |
| Additional functions | | 4.3-inch LCD touch screen/process display LED/sub-process chamber lighting color | |
| Safety device | | Overcurrent leakage protection switch, temperature over-rise prevention, door lock, water leakage sensor | |
| Internal chamber dimensions | | Width 500×Depth 480×Height 480mm | |
| External dimensions | | Width 580×Depth 600×Height 845mm | |
| Weight | | Approx. 87kg | |
| Power supply | | Single-phase AC220V 50/60Hz 14A | |
| Accessories | | Water supply hose, drainage hose, upper cleaning rack removal plug, cleaning agent stainless steel rack, drainage pipe clip, overflow hose, cleaning agent hose | |

Operation screen



Chinese color touch screen operation
Operation details can be confirmed at any time
10 cleaning programs can be called

Washing tank



- ① Upper sliding cleaning rack (optional)
② Lower sliding cleaning rack (optional)

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Made in
Japan

Front LED light and tank color illumination


The current operating status can be understood based on the color illumination indicating each process even from a distance.



Cleaning timetable

Operating conditions: Room temperature: 23±5°C/Water temperature: 23°C/Raw water pressure: 0.1MPa/
Power supply: 200V±5%

Mode **【B】** (rinse, pure water rinse once) Reference time:



Cleaning rack setup example



Sliding cleaning frame

Each pipe ends with a nozzle




Test tube cleaning frame + upper sliding cleaning rack + lower sliding cleaning rack

Options

| No. | Product name | Specifications | Product number |
|-----|--------------------------------|---|----------------|
| ① | Upper sliding cleaning rack | | 291061 |
| ② | Lower sliding cleaning rack | | 291062 |
| ③ | Sliding jet cleaning frame | Do not use sliding cleaning rack Number of nozzles: 36 | 291063 |
| ④ | Test tube cleaning frame | Used with sliding cleaning rack Placement quantity: Φ10.0mm: about 212 pieces×4 frames Φ16.5mm: about 84 pieces×4 frames Φ18.0mm: about 71 pieces×4 frames | 291068 |
| ⑤ | Small pressure rod | For test tube cleaning frame, 4 pcs per set | 291071 |
| ⑥ | Beaker cleaning frame | Used with sliding cleaning rack Placement quantity: 1L beaker, about 8 pieces | 291064 |
| ⑦ | Beaker cleaning frame net | For beaker cleaning frame | 291069 |
| ⑧ | Mesh frame | Used with sliding cleaning rack Width 140×Depth 220×Height 50mm | 291072 |
| ⑨ | Sample vial cleaning frame | 5mL~20mL, diameter 7mm, quantity 96 pieces | OAW56 |
| ⑩ | Sample vial cleaning frame net | Used with sample vial cleaning frame | OAW58 |
| ⑪ | Stand | Width 580×Depth 600×Height 535mm | 291067 |
| ⑫ | Drainage converging pipe | For converging overflow drainage and normal drainage | 291065 |
| ⑬ | Water supply connection unit | IN: R1/2 OUT: G1/2 valve + direct head | 291066 |
| ⑭ | Water Purification Systems | WL220T | 253629 |
| ⑮ | Water inlet assembly | OWH10 (for WL220T) | 253686 |
| ⑯ | Setup pad | OWL50 (for WL220T) Must be used when placed on AWD510 top | 253271 |
| ⑰ | Power cord 4m | OWL52 (for WL220T) | 253273 |
| ⑱ | Ion exchange resin | CPC-P + CPC-E (for WL220T) | 253262 |
| ⑲ | Laboratory cleaning agent | AWL100, general alkaline liquid detergent 2L×4 barrels, alkaline | 291077 |
| ⑳ | Laboratory cleaning agent | AWL200, oil stain alkaline liquid detergent 2L×4 barrels, alkaline | 291078 |
| ㉑ | Laboratory cleaning agent | AWL300, general liquid detergent 2L×4, weak alkaline | 291079 |
| ㉒ | Tank cleaning agent | AWP500, acidic powder detergent 500g×4 boxes, strong acid, for cleaning internal chamber | 291080 |
| ㉓ | Phosphorus-free detergent | General alkaline powder detergent, 8kg | 8190026001 |



① Upper sliding cleaning rack



② Lower sliding cleaning rack



③ Sliding jet cleaning frame



⑲ Laboratory cleaning agent



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Laboratory Washers for Glassware | Fully-Automatic, Large

AW83Z

Made in
Japan

Cleaning time (each process) 1~99 min adjustable

Washing water temp. Supply water temperature~80℃

Fully-Automatic cleaning of large quantities of utensils.



Equipped with a large pressure jet-type fully-automatic washing machine with upper and lower double chambers.

Features

- Cleaning frames for various containers can be placed and cleaned simultaneously.
- Initial wash→detergent wash→rinse→pure water rinse (optional) cleaning mode (3 types) available.
- Jet frame (optional) can clean stubborn stains.
- Use frame fixation, nozzle rotation method to avoid container damage.
- Built-in heater can perform warm water cleaning.
- Liquid detergent is automatically dispensed through a special detergent container. (Powder detergent can be used in portions as needed)

Cleaning room




Specifications

| Model | AW83Z |
|----------------------------------|--|
| Cleaning room | Upper cleaning chamber + lower cleaning chamber 2 chambers Upward and downward two-way pressurized water jet method, rotating nozzles Frame fixed type (replaceable jet frame) |
| Cleaning mode | 3 modes are optional, with optional features in (): ① Pre-wash → Wash → Rinse → (Pure Water Rinse) ② Wash → Rinse → (Pure Water Rinse) ③ Rinse → (Pure Water Rinse) |
| Washing time | Pre-wash: 1~99 min setting Wash: 1~99 min setting Rinse: 1~99 min setting Rinse can be set separately up to 5 times; pure water rinse of 50L can be set 5 times (requires connection to the optional pure water supply device). Washing starts after reaching the specified water volume, and timing begins after reaching the set temperature. |
| Supply water temp. | Room temperature~60℃ |
| Washing water temp. | Supply water temperature~80℃ Pre-wash and wash are set at the same temperature; rinse (pure water rinse) is not heated. |
| Liquid detergent supply method | Pump automatic supply (supply amount adjustable) |
| Liquid detergent tank capacity | 2L |
| Water usage | Pre-wash, wash, rinse: approximately 28L each time |
| Hot water supply method | Heater (built-in: 6KW), or connect to an external hot water supply pipe |
| Supply water pressure | 0.1~0.3MPa |
| Container stand | 2 shelf panels |
| Water supply and drainage method | Water supply: electromagnetic valve switch, water volume controlled by water level gauge (includes overflow prevention) Drainage: pump forced drainage |
| Power supply (50/60Hz) | Three-phase AC200V 30A (30A) |
| Exterior material | Steel plate, surface coated with chemical-resistant paint |
| Interior material | Stainless steel plate |
| External dimensions | Width 860×depth 770×height 1795 mm |
| Internal chamber dimensions | Width 600×depth 630×height 1080 mm (1 chamber) |
| Pump | Washing pump: three-phase AC200V 50/60Hz 355W/560W Drain Pump: Single-phase AC200V 45W |
| Shelf | 550×550 mm maximum load 245N (25kg) |
| Weight | Approx. 220kg |
| Accessories | Water supply pipe (with connector) 2m 1 piece, drain pipe (inner diameter 18mm) 2m 1 piece Special detergent (phosphorus-free) 1kg, measuring spoon (50ml) 1 piece Test tube frame support table 1 piece, tube clamp 1 piece |


Nozzles



Process • Timetable



Dimension diagram (mm)



1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Laboratory Washers for Glassware | Fully-Automatic, Benchtop

AW62

Made in
Japan

Washing time (each process) 0~30min

Water temp. 45~80°C


Benchtop, fully-automatic washer.

A compact, powerful cleaning machine designed for laboratory glassware cleaning needs.


Features

- All processes from washing to rinsing run fully-automatically, with indicator lights for each step.
- The cleaning process and time can be freely selected based on the shape of the container, contamination weight, and cleaning level.
- Cleaning effectiveness is highly influenced by water temperature. This device has a built-in heater to heat the washing water, eliminating the need for boiler piping or water heating systems.
- Compact benchtop type, standard with a rotating table and two-directional high-pressure water spraying mechanism. For pure water rinsing, a pure water supply device is required. The pure water supply device (optional) can serve as device stand and will not occupy additional space.
- If the cleaning liquid does not fully reach the inside of the container and contaminants cannot be cleaned, please replace the spray frame (optional).

Process • Timetable



Dimension diagram (mm)



Specifications

| Type | AW62 |
|------------------------------|--|
| Performance | Washing method |
| | Two-direction pressure water spraying method, fixed nozzle |
| | Frame rotating type (replaceable spray frame) |
| | Cleaning cycle process |
| | Washing (0~30 min) can be set |
| | Rinsing (0~30 min) can be set |
| | Pure water rinse (optional) (connected with pure water supply device) |
| | 20L pure water cleaning |
| | The above processes can be freely selected |
| | Supply water temp. |
| Construction - Specification | Room temp. ~60°C |
| | Washing water temp. |
| | 45~80°C |
| | Hot water supply method |
| | Heater (built-in: 6KW) |
| | Supply water pressure |
| | 0.1~0.3MPa |
| | Container stand |
| | Rotating table (standard) (frame optional) |
| | Water supply and drainage method |
| Accessories | Water supply: Electromagnetic valve switch, water volume controlled by liquid level control switch |
| | Drainage: Natural gravity drainage |
| | Exterior material |
| | Steel plate, surface coated with chemical-resistant paint |
| | Interior material |
| | Stainless steel plate |
| | External dimensions |
| | Width 600×Depth 620×Height 940mm |
| | Internal chamber dimensions |
| | Width 594×depth 572×height 564 mm, effective height: 345 mm |
| | Rotating table |
| | Diameter 550 mm (maximum load bearing 25kg) |
| | Pump |
| | Three-phase AC200V 250W |
| | Door |
| | Up and down open/close type (stops at any position) |
| | Weight |
| | Approx. 90kg |
| | Power supply (50/60Hz) |
| | 3-phase AC200V 30A |
| | Water supply pipe (with connector) 2m 1 piece, drain pipe (inner diameter 25.4mm) 1.5m 1 piece |
| | Special detergent (phosphorus-free) 1kg, measuring spoon (50mL) 1 piece |
| | Resin lid, spray nozzle brush 1 piece |
| | Water supply port component 1 set |

Laboratory Washers for Glassware | Semi-Automatic, Benchtop

AW47

Made in
Japan

Container quantity 450 pieces (test tube 16.5mL) 36 pieces (triangular cone 100mL) 16.5mL)

Washing time 0~60 min can be set

Benchtop, semi-automatic laboratory washer.



Simple and easy to use; it can be configured with 1 personal benchtop washer per room.

Features

- Simply set the time and press the switch for automatic operation.
- Use a two-directional strong pressure spraying method, with rotating nozzles for thorough cleaning (detergent can be used).
- Standard with washing water heater, no need for additional hot water heaters and corresponding piping work.
- Replacing the strong spray frame (optional) allows cleaning of hard-to-clean stains.

Control panel



Options




Spray frame



Test tube frame

Dimension diagram (mm)



Specifications

| Model | | AW47 |
|--------------------------|----------------------------------|--|
| Performance | Washing method | Two-direction pressure spraying method, rotating nozzle (replaceable spray frame) |
| | Supply water temp. | Room temperature~60℃ |
| | Hot water supply method | Heater (built-in: 1KW) |
| | Water supply pressure | 0.1~0.3MPa |
| | Container stand | Stand (standard) (cleaning frame optional) |
| | Water supply and drainage method | Water supply: electromagnetic valve switch Drainage: Natural gravity drainage |
| Structure Specifications | Exterior material | Steel plate, surface coated with chemical-resistant paint |
| | Interior material | Stainless steel plate |
| | External dimensions | Width 450×Depth 490×Height 875mm |
| | Internal chamber dimensions | Width 420×Depth 450×Height 570mm |
| | Pump | AC100V 50/60Hz 200W |
| | Rotating table | Diameter 420 mm |
| | Door | Up and down open/close type |
| | Weight | Approx. 43kg |
| Accessories | Power supply (50/60Hz) | AC100V 15A |
| | | Water supply pipe (with connector) 2m 1 piece |
| | | Drain pipe (inner diameter 25.4mm) 1.5m 1 piece |
| | | Plastic cover 1 piece |
| | | Special detergent (phosphorus-free) 1kg, measuring spoon (50mL) 1 piece Water supply port component 1 set |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers

13 Washers

14 Analysis and
Test Devices

15 Options

Ultrasonic Pipet Washers | Ultrasonic Cleaning, Benchtop

AW31

Made in
Japan

| | | | |
|---------------------|------------------|------------------|-------------------|
| Pipet rack quantity | 264 pieces (1ml) | 216 pieces (5ml) | 136 pieces (10ml) |
|---------------------|------------------|------------------|-------------------|

Powerful dual cleaning mechanism.

Use siphon principle, using ultrasonic waves and detergent for powerful cleaning.

Features

- Built-in high-efficiency cleaning circulation mechanism; the contaminated cleaning liquid inside the pipet is forcibly discharged from the bottom nozzle, fully maximizing the cleaning effect and ensuring no dead angles during cleaning.
- The pipet storage frame in the washing tank maintains a certain rotational motion, allowing sufficient contact between the cleaning item and the ultrasonic waves and spray, achieving uniform cleaning results.
- The spray flow mechanism, ultrasonic generating device, rotary mechanism, and siphon mechanism are efficiently integrated with a rational and compact structural design, making it a high-performance ultrasonic cleaner specifically tailored for pipets.
- By operating and adjusting the water supply valve, the entire process from cleaning to rinsing with clean water can be completed in one washing tank.
- Water-soluble alkaline detergents can be used for cleaning, ensuring very high safety.



Specifications


| Model | | AW31 |
|------------------------------|-----------------------------------|---|
| Pipet quantity | Mother pipet (maximum 500mm) | |
| | 1ml, 264 pieces | |
| | 5ml, 216 pieces | |
| | 10ml, 136 pieces | |
| Water usage | | 1~3L/min |
| Ultrasonic generating device | | High frequency output 55W 28kHz |
| | | Vibrating plate SUS304 stainless steel |
| | | Transducer round ferrite |
| Structure | Washing tank material | Transparent resin |
| | Chamber material | Stainless steel SUS304 |
| | Storage frame dimensions (basket) | Inner diameter Φ 138×height 540mm |
| | | Outer diameter Φ 146×height 570mm |
| | External dimensions | Width 390×Depth 400×Height 844mm |
| | Weight | Body weight approx. 18.5kg, pipe frame approx. 1.0kg |
| Accessories | | Power supply |
| | | AC100V 50/60Hz 2A |
| | | Frame |
| | | Hand pump |
| | | Water supply pipe: Inner diameter 12mm×Length 1.5m |
| | | Drainage pipe: Inner diameter 25mm×Length 1.5m |
| | | Detergent: 1L (low-foam alkaline detergent approx. 1kg) |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Options



Jet cleaning frame (for AW83Z)



Beaker cleaning frame (for AW83Z/62)



Test tube cleaning frame (for AW83Z/62)



Triangular flask cleaning frame (for AW83Z/62)



Jet cleaning frame (for AW62)



Jet cleaning frame (for AW47)



Test tube cleaning frame (for AW47)

Pure water supply unit (for AW62)



Pure water supply unit



Usage examples

Phosphorus-free detergent



For AW83Z

| | Product name | Specifications | Product number |
|------------------------|--|---|----------------|
| Pure water supply unit | External dimensions | Width 420×Depth 600×Height 1,065mm | 291092 |
| | Pump | 355W/560W/(50/60Hz) | |
| | Pure water storage tank | 50L | |
| | Ion exchange column | Mixed type, resin quantity 10L (water collection quantity 1500L, warning light after 30 rinses) | |
| | Jet cleaning frame | 100ml flask, 44 pcs installed | |
| Accessories | Beaker cleaning frame | 50ml beaker, approx. 85 pcs installed | 291081 |
| | Test tube cleaning frame | Φ16.5 test tube, approx. 600 pcs installed | 291082 |
| | Test tube cleaning frame support table | — | 281255 |
| | Triangular flask cleaning frame | 58ml triangular flask, approx. 68 pcs installed | 291083 |
| | Water supply and drainage port change | Specify when placing an order | 281256 |
| Consumables | Dedicated detergent | Phosphorus-free detergent 8kg | 8190026001 |
| | New ion exchange bottle | New product | 000821 |
| | Regenerated ion exchange bottle | Recycling | 000822 |

For AW62

| | Product name | Specifications | Product number |
|-------------|---------------------------------|---|----------------|
| | Pure water supply unit | Built-in 20L pure water bath | 291017 |
| | Jet cleaning frame | 100ml flask, 42 pcs installed | 291086 |
| | Beaker cleaning frame | 50ml beaker, approx. 85 pcs installed | 291081 |
| | Test tube cleaning frame | Φ16.5 test tube, approx. 600 pcs installed | 291082 |
| | Triangular flask cleaning frame | 58ml triangular flask, approx. 68 pcs installed | 291083 |
| Consumables | Dedicated detergent | Phosphorus-free detergent 8kg | 8190026001 |
| | Ion exchange column | Ion exchange resin (3L) | CPCN30010 |

For AW47

| | Product name | Specifications | Product number |
|--|---------------------------------|--|----------------|
| | Jet cleaning frame | 100ml flask, 36 pcs installed | 291090 |
| | Test tube cleaning frame | Φ18.5 test tube, approx. 450 pcs installed | 291091 |
| | Consumables/Dedicated detergent | Phosphorus-free detergent 8kg | 8190026001 |

Analysis and Test Devices

Contents


| | |
|--------------------------------------|-----|
| Thermal Resistance Test System | 226 |
|--------------------------------------|-----|

Thermal Resistance Test System | Thermal Resistance Testing for Power Semiconductor Ceramic Substrates

TE100

Made in
Japan

According to international ISO standards, test the thermal resistance of power semiconductors/semiconductor materials.



The accuracy of sample thermal features (thermal resistance) assessment has been improved!

Features

- Equipment for evaluating the thermal features (thermal resistance) of power device ceramic substrates.
- Tested according to International Organization for Standardization ISO4825-1: 2023-01.
- Capable of testing the heat dissipation features generated based on module structure.
- Capable of testing and evaluating the heat dissipation features of individual substrate materials.

ISO4825-1: 2023-01

Fine ceramics (advanced ceramics, advanced technical ceramics) —Test method for thermal property measurements of metallized ceramic substrates — Part 1:

Fine ceramics (advanced ceramics, advanced technical ceramics) — Test methods for thermal property measurements of metallized ceramic substrates — Part 1: Evaluation of thermal resistance for use in power modules.

Specifications

Main Unit

| Model | | TE100 |
|--|------------------|---------------------|
| Sample size (ISO4825-1: 2023-01 standard) | | 30×30mm |
| Sample load | | 10kg |
| Temp. accuracy | | ≥ 0.01℃ |
| Resistance measurement error | | ±0.1mΩ (70~130Ω) |
| Sampling rate | | 100 times/sec (max) |
| Power supply voltage | | AC100V · 50/60Hz |
| Dimensions | Control unit | W380×D470×H180mm |
| | Measurement unit | W380×D400×H320mm |

※ Exclude monitor, keyboard, and mouse

TEG chip

| Model | | TEG chip |
|-------------------------|--|---------------------|
| Heat generation density | | 1KW/cm ² |
| Maximum input power | | Approximately 250W |
| Temp. rise rate | | 1.4×10 K/sec |
| Dimensions | | W5×D5×H0.35mm |

Cooling water circulation device

| Model | | CFA312C |
|---------------------|--|------------------------------------|
| Circulation method | | External closed system circulation |
| Cooling method | | Air cooling |
| Temp. control range | | -10~60℃ |
| Power supply | | AC220V 8A |
| Dimensions | | W380×D565×H725mm |

1 Sterilizers

2 Granulation
and Spray
Dryers3 Muffle
Furnaces

4 Ovens

5 Incubators

6 Plasma
Equipment7 Water
Purifiers

8 Baths

9 Water
Circulators10 Rotary
Evaporators11 Freeze Dryers
& Cold Traps12 Stirrers &
Shakers


13 Washers

14 Analysis and
Test Devices

15 Options


TEG chip (Consumables)

TEG chips are installed on the samples to be evaluated, such as metallized substrates.



When measuring, open the bottom part to place the sample. (See below diagram)

Measurement method



Standard configuration analysis system (software)

- Simple operation screen consisting of "Setting/Measurement/Result/Help".
- Centralized management of heating of TEG chips and cooling of the CFA312C water circulation device.



Sample testing is available; please contact Yamato Scientific.

Q&A Thermal features evaluation device

Q It can it be used in what kind of market? Also, is the object only metal ceramic substrates?

Mainly aimed at power semiconductor fields such as electric vehicles, charging stations, high-speed trains, etc., contributing to the high thermal conductivity design of semiconductors. Suitable for power semiconductor components such as ceramic substrates, thermal materials, heat sinks, etc.

Q Is sample testing available?

Yes. We can install your sample on a TEG chip and test its thermal resistance; feel free to contact us.

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Contents


| | |
|---------------------------------|-----|
| Stacking Fittings | 230 |
| Stands | 231 |
| Data Loggers | 232 |
| Vacuum Pumps | 233 |
| Air Compressors | 236 |
| Transformers | 236 |
| N ₂ Generators | 237 |

Fitting for Stacking

OD·ODN·ODK


■ OD40C/60C, ODN26C/28C

Dimension diagram (mm)



■ ODK80C/82C/84C

Dimension diagram (mm)




Shelf Plate

Applicable to chamber-type products



Stainless steel punched mesh plate



Stainless steel wire mesh plate

- 1 Sterilizers
- 2 Granulation and Spray Dryers
- 3 Muffle Furnaces
- 4 Ovens
- 5 Incubators
- 6 Plasma Equipment
- 7 Water Purifiers
- 8 Baths
- 9 Water Circulators
- 10 Rotary Evaporators
- 11 Freeze Dryers & Cold Traps
- 12 Stirrers & Shakers
- 13 Washers
- 14 Analysis and Test Devices
- 15 Options


Stands

ON·ONS·OP·OH

ON30C




Dimension diagram (mm)



ON61C




Dimension diagram (mm)



ONS10C/30C/60C




Dimension diagram (mm)



OP42C/62C, OP45C/65C




Dimension diagram (mm)



OH41C/61C



Dimension diagram (mm)



| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

Data Loggers

SR·AL·EX·KR



■ Paper-based data logger (Yokogawa)

| Model | SR10001 | SR10002 | SR10003 | SR10004 | SR10006 |
|----------------------|---|---------|---------|---------|---------|
| Input type | General input: thermocouple, thermistor, DC voltage, DC current | | | | |
| Measurement cycle | Pen type: 125ms/channel; Dot type: 1s/6 points or 2.5s/6 points | | | | |
| Calculation function | Differential calculation between channels, Linear scale, Square root, Deviation correction function | | | | |
| Recording section | Effective recording width: 100mm | | | | |
| Display section | LED 7-segment display (2+5 digits) | | | | |
| Dimensions | W144×H144×D247 | | | | |
| Power supply | 100-240VAC | | | | |



■ Paper-based data logger (CHINO)

| Model | AL4706 | AL4712 | AL4724 |
|---------------------|---|--------|--------|
| Input type | General input: thermocouple, thermistor, DC voltage, DC current | | |
| Measurement cycle | 1 second/6 points, 2 seconds/12 points, 2 seconds/24 points | | |
| Recording method | Color tape dot recording method, 6 colors | | |
| Recording method | Effective recording width: 100mm | | |
| Scheduled recording | Digitally records month, day, time, channel number, value, unit. Interval (h, min) can be set arbitrarily | | |
| Dimensions | W144×H144×D222 | | |
| Power supply | 100-240VAC | | |



■ Paperless data logger (Yokogawa)

| Model | EX102 | EX106 | EX108 | EX110 | EX112 |
|------------------------------|--|-------|-------|-------|-------|
| Input type | General input: thermocouple, thermistor, DC voltage | | | | |
| Display section | 5.7-inch wide viewing angle high-resolution TFT color LCD screen | | | | |
| Minimum measurement cycle | 125ms | | | 1s | |
| Measurement/display accuracy | ±0.05% of RDG (DC Voltage), ±0.15% of RDG (Thermocouple, Thermistor) | | | | |
| Optional accessories | Maximum 2GB CF card, Power measurement recording, Vacuum degree recording, Flow accumulation, F0 value calculation | | | | |
| Dimensions | W144×H144×D184 | | | | |
| Power supply | 100-240VAC | | | | |



■ Paperless data logger (CHINO)

| Model | KR2S2P | KR2S6P |
|------------------------------|--|--------|
| Input type | General input: thermocouple, thermistor, DC voltage | |
| Display section | 5.7-inch wide viewing angle high-resolution TFT color LCD screen | |
| Minimum measurement cycle | 1s/6 points (100ms/4 points) | |
| Measurement/display accuracy | ±0.1%±1 digit (with exceptions) | |
| External storage | CF card (capacity 256MB-8GB), standard configuration 256MB; USB drive (capacity 32M-8GB) | |
| Dimensions | W144×H144×D257 | |
| Power supply | 100-240VAC | |

1 Sterilizers

2 Granulation and Spray Dryers

3 Muffle Furnaces

4 Ovens

5 Incubators

6 Plasma Equipment

7 Water Purifiers

8 Baths

9 Water Circulators

10 Rotary Evaporators

11 Freeze Dryers & Cold Traps

12 Stirrers & Shakers

13 Washers

14 Analysis and Test Devices

15 Options

Vacuum Pumps

GLD·GCD

GLD series



| Model | GLD-N051 | GLD-N137 | GLD-N202 | GLD-N280 |
|-------------------------|---|--------------------------|--------------------------|--------------------------|
| Geometric pumping speed | 60L/min | 165L/min | 240L/min | 336L/min |
| Ultimate pressure | Gas ballast valve closed: 6.7×10^{-1} Pa | | | |
| | Gas ballast valve open: 6.7 Pa | | | |
| Power supply | Single phase AC220V 2.5A | Single phase AC220V 3.5A | Single phase AC220V 4.1A | Single phase AC220V 5.6A |
| Oil used | ULVOIL R-4 | ULVOIL R-7 | | |
| Required oil amount | 500-800mL | 1000mL | 1100mL | 750-1500mL |
| Air inlet | KF25 | | | |
| External dimensions | W150×D396.5×H222.7mm | W170×D461.5×H241.1mm | W170×D489.5×H241.1mm | W181×D496.5×H269mm |
| Weight | 12.5kg | 24.5kg | 26.5kg | 31.5kg |

GCD series



| Model | GCD-051X | GCD-136X | GCD-201X |
|-------------------------|---|--------------------------|--------------------------|
| Geometric pumping speed | 50L/min | 135L/min | 200L/min |
| Ultimate pressure | Gas ballast valve closed: 6.7×10^{-1} Pa | | |
| | Gas ballast valve open: 6.7 Pa | | |
| Power supply | Single phase AC220V 2.4A | Single phase AC220V 3.6A | Single phase AC220V 3.6A |
| Oil used | SO-M | | |
| Required oil amount | 500-800mL | 1000mL | 1100mL |
| Air inlet | KF25 | | |
| External dimensions | W165.5×D419×H222.7mm | W170×D493×H241.1mm | W170×D541.5×H241.1mm |
| Weight | 14.1kg | 25.4kg | 29.4kg |

Recommended model selection table

| Vacuum drying oven model | GLD series | | | | GCD series | | |
|--------------------------|------------|----------|----------|----------|------------|----------|----------|
| | GLD-N051 | GLD-N137 | GLD-N202 | GLD-N280 | GCD-051X | GCD-136X | GCD-201X |
| DP23C | ○ | | | | ○ | | |
| DP33C | ○ | | | | ○ | | |
| DP43C | | ○ | ○ | | | ○ | ○ |
| DP63C | | ○ | ○ | | | ○ | ○ |
| DP83C | | | | ○ | | | ○ |
| DP103C | | | | ○ | | | ○ |
| ADP210C | ○ | | | | ○ | | |
| ADP310C | ○ | | | | ○ | | |

Sterilizers 1

Granulation and Spray Dryers 2

Muffle Furnaces 3

Ovens 4

Incubators 5

Plasma Equipment 6

Water Purifiers 7

Baths 8

Water Circulators 9

Rotary Evaporators 10

Freeze Dryers & Cold Traps 11

Stirrers & Shakers 12

Washers 13

Analysis and Test Devices 14

Options 15

Vacuum Pumps

PK·NeoDry

PK series



| Model | PK250-2 | | PK600-2 | PK1000-2 |
|---|----------------------------------|--|------------------|------------------|
| Geometric pumping speed | 250L/min | | 600L/min | 1000L/min |
| Ultimate pressure | Gas ballast valve closed: 1.0 Pa | | | |
| | Gas ballast valve open: 10 Pa | | | |
| Power supply voltage | Single phase AC200V | | | |
| Air inlet | KF25 | | KF40 | |
| Exhaust vent | KF25 | | | |
| Maximum water vapor processing capacity | 250g/h | | 350g/h | 600g/h |
| External dimensions | W210×D385×H250mm | | W290×D475×H275mm | W315×D530×H275mm |
| Weight | 23kg | | 54kg | 56kg |

NeoDry Series



| Model | NeoDry7E | NeoDry15E | NeoDry30E | NeoDry36E | NeoDry60E |
|---|--------------------------|-----------|-------------------------|-----------|-----------|
| Geometric pumping speed | 110L/min | 250L/min | 500L/min | 600L/min | 1000L/min |
| Ultimate pressure | 1.0Pa | | | | |
| Power supply voltage | Single-phase AC220V~240V | | Three-phase AC220V~240V | | |
| Air inlet | KF25 | | KF40 | | |
| Exhaust vent | KF25 | | | | |
| Maximum water vapor processing capacity | 120g/h | 250g/h | | 350g/h | 600g/h |
| Sound intensity | ≤56dB | | | ≤58dB | ≤60dB |
| Weight | 19kg | 23kg | 25kg | 54kg | 56kg |

Recommended model selection table

| Vacuum drying oven model | PK series | | | NeoDry series | | | | |
|--------------------------|-----------|---------|----------|---------------|-----------|-----------|-----------|-----------|
| | PK250-2 | PK600-2 | PK1000-2 | NeoDry7E | NeoDry15E | NeoDry30E | NeoDry36E | NeoDry60E |
| DP23C | ○ | | | ○ | ○ | | | |
| DP33C | ○ | | | ○ | ○ | | | |
| DP43C | ○ | | | ○ | ○ | | | |
| DP63C | ○ | | | | ○ | | | |
| DP83C | ○ | ○ | | | ○ | ○ | | |
| DP103C | ○ | ○ | | | ○ | ○ | | |
| ADP210C | ○ | | | | ○ | | | |
| ADP310C | ○ | | | | ○ | | | |

Vacuum Pumps

DA·DOP

DA Series



| Model | DA-60D | DA-120S | DA-121DF | DA-241SF |
|-------------------------|---|------------------|--|--------------------------|
| Geometric pumping speed | 60L/min | 120L/min | | 240L/min |
| Ultimate pressure | 3.32×103Pa | 13.3×103Pa | 3.3×103Pa | 16.0×103Pa |
| Power supply | Single phase AC220V 2.4A | | Single phase AC220V 2.3A | Single phase AC220V 2.5A |
| Inlet/outlet diameter | Outer diameter 14×Inner diameter 9 (G3/8) | | Outer diameter 16×Inner diameter 12 (G1/2) | |
| Ambient temp. | 7~40℃ | | 0~40℃ | |
| External dimensions | W156×D358×H238mm | W162×D358×H238mm | W193.5×D411×H241.1mm | W207×D411×H285mm |
| Weight | 19kg | | 26kg | |

DOP Series




| Model | DOP-181SD | DOP-301SB | DOP-400SB | DOP-420SA |
|-------------------------|--------------------------|---|--|---|
| Geometric pumping speed | 180L/min | 300L/min | 400L/min | 420L/min |
| Ultimate pressure | 10.0×103Pa | 8.0×103Pa | 12.0×103Pa | 17.3×103Pa |
| Power supply | Single phase AC220V 2.9A | Three-phase AC200V~230V 2.1A | Three-phase AC200V~230V 2.4A | Three-phase AC200V~230V 3.5A |
| Inlet/outlet diameter | Rc3/8 | Outer diameter 16×Inner diameter 12 (Rc1/2) | Suitable for piping outer diameter Φ16 | Outer diameter 26×Inner diameter 20 (Rc3/4) |
| Ambient temp. | 7~40℃ | 0~40℃ | | |
| External dimensions | W162×D266×H235mm | W315×D443×H231mm | W316×D434×H231mm | W310×D523×H253mm |
| Weight | 12kg | 20kg | 23kg | 33kg |

Recommended model selection table

| Vacuum drying oven model | DA series | | | | DOP series | | | |
|--------------------------|-----------|---------|----------|----------|------------|-----------|-----------|-----------|
| | DA-60D | DA-120S | DA-121DF | DA-241SF | DOP-181SD | DOP-301SB | DOP-400SB | DOP-420SA |
| DP23C | ○ | | | | | | | |
| DP33C | ○ | | | | | | | |
| DP43C | | ○ | ○ | ○ | ○ | ○ | | |
| DP63C | | ○ | ○ | ○ | ○ | ○ | | |
| DP83C | | | | | | ○ | ○ | ○ |
| DP103C | | | | | | ○ | ○ | ○ |
| ADP210C | ○ | | | | | | | |
| ADP310C | ○ | | | | | | | |

- Sterilizers 1
- Granulation and Spray Dryers 2
- Muffle Furnaces 3
- Ovens 4
- Incubators 5
- Plasma Equipment 6
- Water Purifiers 7
- Baths 8
- Water Circulators 9
- Rotary Evaporators 10
- Freeze Dryers & Cold Traps 11
- Stirrers & Shakers 12
- Washers 13
- Analysis and Test Devices 14
- Options 15

Vacuum piping



Air Compressors

TYW-1/2



Small oil-free air compressor

| Model | TYW-1 | TYW-2 |
|--------------------------|----------------------------------|----------|
| Start pressure | 0.45MPa (0.45~0.5MPa adjustable) | |
| Shutdown pressure | 0.75MPa | |
| Maximum working pressure | 1MPa | |
| Exhaust volume | 60L/min | 90L/min |
| Air tank volume | 22L | 38L |
| Air tank size | Φ377×242 | Φ406×340 |
| Sound intensity | ≤65dB | |
| Output power | 550W | 840W |
| Power supply | Single phase AC220V | |

Transformers

TR10·TR20



Voltage conversion transformer (100v)

| Model | TR10-1000 | TR10-3000 | TR10-5000 |
|------------------|--|----------------|-----------|
| Input voltage | Single phase AC220V | | |
| Output voltage | Single-phase AC100V (Japanese power socket) | | |
| Rated power | 1kVA | 3kVA | 5kVA |
| Output interface | Socket (universal) | Terminal block | |
| Safety device | Overcurrent and leakage protection circuit breaker | | |

Voltage conversion transformer (200v)

| Model | TR20-1000 | TR20-3000 | TR20-5000 |
|------------------|--|----------------|-----------|
| Input voltage | Single phase AC220V | | |
| Output voltage | Single phase AC200V | | |
| Rated power | 1kVA | 3kVA | 5kVA |
| Output interface | Socket (universal) | Terminal block | |
| Safety device | Overcurrent and leakage protection circuit breaker | | |

N₂ Generators

NiGen LCMS 40-1




Main features

- It can simultaneously meet the nitrogen and air demands of multiple LC-MS systems.
- Use advanced PSA (Pressure swing adsorption) technology combined with carbon molecular sieve technology, along with unique high-purity nitrogen backflow technology, efficiently removing impurities and hydrocarbons from the air to ensure nitrogen purity.
- Built-in air compressor with proprietary technology, 16,000 h maintenance-free, superior to similar products.
- Patented FastArt technology ensures quick startup, reducing gas generation time.
- Integrated with a 50L nitrogen buffer, reducing compressor startups, extending compressor life, and lowering maintenance costs and noise.
- Special soundproof materials and vibration reduction measures ensure silent operation.
- Equipped with wheels at the base for easy mobility.
- Use a LCD touchscreen menu display, with the CPU connected to the generator via CAN-BUS to monitor operational parameters, alarm prompts, and facilitate operation.

Working principle

- Use advanced PSA (Pressure Swing Adsorption) technology to guide compressed air from the air compressor into the carbon molecular sieve. Oxygen, carbon dioxide, water, and other impurities are removed as they pass through the carbon molecular sieve, allowing only nitrogen to pass through and enter the gas storage tank. After pressure and flow adjustments in the gas storage tank, it can be connected to the LC-MS.



Technical specifications

| Flow | | 40 NL/min |
|-----------------------------|--|-------------|
| Separation technology | PSA, CLAIND efficient carbon molecular sieve columns×2 | |
| Outlet pressure | 0-8 bar/0-115 psi, adjustable | |
| Oxygen purity | Up to 99.9%, adjustable | |
| Dew point | ≤ -40℃ | |
| Air compressor | Built-in, air flow: 135L/min, free of suspended liquids and phthalates | |
| Power requirements | 230Vac±10% 50 Hz 1.6KW | |
| Enclosure protection rating | IP20 | |
| Operating temp. | 5℃-40℃/40℃-100℃ | |
| Environmental humidity | 5%-90% | |
| Maximum altitude | 2000m asl | |
| External dimensions | Height | 130cm/51.1" |
| | Width | 44cm/17.3" |
| | Depth | 110cm/43.3" |
| Weight | 220 kg/489 lbs | |

| | |
|------------------------------|----|
| Sterilizers | 1 |
| Granulation and Spray Dryers | 2 |
| Muffle Furnaces | 3 |
| Ovens | 4 |
| Incubators | 5 |
| Plasma Equipment | 6 |
| Water Purifiers | 7 |
| Baths | 8 |
| Water Circulators | 9 |
| Rotary Evaporators | 10 |
| Freeze Dryers & Cold Traps | 11 |
| Stirrers & Shakers | 12 |
| Washers | 13 |
| Analysis and Test Devices | 14 |
| Options | 15 |

